



Agilent Technologies
R&D and Marketing GmbH & Co. KG
Hewlett-Packard Straße 8
76337 Waldbronn
Deutschland

www.agilent.com

To Whom It May Concern

The purpose of this letter is to provide additional information concerning the injection accuracy of the Agilent 1260 Infinity II Vialsampler (Agilent Product Number: G7129A) and the Agilent 1290 Infinity II Vialsampler (Agilent Product Number: G7129B), not covered in the Product Data sheet (Agilent Publication Number: 5991-6287EN).

On May 10, 2016, Agilent Technologies introduced the new 1260 Infinity II Liquid Chromatography modules and systems, evolving from one of the world's most-trusted LC systems, the Agilent 1100 Series, 1200 Series and 1200 Infinity Series LC.

In addition to the specification, documented in the product data sheet, it is a great pleasure to confirm that the injection accuracy of the Agilent G7129A and G7129B is typically smaller than 1.0 %.

Sincerely yours,

A handwritten signature in blue ink that reads "H. Schulenberg-Schell".

Dr. Helmut Schulenberg-Schell
Director Business Development LPSD
Liquid Phase Analysis Division
Waldbronn, June 23, 2016

A handwritten signature in blue ink that reads "Michael Kraft".

Dr. Michael Kraft
Director Industry Marketing LPSD
Liquid Phase Analysis Division
Waldbronn, June 23, 2016