

IV CANNULA PILOT OVERVIEW

To avoid repeated vein puncture in patients for whom fluids, electrolyte solutions or medication must be given intravenously on a permanent basis, an I.V. cannula is inserted, which can generally be left in place for several days

FEATURES AND BENEFITS

NEEDLE: Siliconised stainless steel back cut needle help clinicians to perform easy and painless venipuncture.

CATHETER: Lubricated with a exclusive process and made of Biocompatible **FEP** material has smooth outer surface and tapered tip for very smooth and painless venipuncture with minimum trauma.

FLASH BACK CHAMBER: It allows quick visualization of flash back of blood and confirmation of catheter placement in vessel.

FLEXIBLE WINGS: The flexible wings provides easy gripping and follows the contours of site allowing more secure fixation.

INJECTION VALVE: It is used for intermittent medication.

NEEDLE COVER: It requires minimum force during the injection of intermittent medication



OPTIONS AVAILABLE

- FEP/PTFE/PUR Catheter, Colored Wings
- Packaging in Coloured Paper
- Medical Kraft Paper/Tyvek

MEDICAL GRADE MATERIALS SAFETY TESTED

- Needle - high quality steel with high chrome-nickel content, manufactured in compliance with AISI standards
- Catheter - **PU, PTFE, FEP**
- Packaging - Blister/ box / carton

QUALITY

- Conforms to ISO 10555 standard

STERILIZATION

- Ethylene oxide

SHELF LIFE

- 5 years

ORDER INFORMATION

CODE NO.	SIZE	COLOR	CATHETER O.D (MM)	CATHETER LENGTH (MM)	WATER FLOW RATE (ML./MIN.)	PACKING INNER/OUTER UNIT
	14G	ORANGE	2.1	270	270	100/1000
	16G	GREY	1.7	180	180	100/1000
	17G	WHITE	1.5	125	125	100/1000
	18G	GREEN	1.3	95	95	100/1000
	20G	PINK	1.1	65	65	100/1000
	22G	BLUE	0.9	36	36	100/1000
	24G	YELLOW	0.7	23	23	100/1000
	26G	VIOLET	0.6	15	15	100/1000



CE 0123 EN ISO 13485



WELLEX
ENGINEERING FOR HEALING



PACKAGING
Tray / Carton



STERILISATION
Ethylene Oxide



QUALITY



SHELF LIFE
5 Years

ORDER INFORMATION CODES FOR ORDERS UPON REQUEST

www.wellexmed.com



IV Cannula
























IV Cannula

Overview

- Intravenous cannulas / Catheters allow the introduction or withdrawal of fluids from the human circulation system. The short flexible and kink-resistance catheter is introduced into a blood vessel over a hollow introducer needle.
- These IV Cannulas are fitted with an injection port that allows the supplementation of a drug into the cannula itself by a syringe without a needle, and prevents back-flow. The color-coded cap is <<flip top>> (Snapcap) which allow a rapid one-finger closure. The flexible wings provide easy and proper fixation with fixation tapes.

Ordering Information

#	Code	Colour	Size	Description
1.	FIV143OR		14G	IV Cannula with wings&Injection port
2.	FIV163GR		16G	IV Cannula with wings&Injection port
3.	FIV183GR		18G	IV Cannula with wings&Injection port
4.	FIV203PK		20G	IV Cannula with wings&Injection port
5.	FIV223BU		22G	IV Cannula with wings&Injection port
6.	FIV243YW		24G	IV Cannula with wings&Injection port
7.	FIV263VO		26G	IV Cannula with wings&Injection port
8.	FIVP14OR		14G	PEN IV Cannula
9.	FIVP16GR		16G	PEN IV Cannula
10.	FIVP18GR		18G	PEN IV Cannula
11.	FIVP20PK		20G	PEN IV Cannula
12.	FIVP22BU		22G	PEN IV Cannula
13.	FIVP24YW		24G	PEN IV Cannula
14.	FIVP26VO		26G	PEN IV Cannula
15.	FIVS14OR		14G	Safety IV Cannula
16.	FIVS16GR		16G	Safety IV Cannula
17.	FIVS18GR		18G	Safety IV Cannula
18.	FIVS20PK		20G	Safety IV Cannula
19.	FIVS22BU		22G	Safety IV Cannula
20.	FIVS24YW		24G	Safety IV Cannula
21.	FIVS26VO		26G	Safety IV Cannula





Add :

Cairo-Ismailia Desert
Road, Ismailia Pub-
licFree Zone 41639
Ismailia City, Egypt.

Tel :

(002) 064 3488080

Fax :

(002) 064 3488070

E-mail :

info@wellexmed.com

Website :

www.wellexmed.com



Add : Cairo-Ismailia Desert Road, Ismailia PublicFree Zone 41639 Ismailia City, Egypt.

Tel : (002) 064 3488080

Fax : (002) 064 3488070

E-mail : info@wellexmed.com

Website : www.wellexmed.com