



Cytology Brush

Used endoscopically to obtain mucosal tissue samples from gastrointestinal tract for microscopic-histopathological examination



CY
not wire guided




CY-DL
wire guided

Art. No.	Guidewire	Working Length (cm)	Working Channel Ø (mm)
CY	0.032"	230	2.8
CY-DL	0.032"	230	2.8

Naso Jejunal Feeding Tube

To feed nutrition through mouth or nose with the tube leading into the jejunal

Tube Polyurethane



Art. No.	Size (Fr/CH)	Working Length (cm)
NJFT-12120	12	120
NJFT-12	12	155
NJFT-14120	12	120
NJFT-14	12	155
NJFT-16120	12	120
NJFT-16	12	155