

Clearing agents are completely miscible with alcohol and paraffin, making them an ideal intermedium in histology and cytology processing. They are used multiple times during tissue processing and section staining; before embedding the tissue in paraffin, for deparaffinization before staining and before covering the stained section using medium for permanent preservation.

Xylene-based medium for use in histology and cytology.

Xylene substitute, aliphatic hydrocarbon-based. Essentially no odor. For use in histology and cytology.

1000 mL	BCN-1L
2500 mL	BCN-2.5L
5000 mL	BCN-5L
10000 mL	BCN-10L

