

Bipolar Forceps

For precise desiccation and coagulation control. All forceps are autoclavable and require a bipolar cord and foot switch for operation. Cold sterilization may damage insulating material.

	Ref. No.
Jeweler's Micro Tip, Straight.....	7-809-1
Length, 4in. (10.2cm.) Tip, 0.4mm. Insulated.	
Requires foot switch and bipolar cord.	
Semkin Micro Tips, with Stops.....	60-1771-001
Length, 5.5in. (14cm.) Tip, 0.7mm. Uninsulated.	
Requires foot switch and bipolar cord.	
Semkin Straight	60-1772-001
Length, 6in. (15.2cm.) Tip, 1.5mm. Insulated.	
Requires foot switch and bipolar cord.	
Gerald Micro Tips	7-809-2
Length, 7in. (18.0cm.) Tip, 1.0mm. Insulated.	
Requires foot switch and bipolar cord.	
Hardy Micro Tips	7-809-3
Bayonet with Stops. Length, 8.25in (21cm). Tip, 0.5mm. Uninsulated	
Requires foot switch and bipolar cord.	
Adson Smooth Tips	7-809-4
Length, 5in. (13cm.) Tip, 1.0mm. Insulated.	
Requires foot switch and bipolar cord.	
Cushing Smooth Tips	7-809-5
Length, 7in. (18.0cm.) Tip, 1.5mm. Insulated.	
Requires foot switch and bipolar cord.	
Cushing Smooth Tips	7-809-6
Bayonet. Length, 7.5in. (19.0cm.) Tip, 0.7mm. Insulated	
Requires foot switch and bipolar cord.	

