
XBit 80

Color Doppler System Datasheet

V1.0



General Information

Dimensions and Weight

- Dimensions of main unit (approx.): 764mm*630mm*1440mm
- Net weight of main unit (approx): 75kg (no probe included)

Electrical Power

- Power supply voltage: Auto adaptable for AC100-240V
- Power supply frequency: 50-60 Hz
- Power consumption: 600 VA
- Stand-by mode
- Battery for option

User Interface

Operation Panel

- Control panel
 - Height adjustable: 15 cm
 - Rotatable with 100° left and right
 - released from the machine body
- Alphanumeric keyboard
- 8 TGC Slides
- Interactive backlit keys
- Integrated speaker
 - Volume adjustable

Touch panel

- 13.3 inch Touch screen
 - customizable layout
 - can show the image real time
 - sensitive to use

Display Screen

- High resolution color LED
- Dimension: standard 23.8 inch
- Resolution: 1920×1080
- Image Area: 800×600
Full Screen: 1120×840
1440×1080
- Brightness and contrast adjustment

System Overview

Applications

- Abdominal (Gynecology & Urology)
- Fetal/OB
- Small Parts
 - Breast
 - Thyroid
 - Scrotum

- Pediatrics
- MSK_Conventional & Superficial
- Cardiac (adult & pediatric)
- Transvaginal

Scanning Method

- Electronic convex
- Electronic linear
- Electronic phased array
- Volume convex

Transducer Types

- Convex probe:
C3-T, C3M-T
- Linear probe:
L7-T, L12-T, L8M-T, L12M-T, L10i-T,
L8M6-T,L8M5-T,L18-T,L7SVA-T,L7R-T
- Trans-vaginal probe:
E7-T, E7MW-T
- Phased array probe:
P2-T, P5-T,P3T-T
- Volume probe:
VC4-T, VE6-T
- Micro convex probe:
MC6-T, MC3-T
- Bi-plane probe:
R7B8-T
- Bi-plane probe:
- TEE probe:
T5-T
- pencil probe:
CW2-T

Image Modes

- B Mode
- B/M mode
- M mode
- 2B Mode
- 4B Mode
- CFM Mode
- 2D Steer
- PD Mode
- DPD Mode
- PW Mode
- B/BC Mode
- Triplex
- Quadplex
- CW Mode
- Free Steering M Mode
- TDI
- Color M Mode
- Curved Panoramic Imaging
- Trapezoidal imaging
- Compound
- SRA
- Elastography
- Stress Echo
- ECG
- Super Needle
- 4D
- Virtual HD
- FHI mode
- AIO

Display Mode

- Quad/Dual display
- Duplex mode
- Triplex mode
- Quadplex mode

Display Annotation

- Hospital name

- Date/Time
- Patient Name and Patient ID
- Gray/Color bar
- Cine guide
- Scanning direction
- Measurement results window
- Transducer type
- Frequency
- Application name
- Menu indication
- Trackball functions indication
- Imaging parameters displayed on the screen

Standard Configuration

- 23.8 inch LED monitor
- 13.3 inch touch screen
- 4 active transducer ports
- 500G integrated hard disk
- DVD-R/W
- USB ports: 8
- TGC
- LGC
- B, 2B, 4B, B/M, B/BC, CFM, PW, Power Doppler/Directional PD, Instant Triplex, Duplex, Quadplex, Trapezoidal, Chroma B&M&PW, Full Screen
- Automatic PW trace and measurement in real time
- Super Image module: FHI, Multiple Compound Imaging, SRA (Speckle Reduction Algorithm), AIO
- Q-image (intelligent image optimization), X-contrast, Q-beam, Q-flow
- Measurement & calculation software packages: General, OB&GYN, Cardiac
- Zoom
 - Zoom navigator
 - Improve the fps
 - High resolution Zoom
 - Real-time Zoom

- Frozen Zoom
- PIP Zoom
- Pan Zoom

Software Options

- 4D module
- Virtual HD/Depth View
- HD Niche/ Smart Volume Slice/SonoCrystal
- 2D steer
- Stress Echo
- Auto EF
- Strain
- Intelligent Doppler
- SonoOB
- SonoContrast
- SonoPW
- SonoColor
- SonoNeedle
- UltraRemote
- SoundFlow
- Virtual Apex
- Static 3D
- MVI Mode
- Elastography
- SWE
- Super Needle
- Curved panoramic
- Color Panoramic
- Curved expanding
- TSS
- HIPPA
- Extended Cardiac Package: ECG Software, Free M, Color M, CW, TDI, IMT
- DICOM 3.0
- HL7
- SonoIMT
- WIFI Function
- Bluetooth

- Biopsy kit: for convex/linear/TV/ Micro-Convex probe respectively
- SonoCarotid
- SonoFollicle
- Remote Service

Hardware Options

- Footswitch
- ECG Lead

Peripherals

- SONY UP-X898MD B&W Video Printer
- SONY UP-D25MD

Imaging Parameters

B Mode

- Gain: 0-255, 5/ step
- Compound: 1, 2,3,4, off
- SRA: off.1.2
- Focus Number: 1-9
- Focus Position: 0-15, depend on the probes
- Full Screen: on/off
- X-contrast: Normal/Enhance/Suppress
- Q-image: 0-4
- Persistence: 0-7
- Density: Low/Mid/High
- 2D Map: 1-20
- Noise Reject: 0-255
- Scan Width: 6-100%
- Image Rotate: 0, 90, 180, 270
- Gamma: 0-8
- Chroma: 1-29, user
- Smooth: 0-7
- Edge enhance: 0-6
- A.power: 0-100%
- Frequency: depend on the probes
- Dynamic: 20-290
- Depth: depend on the probes

- Zoom: 0.05-20, 0.05/ step
- TGC: 8
- Center Line: on/off
- Trapezoidal Mode: on/off
- Biopsy: on/off
- Biopsy Level: 30,45,60
- Super Needle: on/off
- Needle angle: $\pm 30^\circ$
- Curved panoramic imaging: on/off (only for linear probe)
- Elastography: on/off
- 2D steer: $\pm 20^\circ$

M Mode

- Gain: 0-255, 5/ step
- Layout: LR/UD
- Display Format: 1:2, 1:1, 2:1
- Chroma: 1-29, user
- Free Steering M Mode: off, 1, 2, 3
- Color Map: 1-9, user
- 2D map: 1-20
- Dynamic: 20-290, 18/ step
- Speed: 1-8

Color Mode

- Gain: 0-255, 5/ step
- Color Map: 1-20, user
- Color Invert: on/off
- Q-flow: on/off
- Q-beam: on/off
- Persistence: 0-7
- Color Mode: Velocity, Variance, Vel-Var
- Wall Filter: 0-3
- Density: Low/Mid/High
- Wall Thre: 0-15
- Blood Ejection: Smooth, HRes, HRes2, HRes3
- B/BC: on/off
- Frequency: depend on the probes
- Baseline: -3-3

- Scale: depend on the probes
- Steer: -20-20 (only for linear probe)
- PRF: depend on the probes

CPA/DPD Mode

- Gain: 0-255, 5/ step
- Wall Filter: 0-3
- Q-beam: on/off
- Q-flow: on/off
- Wall Thre: 0-15
- Persistence: 0-7
- Frequency: depend on the probes
- PRF: depend on the probes
- Steer: -20-20 (only for linear probe)
- Color Map: 1, user

PW Mode

- Gain: 0-255, 5/ step
- 2D Map: 1-20
- Wall Filter: 0-3
- Spectrum Enhance: 0-3
- Dynamic Range: 46-67, 3/ step
- Invert: on/off
- Display format: 1:2, 1:1, 2:1
- Triplex: on/off
- Quadplex: on/off
- Auto Cal Parameter: on/off
- DTrace Smooth: 0-3
- Threshold: 1-5
- DVmean: on/off
- DVmax: on/off
- Trace area: above, below, all
- Audio: 0-100%
- Color Map: 1-29, user
- QuickAngle: on/off
- Auto Cal: on/off
- Freq: depend on the probes
- Baseline: -5-5
- PRF: depend on the probes

- Steer: -30-30 (only for linear probe)
- Speed: 1-8

CW Mode

- Gain: 0-255, 5/ step
- 2D Map: 1-20
- Spectrum Enhance: 0-3
- Dynamic: 46-67, 3/ step
- Audio: 0-100%
- Wall Filter: 0-3
- Color Map: 1-9, user
- QuickAngle: on/off
- Baseline: -5-5
- PRF: depend on the probes
- Speed: depend on the probes

SonoContrast

- SonoContrast provides exceptional Contrast agent detecting capability, not only extracts second harmonic, but also non-linear fundamental signals
- Available for convex, linear, phased array and endocavitary probe
- Available for abdomen, GYN, URO, Thyroid, Breast,
- Analysis package
- mechanical index
- Timer1: on/off
- Timer2: on/off
- Retro capture and Pro capture storage
- Comparative analysis on complex curves
- Dual live: side by side displays tissue image and contrast image
- Mix: mix contrast image with tissue image
- Mix map: 7 types
- Persistence: 8 steps
- Dynamic range: 20-290
- 2D map: 20 types
- Chroma: 30 types
- Supports U/D Flip and L/R Flip
- Rotation: 90 degrees/ step

- transpose position of contrast and tissue image
- Line density: low/ Mid /high

Triplex Mode

- B+C+D
- Available on all probes

Quadplex Mode

- B+C+D+auto trace
- DTrace Calc Parameters: Vs, Vd, TAMAX, VTI, Time, RI, PI, S/D, HR

Technology and Function

Fusion Harmonic Imaging

- Available on all probes
- FHI key ON/OFF
- Second active multi-frequency

Trapezoidal

- Available for linear probe
- combined with compound algorithm space

SonoOB

- Automatic measurement: BPD, HC, AC, FL, NT ,OFD
- Efficiency and accuracy

HIPPA

- Password to get into the system
- User define

Curved panoramic

- Real time
- Support measurement
- Erasable design
- Color map: 30 types

Elastography

Available on linear, convex, endocavitary

- Dual images simultaneous
- Modifiable ROI
- Support strain ratio measurement
- Real-time display of pressure column
- Quantitative comparison

Dynamic focusing

- Wider the focus area
provide image more detail and higher resolution

Post Processing for raw data

- Support measurement

- Adjust the gain, TGC, 2D map, chroma, dynamic range, invert etc.

Stress Echo

- Available on phased array probes
- 36 factory protocols
- User-defined protocols
- Analysis system: wall motion scoring
- Professional report

Cineloop

- Cine loops: up to 2000 frames
- Perspective: Save cine in real time
- Retrospective: Save cine loops in frozen mode
- Support 2D, M, PW, CFM, CPA, DPD, CW, Color M, Free Steering M
- Simultaneous and independent review in duplex mode
- Cineloop auto/manual
- Variable cine playback speed
- User-define start and end frame of cine storage
- User-define start and end frame of cine review
- storage in hard disk and display in real-time modes
- Slide show: slide show function

Storage

- 500GB integrated SSD
- DVD-R/W driver
- USB ports
- Still images storage format: IMAG
- Still images export format: BMP, JPG, DCM, PNG, TIFF
- Cine loops storage format: CINE
- Cine loops export format: AVI
- Fast storage setting
- System suitable to avoid the loss of data / images

EasyView

- Image review Layout: 1×1, 2×2
- Image management

Exam Review

- Search Exam
- Exam review: patient view, study view
- Exam management
 - Delete selected exam
 - Export selected exam
 - Backup selected exam
 - Recover from the backup exam
 - Selected all
 - Expand all
 - Collapse all
 - Edit selected Exam
 - Review selected Exam
 - Continue selected Exam

Connectivity

- Ethernet work connection
- USB for USB Device
- DICOM support(option)
 - Verify
 - Print
 - Store
 - Worklist
 - Structure report
 - MPPS
 - Query/retrieve

Measurement and Calculations

General Measurement Package

- Software packages for various specific clinical use
- Comprehensive analysis methods
- Clinical analysis reports

General measurement package

- B mode General measurement
 - Distance
 - Length__Area (Ellipse)
 - Length__Area (Trace)
 - Volume (1 Distance)
 - Volume (2 Distance)

Volume (3 Distance)	MDistance
Volume (1 Ellipse)	MTime
Volume (2 Ellipse)	Velocity
Volume (1 Distance 1 Ellipse)	HR
Ratio	HR Manual
Angle	● PW mode GYN measurement
Strain Ratio	Umb A
HR Manual	MCA
SonoColor	Uterine Artery
SWE	Fetal AO
● M mode normal measurement	HR Manual
MDistance	● B mode OB measurement
MTime	Distance
Velocity	GS
HR	CRL
HR Manual	BPD
● PW mode Normal measurement	Auto BPD
Velocity	AC(Ellipse)
Distance	Auto AC
Peak	HC
Auto Trace	Auto HC
Manual Trace	FL
StD%	Auto FL
StA%	Humerus
Area	OFD
ICA/CCA	Auto OFD
HR	NT
Volume Flow	Free NT
HR Manual	Fetal Biometry
Clinical Analysis Packages	Fetal Long Bones
● B mode GYN measurement	Fetal Cranium
Distance	OB Others
UT	Z Score
Cervix Vol.	AFI
ENDO	Ductus Venosus
OV Volume	CX_L
FO_D	Aorta
FO Auto	Descending Aorta
Uterine Artery	MCA
HR Manual	Umb A
Strain Ratio	Uterine Artery
● M mode GYN measurement	Pulmonary Artery
	Fetal Select
	HR Manual

Strain Ratio	Velocity
● M mode OB measurement	HR
MDistance	HR Manual
MTime	● PW mode Vessel measurement
Velocity	CCA
HR	ICA
HR Manual	ECA
FHR	Vertebral A
● PW mode OB measurement	INT IIL
Umb A	EXT IL
Aorta	ILIAC
Descending Aorta	CFA
Left Uterine Artery	ProFun
Right Uterine Artery	LTCIR
Pulmonary Artery	SFA
MCA	Pop A
FHR	ATA
HR Manual	PTA
Duct Venosus	PERON
● B mode Vascular measurement	DRPED
IMT (Auto)	HR
IMT Mean	Volume Flow
SonoColor	HR Manual
CCA	● B mode URO measurement
ICA	Distance
ECA	Void Vol.
Vertebral A	Prostate Vol.
EXT IL	Kidney Volume
INT IIL	T-Zone Vol.
ILIAC	Bladder Vol.
CFA	StA%
ProFun	StD%
LTCIR	Vessel Area
SFA	Vessel Dis
Pop A	HR Manual
ATA	Strain Ratio
PTA	● M mode URO measurement
PERON	MDistance
DRPED	MTime
HR Manual	Velocity
Strain Ratio	HR
● M mode Vessel measurement	HR Manual
MDistance	● PW mode URO measurement
MTime	Velocity

Acceleration	StA%
Distance	Area
Peak	ICA/CCA
Auto Trace	HR
Manual Trace	Volume Flow
StD%	HR Manual
StA%	● B mode Pediatrics measurement
Area	HIP
ICA/CCA	Vol(3Dis)
HR	HR Manual
Volume Flow	Strain Ratio
HR Manual	● M mode Pediatrics measurement
● B mode Small Parts measurement	MDistance
Distance	MTime
Length__Area (Ellipse)	Velocity
Length__Area (Trace)	HR
Volume (1 Distance)	HR Manual
Volume (2 Distance)	● PW mode Pediatrics measurement
Volume (3 Distance)	Velocity
Volume (1 Ellipse)	Distance
Volume (2 Ellipse)	Peak
Volume (1 Distance 1 Ellipse)	Auto Trace
Ratio	Manual Trace
Angle	StD%
Strain Ratio	StA%
Breast	Area
Auto Breast	ICA/CCA
Thyroid	HR
Auto Thyroid	Volume Flow
HR Manual	HR Manual
SWE	● B mode Carotid measurement
● M mode Small Parts measurement	Subclavian A
MDistance	CCA
MTime	Bulb
Velocity	ICA
HR	ECA
HR Manual	Vertebral A
● PW mode Small Parts measurement	General Measurement
Velocity	Strain Ratio
Distance	HR Manual
Peak	● M mode Carotid measurement
Auto Trace	MDistance
Manual Trace	MTime
StD%	Velocity

HR	Aortic Valve
HR Manual	Tricuspid Valve
● PW mode Carotid measurement	Pulmonary Valve
Subclavian A	RV/LV
CCA	LV Mass
Bulb	TAPSE
ICA	Vp
ECA	HR Manual
Vertebral A	● PW mode Cardiac measurement
General Measurement	Velocity
ICA/CCA	Acceleration
HR	Time
Volume Flow	Slope
HR Manual	HR
● B mode Cardiac measurement	ED/PS
Auto EF	Mitral Valve
Teichholz	Aortic
Simpson SP	Tricuspid Valve
Simpson Biplane	Pulmonary Valve
Modify Simpson	Pulmonary Vein
Cube	PISA
Bullet Volume	Qp/Qs
Gibson	Tei Index
Mitral Valve	TDI
Aortic Valve	HR Manual
Pulmonary Valve	● B mode Abdomen measurement
Tricuspid Valve	CBD
LVOT	GB Wall
RVOT	Liver Length
PISA	Artery
LV Mass	Spleen
Qp/Qs	Renal Vol.
RV/LV	GB Volume
IVC	Iliac
RA/LA	HR Manual
AO/LA	Strain Ratio
HR Manual	● M mode Abdomen measurement
● M mode Cardiac measurement	MDistance
MDistance	MTime
MTime	Velocity
Slope	HR
HR	HR Manual
Left Ventricle	● PW mode Abdomen measurement
Mitral Valve	Velocity

- Acceleration
- Distance
- Peak
- Auto Trace
- Manual Trace
- StD%
- StA%
- Area
- ICA/CCA
- HR
- RAR
- Volume Flow
- HR Manual
- B mode TCD measurement
 - ICA
 - CS
 - MCA
 - ACA
 - PCA
 - ACOA
 - PCOA
 - OA
 - Vertebral A
 - BA
 - PICA
 - HR Manual
- PW mode TCD measurement
 - ICA
 - CS
 - MCA
 - ACA
 - PCA
 - ACOA
 - PCOA
 - OA
 - Vertebral A
 - BA
 - PICA
 - HR Manual

SYSTEM SETUP

By using system setup, users could

- Customize hospital information
- Customize language

- Customize fast storage time
- Customize color map
- Customize functions to Footswitch, P1 key, Print key
- Customize functions to alphanumeric 0~9
- Customize PC and Video Print
- Option
- Customize Measure
- Customize Comment library
- Customize Report

User Define Functions

- By user-define function, users could customize user-define preset, including
 - Applications name, Presets name
 - Applications exam type
 - Imaging parameters

Multi-language Display Interface

- English
- Chinese
- Other languages

Note: other languages for detailed, please contact CHISON.

Inputs and outputs

- AC Power In: 1
- AC power Out: 1
- Power Button: 1
- USB Port: 8个
- Ethernet: 1
- Remote Control: 1
- S-Video Out: 1
- Audio: L,R
- DVI: 1
- VGA Out: 1
- Video Out: 1
- Footswitch Port: 2
- Ground pole: 1

Operating conditions

- Ambient temperature: 10°C to 40°C
- Relative humidity: 30% to 75% (no condensation)
- Atmospheric pressure: 700 hPa to 1060 hPa

Storage and Transport conditions

- Ambient temperature: -5°C to 40°C
- Relative humidity: ≤80% (no condensation)
- Atmospheric pressure: 700 hPa to 1060 hPa

Not all features or specifications described in this document may be available.

CHISON Medical Technologies Co., Ltd. reserves the right to make changes in specifications and features shown herein, or discontinues the product at any time without notice or obligation. For the most current information, Contact CHISON Representative.