

Save \$650 & get a lab-themed hidden object puzzle Save and claim >



Home > Shop All Products > Tubes > PCR Tubes > MicroAmp™ Optical 8-Tube Strip,...



[Certificates](#)


Applied Biosystems™

MicroAmp™ Optical 8-Tube Strip, 0.2 mL

The Applied Biosystems™ MicroAmp™ Optical 8-Tube Strip is optimally designed to fit snugly in our Applied Biosystems™ Real-Time PCR [Read more](#)

Have Questions? [Contact Us](#)

Catalog Number	Quantity
4316567	125 Strips

alog number 4316567 

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)

ThermoFisher
SCIENTIFIC



Quantity: 125 Strips

[Product Overview](#)

[Documents](#)

[FAQ](#)

The Applied Biosystems™ MicroAmp™ Optical 8-Tube Strip is optimally designed to fit snugly in our Applied Biosystems™ Real-Time PCR Systems or Thermal Cyclers, maximizing thermal conductivity and providing precise control of cycling temperatures. Requires a tray or retainer tray specific to your instrument.

Note: See user's manual or package insert for limited label license, and trademark information.

Specifications

Certifications/Compliance	Certified DNA/Rnase Free
For Use With (Equipment)	Thermal Cyclers or Real-Time PCR Systems
Format	(125) 8-Tube Strips
Product Line	MicroAmp™
Quantity	125 Strips
Reaction Speed	Standard
Shipping Condition	Room Temperature
Volume (Metric) Working	0.2 mL
Closure Type	Not Included
Color	Optical
Material	Polypropylene
Product Type	Optical 8-Tube Strip
Volume (Metric)	0.25 mL
Unit Size	Each

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



Storage: Room temperature

Documents & Downloads

Certificates

Lot #	Certificate Type	Date	Catalog Number(s)
I52F24Q223	Certificate of Analysis	May 27, 2025	4316567
I52F24Q113	Certificate of Analysis	May 27, 2025	4316567
I52F24Q321	Certificate of Analysis	May 27, 2025	4316567
I52F24Q122	Certificate of Analysis	May 27, 2025	4316567
I52F24Q123	Certificate of Analysis	May 27, 2025	4316567

5 results displayed, search above for a specific certificate

[Request a Certificate](#)

Product Information

Manuals



- [User Bulletin: Understanding the Differences in PCR Plastics Format](#)
- [User Bulletin: How to Use MicroAmp Reaction Plates, Tube Strips, and Tubes](#)

Save \$650 & get a lab-themed hidden object puzzle [Save and claim >](#)



Brochures



 [Technical reference: MicroAmp plastics instrument compatibility chart](#)

Frequently asked questions (FAQs)

Can I run strip tubes on the ViiA 7 instrument?



What is the difference between the MicroAmp Optical Tubes (P/N801-0933) and the ABI PRISM Optical Tubes (P/N 4316567)? Is there any difference in the quality of the plastic used for the two different types?



What does the SeqStudio Starter Kit (PN A35000) contain?

