

WANCARE

ULTRASOUND ECG gel

USAGE

Ensures clearer and uninterrupted transmission of ultrasonographic waves by disconnecting the air between the skin and the prob during any and all Ultrasonography, Doppler ECG or effort test.

PROPERTIES

pH balanced between 5,00 - 7,00.

No toxic effect. Not irritate the skin.

Water soluble. Easy to clean.

Not contain formaldehyde, stabilizers, oily substances and salt. Odorless. No damage and corrosion to the ECG electrode.



COMPOSITION

Deionized water, Carbormer, TriEthanol Amine, MonoPropylene Glycol, Benzyl Alcohol, Methylchloroisothiazolinone

STORAGE CONDITIONS

Please store under 35°C keep away from direct sunlight, store in a cool and dry place.

SIZE	REF NO
250ML	KAF G31-2
500ML	KAF G31-1
1LT	KAF G31
5LT	KAF G31-3

KAFGRUP
SAĞLIK HİZMETLERİ

Atakent Mh. 221. Sk. No:3A Rota Office A Blok D.83
K.Çekmece/İSTANBUL
T:+90 212 471 4 200(pbx) F:+90 212 471 4 201
www.kafgrup.com / info@kafgrup.com

KAF GRUPOĞRETİM HİZMETLERİ
İNŞAAT SANAYİ VE TİCARET LTD. ŞTİ.
Atakent Mh. 221. Sk. No:3A Rota Office
A Blok D.83 K.Çekmece / İSTANBUL
Tel: 0212 471 42 00 / fax: 0212 471 42 01
Halkalı V.B.: 486 053 3864