

## Technical Data Sheet


**Model No.:N7603**  
**BUDI-DI: 697442484L30QE**

- **Product Name:** C-Section Pack
- **Classification:** Class Is, according to Annex IX, MDD 93/42/EEC
- **Description:** C-section pack is mainly made of AS-Only SMMS, by cutting, punching and processing with certain protective characteristics, including 7 accessories: table cover, reinforced surgical gown, C-section drape, baby drape, etc.
- **Intended use:** C-section pack is used for clinical surgery in the obstetrics department in medical institutions.
- **Material:** 45gsm AS-Only SMMS
- **Approvals:** CE, ISO 13485, EN13795-1
- **Manufacturer:** Leboo Healthcare Products Limited
- **Country of origin:** China
- **Specification:** Detail see specification as below.

Component	Description	Color	Quantity
Table Cover	150x200cm (75x200cm)	Blue	1
SMMS Reinforced Surgical Gown	L 125x145cm	Blue	2
C-Section Drape	200x300cm, with: Fenestration 30x30cm, adhesive,liquid collection pouch 75x75cm, tube holders	Blue	1
BABY Drape	75x90cm	/	1
Umbilical clip	/	/	1
Hand Towel	30x40cm	White	2
Medical sheet	100x100cm	Blue	1

### ● Packaging:

- (1) Carton Requirements: Five layers corrugated
- (2) Packing Quantity: 1pc/header pouch, 6pcs/ctn  
5 Layers Carton(Paper)



### ● Storage:

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5°C to +45°C and with relative humidity below 80%.

Design:

Check:

Approve:

Date: 2022-9-13

#### ● Shelf Life:

Shelf life is 36 months from the date of manufacture when stored as stated above.

#### ● Cautions and disposal:

- (1) Do not contact with heavy oils, sparks or flames;
- (2) Do not use in the environment with conditions of excessive heat;
- (3) Replace if damaged, heavily contaminated or in accordance with local practice;
- (4) Hand and dispose of contaminated products with care and in accordance with national regulations.

CE 0197 CE Symbol



Medical Device



Keep away from sunlight



Do not reuse



Keep dry



Don't use when packing damaged  
and consult instructions for use



Temperature limit



EO Sterile



H.O.: A-1506, Lar Valley International, No.168 Guang'anmen Wai St., Xicheng District, Beijing 100055, China  
T: +86 10 63252899 F: +86 10 63252899 E: leboo@leboochina.com.cn W: www.leboo.com.cn

Design:

Check:

Approve:

Date: 2022-9-13