


Patient Monitor

With people in mind

P Series P22 / P18 / P15 / P12

Powerful. Practical. Pioneer.





Winner of the world's leading design award

P Series won the Red Dot Design Award: Product Design 2021. This international award confirms global recognition of the superiority of the P Series design and calls attention to the Biolight vision of continuous development of medical devices for the benefit of human health.


















NMT (Neuromuscular transmission)

NMT helps in determining the depth of neuromuscular blockade.



rSO₂ (regional oxygen saturation)

A non-invasive method of directly measuring cerebral perfusion




RM (respiratory mechanics)

Help understand the ventilator / anesthesia machine operation and patient's respiratory function




BIS (bispectral index)

Used to determine the depth of sedation or hypnosis.



EEG (electroencephalograghy)


Used to determine changes in brain activity. Findings assist in diagnosing brain disorder.



C.O. (cardio output)

A gold standard for monitoring the critically ill patient to ensure tissue oxygenation.

Monitoring parameters of all vital organs



P1. Patient-First.

Mobile Care. Continuous Care.





Multiple use and flexibility

(stand-alone use or and as part of Host Monitor P12 – P22)

Monitoring also for intra-hospital transfer

Just slide the P1 transport monitor into the patient monitor, the P Series can implement the seemless transfer from operating room, ICU, to the recovery room.



P1



Screen (inch)	5.5 inch touch screen
Weight	< 0.9 kg
Battery	8 hours continuous operation
Resolution	1280*720
Built in module	CO ₂

Expanded module Masimo, Nellcor SpO₂ Suntech NIBP Sidestream/Mainstream CO₂







EWS (Early Warning Score)

Helps identify early signs of deterioration and quickly determine the degree of illness.

GCS (Glasgow Coma Scale)

Provides a practical method for assessment of conscious level in response to defined stimuli.






Sepsis View

Helps to score and guide the treatment of sepsis based on the latest SSC guidelines.

BOA (Balance of Anesthesia)

Focus on relevant parameters over the course of the surgery.








Advanced Seamless Network Connection


Anyview P Series integrates the hospital network to achieve efficient management of patient data through high-end connection technology.

It connects to other bedside medical devices such as anesthesia machine through Biolink module. Via central monitoring station, the patient data and other beside devices data are shared seamlessly with the hospital network, and support the clincian the diagnosis.



Unique Drip Monitoring





Accessories P Series

Module Rack Extension of module capacity



Transfer HolderFor intra-hospital transfer



LÖWENSTEIN medical

Sales + Service

Löwenstein Medical Arzbacher Straße 80 56130 Bad Ems, Germany T. +49 2603 9600-0 F. +49 2603 9600-50 info@loewensteinmedical.com loewensteinmedical.com









Guangdong Biolight Meditech Co.,Ltd.

NO.2 Innovation First Road,
Technical Innovation Coast,
Hi-tech Zone, Zhuhai, P.R. China.
Postcode: 519085
T. +86-756-3399935
F. +86-756-3399911
overseas@blt.com.cn
blt.com.cn



