

Trocar pyramidal sharp



Obturator round



Trocar Sleeve / Sheath



Arthroscope Sheath and Obturator

- 1 stopcocks fixed
- Working length of Sheath: 130mm
- $\varnothing$ 6mm diameter, Highflow, Shaft Trocar
- [For Endoscopes: $\varnothing$ 4mm, 175 length, 30° direction of view]
- Obturator diameter $\varnothing$ 4mm, length 175mm



Straight Arthroscopic Forceps (Punch)

Model: 542-701C

- Tip profile: 2,28mm
- Bite Width: 1,59mm
- Tip width: 3,25mm
- Working Length: 13cm



Angled Arthroscopic Forceps (Punch)

Model: 542-007C

- Direction: Right
- Tip profile: 2,52mm
- Bite Width: 3,17mm
- Tip width: 5,05mm
- Working Length: 13cm



Angled Arthroscopic Forceps (Punch)

Model: 542-006C

- Direction: Left
- Tip profile: 2,52mm
- Bite Width: 3,17mm
- Tip width: 5,05mm
- Working Length: 13cm





Oval Arthroscopic Forceps

Model: 542-031C

- Tip profile: 1,65mm
- Bite Width: 4,39mm
- Tip width: 6,57mm
- Working Length: 13cm



Oval Arthroscopic Forceps

Model: 542-031C

- Tip profile: 1,65mm
- Bite Width: 4,39mm
- Tip width: 6,57mm
- Working Length: 13cm



Straight Arthroscopic Forceps

Model: 542-001C

- Tip profile: 2,52mm
- Bite Width: 3,17mm
- Tip width: 5,05mm
- Working Length: 13cm



Arthroscopic Punch

Model: 542-8864C


- Retrgrade right-cutting
- Jaw Section (H x W): 4,0 x 5,7mm
- Cut width: 2,3mm
- Working Length: 13cm



Straight Arthroscopic Scissors

Model: 542-061C

- Tip profile: 2,71mm
- Bite Width: 0,76mm
- Tip width: 2,75mm
- Working Length: 13cm



Microfactor Instruments

Model: 577-001

- $\varnothing$ 4,0mm
- Microfactor 35°
- Length: 160mm



Arthroscopic Rasp, convex

Model: 572-047

- convex
- rough tip
- length: 140mm



Arthroscopic Rasp, concave

Model: 572-048


- concave
- rough tip
- length: 140mm



Arthroscopic Instruments, for palpation

Model: 572-028

- hook probe 5mm graduated
- length: 160mm



PROWATM
MEDICAL INSTRUMENTS



Inflow Irrigation Cannula

Model: 515-521

- $\varnothing$ 8mm,
- Capacity 7mm
- Length: 168mm

