

Metacross[®]
PTA Balloon Dilatation Catheter

Crosstella[®]
PTA Balloon Dilatation Catheter

Croserio[®]
PTA Balloon Dilatation Catheter



MAXIMIZE
YOUR OPTIONS
FOR **RX AND OTW**
PERFORMANCE



Digitally signed by Grabazei Alexandru
Date: 2020.04.09 14:33:58 EEST
Reason: MoldSign Signature
Location: Moldova

TERUMO
INTERVENTIONAL
SYSTEMS

CHOOSE FROM A FULL LINE OF RX AND OTW OPTIONS

Designed for optimal performance

METACROSS® PTA Balloon Dilatation Catheter (0.035")

CROSSTELLA® PTA Balloon Dilatation Catheter (0.018")

CROSPERIO® PTA Balloon Dilatation Catheter (0.014")

Providing end-to-end performance with less effort and more options, our rapid exchange (RX) and over-the-wire (OTW) PTA balloon dilatation catheters let you choose the optimal solution for you and your patients.

- 3 guidewire sizes available (0.014", 0.018", and 0.035")¹⁻⁶
- Smaller balloon profiles help ease crossability^{7,*}
 - Supports re-crossability within the same procedure¹⁻⁶
- Well-designed shaft and stiff tapered core wire engineered for maximizing pushability[†]
 - Effective transition of force from proximal to distal end



Catheter	Guidewire Compatibility (in)	Usable Shaft Length (cm)	
Metacross RX	0.035	135	
Metacross OTW	0.035	65/90/135	
Crostella RX	0.018	150	Pro
Crostella OTW	0.018	150	
Crosperio RX	0.014	150	
Crosperio OTW	0.014	65/150	

* Diameter range: METACROSS® RX/OTW (3-12 mm); CROSSTELLA® RX/OTW (2-6 mm); CROSPERIO® RX/OTW (1.5-4 mm).


† RX design feature.

‡ METACROSS® RX/OTW PTA Balloon Catheter illustrations shown.

PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**




CHOOSE FROM A FULL LINE OF RX AND OTW OPTIONS

Designed for optimal performance

- METACROSS®** PTA Balloon Dilatation Catheter (0.035")
- CROSSTELLA®** PTA Balloon Dilatation Catheter (0.018")
- CROSPERIO®** PTA Balloon Dilatation Catheter (0.014")

Providing end-to-end performance with less effort and more options, our rapid exchange (RX) and over-the-wire (OTW) PTA balloon dilatation catheters let you choose the optimal solution for you and your patients.

- 3 guidewire sizes available (0.014", 0.018", and 0.035")¹⁻⁶
- Smaller balloon profiles help ease crossability^{7,*}
 - Supports re-crossability within the same procedure¹⁻⁶
- Well-designed shaft and stiff tapered core wire engineered for maximizing pushability[†]
 - Effective transition of force from proximal to distal end



RAPID EXCHANGE: BETTER HANDLING BY DESIGN

Achieve optimal results with Rapid Exchange

- Enables faster, more convenient exchange of devices
- Supports the use of shorter guidewires
- Encourages procedural efficiency, reducing fluoroscopy time and radiation dose
- Optimizes control, allowing a single operator to mitigate lost wire position


Metacross® RX
PTA Balloon Dilatation Catheter

ENGINEERED FOR PERFORMANCE

Excellent crossability with innovative catheter design⁷

- Designed to provide best-in-class pushability⁷
- Larger diameter proximal shaft⁷
- Smallest balloon profile among tested semi-compliant balloons⁷

METACROSS® RX has excellent crossability compared to select competitors



Catheter	Guidewire Compatibility (in)	Usable Shaft Length (cm)	Shaft Diameter		Proximal Tapered	Proximal Marker (mm)	Center (mm)	Distal Marker (mm)	Distal Tapered (mm)	Distal (mm)
			Proximal	Distal						
Metacross® RX	0.035	135	Proximal 3.9 Fr (1.30 mm)	Distal 5.7 Fr (1.90 mm)	1.7	1.7	1.6	1.7	1.4	1.02
Metacross® OTW	0.035	65/90/135	5 Fr compatible: 5.4 Fr (1.80 mm) 6/7 Fr compatible: 5.7 Fr (1.90 mm)		1.7	1.7	1.6	1.7	1.4	1.07
Crostella® RX	0.018	150	Proximal 3.4 Fr (1.12 mm) or 4.1 Fr (1.35 mm)	Distal 4 Fr (1.32 mm)	1.28	1.32	1.28	1.36	1.29	0.62
Crostella® OTW	0.018	150	4 Fr compatible: 3.9 Fr (1.30 mm) 5 Fr compatible: 4.3 Fr (1.43 mm)		1.19	1.30	1.26	1.30	1.01	0.58
Crosperio® RX	0.014	150	Proximal 3.5 Fr (1.16 mm)	Distal 3.5 Fr (1.16 mm)	0.85	0.85	0.82	0.86	0.83	0.47
Crosperio® OTW	0.014	65/150	3.8 Fr (1.28 mm)		0.90	0.87	0.81	0.83	0.76	0.46

* Diameter range: METACROSS® RX/OTW (3-12 mm); CROSSTELLA® RX/OTW (2-6 mm); CROSPERIO® RX/OTW (1.5-4 mm).
† RX design feature.
‡ METACROSS® RX/OTW PTA Balloon Catheter illustrations shown.

METACROSS® RX

CROSTELLA® RX

CROSPERIO® RX

See inserts for ordering information.

Metacross®

PTA Balloon Dilatation Catheter

Crosstella®

PTA Balloon Dilatation Catheter

Crosperio®

PTA Balloon Dilatation Catheter



FIND OUT MORE



Phone: 800.862.4143



terumo.com



Fax: 800.411.5870

Indications

CROSPERIO® RX PTA Dilatation Balloon Catheter, CROSPERIO® OTW PTA Dilatation Balloon Catheter, CROSTELLA® RX PTA Dilatation Balloon Catheter, CROSTELLA® OTW PTA Dilatation Balloon Catheter, METACROSS® RX PTA Dilatation Balloon Catheter, and METACROSS® OTW PTA Dilatation Balloon Catheter are intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native and synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

References:

1. METACROSS® RX IFU, January 2016. **2.** METACROSS® OTW IFU, March 2016.
3. CROSTELLA® RX IFU, February 2016. **4.** CROSTELLA® OTW IFU, May 2016.
5. CROSPERIO® RX IFU, March 2016. **6.** CROSPERIO® OTW IFU, March 2016. **7.** Data on file.

©2018 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-00605

TERUMO
INTERVENTIONAL
SYSTEMS

Metacross® RX

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

METACROSS® RX 0.035" PTA Balloon Dilatation Catheter – 135 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
3	BD-L30020CR	BD-L30040CR	BD-L30060CR	BD-L30080CR	BD-L30100CR	BD-L30120CR	BD-L30150CR	BD-L30200CR
4	BD-L40020CR	BD-L40040CR	BD-L40060CR	BD-L40080CR	BD-L40100CR	BD-L40120CR	BD-L40150CR	BD-L40200CR
5	BD-L50020CR	BD-L50040CR	BD-L50060CR	BD-L50080CR	BD-L50100CR	BD-L50120CR	BD-L50150CR	BD-L50200CR
6	BD-L60020CR	BD-L60040CR	BD-L60060CR	BD-L60080CR	BD-L60100CR	BD-L60120CR	BD-L60150CR	BD-L60200CR
7	BD-L70020CR	BD-L70040CR	BD-L70060CR	BD-L70080CR	BD-L70100CR	BD-L70120CR	BD-L70150CR	BD-L70200CR
8	BD-L80020CR	BD-L80040CR	BD-L80060CR	BD-L80080CR	—	—	—	—
9	BD-L90020CR	BD-L90040CR	BD-L90060CR	BD-L90080CR	—	—	—	—
10	BD-LX0020CR	BD-LX0040CR	BD-LX0060CR	BD-LX0080CR	—	—	—	—
12	BD-LX2020CR	BD-LX2040CR	—	—	—	—	—	—

6F

7F



Metacross® OTW

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 65 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120		150	200
3	BD-L30020SW	BD-L30040SW	BD-L30060SW	BD-L30080SW	BD-L30100SW	BD-L30120SW	BD-L30150SW	BD-L30200SW	5Fr
4	BD-L40020SW	BD-L40040SW	BD-L40060SW	BD-L40080SW	BD-L40100SW	BD-L40120SW	BD-L40150SW	BD-L40200SW	
5	BD-L50020SW	BD-L50040SW	BD-L50060SW	BD-L50080SW	BD-L50100SW	BD-L50120SW	BD-L50150SW	BD-L50200SW	
6	BD-L60020SW	BD-L60040SW	BD-L60060SW	BD-L60080SW	BD-L60100SW	BD-L60120SW	BD-L60150SW	BD-L60200SW	
7	BD-L70020SW	BD-L70040SW	BD-L70060SW	BD-L70080SW	BD-L70100SW	BD-L70120SW	BD-L70150SW	BD-L70200SW	
8	BD-L80020SW	BD-L80040SW	BD-L80060SW	BD-L80080SW	—	—	—	—	6Fr
9	BD-L90020SW	BD-L90040SW	BD-L90060SW	BD-L90080SW	—	—	—	—	
10	BD-LX0020SW	BD-LX0040SW	BD-LX0060SW	BD-LX0080SW	—	—	—	—	7Fr
12	BD-LX2020SW	BD-LX2040SW	—	—	—	—	—	—	

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 90 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120		150	200
3	BD-L30020MW	BD-L30040MW	BD-L30060MW	BD-L30080MW	BD-L30100MW	BD-L30120MW	BD-L30150MW	BD-L30200MW	5Fr
4	BD-L40020MW	BD-L40040MW	BD-L40060MW	BD-L40080MW	BD-L40100MW	BD-L40120MW	BD-L40150MW	BD-L40200MW	
5	BD-L50020MW	BD-L50040MW	BD-L50060MW	BD-L50080MW	BD-L50100MW	BD-L50120MW	BD-L50150MW	BD-L50200MW	
6	BD-L60020MW	BD-L60040MW	BD-L60060MW	BD-L60080MW	BD-L60100MW	BD-L60120MW	BD-L60150MW	BD-L60200MW	
7	BD-L70020MW	BD-L70040MW	BD-L70060MW	BD-L70080MW	BD-L70100MW	BD-L70120MW	BD-L70150MW	BD-L70200MW	
8	BD-L80020MW	BD-L80040MW	BD-L80060MW	BD-L80080MW	—	—	—	—	6Fr
9	BD-L90020MW	BD-L90040MW	BD-L90060MW	BD-L90080MW	—	—	—	—	
10	BD-LX0020MW	BD-LX0040MW	BD-LX0060MW	BD-LX0080MW	—	—	—	—	7Fr
12	BD-LX2020MW	BD-LX2040MW	—	—	—	—	—	—	

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 135 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120		150	200
3	BD-L30020CW	BD-L30040CW	BD-L30060CW	BD-L30080CW	BD-L30100CW	BD-L30120CW	BD-L30150CW	BD-L30200CW	5Fr
4	BD-L40020CW	BD-L40040CW	BD-L40060CW	BD-L40080CW	BD-L40100CW	BD-L40120CW	BD-L40150CW	BD-L40200CW	
5	BD-L50020CW	BD-L50040CW	BD-L50060CW	BD-L50080CW	BD-L50100CW	BD-L50120CW	BD-L50150CW	BD-L50200CW	
6	BD-L60020CW	BD-L60040CW	BD-L60060CW	BD-L60080CW	BD-L60100CW	BD-L60120CW	BD-L60150CW	BD-L60200CW	
7	BD-L70020CW	BD-L70040CW	BD-L70060CW	BD-L70080CW	BD-L70100CW	BD-L70120CW	BD-L70150CW	BD-L70200CW	
8	BD-L80020CW	BD-L80040CW	BD-L80060CW	BD-L80080CW	—	—	—	—	6Fr
9	BD-L90020CW	BD-L90040CW	BD-L90060CW	BD-L90080CW	—	—	—	—	
10	BD-LX0020CW	BD-LX0040CW	BD-LX0060CW	BD-LX0080CW	—	—	—	—	7Fr
12	BD-LX2020CW	BD-LX2040CW	—	—	—	—	—	—	

FIND OUT MORE



Phone: 800.862.4143



terumo.com



Fax: 800.411.5870

