

**PAGINA 0**

**CREATA PENTRU A  
PUTEA SEMNA  
ELECTRONIC  
DOCUMENTUL**



RELATED PRODUCTS

- Chronic Total Occlusion Devices
- Cardiovascular Products

REQUEST MORE INFORMATION

CONTACT US



OVERVIEW

Target a vessel's true lumen from the subintimal space by using the Enteer™ re-entry catheter to cross chronic total occlusions (CTOs). Leverage the self-orienting balloon that, when combined with an Enteer™ guidewire, enables re-entry. With two balloon sizes and three guidewire options, the Enteer™ re-entry system can be used to access and cross lesions above or below the knee, without any capital equipment investment.

INDICATIONS

The Enteer™ Re-entry Guidewire is intended to facilitate placement of balloon dilatation catheters or other intravascular devices during percutaneous transluminal angioplasty (PTA). The Enteer™ Guidewire is not to be used in cerebral blood vessels.

The Enteer™ Re-entry Catheter is indicated for directing, steering, controlling and supporting a guidewire in order to access discrete regions of the peripheral vasculature.

- case and minimises guidewire exchanges.
- With three offerings, the 28° Enteer™ guidewires are designed to find the Enteer™ catheter exit port and re-enter true lumen.
- When inflated, the flat shaped balloon self-oriens toward the true lumen within the subintimal space.



180° and offset exit ports on each side of device allow the Enteer™ guidewire to re-enter the correct port into the true lumen.



MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device.

MODEL SPECIFICATIONS

RE-ENTRY BALLOON CATHETER

Reference Number	Balloon Size (W × H × L mm)	Working Length (cm)	Guidewire Compatibility (in)	Crossing Profile (max in)
ENB-375-20-135	ATK 3.75 × 1.5 × 20	135	≤ 0.018	0.066
ENB-275-20-150	BTK 2.75 × 1.0 × 20	150	≤ 0.018	0.066

GUIDEWIRE

Reference Number	Description	Diameter (in)	Length (cm)	Tip Reach (mm)
ENW-FX-014-300	Flexible	0.014	300	1.5
ENW-SD-014-300	Standard	0.014	300	1.5
ENW-SF-014-300	Stiff	0.014	300	2.5

See product catalog for complete, detailed product information.

HEALTHCARE PROFESSIONALS

- Overview
- Products
- Therapies & Procedures
- Integrated Health Solutions

PATIENTS

- Overview

OUR COMPANY

- Who We Are
- Key Facts
- Leadership
- History
- Mission
- Locations
- Governance
- Careers

OUR IMPACT

- Our Impact
- Healthcare access
- Inclusion, Diversity & Equity
- Protecting our planet
- ESG Report