

PhysioGo 200A

Ultrasound therapy



Features

| product code | A-US-AST-PHG200A |
|--|------------------|
| color display with touch panel | 7" |
| independent treatment channels | 1 |
| manual mode | \checkmark |
| disease entities selected by name or medical field | \checkmark |
| preset treatment programs database | \checkmark |
| user-defined programs database | \checkmark |
| favorite programs | \checkmark |
| possibility of program names edition | \checkmark |
| encyclopedia describing the treatment methodology | \checkmark |
| statistics of performed treatment procedures | \checkmark |
| buzzer sound volume regulation | \checkmark |

Ultrasound therapy

| waterproof ultrasound heads | \checkmark |
|---|--------------|
| continuous / pulse emission | \checkmark |
| ultrasound head contact control (effective treatment time measured) | \checkmark |
| head sensitivity calibration according to the needs | \checkmark |

Preset treatment programs

| built-in treatment programs for ultrasound therapy |
|--|
| user configurable programs |
| favorite programs |

Ultrasound therapy technical parameters

| operating frequency | 1 & 3,5 MHz |
|-----------------------------------|-----------------------|
| effective radiation area | 1 cm², 4 cm² |
| maximum ultrasound wave intensity | 2/3 W/cm ² |
| frequency in pulse mode | 16 Hz, 48 Hz, 100 Hz |
| duty factor in pulse mode | 5 - 75 %, step 5% |
| treatment timer | 30 s - 30 minutes |

General technical parameters

dimensions

34 x 28 x 11-16 cm

58 50 \checkmark



| device weight | 6 kg |
|---------------------------------|---------------------------------|
| battery type | Li-lon |
| battery capacity | 2250 mAh |
| power supply, power consumption | 230 V, 50/60 Hz, 75 W, 90 VA |

