


Figure 3: Patient unit, Inspiratory section



Servo-i Ventilator Spare Parts List

Table 3: Patient unit, Inspiratory section

ID	Part no.	Description	Comments
60	6671137	Gas module O2, Type III	
60A	6693792	Nozzle unit	
60B	6039528	Rubber seal for bacteria filter	
60C	6169416	Bacteria Filter for Gas Module	
60D	6169358	Non-return valve, 2 pcs	
60E	6532654	Module adapter O2, with O-rings	Fits only modules with 2 screws in end plate
61	6671135	Gas module Air, Type III	
61A	6693792	Nozzle unit	
61B	6039528	Rubber seal for bacteria filter	
61C	6169416	Bacteria Filter for Gas Module	
61D	6169358	Non-return valve, 2 pcs	
61E	6532647	Module adapter AIR, with O-rings	Fits only modules wit 2 screws in the end plate
62	6488113	Module bracket	
63	6670330	Pull magnet incl. bracket	
64	9303892	Filter for O2 cell, 25 pcs	
65	6532662	Inspiratory pipe	
66	6532977	Tube	
67	6640044	O2 cell	
68	6448612	Connector muff	
69	6487958	Cable for O2 cell	
70	6888750	O2 cell holder	

LAP SUCTION IRRIGATION PUMP 2216

SUCTION IRRIGATION TUBE SETS

Irrigation			Product no.
Accessories – reusable	Tube set with piercing spike	with Luer-Lock connection, autoclavable, incl. 10 spare membranes, reusable for 20 reprocessing cycles	8171223
	Irrigation tube set with piercing spike, L 3 m	(PACK= 10 PCS), with Luer-Lock connection, sterile, disposable product	4171223
Accessories – disposable	Tube set with CARE-LOCK, L 3 m	(PACK= 10 PCS), sterile, disposable product	4171224
	Suction irrigation tube set	(PACK= 10 PCS), with piercing spike, sterile disposable product, PVC	4171225
	Day irrigation tube set	(PACK= 10 PCS), sterile, disposable product, PVC	4171226
	Patient irrigation tube set, sterile	(PACK= 10 PCS), with Luer-Lock connection, sterile, disposable product, PVC	4171227
Note: Day tube and patient tube can only be used together. A new patient tube is required for each patient; the day tube can be used for up to 20 patients			

Suction			Product no.
Accessories – reusable	Tube set for vacuum, incl. filter	with integrated hygiene filter, PACK= 10 PCS, service life: 30 days (connected to irrigation pump)	2206207
	Irrigation tube	L 1.2 m, with Luer-Lock, for suction irrigation handle 8385.901 and tube set 8171223 or 8170.223	817123
	Suction tube L 3 m	for suction irrigation handle 8385.901	8171402
	Secretion container, 2 l	2 liters, with lid	8170.655
	Secretion container, 2 l	2 liters, no lid	8170.6551
Accessories – disposable	Secretion container, 2.5 l	2.5 liters (PACK= 2 PCS)	41706566
	Secretion container bndl consisting of:	Secretion container 2.5 l (41706566) and secretion bag (PACK= 90 PCS) (41706564)	4170657
	Tube set for suction (2 connections)	sterile, disposable product, (PACK= 10 PCS) for connecting instruments to collecting container	T0503-01
	Angled 90° plug-in coupling	(PACK= 10 PCS) For connecting vacuum tube 2206207 to disposable container systems	0771019

Lotul : 212

Notes Help

2 Spare parts

Fit To Window

Pictures



Parts List

IR transmitter without charging unit

GETINGE

115091CE

IR transmitter without charging unit

ID	Part-No	Description
1	2270294	Battery pack 10,8V/ 550 mAh
	11360499	Tallow-drop screw DIN7985-M3x20-A2-50-H
	98019444	Escutcheon
	115091CE	IR transmitter without charging unit
100		Wiring diagram
110		Circuit diagram
120		Hydraulic diagram
210		service manual customer
240		recycling manual