


DOUBLE BLOOD BAG



This system is used for separation of two components from whole blood. This double system includes one primary bag with anticoagulant CPDA-1 Solutions USP and one empty satellite bag.

Sizes and Specification:



Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
DC-250W	250ml	1x200ml	16G	CPDA-1	35ml	5/80
DC-350W	350ml	1x300ml	16G	CPDA-1	49ml	5/80
DC-450W	450ml	1x300ml	16G	CPDA-1	63ml	5/80
DC-500W	500ml	1x300ml	16G	CPDA-1	70ml	5/80

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
DC-250E	250ml	1x200ml	16G	CPDA-1	35ml	5/80
DC-350E	350ml	1x300ml	16G	CPDA-1	49ml	5/80
DC-450E	450ml	1x300ml	16G	CPDA-1	63ml	5/80
DC-500E	500ml	1x300ml	16G	CPDA-1	70ml	5/80