



## Leptospira Medium Base, Korthof, Modified

M457

Leptospira Medium is used for isolation, cultivation and maintenance of *Leptospira* species.

### Composition\*\*

| Ingredients                      | Gms / Litre |
|----------------------------------|-------------|
| Peptic digest of animal tissue   | 0.800       |
| Sodium chloride                  | 1.400       |
| Sodium bicarbonate               | 0.020       |
| Potassium chloride               | 0.040       |
| Calcium chloride                 | 0.040       |
| Monopotassium hydrogen phosphate | 0.240       |
| Disodium hydrogen phosphate      | 0.880       |
| Final pH ( at 25°C)              | 7.2±0.2     |

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

#### 1)Preparation of Base:

Suspend 3.42 grams of M457 in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Distribute in 100 ml amounts in flasks. Sterilize by autoclaving at 115°C for 15 minutes. Cool to 55°C.

#### 2) Preparation of Haemoglobin Solution:

To the rabbit blood clot, after removing serum, add equal volume of distilled water. Freeze and thaw repeatedly to haemolyse the corpuscles. Sterilize by Seitz or millipore filtration.

#### 3) Complete Medium:

To 100 ml sterile base, add sterile 8 ml inactivated blood serum and 0.8 ml sterile haemoglobin solution. Mix thoroughly. Distribute if desired in 2-3 ml amount in sterile screw capped Bijou bottles/tubes. Test for sterility by incubating at 37°C.

### Principle And Interpretation

Leptospirosis is an acute febrile disease caused by members of the genus *Leptospira* (1,2). Direct culture of blood is the most reliable way to detect *Leptospira* during the first week of illness. After the first week of illness and for several months thereafter, leptospire may be isolated by direct culture of undiluted urine specimens. By autopsy, leptospire may be isolated from kidney and liver tissues as well as from blood and urine. Leptospira Medium Base, Korthof, Modified is formulated as described by Korthof (3, 4) for cultivation and maintenance of *Leptospira* species.

Peptic digest of animal tissue provide amino acids and other nitrogenous substances to support bacterial growth. Haemoglobin solution and inactivated blood serum provide additional sources of nutrients to the Leptospire. The salts supply essential nutrients for the growth of the organisms. Phosphates form buffering system while sodium chloride maintains osmotic equilibrium and also provides essential ions.

All cultures are incubated at room temperature in the dark for up to 6 weeks. The organisms grow below the surface. Material collected from a few centimeters below the surface of broth cultures should be examined weekly for the presence of growth using a direct wet preparation under dark field illumination. Leptospire will exhibit corkscrew like motility (1).

Examine the tubes for growth every 5-7 days. Growth occurs as a ringed area (disc) 1-3 cm below the surface of the medium. The absence of a ringed area of growth doesn't necessarily mean leptospire are not present. Remove a small amount of growth from the disc area and examine microscopically (gram stain is not satisfactory). Microcolonies can be fixed with methanol and stained with Giemsa's stain to show rod forms (3).

### Quality Control

**Appearance**

Off-white to yellow homogeneous free flowing powder

**Colour and Clarity of prepared medium**

Yellowish brown coloured, clear to slightly opalescent solution after addition of serum and haemoglobin

**Reaction**

Reaction of 0.342% w/v aqueous solution at 25°C. pH : 7.2±0.2

**pH**

7.00-7.40

**Cultural Response**

M457: Cultural characteristics observed with added inactivated blood serum and sterile haemoglobin solution, after an incubation at 30°C for upto 2-7days.

**Organism****Growth**

*Leptospira interrogans* luxuriant

*sero.grippotyphosa*

*Leptospira interrogans sero.* luxuriant

*Australis*

*Leptospira interrogans sero.* luxuriant

*Canicola*

**Storage and Shelf Life**

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

**Reference**

1. Forbes B. A., Sahm A. S., and Weissfeld D. F., Bailey & Scotts Diagnostic Microbiology, 10th Ed., 1998, Mosby, Inc., St. Louis, Mo.
2. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
3. Korthof G., 1932, Zentralbl. Bakteriol. Parasitenkd. Infektionskr. Hyg. Abt. I. Orig., 125:429.
4. Rechigl M. Jr. (Ed.), 1978, Handbook Series in Nutrition and Food, Vol. III, CRC Press.

Revision : 1 / 2011

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.