

EC Declaration of Conformity

For the following equipment:

U2 Total Knee System

(Product list is shown as following pages)

UMDNS code: 16-096 Prosthesis, Joint, Knee, Total

GMDN code: 33664 Cruciate-retaining total knee prosthesis

33665 Posterior-stabilized total knee prosthesis

(Model, Designation)

is herewith confirmed to comply with the requirements set out in the Council Directive on the harmonization of the Laws of the Member States concerning **Medical Device Directive 93/42/EEC As Amended by 2007/47/EC and 2005/50/EC** with the compliance of conformity assessment Annex II to be certified by British Standards Institution (notify body number – 2797).

For the evaluation regarding the Class III product safety aspects, the following **harmonized standards** are applied:

EN ISO 13485, EN ISO 14630, EN ISO 21534, EN ISO 21536, EN ISO 14971, EN ISO 11137-1, EN ISO 11137-2, EN ISO 11737-1, EN ISO 11737-2, EN ISO 11135-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 10993-1, EN ISO 10993-3, EN ISO 10993-5, EN ISO 10993-6, EN ISO 10993-7, EN ISO 10993-11, EN 1041, EN ISO 15223-1, EN 62366-1

The following European Authorized Representative is stated to the declaration:

mdi Europa GmbH / Langenhagener Straße 71, 30855 Langenhagen, Germany

(Company Name/ Address)

This declaration is issued under the sole responsibility of United Orthopedic Corporation. The following person is responsible for the compliance of declaration:

United Orthopedic Corporation/No. 57, Park Ave. 2, Science Park, Hsinchu City 30075, Taiwan

(Manufacturer Name/ Address)

Jason Lin / President

(Name/ Position)


(Legal Signature)

2020.09/15

(Date)

Catalog No.	Description
2401-1010	UKNEE Patella, in set, small
2401-1020	UKNEE Patella, in set, medium
2401-1030	UKNEE Patella, in set, large
2401-1040	UKNEE Patella, in set, X-large
2401-3010	UKNEE Patella, on set, small
2401-3020	UKNEE Patella, on set, medium
2401-3030	UKNEE Patella, on set, large
2103-1110	Femoral component, CR, porous coated, #1, left
2103-1120	Femoral component, CR, porous coated, #2, left
2103-1130	Femoral component, CR, porous coated, #3, left
2103-1140	Femoral component, CR, porous coated, #4, left
2103-1150	Femoral component, CR, porous coated, #5, left
2103-1160	Femoral component, CR, porous coated, #6, left
2103-1170	Femoral component, CR, porous coated, #7, left
2103-1115	Femoral component, CR, porous coated, #1.5, left
2103-1125	Femoral component, CR, porous coated, #2.5, left
2103-1135	Femoral component, CR, porous coated, #3.5, left
2103-1145	Femoral component, CR, porous coated, #4.5, left
2103-1155	Femoral component, CR, porous coated, #5.5, left
2103-1165	Femoral component, CR, porous coated, #6.5, left
2103-1210	Femoral component, CR, porous coated, #1, right
2103-1220	Femoral component, CR, porous coated, #2, right
2103-1230	Femoral component, CR, porous coated, #3, right
2103-1240	Femoral component, CR, porous coated, #4, right
2103-1250	Femoral component, CR, porous coated, #5, right
2103-1260	Femoral component, CR, porous coated, #6, right
2103-1270	Femoral component, CR, porous coated, #7, right
2103-1215	Femoral component, CR, porous coated, #1.5, right
2103-1225	Femoral component, CR, porous coated, #2.5, right
2103-1235	Femoral component, CR, porous coated, #3.5, right
2103-1245	Femoral component, CR, porous coated, #4.5, right
2103-1255	Femoral component, CR, porous coated, #5.5, right
2103-1265	Femoral component, CR, porous coated, #6.5, right
2103-1310	Femoral component, CR, cemented, #1, left
2103-1320	Femoral component, CR, cemented, #2, left
2103-1330	Femoral component, CR, cemented, #3, left
2103-1340	Femoral component, CR, cemented, #4, left
2103-1350	Femoral component, CR, cemented, #5, left

Catalog No.	Description
2103-1360	Femoral component, CR, cemented, #6, left
2103-1370	Femoral component, CR, cemented, #7, left
2103-1315	Femoral component, CR, cemented, #1.5, left
2103-1325	Femoral component, CR, cemented, #2.5, left
2103-1335	Femoral component, CR, cemented, #3.5, left
2103-1345	Femoral component, CR, cemented, #4.5, left
2103-1355	Femoral component, CR, cemented, #5.5, left
2103-1365	Femoral component, CR, cemented, #6.5, left
2103-1410	Femoral component, CR, cemented, #1, right
2103-1420	Femoral component, CR, cemented, #2, right
2103-1430	Femoral component, CR, cemented, #3, right
2103-1440	Femoral component, CR, cemented, #4, right
2103-1450	Femoral component, CR, cemented, #5, right
2103-1460	Femoral component, CR, cemented, #6, right
2103-1470	Femoral component, CR, cemented, #7, right
2103-1415	Femoral component, CR, cemented, #1.5, right
2103-1425	Femoral component, CR, cemented, #2.5, right
2103-1435	Femoral component, CR, cemented, #3.5, right
2103-1445	Femoral component, CR, cemented, #4.5, right
2103-1455	Femoral component, CR, cemented, #5.5, right
2103-1465	Femoral component, CR, cemented, #6.5, right
2103-3110	Femoral component, PS, , #1, left
2103-3120	Femoral component, PS, , #2, left
2103-3130	Femoral component, PS, , #3, left
2103-3140	Femoral component, PS, , #4, left
2103-3150	Femoral component, PS, , #5, left
2103-3160	Femoral component, PS, , #6, left
2103-3170	Femoral component, PS, , #7, left
2103-3115	Femoral component, PS, #1.5, left
2103-3125	Femoral component, PS, #2.5, left
2103-3135	Femoral component, PS, #3.5, left
2103-3145	Femoral component, PS, #4.5, left
2103-3155	Femoral component, PS, #5.5, left
2103-3165	Femoral component, PS, #6.5, left
2103-3210	Femoral component, PS, , #1, right
2103-3220	Femoral component, PS, , #2, right
2103-3230	Femoral component, PS, , #3, right
2103-3240	Femoral component, PS, , #4, right
2103-3250	Femoral component, PS, , #5, right

Catalog No.	Description
2103-3260	Femoral component, PS, , #6, right
2103-3270	Femoral component, PS, , #7, right
2103-3215	Femoral component, PS, #1.5, right
2103-3225	Femoral component, PS, #2.5, right
2103-3235	Femoral component, PS, #3.5, right
2103-3245	Femoral component, PS, #4.5, right
2103-3255	Femoral component, PS, #5.5, right
2103-3265	Femoral component, PS, #6.5, right
2203-3010	Tibial baseplate, cemented, #1
2203-3020	Tibial baseplate, cemented, #2
2203-3030	Tibial baseplate, cemented, #3
2203-3040	Tibial baseplate, cemented, #4
2203-3050	Tibial baseplate, cemented, #5
2203-3060	Tibial baseplate, cemented, #6
2203-3200	Tibial baseplate, CMA, #0
2203-3210	Tibial baseplate, CMA, #1
2203-3220	Tibial baseplate, CMA, #2
2203-3230	Tibial baseplate, CMA, #3
2203-3240	Tibial baseplate, CMA, #4
2203-3250	Tibial baseplate, CMA, #5
2203-3260	Tibial baseplate, CMA, #6
2203-3270	Tibial baseplate, CMA, #7
2303-1201	Tibial insert, CR, plained, #0, 9mm thick
2303-1202	Tibial insert, CR, plained, #0, 11mm thick
2303-1203	Tibial insert, CR, plained, #0, 13mm thick
2303-1204	Tibial insert, CR, plained, #0, 15mm thick
2303-1205	Tibial insert, CR, plained, #0, 18mm thick
2303-1200	Tibial insert, CR, plained, #0, 17mm thick
2303-1206	Tibial insert, CR, plained, #0, 10mm thick
2303-1207	Tibial insert, CR, plained, #0, 12mm thick
2303-1208	Tibial insert, CR, plained, #0, 14mm thick
2303-1209	Tibial insert, CR, plained, #0, 16mm thick
2303-1211	Tibial insert, CR, plained, #1, 9mm thick
2303-1212	Tibial insert, CR, plained, #1, 11mm thick
2303-1213	Tibial insert, CR, plained, #1, 13mm thick
2303-1214	Tibial insert, CR, plained, #1, 15mm thick
2303-1215	Tibial insert, CR, plained, #1, 18mm thick
2303-1210	Tibial insert, CR, plained, #1, 17mm thick
2303-1216	Tibial insert, CR, plained, #1, 10mm thick

Catalog No.	Description
2303-1217	Tibial insert, CR, plained, #1, 12mm thick
2303-1218	Tibial insert, CR, plained, #1, 14mm thick
2303-1219	Tibial insert, CR, plained, #1, 16mm thick
2303-1221	Tibial insert, CR, plained, #2, 9mm thick
2303-1222	Tibial insert, CR, plained, #2, 11mm thick
2303-1223	Tibial insert, CR, plained, #2, 13mm thick
2303-1224	Tibial insert, CR, plained, #2, 15mm thick
2303-1225	Tibial insert, CR, plained, #2, 18mm thick
2303-1220	Tibial insert, CR, plained, #2, 17mm thick
2303-1226	Tibial insert, CR, plained, #2, 10mm thick
2303-1227	Tibial insert, CR, plained, #2, 12mm thick
2303-1228	Tibial insert, CR, plained, #2, 14mm thick
2303-1229	Tibial insert, CR, plained, #2, 16mm thick
2303-1231	Tibial insert, CR, plained, #3, 9mm thick
2303-1232	Tibial insert, CR, plained, #3, 11mm thick
2303-1233	Tibial insert, CR, plained, #3, 13mm thick
2303-1234	Tibial insert, CR, plained, #3, 15mm thick
2303-1235	Tibial insert, CR, plained, #3, 18mm thick
2303-1230	Tibial insert, CR, plained, #3, 17mm thick
2303-1236	Tibial insert, CR, plained, #3, 10mm thick
2303-1237	Tibial insert, CR, plained, #3, 12mm thick
2303-1238	Tibial insert, CR, plained, #3, 14mm thick
2303-1239	Tibial insert, CR, plained, #3, 16mm thick
2303-1241	Tibial insert, CR, plained, #4, 9mm thick
2303-1242	Tibial insert, CR, plained, #4, 11mm thick
2303-1243	Tibial insert, CR, plained, #4, 13mm thick
2303-1244	Tibial insert, CR, plained, #4, 15mm thick
2303-1245	Tibial insert, CR, plained, #4, 18mm thick
2303-1240	Tibial insert, CR, plained, #4, 17mm thick
2303-1246	Tibial insert, CR, plained, #4, 10mm thick
2303-1247	Tibial insert, CR, plained, #4, 12mm thick
2303-1248	Tibial insert, CR, plained, #4, 14mm thick
2303-1249	Tibial insert, CR, plained, #4, 16mm thick
2303-1251	Tibial insert, CR, plained, #5, 9mm thick
2303-1252	Tibial insert, CR, plained, #5, 11mm thick
2303-1253	Tibial insert, CR, plained, #5, 13mm thick
2303-1254	Tibial insert, CR, plained, #5, 15mm thick
2303-1255	Tibial insert, CR, plained, #5, 18mm thick
2303-1250	Tibial insert, CR, plained, #5, 17mm thick

Catalog No.	Description
2303-1256	Tibial insert, CR, plained, #5, 10mm thick
2303-1257	Tibial insert, CR, plained, #5, 12mm thick
2303-1258	Tibial insert, CR, plained, #5, 14mm thick
2303-1259	Tibial insert, CR, plained, #5, 16mm thick
2303-1261	Tibial insert, CR, plained, #6, 9mm thick
2303-1262	Tibial insert, CR, plained, #6, 11mm thick
2303-1263	Tibial insert, CR, plained, #6, 13mm thick
2303-1264	Tibial insert, CR, plained, #6, 15mm thick
2303-1265	Tibial insert, CR, plained, #6, 18mm thick
2303-1260	Tibial insert, CR, plained, #6, 17mm thick
2303-1266	Tibial insert, CR, plained, #6, 10mm thick
2303-1267	Tibial insert, CR, plained, #6, 12mm thick
2303-1268	Tibial insert, CR, plained, #6, 14mm thick
2303-1269	Tibial insert, CR, plained, #6, 16mm thick
2303-1271	Tiaial insert, CR, plained, #7, 9mm thick
2303-1272	Tiaial insert, CR, plained, #7, 11mm thick
2303-1273	Tiaial insert, CR, plained, #7, 13mm thick
2303-1274	Tiaial insert, CR, plained, #7, 15mm thick
2303-1275	Tiaial insert, CR, plained, #7, 18mm thick
2303-1270	Tibial insert, CR, plained, #7, 17mm thick
2303-1276	Tibial insert, CR, plained, #7, 10mm thick
2303-1277	Tibial insert, CR, plained, #7, 12mm thick
2303-1278	Tibial insert, CR, plained, #7, 14mm thick
2303-1279	Tibial insert, CR, plained, #7, 16mm thick
2303-1601	XPE Tibial insert, CR, #0, 9mm thick
2303-1602	XPE Tibial insert, CR, #0, 11mm thick
2303-1603	XPE Tibial insert, CR, #0, 13mm thick
2303-1604	XPE Tibial insert, CR, #0, 15mm thick
2303-1605	XPE Tibial insert, CR, #0, 18mm thick
2303-1606	XPE Tibial insert, CR, #0, 10mm thick
2303-1607	XPE Tibial insert, CR, #0, 12mm thick
2303-1608	XPE Tibial insert, CR, #0, 14mm thick
2303-1609	XPE Tibial insert, CR, #0, 16mm thick
2303-1600	XPE Tibial insert, CR, #0, 17mm thick
2303-1611	XPE Tibial Insert, CR, #1, 9mm thick
2303-1612	XPE Tibial Insert, CR, #1, 11mm thick
2303-1613	XPE Tibial Insert, CR, #1, 13mm thick
2303-1614	XPE Tibial Insert, CR, #1, 15mm thick
2303-1615	XPE Tibial Insert, CR, #1, 18mm thick

Catalog No.	Description
2303-1616	XPE Tibial insert, CR, #1, 10mm thick
2303-1617	XPE Tibial insert, CR, #1, 12mm thick
2303-1618	XPE Tibial insert, CR, #1, 14mm thick
2303-1619	XPE Tibial insert, CR, #1, 16mm thick
2303-1610	XPE Tibial insert, CR, #1, 17mm thick
2303-1621	XPE Tibial Insert, CR, #2, 9mm thick
2303-1622	XPE Tibial Insert, CR, #2, 11mm thick
2303-1623	XPE Tibial Insert, CR, #2, 13mm thick
2303-1624	XPE Tibial Insert, CR, #2, 15mm thick
2303-1625	XPE Tibial Insert, CR, #2, 18mm thick
2303-1626	XPE Tibial insert, CR, #2, 10mm thick
2303-1627	XPE Tibial insert, CR, #2, 12mm thick
2303-1628	XPE Tibial insert, CR, #2, 14mm thick
2303-1629	XPE Tibial insert, CR, #2, 16mm thick
2303-1620	XPE Tibial insert, CR, #2, 17mm thick
2303-1631	XPE Tibial Insert, CR, #3, 9mm thick
2303-1632	XPE Tibial Insert, CR, #3, 11mm thick
2303-1633	XPE Tibial Insert, CR, #3, 13mm thick
2303-1634	XPE Tibial Insert, CR, #3, 15mm thick
2303-1635	XPE Tibial Insert, CR, #3, 18mm thick
2303-1636	XPE Tibial insert, CR, #3, 10mm thick
2303-1637	XPE Tibial insert, CR, #3, 12mm thick
2303-1638	XPE Tibial insert, CR, #3, 14mm thick
2303-1639	XPE Tibial insert, CR, #3, 16mm thick
2303-1630	XPE Tibial insert, CR, #3, 17mm thick
2303-1641	XPE Tibial Insert, CR, #4, 9mm thick
2303-1642	XPE Tibial Insert, CR, #4, 11mm thick
2303-1643	XPE Tibial Insert, CR, #4, 13mm thick
2303-1644	XPE Tibial Insert, CR, #4, 15mm thick
2303-1645	XPE Tibial Insert, CR, #4, 18mm thick
2303-1646	XPE Tibial insert, CR, #4, 10mm thick
2303-1647	XPE Tibial insert, CR, #4, 12mm thick
2303-1648	XPE Tibial insert, CR, #4, 14mm thick
2303-1649	XPE Tibial insert, CR, #4, 16mm thick
2303-1640	XPE Tibial insert, CR, #4, 17mm thick
2303-1651	XPE Tibial Insert, CR, #5, 9mm thick
2303-1652	XPE Tibial Insert, CR, #5, 11mm thick
2303-1653	XPE Tibial Insert, CR, #5, 13mm thick
2303-1654	XPE Tibial Insert, CR, #5, 15mm thick

Catalog No.	Description
2303-1655	XPE Tibial Insert, CR, #5, 18mm thick
2303-1656	XPE Tibial insert, CR, #5, 10mm thick
2303-1657	XPE Tibial insert, CR, #5, 12mm thick
2303-1658	XPE Tibial insert, CR, #5, 14mm thick
2303-1659	XPE Tibial insert, CR, #5, 16mm thick
2303-1650	XPE Tibial insert, CR, #5, 17mm thick
2303-1661	XPE Tibial Insert, CR, #6, 9mm thick
2303-1662	XPE Tibial Insert, CR, #6, 11mm thick
2303-1663	XPE Tibial Insert, CR, #6, 13mm thick
2303-1664	XPE Tibial Insert, CR, #6, 15mm thick
2303-1665	XPE Tibial Insert, CR, #6, 18mm thick
2303-1666	XPE Tibial insert, CR, #6, 10mm thick
2303-1667	XPE Tibial insert, CR, #6, 12mm thick
2303-1668	XPE Tibial insert, CR, #6, 14mm thick
2303-1669	XPE Tibial insert, CR, #6, 16mm thick
2303-1660	XPE Tibial insert, CR, #6, 17mm thick
2303-1671	XPE Tibial insert, CR, #7, 9mm thick
2303-1672	XPE Tibial insert, CR, #7, 11mm thick
2303-1673	XPE Tibial insert, CR, #7, 13mm thick
2303-1674	XPE Tibial insert, CR, #7, 15mm thick
2303-1675	XPE Tibial insert, CR, #7, 18mm thick
2303-1676	XPE Tibial insert, CR, #7, 10mm thick
2303-1677	XPE Tibial insert, CR, #7, 12mm thick
2303-1678	XPE Tibial insert, CR, #7, 14mm thick
2303-1679	XPE Tibial insert, CR, #7, 16mm thick
2303-1670	XPE Tibial insert, CR, #7, 17mm thick
2303-3001	Tibial insert, PS, #0, 9mm thick
2303-3002	Tibial insert, PS, #0, 11mm thick
2303-3003	Tibial insert, PS, #0, 13mm thick
2303-3004	Tibial insert, PS, #0, 15mm thick
2303-3005	Tibial insert, PS, #0, 18mm thick
2303-3006	Tibial insert, PS, #0, 10mm thick
2303-3007	Tibial insert, PS, #0, 12mm thick
2303-3008	Tibial insert, PS, #0, 14mm thick
2303-3009	Tibial insert, PS, #0, 16mm thick
2303-3000	Tibial insert, PS, #0, 17mm thick
2303-3011	Tibial insert, PS, #1, 9mm thick
2303-3012	Tibial insert, PS, #1, 11mm thick

Catalog No.	Description
2303-3013	Tibial insert, PS, #1, 13mm thick
2303-3014	Tibial insert, PS, #1, 15mm thick
2303-3015	Tibial insert, PS, #1, 18mm thick
2303-3016	Tibial insert, PS, #1, 10mm thick
2303-3017	Tibial insert, PS, #1, 12mm thick
2303-3018	Tibial insert, PS, #1, 14mm thick
2303-3019	Tibial insert, PS, #1, 16mm thick
2303-3010	Tibial insert, PS, #1, 17mm thick
2303-3021	Tibial insert, PS, #2, 9mm thick
2303-3022	Tibial insert, PS, #2, 11mm thick
2303-3023	Tibial insert, PS, #2, 13mm thick
2303-3024	Tibial insert, PS, #2, 15mm thick
2303-3025	Tibial insert, PS, #2, 18mm thick
2303-3026	Tibial insert, PS, #2, 10mm thick
2303-3027	Tibial insert, PS, #2, 12mm thick
2303-3028	Tibial insert, PS, #2, 14mm thick
2303-3029	Tibial insert, PS, #2, 16mm thick
2303-3020	Tibial insert, PS, #2, 17mm thick
2303-3031	Tibial insert, PS, #3, 9mm thick
2303-3032	Tibial insert, PS, #3, 11mm thick
2303-3033	Tibial insert, PS, #3, 13mm thick
2303-3034	Tibial insert, PS, #3, 15mm thick
2303-3035	Tibial insert, PS, #3, 18mm thick
2303-3036	Tibial insert, PS, #3, 10mm thick
2303-3037	Tibial insert, PS, #3, 12mm thick
2303-3038	Tibial insert, PS, #3, 14mm thick
2303-3039	Tibial insert, PS, #3, 16mm thick
2303-3030	Tibial insert, PS, #3, 17mm thick
2303-3041	Tibial insert, PS, #4, 9mm thick
2303-3042	Tibial insert, PS, #4, 11mm thick
2303-3043	Tibial insert, PS, #4, 13mm thick
2303-3044	Tibial insert, PS, #4, 15mm thick
2303-3045	Tibial insert, PS, #4, 18mm thick
2303-3046	Tibial insert, PS, #4, 10mm thick
2303-3047	Tibial insert, PS, #4, 12mm thick
2303-3048	Tibial insert, PS, #4, 14mm thick
2303-3049	Tibial insert, PS, #4, 16mm thick
2303-3040	Tibial insert, PS, #4, 17mm thick
2303-3051	Tibial insert, PS, #5, 9mm thick

Catalog No.	Description
2303-3052	Tibial insert, PS, #5, 11mm thick
2303-3053	Tibial insert, PS, #5, 13mm thick
2303-3054	Tibial insert, PS, #5, 15mm thick
2303-3055	Tibial insert, PS, #5, 18mm thick
2303-3056	Tibial insert, PS, #5, 10mm thick
2303-3057	Tibial insert, PS, #5, 12mm thick
2303-3058	Tibial insert, PS, #5, 14mm thick
2303-3059	Tibial insert, PS, #5, 16mm thick
2303-3050	Tibial insert, PS, #5, 17mm thick
2303-3061	Tibial insert, PS, #6, 9mm thick
2303-3062	Tibial insert, PS, #6, 11mm thick
2303-3063	Tibial insert, PS, #6, 13mm thick
2303-3064	Tibial insert, PS, #6, 15mm thick
2303-3065	Tibial insert, PS, #6, 18mm thick
2303-3066	Tibial insert, PS, #6, 10mm thick
2303-3067	Tibial insert, PS, #6, 12mm thick
2303-3068	Tibial insert, PS, #6, 14mm thick
2303-3069	Tibial insert, PS, #6, 16mm thick
2303-3060	Tibial insert, PS, #6, 17mm thick
2303-3071	Tibial insert, PS, #7, 9mm thick
2303-3072	Tibial insert, PS, #7, 11mm thick
2303-3073	Tibial insert, PS, #7, 13mm thick
2303-3074	Tibial insert, PS, #7, 15mm thick
2303-3075	Tibial insert, PS, #7, 18mm thick
2303-3076	Tibial insert, PS, #7, 10mm thick
2303-3077	Tibial insert, PS, #7, 12mm thick
2303-3078	Tibial insert, PS, #7, 14mm thick
2303-3079	Tibial insert, PS, #7, 16mm thick
2303-3070	Tibial insert, PS, #7, 17mm thick
2303-3601	XPE Tibial insert, PS, #0, 9mm thick
2303-3602	XPE Tibial insert, PS, #0, 11mm thick
2303-3603	XPE Tibial insert, PS, #0, 13mm thick
2303-3604	XPE Tibial insert, PS, #0, 15mm thick
2303-3605	XPE Tibial insert, PS, #0, 18mm thick
2303-3606	XPE Tibial insert, PS, #0, 10mm thick
2303-3607	XPE Tibial insert, PS, #0, 12mm thick
2303-3608	XPE Tibial insert, PS, #0, 14mm thick
2303-3609	XPE Tibial insert, PS, #0, 16mm thick
2303-3600	XPE Tibial insert, PS, #0, 17mm thick

Catalog No.	Description
2303-3611	XPE Tibial Insert, PS, #1, 9mm thick
2303-3612	XPE Tibial Insert, PS, #1, 11mm thick
2303-3613	XPE Tibial Insert, PS, #1, 13mm thick
2303-3614	XPE Tibial Insert, PS, #1, 15mm thick
2303-3615	XPE Tibial Insert, PS, #1, 18mm thick
2303-3616	XPE Tibial insert, PS, #1, 10mm thick
2303-3617	XPE Tibial insert, PS, #1, 12mm thick
2303-3618	XPE Tibial insert, PS, #1, 14mm thick
2303-3619	XPE Tibial insert, PS, #1, 16mm thick
2303-3610	XPE Tibial insert, PS, #1, 17mm thick
2303-3621	XPE Tibial Insert, PS, #2, 9mm thick
2303-3622	XPE Tibial Insert, PS, #2, 11mm thick
2303-3623	XPE Tibial Insert, PS, #2, 13mm thick
2303-3624	XPE Tibial Insert, PS, #2, 15mm thick
2303-3625	XPE Tibial Insert, PS, #2, 18mm thick
2303-3626	XPE Tibial insert, PS, #2, 10mm thick
2303-3627	XPE Tibial insert, PS, #2, 12mm thick
2303-3628	XPE Tibial insert, PS, #2, 14mm thick
2303-3629	XPE Tibial insert, PS, #2, 16mm thick
2303-3620	XPE Tibial insert, PS, #2, 17mm thick
2303-3631	XPE Tibial Insert, PS, #3, 9mm thick
2303-3632	XPE Tibial Insert, PS, #3, 11mm thick
2303-3633	XPE Tibial Insert, PS, #3, 13mm thick
2303-3634	XPE Tibial Insert, PS, #3, 15mm thick
2303-3635	XPE Tibial Insert, PS, #3, 18mm thick
2303-3636	XPE Tibial insert, PS, #3, 10mm thick
2303-3637	XPE Tibial insert, PS, #3, 12mm thick
2303-3638	XPE Tibial insert, PS, #3, 14mm thick
2303-3639	XPE Tibial insert, PS, #3, 16mm thick
2303-3630	XPE Tibial insert, PS, #3, 17mm thick
2303-3641	XPE Tibial Insert, PS, #4, 9mm thick
2303-3642	XPE Tibial Insert, PS, #4, 11mm thick
2303-3643	XPE Tibial Insert, PS, #4, 13mm thick
2303-3644	XPE Tibial Insert, PS, #4, 15mm thick
2303-3645	XPE Tibial Insert, PS, #4, 18mm thick
2303-3646	XPE Tibial insert, PS, #4, 10mm thick
2303-3647	XPE Tibial insert, PS, #4, 12mm thick
2303-3648	XPE Tibial insert, PS, #4, 14mm thick
2303-3649	XPE Tibial insert, PS, #4, 16mm thick

Catalog No.	Description
2303-3640	XPE Tibial insert, PS, #4, 17mm thick
2303-3651	XPE Tibial Insert, PS, #5, 9mm thick
2303-3652	XPE Tibial Insert, PS, #5, 11mm thick
2303-3653	XPE Tibial Insert, PS, #5, 13mm thick
2303-3654	XPE Tibial Insert, PS, #5, 15mm thick
2303-3655	XPE Tibial Insert, PS, #5, 18mm thick
2303-3656	XPE Tibial insert, PS, #5, 10mm thick
2303-3657	XPE Tibial insert, PS, #5, 12mm thick
2303-3658	XPE Tibial insert, PS, #5, 14mm thick
2303-3659	XPE Tibial insert, PS, #5, 16mm thick
2303-3650	XPE Tibial insert, PS, #5, 17mm thick
2303-3661	XPE Tibial Insert, PS, #6, 9mm thick
2303-3662	XPE Tibial Insert, PS, #6, 11mm thick
2303-3663	XPE Tibial Insert, PS, #6, 13mm thick
2303-3664	XPE Tibial Insert, PS, #6, 15mm thick
2303-3665	XPE Tibial Insert, PS, #6, 18mm thick
2303-3666	XPE Tibial insert, PS, #6, 10mm thick
2303-3667	XPE Tibial insert, PS, #6, 12mm thick
2303-3668	XPE Tibial insert, PS, #6, 14mm thick
2303-3669	XPE Tibial insert, PS, #6, 16mm thick
2303-3660	XPE Tibial insert, PS, #6, 17mm thick
2303-3671	XPE Tibial insert, PS, #7, 9mm thick
2303-3672	XPE Tibial insert, PS, #7, 11mm thick
2303-3673	XPE Tibial insert, PS, #7, 13mm thick
2303-3674	XPE Tibial insert, PS, #7, 15mm thick
2303-3675	XPE Tibial insert, PS, #7, 18mm thick
2303-3676	XPE Tibial insert, PS, #7, 10mm thick
2303-3677	XPE Tibial insert, PS, #7, 12mm thick
2303-3678	XPE Tibial insert, PS, #7, 14mm thick
2303-3679	XPE Tibial insert, PS, #7, 16mm thick
2303-3670	XPE Tibial insert, PS, #7, 17mm thick
2303-1401	XPE Tibial insert, UC, #0, 9mm thick
2303-1402	XPE Tibial insert, UC, #0, 11mm thick
2303-1403	XPE Tibial insert, UC, #0, 13mm thick
2303-1404	XPE Tibial insert, UC, #0, 15mm thick
2303-1405	XPE Tibial insert, UC, #0, 18mm thick
2303-1411	XPE Tibial insert, UC, #1, 9mm thick
2303-1412	XPE Tibial insert, UC, #1, 11mm thick
2303-1413	XPE Tibial insert, UC, #1, 13mm thick

Catalog No.	Description
2303-1414	XPE Tibial insert, UC, #1, 15mm thick
2303-1415	XPE Tibial insert, UC, #1, 18mm thick
2303-1421	XPE Tibial insert, UC, #2, 9mm thick
2303-1422	XPE Tibial insert, UC, #2, 11mm thick
2303-1423	XPE Tibial insert, UC, #2, 13mm thick
2303-1424	XPE Tibial insert, UC, #2, 15mm thick
2303-1425	XPE Tibial insert, UC, #2, 18mm thick
2303-1431	XPE Tibial insert, UC, #3, 9mm thick
2303-1432	XPE Tibial insert, UC, #3, 11mm thick
2303-1433	XPE Tibial insert, UC, #3, 13mm thick
2303-1434	XPE Tibial insert, UC, #3, 15mm thick
2303-1435	XPE Tibial insert, UC, #3, 18mm thick
2303-1441	XPE Tibial insert, UC, #4, 9mm thick
2303-1442	XPE Tibial insert, UC, #4, 11mm thick
2303-1443	XPE Tibial insert, UC, #4, 13mm thick
2303-1444	XPE Tibial insert, UC, #4, 15mm thick
2303-1445	XPE Tibial insert, UC, #4, 18mm thick
2303-1451	XPE Tibial insert, UC, #5, 9mm thick
2303-1452	XPE Tibial insert, UC, #5, 11mm thick
2303-1453	XPE Tibial insert, UC, #5, 13mm thick
2303-1454	XPE Tibial insert, UC, #5, 15mm thick
2303-1455	XPE Tibial insert, UC, #5, 18mm thick
2303-1461	XPE Tibial insert, UC, #6, 9mm thick
2303-1462	XPE Tibial insert, UC, #6, 11mm thick
2303-1463	XPE Tibial insert, UC, #6, 13mm thick
2303-1464	XPE Tibial insert, UC, #6, 15mm thick
2303-1465	XPE Tibial insert, UC, #6, 18mm thick
2303-1471	XPE Tibial insert, UC, #7, 9mm thick
2303-1472	XPE Tibial insert, UC, #7, 11mm thick
2303-1473	XPE Tibial insert, UC, #7, 13mm thick
2303-1474	XPE Tibial insert, UC, #7, 15mm thick
2303-1475	XPE Tibial insert, UC, #7, 18mm thick
2303-1406	XPE Tibial insert, UC, #0, 10mm thick
2303-1407	XPE Tibial insert, UC, #0, 12mm thick
2303-1408	XPE Tibial insert, UC, #0, 14mm thick
2303-1409	XPE Tibial insert, UC, #0, 16mm thick
2303-1400	XPE Tibial insert, UC, #0, 17mm thick
2303-1416	XPE Tibial insert, UC, #1, 10mm thick
2303-1417	XPE Tibial insert, UC, #1, 12mm thick

Catalog No.	Description
2303-1418	XPE Tibial insert, UC, #1, 14mm thick
2303-1419	XPE Tibial insert, UC, #1, 16mm thick
2303-1410	XPE Tibial insert, UC, #1, 17mm thick
2303-1426	XPE Tibial insert, UC, #2, 10mm thick
2303-1427	XPE Tibial insert, UC, #2, 12mm thick
2303-1428	XPE Tibial insert, UC, #2, 14mm thick
2303-1429	XPE Tibial insert, UC, #2, 16mm thick
2303-1420	XPE Tibial insert, UC, #2, 17mm thick
2303-1436	XPE Tibial insert, UC, #3, 10mm thick
2303-1437	XPE Tibial insert, UC, #3, 12mm thick
2303-1438	XPE Tibial insert, UC, #3, 14mm thick
2303-1439	XPE Tibial insert, UC, #3, 16mm thick
2303-1430	XPE Tibial insert, UC, #3, 17mm thick
2303-1446	XPE Tibial insert, UC, #4, 10mm thick
2303-1447	XPE Tibial insert, UC, #4, 12mm thick
2303-1448	XPE Tibial insert, UC, #4, 14mm thick
2303-1449	XPE Tibial insert, UC, #4, 16mm thick
2303-1440	XPE Tibial insert, UC, #4, 17mm thick
2303-1456	XPE Tibial insert, UC, #5, 10mm thick
2303-1457	XPE Tibial insert, UC, #5, 12mm thick
2303-1458	XPE Tibial insert, UC, #5, 14mm thick
2303-1459	XPE Tibial insert, UC, #5, 16mm thick
2303-1450	XPE Tibial insert, UC, #5, 17mm thick
2303-1466	XPE Tibial insert, UC, #6, 10mm thick
2303-1467	XPE Tibial insert, UC, #6, 12mm thick
2303-1468	XPE Tibial insert, UC, #6, 14mm thick
2303-1469	XPE Tibial insert, UC, #6, 16mm thick
2303-1460	XPE Tibial insert, UC, #6, 17mm thick
2303-1476	XPE Tibial insert, UC, #7, 10mm thick
2303-1477	XPE Tibial insert, UC, #7, 12mm thick
2303-1478	XPE Tibial insert, UC, #7, 14mm thick
2303-1479	XPE Tibial insert, UC, #7, 16mm thick
2303-1470	XPE Tibial insert, UC, #7, 17mm thick
2403-1010	Patella, On set, X-Small, 26 mm Dia. 7mm thick
2403-1020	Patella, On set, Small, 29 mm Dia. 8mm thick
2403-1030	Patella, On set, Medium, 32 mm Dia. 8.5mm thick
2403-1040	Patella, On set, Large, 35 mm Dia. 9mm thick
2403-1050	Patella, On set, X-Large, 38 mm Dia. 9.5mm thick
2403-1060	Patella, On set, XX-Large, 41 mm Dia. 10mm thick

Catalog No.	Description
2403-1070	Patella, On set, E-Large, 44 mm Dia. 10.5mm thick
2403-3210	XPE Patella, On set, X-Small, 26 mm Dia. 7mm thick
2403-3220	XPE Patella, On set, Small, 29 mm Dia. 8mm thick
2403-3230	XPE Patella, On set, Medium, 32 mm Dia. 8.5mm thick
2403-3240	XPE Patella, On set, Large, 35 mm Dia. 9mm thick
2403-3250	XPE Patella, On set, X-Large, 38 mm Dia. 9.5mm thick
2403-3260	XPE Patella, On set, XX-Large, 41 mm Dia. 10mm thick
2403-3270	XPE Patella, On set, E-Large, 44 mm Dia. 10.5mm thick
2403-3010	XPE Patella, Small, 22mm Dia. 8mm thick
2403-3020	XPE Patella, Medium, 25mm Dia. 10mm thick
2403-3030	XPE Patella, Large, 28mm Dia. 10mm thick
2403-3040	XPE Patella, X-Large, 32mm Dia. 10mm thick
2203-1001	All poly tibial component, CR, #0, 9mm thick
2203-1002	All poly tibial component, CR, #0, 11mm thick
2203-1003	All poly tibial component, CR, #0, 13mm thick
2203-1004	All poly tibial component, CR, #0, 15mm thick
2203-1005	All poly tibial component, CR, #0, 18mm thick
2203-1011	All poly tibial component, CR, #1, 9mm thick
2203-1012	All poly tibial component, CR, #1, 11mm thick
2203-1013	All poly tibial component, CR, #1, 13mm thick
2203-1014	All poly tibial component, CR, #1, 15mm thick
2203-1015	All poly tibial component, CR, #1, 18mm thick
2203-1016	All poly tibial component, CR, #1, 10mm thick
2203-1017	All poly tibial component, CR, #1, 12mm thick
2203-1018	All poly tibial component, CR, #1, 14mm thick
2203-1019	All poly tibial component, CR, #1, 16mm thick
2203-1010	All poly tibial component, CR, #1, 17mm thick
2203-1021	All poly tibial component, CR, #2, 9mm thick
2203-1022	All poly tibial component, CR, #2, 11mm thick
2203-1023	All poly tibial component, CR, #2, 13mm thick
2203-1024	All poly tibial component, CR, #2, 15mm thick
2203-1025	All poly tibial component, CR, #2, 18mm thick
2203-1026	All poly tibial component, CR, #2, 10mm thick
2203-1027	All poly tibial component, CR, #2, 12mm thick
2203-1028	All poly tibial component, CR, #2, 14mm thick
2203-1029	All poly tibial component, CR, #2, 16mm thick
2203-1020	All poly tibial component, CR, #2, 17mm thick
2203-1031	All poly tibial component, CR, #3, 9mm thick
2203-1032	All poly tibial component, CR, #3, 11mm thick

Catalog No.	Description
2203-1033	All poly tibial component, CR, #3, 13mm thick
2203-1034	All poly tibial component, CR, #3, 15mm thick
2203-1035	All poly tibial component, CR, #3, 18mm thick
2203-1036	All poly tibial component, CR, #3, 10mm thick
2203-1037	All poly tibial component, CR, #3, 12mm thick
2203-1038	All poly tibial component, CR, #3, 14mm thick
2203-1039	All poly tibial component, CR, #3, 16mm thick
2203-1030	All poly tibial component, CR, #3, 17mm thick
2203-1041	All poly tibial component, CR, #4, 9mm thick
2203-1042	All poly tibial component, CR, #4, 11mm thick
2203-1043	All poly tibial component, CR, #4, 13mm thick
2203-1044	All poly tibial component, CR, #4, 15mm thick
2203-1045	All poly tibial component, CR, #4, 18mm thick
2203-1046	All poly tibial component, CR, #4, 10mm thick
2203-1047	All poly tibial component, CR, #4, 12mm thick
2203-1048	All poly tibial component, CR, #4, 14mm thick
2203-1049	All poly tibial component, CR, #4, 16mm thick
2203-1040	All poly tibial component, CR, #4, 17mm thick
2203-1051	All poly tibial component, CR, #5, 9mm thick
2203-1052	All poly tibial component, CR, #5, 11mm thick
2203-1053	All poly tibial component, CR, #5, 13mm thick
2203-1054	All poly tibial component, CR, #5, 15mm thick
2203-1055	All poly tibial component, CR, #5, 18mm thick
2203-1056	All poly tibial component, CR, #5, 10mm thick
2203-1057	All poly tibial component, CR, #5, 12mm thick
2203-1058	All poly tibial component, CR, #5, 14mm thick
2203-1059	All poly tibial component, CR, #5, 16mm thick
2203-1050	All poly tibial component, CR, #5, 17mm thick
2203-1061	All poly tibial component, CR, #6, 9mm thick
2203-1062	All poly tibial component, CR, #6, 11mm thick
2203-1063	All poly tibial component, CR, #6, 13mm thick
2203-1064	All poly tibial component, CR, #6, 15mm thick
2203-1065	All poly tibial component, CR, #6, 18mm thick
2203-1066	All poly tibial component, CR, #6, 10mm thick
2203-1067	All poly tibial component, CR, #6, 12mm thick
2203-1068	All poly tibial component, CR, #6, 14mm thick
2203-1069	All poly tibial component, CR, #6, 16mm thick
2203-1060	All poly tibial component, CR, #6, 17mm thick
2203-1071	All poly tibial component, CR, #7, 9mm thick

Catalog No.	Description
2203-1072	All poly tibial component, CR, #7, 11mm thick
2203-1073	All poly tibial component, CR, #7, 13mm thick
2203-1074	All poly tibial component, CR, #7, 15mm thick
2203-1075	All poly tibial component, CR, #7, 18mm thick
2203-1076	All poly tibial component, CR, #7, 10mm thick
2203-1077	All poly tibial component, CR, #7, 12mm thick
2203-1078	All poly tibial component, CR, #7, 14mm thick
2203-1079	All poly tibial component, CR, #7, 16mm thick
2203-1070	All poly tibial component, CR, #7, 17mm thick
2203-1201	All poly tibial component, PS, #0, 9mm thick
2203-1202	All poly tibial component, PS, #0, 11mm thick
2203-1203	All poly tibial component, PS, #0, 13mm thick
2203-1204	All poly tibial component, PS, #0, 15mm thick
2203-1205	All poly tibial component, PS, #0, 18mm thick
2203-1211	All poly tibial component, PS, #1, 9mm thick
2203-1212	All poly tibial component, PS, #1, 11mm thick
2203-1213	All poly tibial component, PS, #1, 13mm thick
2203-1214	All poly tibial component, PS, #1, 15mm thick
2203-1215	All poly tibial component, PS, #1, 18mm thick
2203-1216	All poly tibial component, PS, #1, 10mm thick
2203-1217	All poly tibial component, PS, #1, 12mm thick
2203-1218	All poly tibial component, PS, #1, 14mm thick
2203-1219	All poly tibial component, PS, #1, 16mm thick
2203-1210	All poly tibial component, PS, #1, 17mm thick
2203-1221	All poly tibial component, PS, #2, 9mm thick
2203-1222	All poly tibial component, PS, #2, 11mm thick
2203-1223	All poly tibial component, PS, #2, 13mm thick
2203-1224	All poly tibial component, PS, #2, 15mm thick
2203-1225	All poly tibial component, PS, #2, 18mm thick
2203-1226	All poly tibial component, PS, #2, 10mm thick
2203-1227	All poly tibial component, PS, #2, 12mm thick
2203-1228	All poly tibial component, PS, #2, 14mm thick
2203-1229	All poly tibial component, PS, #2, 16mm thick
2203-1220	All poly tibial component, PS, #2, 17mm thick
2203-1231	All poly tibial component, PS, #3, 9mm thick
2203-1232	All poly tibial component, PS, #3, 11mm thick
2203-1233	All poly tibial component, PS, #3, 13mm thick
2203-1234	All poly tibial component, PS, #3, 15mm thick
2203-1235	All poly tibial component, PS, #3, 18mm thick

Catalog No.	Description
2203-1236	All poly tibial component, PS, #3, 10mm thick
2203-1237	All poly tibial component, PS, #3, 12mm thick
2203-1238	All poly tibial component, PS, #3, 14mm thick
2203-1239	All poly tibial component, PS, #3, 16mm thick
2203-1230	All poly tibial component, PS, #3, 17mm thick
2203-1241	All poly tibial component, PS, #4, 9mm thick
2203-1242	All poly tibial component, PS, #4, 11mm thick
2203-1243	All poly tibial component, PS, #4, 13mm thick
2203-1244	All poly tibial component, PS, #4, 15mm thick
2203-1245	All poly tibial component, PS, #4, 18mm thick
2203-1246	All poly tibial component, PS, #4, 10mm thick
2203-1247	All poly tibial component, PS, #4, 12mm thick
2203-1248	All poly tibial component, PS, #4, 14mm thick
2203-1249	All poly tibial component, PS, #4, 16mm thick
2203-1240	All poly tibial component, PS, #4, 17mm thick
2203-1251	All poly tibial component, PS, #5, 9mm thick
2203-1252	All poly tibial component, PS, #5, 11mm thick
2203-1253	All poly tibial component, PS, #5, 13mm thick
2203-1254	All poly tibial component, PS, #5, 15mm thick
2203-1255	All poly tibial component, PS, #5, 18mm thick
2203-1256	All poly tibial component, PS, #5, 10mm thick
2203-1257	All poly tibial component, PS, #5, 12mm thick
2203-1258	All poly tibial component, PS, #5, 14mm thick
2203-1259	All poly tibial component, PS, #5, 16mm thick
2203-1250	All poly tibial component, PS, #5, 17mm thick
2203-1261	All poly tibial component, PS, #6, 9mm thick
2203-1262	All poly tibial component, PS, #6, 11mm thick
2203-1263	All poly tibial component, PS, #6, 13mm thick
2203-1264	All poly tibial component, PS, #6, 15mm thick
2203-1265	All poly tibial component, PS, #6, 18mm thick
2203-1266	All poly tibial component, PS, #6, 10mm thick
2203-1267	All poly tibial component, PS, #6, 12mm thick
2203-1268	All poly tibial component, PS, #6, 14mm thick
2203-1269	All poly tibial component, PS, #6, 16mm thick
2203-1260	All poly tibial component, PS, #6, 17mm thick
2203-1271	All poly tibial component, PS, #7, 9mm thick
2203-1272	All poly tibial component, PS, #7, 11mm thick
2203-1273	All poly tibial component, PS, #7, 13mm thick
2203-1274	All poly tibial component, PS, #7, 15mm thick

Catalog No.	Description
2203-1275	All poly tibial component, PS, #7, 18mm thick
2203-1276	All poly tibial component, PS, #7, 10mm thick
2203-1277	All poly tibial component, PS, #7, 12mm thick
2203-1278	All poly tibial component, PS, #7, 14mm thick
2203-1279	All poly tibial component, PS, #7, 16mm thick
2203-1270	All poly tibial component, PS, #7, 17mm thick
2203-1401	All poly tibial component, UC, #0, 9mm thick
2203-1402	All poly tibial component, UC, #0, 11mm thick
2203-1403	All poly tibial component, UC, #0, 13mm thick
2203-1404	All poly tibial component, UC, #0, 15mm thick
2203-1405	All poly tibial component, UC, #0, 18mm thick
2203-1411	All poly tibial component, UC, #1, 9mm thick
2203-1412	All poly tibial component, UC, #1, 11mm thick
2203-1413	All poly tibial component, UC, #1, 13mm thick
2203-1414	All poly tibial component, UC, #1, 15mm thick
2203-1415	All poly tibial component, UC, #1, 18mm thick
2203-1416	All poly tibial component, UC, #1, 10mm thick
2203-1417	All poly tibial component, UC, #1, 12mm thick
2203-1418	All poly tibial component, UC, #1, 14mm thick
2203-1419	All poly tibial component, UC, #1, 16mm thick
2203-1410	All poly tibial component, UC, #1, 17mm thick
2203-1421	All poly tibial component, UC, #2, 9mm thick
2203-1422	All poly tibial component, UC, #2, 11mm thick
2203-1423	All poly tibial component, UC, #2, 13mm thick
2203-1424	All poly tibial component, UC, #2, 15mm thick
2203-1425	All poly tibial component, UC, #2, 18mm thick
2203-1426	All poly tibial component, UC, #2, 10mm thick
2203-1427	All poly tibial component, UC, #2, 12mm thick
2203-1428	All poly tibial component, UC, #2, 14mm thick
2203-1429	All poly tibial component, UC, #2, 16mm thick
2203-1420	All poly tibial component, UC, #2, 17mm thick
2203-1431	All poly tibial component, UC, #3, 9mm thick
2203-1432	All poly tibial component, UC, #3, 11mm thick
2203-1433	All poly tibial component, UC, #3, 13mm thick
2203-1434	All poly tibial component, UC, #3, 15mm thick
2203-1435	All poly tibial component, UC, #3, 18mm thick
2203-1436	All poly tibial component, UC, #3, 10mm thick
2203-1437	All poly tibial component, UC, #3, 12mm thick
2203-1438	All poly tibial component, UC, #3, 14mm thick

Catalog No.	Description
2203-1439	All poly tibial component, UC, #3, 16mm thick
2203-1430	All poly tibial component, UC, #3, 17mm thick
2203-1441	All poly tibial component, UC, #4, 9mm thick
2203-1442	All poly tibial component, UC, #4, 11mm thick
2203-1443	All poly tibial component, UC, #4, 13mm thick
2203-1444	All poly tibial component, UC, #4, 15mm thick
2203-1445	All poly tibial component, UC, #4, 18mm thick
2203-1446	All poly tibial component, UC, #4, 10mm thick
2203-1447	All poly tibial component, UC, #4, 12mm thick
2203-1448	All poly tibial component, UC, #4, 14mm thick
2203-1449	All poly tibial component, UC, #4, 16mm thick
2203-1440	All poly tibial component, UC, #4, 17mm thick
2203-1451	All poly tibial component, UC, #5, 9mm thick
2203-1452	All poly tibial component, UC, #5, 11mm thick
2203-1453	All poly tibial component, UC, #5, 13mm thick
2203-1454	All poly tibial component, UC, #5, 15mm thick
2203-1455	All poly tibial component, UC, #5, 18mm thick
2203-1456	All poly tibial component, UC, #5, 10mm thick
2203-1457	All poly tibial component, UC, #5, 12mm thick
2203-1458	All poly tibial component, UC, #5, 14mm thick
2203-1459	All poly tibial component, UC, #5, 16mm thick
2203-1450	All poly tibial component, UC, #5, 17mm thick
2203-1461	All poly tibial component, UC, #6, 9mm thick
2203-1462	All poly tibial component, UC, #6, 11mm thick
2203-1463	All poly tibial component, UC, #6, 13mm thick
2203-1464	All poly tibial component, UC, #6, 15mm thick
2203-1465	All poly tibial component, UC, #6, 18mm thick
2203-1466	All poly tibial component, UC, #6, 10mm thick
2203-1467	All poly tibial component, UC, #6, 12mm thick
2203-1468	All poly tibial component, UC, #6, 14mm thick
2203-1469	All poly tibial component, UC, #6, 16mm thick
2203-1460	All poly tibial component, UC, #6, 17mm thick
2203-1471	All poly tibial component, UC, #7, 9mm thick
2203-1472	All poly tibial component, UC, #7, 11mm thick
2203-1473	All poly tibial component, UC, #7, 13mm thick
2203-1474	All poly tibial component, UC, #7, 15mm thick
2203-1475	All poly tibial component, UC, #7, 18mm thick
2203-1476	All poly tibial component, UC, #7, 10mm thick
2203-1477	All poly tibial component, UC, #7, 12mm thick

Catalog No.	Description
2203-1478	All poly tibial component, UC, #7, 14mm thick
2203-1479	All poly tibial component, UC, #7, 16mm thick
2203-1470	All poly tibial component, UC, #7, 17mm thick
2303-1800	E-XPE Tibial insert, CR, #0, 17mm
2303-1801	E-XPE Tibial insert, CR, #0, 9mm
2303-1802	E-XPE Tibial insert, CR, #0, 11mm
2303-1803	E-XPE Tibial insert, CR, #0, 13mm
2303-1804	E-XPE Tibial insert, CR, #0, 15mm
2303-1805	E-XPE Tibial insert, CR, #0, 18mm
2303-1806	E-XPE Tibial insert, CR, #0, 10mm
2303-1807	E-XPE Tibial insert, CR, #0, 12mm
2303-1808	E-XPE Tibial insert, CR, #0, 14mm
2303-1809	E-XPE Tibial insert, CR, #0, 16mm
2303-1810	E-XPE Tibial insert, CR, #1, 17mm
2303-1811	E-XPE Tibial insert, CR, #1, 9mm
2303-1812	E-XPE Tibial insert, CR, #1, 11mm
2303-1813	E-XPE Tibial insert, CR, #1, 13mm
2303-1814	E-XPE Tibial insert, CR, #1, 15mm
2303-1815	E-XPE Tibial insert, CR, #1, 18mm
2303-1816	E-XPE Tibial insert, CR, #1, 10mm
2303-1817	E-XPE Tibial insert, CR, #1, 12mm
2303-1818	E-XPE Tibial insert, CR, #1, 14mm
2303-1819	E-XPE Tibial insert, CR, #1, 16mm
2303-1820	E-XPE Tibial insert, CR, #2, 17mm
2303-1821	E-XPE Tibial insert, CR, #2, 9mm
2303-1822	E-XPE Tibial insert, CR, #2, 11mm
2303-1823	E-XPE Tibial insert, CR, #2, 13mm
2303-1824	E-XPE Tibial insert, CR, #2, 15mm
2303-1825	E-XPE Tibial insert, CR, #2, 18mm
2303-1826	E-XPE Tibial insert, CR, #2, 10mm
2303-1827	E-XPE Tibial insert, CR, #2, 12mm
2303-1828	E-XPE Tibial insert, CR, #2, 14mm
2303-1829	E-XPE Tibial insert, CR, #2, 16mm
2303-1830	E-XPE Tibial insert, CR, #3, 17mm
2303-1831	E-XPE Tibial insert, CR, #3, 9mm
2303-1832	E-XPE Tibial insert, CR, #3, 11mm
2303-1833	E-XPE Tibial insert, CR, #3, 13mm
2303-1834	E-XPE Tibial insert, CR, #3, 15mm
2303-1835	E-XPE Tibial insert, CR, #3, 18mm

Catalog No.	Description
2303-1836	E-XPE Tibial insert, CR, #3, 10mm
2303-1837	E-XPE Tibial insert, CR, #3, 12mm
2303-1838	E-XPE Tibial insert, CR, #3, 14mm
2303-1839	E-XPE Tibial insert, CR, #3, 16mm
2303-1840	E-XPE Tibial insert, CR, #4, 17mm
2303-1841	E-XPE Tibial insert, CR, #4, 9mm
2303-1842	E-XPE Tibial insert, CR, #4, 11mm
2303-1843	E-XPE Tibial insert, CR, #4, 13mm
2303-1844	E-XPE Tibial insert, CR, #4, 15mm
2303-1845	E-XPE Tibial insert, CR, #4, 18mm
2303-1846	E-XPE Tibial insert, CR, #4, 10mm
2303-1847	E-XPE Tibial insert, CR, #4, 12mm
2303-1848	E-XPE Tibial insert, CR, #4, 14mm
2303-1849	E-XPE Tibial insert, CR, #4, 16mm
2303-1850	E-XPE Tibial insert, CR, #5, 17mm
2303-1851	E-XPE Tibial insert, CR, #5, 9mm
2303-1852	E-XPE Tibial insert, CR, #5, 11mm
2303-1853	E-XPE Tibial insert, CR, #5, 13mm
2303-1854	E-XPE Tibial insert, CR, #5, 15mm
2303-1855	E-XPE Tibial insert, CR, #5, 18mm
2303-1856	E-XPE Tibial insert, CR, #5, 10mm
2303-1857	E-XPE Tibial insert, CR, #5, 12mm
2303-1858	E-XPE Tibial insert, CR, #5, 14mm
2303-1859	E-XPE Tibial insert, CR, #5, 16mm
2303-1860	E-XPE Tibial insert, CR, #6, 17mm
2303-1861	E-XPE Tibial insert, CR, #6, 9mm
2303-1862	E-XPE Tibial insert, CR, #6, 11mm
2303-1863	E-XPE Tibial insert, CR, #6, 13mm
2303-1864	E-XPE Tibial insert, CR, #6, 15mm
2303-1865	E-XPE Tibial insert, CR, #6, 18mm
2303-1866	E-XPE Tibial insert, CR, #6, 10mm
2303-1867	E-XPE Tibial insert, CR, #6, 12mm
2303-1868	E-XPE Tibial insert, CR, #6, 14mm
2303-1869	E-XPE Tibial insert, CR, #6, 16mm
2303-1870	E-XPE Tibial insert, CR, #7, 17mm
2303-1871	E-XPE Tibial insert, CR, #7, 9mm
2303-1872	E-XPE Tibial insert, CR, #7, 11mm
2303-1873	E-XPE Tibial insert, CR, #7, 13mm
2303-1874	E-XPE Tibial insert, CR, #7, 15mm

Catalog No.	Description
2303-1875	E-XPE Tibial insert, CR, #7, 18mm
2303-1876	E-XPE Tibial insert, CR, #7, 10mm
2303-1877	E-XPE Tibial insert, CR, #7, 12mm
2303-1878	E-XPE Tibial insert, CR, #7, 14mm
2303-1879	E-XPE Tibial insert, CR, #7, 16mm
2303-3800	E-XPE Tibial insert, PS, #0, 17mm
2303-3801	E-XPE Tibial insert, PS, #0, 9mm
2303-3802	E-XPE Tibial insert, PS, #0, 11mm
2303-3803	E-XPE Tibial insert, PS, #0, 13mm
2303-3804	E-XPE Tibial insert, PS, #0, 15mm
2303-3806	E-XPE Tibial insert, PS, #0, 10mm
2303-3807	E-XPE Tibial insert, PS, #0, 12mm
2303-3808	E-XPE Tibial insert, PS, #0, 14mm
2303-3809	E-XPE Tibial insert, PS, #0, 16mm
2303-3810	E-XPE Tibial insert, PS, #1, 17mm
2303-3811	E-XPE Tibial insert, PS, #1, 9mm
2303-3812	E-XPE Tibial insert, PS, #1, 11mm
2303-3813	E-XPE Tibial insert, PS, #1, 13mm
2303-3814	E-XPE Tibial insert, PS, #1, 15mm
2303-3815	E-XPE Tibial insert, PS, #1, 18mm
2303-3816	E-XPE Tibial insert, PS, #1, 10mm
2303-3817	E-XPE Tibial insert, PS, #1, 12mm
2303-3818	E-XPE Tibial insert, PS, #1, 14mm
2303-3819	E-XPE Tibial insert, PS, #1, 16mm
2303-3820	E-XPE Tibial insert, PS, #2, 17mm
2303-3821	E-XPE Tibial insert, PS, #2, 9mm
2303-3822	E-XPE Tibial insert, PS, #2, 11mm
2303-3823	E-XPE Tibial insert, PS, #2, 13mm
2303-3824	E-XPE Tibial insert, PS, #2, 15mm
2303-3825	E-XPE Tibial insert, PS, #2, 18mm
2303-3826	E-XPE Tibial insert, PS, #2, 10mm
2303-3827	E-XPE Tibial insert, PS, #2, 12mm
2303-3828	E-XPE Tibial insert, PS, #2, 14mm
2303-3829	E-XPE Tibial insert, PS, #2, 16mm
2303-3830	E-XPE Tibial insert, PS, #3, 17mm
2303-3831	E-XPE Tibial insert, PS, #3, 9mm
2303-3832	E-XPE Tibial insert, PS, #3, 11mm
2303-3833	E-XPE Tibial insert, PS, #3, 13mm
2303-3834	E-XPE Tibial insert, PS, #3, 15mm

Catalog No.	Description
2303-3835	E-XPE Tibial insert, PS, #3, 18mm
2303-3836	E-XPE Tibial insert, PS, #3, 10mm
2303-3837	E-XPE Tibial insert, PS, #3, 12mm
2303-3838	E-XPE Tibial insert, PS, #3, 14mm
2303-3839	E-XPE Tibial insert, PS, #3, 16mm
2303-3840	E-XPE Tibial insert, PS, #4, 17mm
2303-3841	E-XPE Tibial insert, PS, #4, 9mm
2303-3842	E-XPE Tibial insert, PS, #4, 11mm
2303-3843	E-XPE Tibial insert, PS, #4, 13mm
2303-3844	E-XPE Tibial insert, PS, #4, 15mm
2303-3845	E-XPE Tibial insert, PS, #4, 18mm
2303-3846	E-XPE Tibial insert, PS, #4, 10mm
2303-3847	E-XPE Tibial insert, PS, #4, 12mm
2303-3848	E-XPE Tibial insert, PS, #4, 14mm
2303-3849	E-XPE Tibial insert, PS, #4, 16mm
2303-3850	E-XPE Tibial insert, PS, #5, 17mm
2303-3851	E-XPE Tibial insert, PS, #5, 9mm
2303-3852	E-XPE Tibial insert, PS, #5, 11mm
2303-3853	E-XPE Tibial insert, PS, #5, 13mm
2303-3854	E-XPE Tibial insert, PS, #5, 15mm
2303-3855	E-XPE Tibial insert, PS, #5, 18mm
2303-3856	E-XPE Tibial insert, PS, #5, 10mm
2303-3857	E-XPE Tibial insert, PS, #5, 12mm
2303-3858	E-XPE Tibial insert, PS, #5, 14mm
2303-3859	E-XPE Tibial insert, PS, #5, 16mm
2303-3860	E-XPE Tibial insert, PS, #6, 17mm
2303-3861	E-XPE Tibial insert, PS, #6, 9mm
2303-3862	E-XPE Tibial insert, PS, #6, 11mm
2303-3863	E-XPE Tibial insert, PS, #6, 13mm
2303-3864	E-XPE Tibial insert, PS, #6, 15mm
2303-3865	E-XPE Tibial insert, PS, #6, 18mm
2303-3866	E-XPE Tibial insert, PS, #6, 10mm
2303-3867	E-XPE Tibial insert, PS, #6, 12mm
2303-3868	E-XPE Tibial insert, PS, #6, 14mm
2303-3869	E-XPE Tibial insert, PS, #6, 16mm
2303-3870	E-XPE Tibial insert, PS, #7, 17mm
2303-3871	E-XPE Tibial insert, PS, #7, 9mm
2303-3872	E-XPE Tibial insert, PS, #7, 11mm
2303-3873	E-XPE Tibial insert, PS, #7, 13mm

Catalog No.	Description
2303-3874	E-XPE Tibial insert, PS, #7, 15mm
2303-3875	E-XPE Tibial insert, PS, #7, 18mm
2303-3876	E-XPE Tibial insert, PS, #7, 10mm
2303-3877	E-XPE Tibial insert, PS, #7, 12mm
2303-3878	E-XPE Tibial insert, PS, #7, 14mm
2303-3879	E-XPE Tibial insert, PS, #7, 16mm
2303-1700	E-XPE Tibial insert, UC, #0, 17mm
2303-1701	E-XPE Tibial insert, UC, #0, 9mm
2303-1702	E-XPE Tibial insert, UC, #0, 11mm
2303-1703	E-XPE Tibial insert, UC, #0, 13mm
2303-1704	E-XPE Tibial insert, UC, #0, 15mm
2303-1705	E-XPE Tibial insert, UC, #0, 18mm
2303-1706	E-XPE Tibial insert, UC, #0, 10mm
2303-1707	E-XPE Tibial insert, UC, #0, 12mm
2303-1708	E-XPE Tibial insert, UC, #0, 14mm
2303-1709	E-XPE Tibial insert, UC, #0, 16mm
2303-1710	E-XPE Tibial insert, UC, #1, 17mm
2303-1711	E-XPE Tibial insert, UC, #1, 9mm
2303-1712	E-XPE Tibial insert, UC, #1, 11mm
2303-1713	E-XPE Tibial insert, UC, #1, 13mm
2303-1714	E-XPE Tibial insert, UC, #1, 15mm
2303-1715	E-XPE Tibial insert, UC, #1, 18mm
2303-1716	E-XPE Tibial insert, UC, #1, 10mm
2303-1717	E-XPE Tibial insert, UC, #1, 12mm
2303-1718	E-XPE Tibial insert, UC, #1, 14mm
2303-1719	E-XPE Tibial insert, UC, #1, 16mm
2303-1720	E-XPE Tibial insert, UC, #2, 17mm
2303-1721	E-XPE Tibial insert, UC, #2, 9mm
2303-1722	E-XPE Tibial insert, UC, #2, 11mm
2303-1723	E-XPE Tibial insert, UC, #2, 13mm
2303-1724	E-XPE Tibial insert, UC, #2, 15mm
2303-1725	E-XPE Tibial insert, UC, #2, 18mm
2303-1726	E-XPE Tibial insert, UC, #2, 10mm
2303-1727	E-XPE Tibial insert, UC, #2, 12mm
2303-1728	E-XPE Tibial insert, UC, #2, 14mm
2303-1729	E-XPE Tibial insert, UC, #2, 16mm
2303-1730	E-XPE Tibial insert, UC, #3, 17mm
2303-1731	E-XPE Tibial insert, UC, #3, 9mm
2303-1732	E-XPE Tibial insert, UC, #3, 11mm

Catalog No.	Description
2303-1733	E-XPE Tibial insert, UC, #3, 13mm
2303-1734	E-XPE Tibial insert, UC, #3, 15mm
2303-1735	E-XPE Tibial insert, UC, #3, 18mm
2303-1736	E-XPE Tibial insert, UC, #3, 10mm
2303-1737	E-XPE Tibial insert, UC, #3, 12mm
2303-1738	E-XPE Tibial insert, UC, #3, 14mm
2303-1739	E-XPE Tibial insert, UC, #3, 16mm
2303-1740	E-XPE Tibial insert, UC, #4, 17mm
2303-1741	E-XPE Tibial insert, UC, #4, 9mm
2303-1742	E-XPE Tibial insert, UC, #4, 11mm
2303-1743	E-XPE Tibial insert, UC, #4, 13mm
2303-1744	E-XPE Tibial insert, UC, #4, 15mm
2303-1745	E-XPE Tibial insert, UC, #4, 18mm
2303-1746	E-XPE Tibial insert, UC, #4, 10mm
2303-1747	E-XPE Tibial insert, UC, #4, 12mm
2303-1748	E-XPE Tibial insert, UC, #4, 14mm
2303-1749	E-XPE Tibial insert, UC, #4, 16mm
2303-1750	E-XPE Tibial insert, UC, #5, 17mm
2303-1751	E-XPE Tibial insert, UC, #5, 9mm
2303-1752	E-XPE Tibial insert, UC, #5, 11mm
2303-1753	E-XPE Tibial insert, UC, #5, 13mm
2303-1754	E-XPE Tibial insert, UC, #5, 15mm
2303-1755	E-XPE Tibial insert, UC, #5, 18mm
2303-1756	E-XPE Tibial insert, UC, #5, 10mm
2303-1757	E-XPE Tibial insert, UC, #5, 12mm
2303-1758	E-XPE Tibial insert, UC, #5, 14mm
2303-1759	E-XPE Tibial insert, UC, #5, 16mm
2303-1760	E-XPE Tibial insert, UC, #6, 17mm
2303-1761	E-XPE Tibial insert, UC, #6, 9mm
2303-1762	E-XPE Tibial insert, UC, #6, 11mm
2303-1763	E-XPE Tibial insert, UC, #6, 13mm
2303-1764	E-XPE Tibial insert, UC, #6, 15mm
2303-1765	E-XPE Tibial insert, UC, #6, 18mm
2303-1766	E-XPE Tibial insert, UC, #6, 10mm
2303-1767	E-XPE Tibial insert, UC, #6, 12mm
2303-1768	E-XPE Tibial insert, UC, #6, 14mm
2303-1769	E-XPE Tibial insert, UC, #6, 16mm
2303-1770	E-XPE Tibial insert, UC, #7, 17mm
2303-1771	E-XPE Tibial insert, UC, #7, 9mm

Catalog No.	Description
2303-1772	E-XPE Tibial insert, UC, #7, 11mm
2303-1773	E-XPE Tibial insert, UC, #7, 13mm
2303-1774	E-XPE Tibial insert, UC, #7, 15mm
2303-1775	E-XPE Tibial insert, UC, #7, 18mm
2303-1776	E-XPE Tibial insert, UC, #7, 10mm
2303-1777	E-XPE Tibial insert, UC, #7, 12mm
2303-1778	E-XPE Tibial insert, UC, #7, 14mm
2303-1779	E-XPE Tibial insert, UC, #7, 16mm
2403-5210	E-XPE Patellar, on set, 3 pegs, X-small, Ø26 mm
2403-5220	E-XPE Patellar, on set, 3 pegs, small, Ø29 mm
2403-5230	E-XPE Patellar, on set, 3 pegs, medium, Ø32 mm
2403-5240	E-XPE Patellar, on set, 3 pegs, large, Ø35 mm
2403-5250	E-XPE Patellar, on set, 3 pegs, X-large, Ø38 mm
2403-5260	E-XPE Patellar, on set, 3 pegs, XX-large, Ø41 mm
2403-5270	E-XPE Patellar, on set, 3 pegs, E-large, Ø44 mm
2403-5010	E-XPE Patella, small, Ø22mm
2403-5020	E-XPE Patella, medium Ø25mm
2403-5030	E-XPE Patella, large Ø28mm
2403-5040	E-XPE Patella, X-large Ø32mm
2203-3400	Tibial baseplate, TPS+, #0
2203-3410	Tibial baseplate, TPS+, #1
2203-3420	Tibial baseplate, TPS+, #2
2203-3430	Tibial baseplate, TPS+, #3
2203-3440	Tibial baseplate, TPS+, #4
2203-3450	Tibial baseplate, TPS+, #5
2203-3460	Tibial baseplate, TPS+, #6
2203-3470	Tibial baseplate, TPS+, #7

-End of List-