

LEADING THE WAY IN GUIDE WIRE TECHNOLOGY

BIFURCATIONS

TRACKING ACROSS
steep angulated takeoffs

TORTUOUS VESSELS

FULL POLYMER COVER
for smooth tracking

CALCIFIED LESIONS

POLYMER COVER AND
HYDROPHILIC COATING
for outstanding lubricity

CTOs

NAVIGATING LESIONS
WITH PILOT 50

PENETRATING LESIONS
WITH PILOT 200

ABBOTT GUIDE WIRES
— SPECIALTY —

HI-TORQUE

PILOT GUIDE WIRES

The right tool for a
“soft” approach



HI-TORQUE PILOT ORDERING INFORMATION

HI-TORQUE PILOT 50	COVER	TIP LOAD	COATING	CORE	TIP STYLE
	FULL POLYMER	1.5g	HYDROPHILIC	DURASTEEL STAINLESS STEEL	CORE-TO-TIP


PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	PENETRATION POWER (kg/in ²)
HI-TORQUE PILOT 50	1010480-H	190	Straight	3	10
HI-TORQUE PILOT 50	1010480-HJ	190	J	3	10
HI-TORQUE PILOT 50	1010483-H	300	Straight	3	10
HI-TORQUE PILOT 50	1010483-HJ	300	J	3	10

HI-TORQUE PILOT 150	COVER	TIP LOAD	COATING	CORE	TIP STYLE
	FULL POLYMER	2.7g	HYDROPHILIC	DURASTEEL STAINLESS STEEL	CORE-TO-TIP

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	PENETRATION POWER (kg/in ²)
HI-TORQUE PILOT 150	1010481-H	190	Straight	3	20
HI-TORQUE PILOT 150	1010481-HJ	190	J	3	20
HI-TORQUE PILOT 150	1010484-H	300	Straight	3	20
HI-TORQUE PILOT 150	1010484-HJ	300	J	3	20

HI-TORQUE PILOT 200	COVER	TIP LOAD	COATING	CORE	TIP STYLE
	FULL POLYMER	4.1g	HYDROPHILIC	DURASTEEL STAINLESS STEEL	CORE-TO-TIP

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	PENETRATION POWER (kg/in ²)
HI-TORQUE PILOT 200	1010482-H	190	Straight	3	30
HI-TORQUE PILOT 200	1010482-HJ	190	J	3	30
HI-TORQUE PILOT 200	1010485-H	300	Straight	3	30
HI-TORQUE PILOT 200	1010485-HJ	300	J	3	30



Caution: These products are intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available), at eifu.abbottvascular.com or at Manuals.sjm.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. **Information contained herein for distribution outside the U.S. only.** Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Tests performed by and data on file at Abbott. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

HI-TORQUE PILOT, HI-TORQUE PILOT 50, HI-TORQUE PILOT 150, and HI-TORQUE PILOT 200 are trademarks of the Abbott Group of Companies.

www.Vascular.Abbott

©2018 Abbott. All rights reserved. AP2945211-OUS Rev. A



Document# Rev
AP2945211-OUS A

AP2945211- OUS A

AP2945211-OUS
Guide Wire Brochure: Pilot
A
Effective Date: 3/2/2018

Document #: AP2945211-OUS
Current Revision: A
Doc Owner: Howell, Laura
Last Review: 3/2/2018
Classification:
Superceded By Doc. # (Supercede only):
Supporting Documentation Attached? Yes ☐ No ☒
Site(s) where doc used: CA
N/A ☒
Material Disposition: On Order: Rec'g Insp: Stock: FU: FS:
Memo:
Material Division:
Associated Products (product name(s) only): *** NONE *** (N)

Title: Guide Wire Brochure: Pilot
Future Revision:
Change Type: New Release
CO Number: 419146
Next Review: 3/1/2021
Retention:

UOM: Alt. UOM: Qty: 0

Reference Documents

Document	Location
----------	----------

Approval Signatures

Level	Role / Function	Name	Site	Signature	Date	Comments
1	Doc Control/Labeling	Howell, Laura	CA	Howell, Laura	2/20/2018 12:08:27 PM	
2	Clinical Engineer	Kaveney, Kathy A	CA	Kaveney, Kathy A	2/20/2018 1:48:30 PM	
2	Customer Service	Martinez, Catherine L	CA	Martinez, Catherine L	2/28/2018 8:01:52 PM	
2	Doc Owner (AP2945211-OUS)	Howell, Laura	CA	Howell, Laura	2/26/2018 8:30:00 AM	
2	Legal IP	Devore, Ronald D	CA	Devore, Ronald D	3/1/2018 11:01:14 PM	
2	Marketing Com	Navarro, Tricia M	CA	Navarro, Tricia M	2/20/2018 12:28:26 PM	
2	Originator	Howell, Laura	CA	Howell, Laura	2/26/2018 8:31:28 AM	
2	R&D	Durcan, Jonathan P	CA	Durcan, Jonathan P	3/1/2018 7:12:30 PM	
2	RA (A & P)	McFerran, Aurora B	CA	McFerran, Aurora B	2/22/2018 11:53:07 AM	

Revision History

Revision	Date	Description of Change	CO Number	Originator
A	3/2/2018 8:11:55 PM	New Document Purpose : to provide information regarding Guide Wires Geography: OUS Intended Audience: Physicians and Healthcare professionals Method of distribution: Electronic/presentation only & leave behind	419146	Howell, Laura