

PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**



HEALTHCARE PROFESSIONALS
Sentrant Introducer Sheath
 EVAR and TEVAR Procedures

REQUEST MORE INFORMATION

[CONTACT US](#)



AORTIC CATALOGUE

Choose from a complete portfolio of market-leading aortic vascular therapies to treat aneurysm disease.

[DOWNLOAD CATALOGUE](#)

RELATED PRODUCTS

- Aortic Stent Graft Products
- Cardiovascular Products

OVERVIEW

Achieve improved seal in both straightforward and complex anatomies with Sentrant™ introducer sheath.*

INDICATIONS

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic Europe website.

The Medtronic Sentrant™ Introducer Sheaths with Hydrophilic Coating are intended to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature and to minimise blood loss associated with such insertions.

[Indications, Safety, and Warnings](#)



SENTRANT INTRODUCER SHEATH FEATURES

The Sentrant introducer sheath is the choice for superior hemostasis.*

- Smooth dilator to sheath transition
- Radiopaque marker band
- Reinforced coiled tubing
- Two working lengths: 28 cm and 64 cm
- Hydrophilic coating
- Locking mechanism on dilator handle
- Guidewire diameter: 0.035" / 0.89 mm


Compatible with Endurant™ and Valiant Navion™ stent graft delivery systems

ENGINEERED TO DELIVER PROCEDURAL CONFIDENCE

- EnsureSeal Technology delivers superior leak resistance versus competitors*
- Coil-reinforced tubing for added stability and kink resistance
- Maintains lubricity after multiple insertions
- Radiopaque dilator shaft and sheath tip for accurate visualization and guidance
- 64 cm configuration available for broader anatomical applicability
- Compatible with a wide range of endovascular portfolios

SUPERIOR LEAK RESISTANCE VERSUS COOK CHECK-FLO PERFORMER™ INTRODUCER SHEATH AND GORE DRYSEAL™ FLEX INTRODUCER SHEATH*

Comparison of leak resistance for the Medtronic Sentrant introducer sheath, Cook Check Flo Performer introducer, and Gore DrySeal Flex sheath.* The Sentrant introducer sheath had zero leakage on 0.035" / 0.89 mm guidewire.



MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the Medtronic Manual Library or in the product labeling supplied with each device.

MODEL SPECIFICATIONS

ORDERING INFORMATION

Model Number	Sheath Size (Fr)	Working Length (cm)
SENSH1228W	12	28
SENSH1428W	14	28
SENSH1628W	16	28
SENSH1828W	18	28
SENSH2028W	20	28
SENSH2228W	22	28
SENSH2428W	24	28
SENSH2628W	26	28
SENSH1264W	12	64
SENSH1464W	14	64
SENSH1664W	16	64
SENSH1864W	18	64
SENSH2064W	20	64
SENSH2264W	22	64
SENSH2464W	24	64
SENSH2664W	26	64

* Leak Resistance Bench Test Data on file at Medtronic. Test data not indicative of clinical performance. Bench Test compared Cook Check-Flo Performer™ Introducer Sheath 16 Fr and Gore DrySeal™ Flex 16 Fr to Sentrant™ 16 Fr. Graph shows average leakage over an extrapolated 2.5 hour procedure.

*** Third party brands are trademarks of their respective owners.

Disclaimer: 64 cm working length sheath may not be available in your region or country. Please consult the approved indications for use. Content on specific Medtronic products is not intended for users in markets that do not have authorisation for use.

HEALTHCARE PROFESSIONALS

- Overview
- Products
- Therapies & Procedures
- Integrated Health Solutions

PATIENTS

Overview

OUR COMPANY

- Who We Are
- Key Facts
- Leadership
- History
- Mission
- Locations
- Governance
- Careers

OUR IMPACT

- Our Impact
- Healthcare access
- Inclusion, Diversity & Equity
- Protecting our planet
- ESG Report