Detro® San AF Wipes Alcohol Free

USAGE AREAS:

- It is used for rapid disinfection of incubators, dental aerators and medical devices used in ENT examinations.
- It is compatible with medical devices/equipments made of polystirene, stainless steel, pvc, acrylic material, polycarbon, ceramic, ABS, silicone, HDPE, rubber, and glass

PRODUCT PROPERTIES:

- Ready to use product.
- · Aldehyde, phenol and alcohol free.
- Wide spectrum of microbiological activity.
- Detro® San AF provides less product consumption thanks to the spray head which can produce foam.

93/42 EEC Medical device directive class IIa.

COMPOSITION:

• Bis (3-aminopropyl) dodecylamine, didecyldimethylammonium chloride, benzalconium chloride, auxiliary substances.

MICROBIOLOGICAL ACTIVITY:

• Bactericidal, fungicidal, virucidal, tuberculocidal in 1 minute.















• Sprayed onto the medical device to cover al the surface. The contact time is completed and the surface is wiped.

• A sufficient amount of wipes is taken and applied directly to the surface to be wiped. Wait for 1 minute for effect.

QUANTITY	PACKING TYPE	BARCODE
750 ML	9 x BOTTLE	8680097102461
3	9 x BOTTLE	8680097100757
5L	4 x DRUM	8680097101167

QUANTITY PACKING TYPE BARCODE

150
pcs.

CYLINDER BOX

8680097100764

PACKAGE: