

## orthoBAND NORMAL plaster of Paris bandage

TD-04-I.2.c-1.2

edition I

date: 2021-04-22

amendment of: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

**CLASSIFICATION** 

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

## **INTENDED PURPOSE**

 disposable device intended for use as a full or partial dressing for stiffening/immobilisation in orthopaedic, ambulatory and surgical procedures

## **FEATURES**

- gypsum substance carrier made of 100 % cotton gauze
- min. 94% white natural gypsum, carrier sprayed twice with gypsum mass - carrier material covering both sides
- carrier is wound on a "X" shaped plastic core
- prior to application requires soaking in a container of water at the recommended temperature of 20°C for approximately 5 seconds
- initial setting time: 5-6 minutes
- complete stiffening of the plaster of Paris bandage takes place within 24 hours
- available with a straight edge (6 cm, 8 cm) or serrated edge (10 cm, 12 cm, 14 cm, 15 cm, 20 cm)
- does not absorb X-rays
- does not cause allergic reaction
- white colour
- single-use
- non-sterile
- shelf-life: 5 years



CATALOGUE NUMBER	DIMENSIONS	UNIT PACKAGING (SALES PACKAGING) FOIL BAG	BULK PACKAGING CARTON
OB06300-N	6 cm x 3 m	1 pc.	48 x 1 pc.
OB08300-N	8 cm x 3 m	1 pc.	48 x 1 pc.
OB10300-N	10 cm x 3 m	1 pc.	48 x 1 pc.
OB12300-N	12 cm x 3 m	1 pc.	48 x 1 pc.
OB14300-N	14 cm x 3 m	1 pc.	36 x 1 pc.
OB15300-N	15 cm x 3 m	1 pc.	36 x 1 pc.
OB20300-N	20 cm x 3 m	1 pc.	24 x 1 pc.