



Vascular and Cardiothoracic Surgery Solutions

Over 50 years of technical and clinical
experience in vascular graft technology

This document is intended to provide information to an
international audience.

Hemacarotid Patch Knitted

Knitted Patches



Material*:

- Collagen-coated polyester patch
- Reverse lock-knit knitting technique
- Magnetic Resonance Safe
- Water permeability : < 5ml/cm²/min

Hemacarotid Patch Knitted

Dimensions	Reference
6 mm x 75 mm	HEK06/75CP (1)
8 mm x 75 mm	HEK08/75CP (1)
10 mm x 75 mm	HEK10/75CP (1)
12 mm x 75 mm	HEK12/75CP (1)
14 mm x 75 mm	HEK14/75CP (1)