

工鼻/HME



BY-1

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
9.2(@Vt500mL)
- 阻力
Resistance(pa)
140 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
/
- 病毒过滤率
Viral filtration efficiency*(%)
/



BY-1

- 推荐潮气量
Recommended tidal volume(ml)
75-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
11.2(@Vt500mL)
- 阻力
Resistance(pa)
30 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
/
- 病毒过滤率
Viral filtration efficiency*(%)
/



BY-RGB-P-3

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
8.8(@Vt500mL)
- 阻力
Resistance(pa)
65 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
/
- 病毒过滤率
Viral filtration efficiency*(%)
/

合式人工鼻/HME Filter



BY-GL-2

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
10.0(@Vt500mL)
- 阻力
Resistance(pa)
105 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
99.99
- 病毒过滤率
Viral filtration efficiency*(%)
99.99



BY-GL-3

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
10.9(@Vt500mL)
- 阻力
Resistance(pa)
105 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
99.99
- 病毒过滤率
Viral filtration efficiency*(%)
99.99



BY-GL-4

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
9.9(@Vt500mL)
- 阻力
Resistance(pa)
110 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
99.99
- 病毒过滤率
Viral filtration efficiency*(%)
99.99



BY-GL-5

- 推荐潮气量
Recommended tidal volume(ml)
150-1500
- 保湿效果
Moisture loss*(mg H₂O/L air)
9.2(@Vt500mL)
- 阻力
Resistance(pa)
100 at 30L/min
- 细菌过滤率
Bacterial filtration efficiency*(%)
99.99
- 病毒过滤率
Viral filtration efficiency*(%)
99.99