



KANSUK ACD-A ANTİKOAGÜLAN SOLÜSYONU

KANSUK ACD-A ANTICOAGULANT SOLUTION



Endikasyonu;

20000006525 ACD-A: ACD-A solüsyonu aferez prosedüründe kullanılan 500 ml antikoagulan solüsyonudur ve vücut dışında toplanan kanın koagülasyonunu önler.

Indication;

20000006525 ACD-A : *The ACD-A solution is an 500 ml anticoagulant solution used in apheresis procedures and it prevents coagulation of the collected blood outside the body.*

Dikkat: Aferez cihazları ile birlikte kullanım içindir. Bozulmaya ilişkin gözle görünür bir kötüleşme varsa kullanmayınız. Kullanmadan önce Kullanma Kılavuzunu okuyunuz. Direkt damardan vermeyiniz. Delmeyiniz. Temiz, serin ve kuru yerde saklayınız. Işıktan sakınınız. APİROJEN, NONTOKSİK

Caution: *This is for use with apheresis devices. Do not use if there is a visible worsening related to deterioration. Read user instructions before use. Do not transfuse directly. Do not pierce. Keep in clean, cool and dry place. Avoid from the light. NON-PYROGENIC, NONTOKSİK*

AYRINTILI BİLGİ İÇİN FİRMAMIZI ARAYINIZ / PLEASE CONTACT US FOR FURTHER INFORMATION

CE1984

Kansuk Laboratuvarı Sanayi ve Ticaret A.Ş.
Yassören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Türkiye/Turkey
Tel / Phone: +90 212 5921576 Faks / Fax: +90 212 5803772
Web: www.kansuk.com E-mail: kansuk@kansuk.com

