


AUSTROMAT i


DEKEMA
DER OFEN.



AUSTROMAT™/ HIGHLIGHTS

- Capacitive multi-touch Panel
- 10,4" high resolution Display
- Revolutionary, easy to use Pad Operation
- Soft touch Lift Movement
- Internet Database Connection

- Remote Control via Network
- Integrated Tele Service
- Quality Management
- Software Update & Backup via Internet
- Aluminum Housing in Black or Silver




AUSTROMAT™ 624i™

- Firing Furnace up to 1200°C
- Fire & Glaze with and without Vacuum
- Crystallize & Infiltrate with and without Vacuum



AUSTROMAT™ 664i™

- High-Temperature Furnace up to 1560°C
- Sinter, Glaze & Crystallize under Atmosphere
- Chairside Speed Mode
- Space for up to six unit bridges or up to 20 units in one level (Ø 65mm)
- Optional Accessories to sinter in two Levels



AUSTROMAT™ 654i™

- Combi Furnace up to 1200°C
- Press, Fire & Glaze with and without Vacuum
- Crystallize & Infiltrate with and without Vacuum
- Press several Plungers/Colors in one go with the optional trixpress™



AUSTROMAT™ 674i™

- High-Temperature Furnace up to 1560°C
- Sinter, Glaze & Crystallize under Atmosphere
- Chairside Speed Mode
- Space for all common bridge sizes or up to 40 units in one level (Ø 100mm)
- Optional Accessories to sinter in two Levels

 368 °C	# 00000001249		17/03/02	01:24:12
	Chairside Program			


Idle

Handle

Fire

Preheat 

980 °C

Load 

Start 

1500 °C

Stop 

Object temperature: 368 °C

Duration: 00:38

Target time: 02:04:21



New



Copy



Cut



Paste



Delete



Setup



ESC



innovation

meets tradition and vision



inspiration

meets art and technology

