

Dragonfly™ OPTIS™

Imaging Catheter



PRODUCT FEATURES


- OCT imaging catheter for OPTIS™ Next Imaging Systems, OPTIS™ Imaging Systems, ILUMIEN™ OPTIS™ and ILUMIEN™ PCI Optimization Systems
- Proprietary sheath design provides continuous calibration throughout pullback
- Rapid exchange tip for easier guidewire loading
- Radiopaque markers at the distal tip, imaging lens and 50 mm proximal to the lens aid the physician in aligning the OCT pullback to region of interest
- Shaft markers at 90 cm and 100 cm notify the physician when the distal tip is exiting the guide catheter
- 2.7 F low-profile tip
- Hydrophilic coating
- Compatible with standard 0.014" guide wires
- Compatible with D.2 and E.2 software or higher

SPECIFICATIONS

Usable length: 135 cm

Outer diameter: 2.7 F (distal)

Wire lumen: 0.014"



Tests performed by and data on file at Abbott.

Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

©2021 Abbott. All rights reserved. MAT-2104452 v1.0



ORDERING INFORMATION

Reorder Number	Description
C408645	Dragonfly™ OPTIS™ Imaging Catheter
C408641	Dragonfly™ OPTIS™ Imaging Catheter and Sterile DOC Cover
C408646	Dragonfly™ OPTIS™ Imaging Catheter, Sterile DOC Cover and 1 Sterile 3 ml Syringe
200-700-00	Sterile DOC Cover
C408655	Sterile 3 ml Syringe

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott. Data on file at Abbott.

Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2104452 v1.0

