

# Angiographic Syringes

www.scw-medical.com



Single



DSA



Dual

| Product No. |             | Volume(ml) | Description                          | Models        |               |
|-------------|-------------|------------|--------------------------------------|---------------|---------------|
| Single      | Dual        |            |                                      | Single        | Dual          |
| 42.16.20000 | 42.16.20001 | 200        | Used for CT: CT 9000 & CT 9000 ADV*  | CT-200-B1     | CT-200/200-B1 |
| 42.16.20002 |             | 150        | Used for DSA: Angiomat 6000          | DSA-150-B1    |               |
| 42.16.20003 |             | 150        | Used for CT & DSA: Angiomat Illumena | CT/DSA-150-B1 |               |
| 42.16.20004 |             | 100        | Used for DSA: Angiomat 3000          | DSA-100-B1    |               |
| 42.16.20005 | 42.16.20006 | 60         | Used for MRI: Optistar               | MRI-60-B1     | MRI-60/60-B1  |

| Product No. |             | Volume(ml) | Description          | Models    |               |
|-------------|-------------|------------|----------------------|-----------|---------------|
| Single      | Dual        |            |                      | Single    | Dual          |
| 42.16.50000 | 42.16.50001 | 200        | Used for CT: Empower | CT-200-EM | CT-200/200-EM |