



## Glass staining dishes

### General features

Glass staining dishes to use for special stains.

They resist to each stain solution and solvent, but they have not to be used with microwave oven.



Distributor: Bio-Optica Milano S.p.A.

Code	Description	Capacity	Size (mm)	Liquid volume
10-1410	Hellendhal dish	8 slides	50 x 30 x 90	70 ml
10-1610	Schefferdecker dish	8 slides	55 x 50 x 60	80 ml
10-1710	Coplin dish	6 slides	30 x 30 x 70	50 ml
10-2000	Complete staining set	10 slides	81 x 101 x 70	90 ml
10-1820	Spare dish and rack for 10-2000 set	10 slides	Dish 81 x 101 x 70 Rack 68 x 87 x 50	90 ml
10-1910	Spare rack for 10-2000 set	10 slides	68 x 87 x 50	-
10-2010	Spare tweezer for 10-2000 set	-	-	-

### Technical features

For this kind of product, special precautions for handling and storage are not necessary.

Validity date: not applicable.

The products are realized for generic lab uses, ideal for the stain of microscope slides.

The product don't request CE mark.