

## 7 Specifications

Item		Description
DICOM communication	LI output	DICOM Print support
	Network transmission (optional)	DICOM Storage support Modality: RF/XA/CR/DX
	Media output	CD-R, DVD-R (DICOM format)
	MWM (optional)	Receives examination information from server.
	MPPS (optional)	Transmits examination information to server.
	RDSR (option)	Transmits radiation dose information in the report format to server.
	Bar code reader (optional)	Inputs patient information from bar code reader and performs search with MWM (option).
Maintenance	Remote maintenance	Provides remote system maintenance work (communication link equipment is required separately).
Dimension	Control cabinet	W700×H1900×D500 mm
	Operation cabinet	W160×H575×D470 mm
Mass	Control cabinet	320 kg
	Operation cabinet	35 kg
Rated power supply voltage		Single phase 200/220/230/240 V, 50/60 Hz, 7.0 kVA