

Coronary Peripheral

EIFU

Microcatheters



PERIPHERAL



MICROCATHETERS



NEUROVASCULAR



ANGIOGRAPHIC



ALL PRODUCTS

Neurovascular

Angiographic

CHRONIC TOTAL OCCLUSION GUIDE WIRES

Full Product Listing

HI Gaia® Series

Products

choice for the need of high penetration while maintaining maneuverability. Chronic Total Occlusion indicated. Guide wire with ACT ONE low tip load to improve deflection and torqueability.

Technologies

Training & Education





1.7gf

ASAHI Gaia First

3.5gf

ASAHI Gaia Second

4.5gf

ASAHI Gaia Third

TIP RADIOPACITY


15 cm

COATING

HYDROPHILIC* WITH UNCOATED DISTAL TIP

Hydrophilic 40 cm – uncoated distal tip

TIP SHAPE



Mini Pre-shape




ASAHI Gaia First TUBE CLIP



ASAHI Gaia Second TUBE CLIP



ASAHI Gaia Third TUBE CLIP



ACT ONE >



Ordering Information

To order, please call us toll-free at 855-286-9473, or email us at customersupport@asahi-intecc.com. To see a complete listing of all available ASAHI products, go to our [Full Product Listing](#) page.

PRODUCT	CATALOG NO.	TIP SHAPE	LENGTH (CM)	DIAMETER (INCH(MM))	RADIOPAQUE SEGMENT (CM)	SPRING COIL LENGTH (CM)	COATING
ASAHI Gaia First	AHW14R007P	Pre-Shape	190	0.010/0.014 (0.26/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip
ASAHI Gaia First, 300 cm	AHW14R307P	Pre-Shape	300	0.010/0.014 (0.26/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip
ASAHI Gaia Second	AHW14R008P	Pre-Shape	190	0.011/0.014 (0.28/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip
ASAHI Gaia Second, 300 cm	AHW14R308P	Pre-Shape	300	0.011/0.014 (0.28/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip

PRODUCT	CATALOG NO.	TIP SHAPE	LENGTH (CM)	DIAMETER (INCH(MM))	RADIOPAQUE SEGMENT (CM)	SPRING COIL LENGTH (CM)	COATING
ASAHI Gaia Third	AHW14R011P	Pre- Shape	190	0.012/0.014 (0.30/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip
ASAHI Gaia Third, 300 cm	AHW14R311P	Pre- Shape	300	0.012/0.014 (0.30/0.36)	15	15	Hydrophilic* with Uncoated Distal Tip

*Coated with SLIP-COAT® coating

Indications for Use

ASAHI PTCA Guide Wires are intended to facilitate the placement of balloon dilation catheters during percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA). The ASAHI PTCA Guide Wires are not to be used in the neurovasculature.

[< Back to Coronary](#)

[BACK TO TOP ^](#)



[PRODUCT LISTING >](#)

ASAHI FULL PRODUCT LISTING

Select from our comprehensive portfolio of purpose-driven tools to master all your interventional cases.

FOUR CORE TECHNOLOGIES
FULL PRODUCT LISTING
CUSTOMER SERVICE
QUALITY POLICY

DOWNLOAD APP



CONTACT ASAHI INTECC

Toll-free: 855-286-9473
customersupport@asahi-intecc.com

[CONTACT FORM >](#)

Copyright © 2024 ASAHI INTECC USA, INC. All Rights Reserved | [PRIVACY POLICY](#) | [COOKIE POLICY](#)
Please read the terms of use for further details. Unless otherwise specified, all product names appearing in this website are trademarks or registered trademarks owned by or licensed to ASAHI. All other registered trademarks or trademarks are the property of their respective owners. No use of any ASAHI trademark, trade name, or trade dress in this website may be made without the prior written authorization of ASAHI, except to identify the product or services of the company.