


Vitreoretinal Cannula

For Posterior Segment (retinal) Surgical Procedures

OC 251 BACK FLUSH FLUTE HANDLE
WITH SILICONE TIP CANNULA OC 138

OC 251-20G OC 251-23G OC 251-25G



OC 142 BACK FLUSH FLUTE HANDLE WITH ACTIVE
& PASSIVE ASPIRATION




OC 257 BACK FLUSH FLUTE HANDLE WITH ACTIVE
ASPIRATION & PASSIVE ASPIRATION
WITH SILICONE TIP CANNULA
OC 138



OC 257-20G OC 257-23G OC 257-25G




OC 176 PFCL I/A CANNULA
(PERFLUOROCARBON CO-AXIAL I/A CANNULA)
23G x 37mm INSIDE 30G OPEN END IRRIGATION PORT
& OUTER ASPIRATION PORT (0.4MM)
OC 176-23G OC 176-25G (INSIDE 32G)




OC 177 PFCL I/A CANNULA
(PERFLUOROCARBON I/A CANNULA)
25G x 37mm OPEN END IRRIGATION PORT
WITH 25G ASPIRATION PORT
OC 177-25G



OC 178 VITREORETINAL TAPERED NEEDLE
20G Tapers to 30G Straight tip
4mm / 5mm Straight
OC 178-20G



OC 179 VITREORETINAL TAPERED NEEDLE
20G Taper to 30G Straight tip
4mm / 5mm Angled 45°
OC 179-20G



OC 201 VITREORETINAL TAPERED NEEDLE
20G Tapers to 30G Curved tip
4mm / 5mm Straight
OC 201-20G

