



Technical Data

****L.mono Enrichment Supplement I**

FD214

For selective differentiation of *Listeria monocytogenes* from other *Listeria* species, as per ISO Committee.

Composition

Per vial sufficient for 500 ml medium

Ingredients

L – phosphatidylinositol

Distilled water

Concentration

1g

25ml

Directions:

Thaw the contents of 1 vial of L. mono Enrichment Supplement I at room temperature. Aseptically add the sterile contents to 460 ml of sterile, molten, cooled (45-50°C) L. mono Differential Agar Base [M1540](#) / L. mono Differential HiVeg™ Agar Base [MV1540](#) along with sterile rehydrated contents of 1 vial each of L. mono Selective Supplement I [FD212](#) and L. mono Selective Supplement II [FD213](#) . Mix well and pour into sterile petri plates.

Storage and Shelf Life

On receipt product should be stored at -20°C. Use before the expiry date on the label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.