


**BALTON**<sup>®</sup>

[www.balton.pl](http://www.balton.pl)



**Anaesthesiology**

## EPIDURAL ANAESTHESIA SET SMALL

### Set components:

- Tuohy needle
- epidural catheter with hub
- flat epidural filter 0.2  $\mu\text{m}$  with self-adhesive element
- sleeve



### Optional component:



system  
for catheter fixation



10 ml low resistance syringe

The epidural set allows you to administer anesthetics continuously or repeatedly.

Reference:	Catheter size [G]	Tuohy needle size [G]	Tuohy needle length [mm]
ZZOM19G	22 (pediatric)	19	50
<b>ZZOM18G</b>	20	18	80
<b>ZZOM17G</b>	19	17	80
ZZOM16G	18	16	80

When ordering use Ref. Number, e.g.: ZZOM18G.

NOTE: The set can be purchased together with a low-resistance syringe.

While ordering, add the letters "SN" to the reference number – for example: ZZOM18GSN.

Special order can contain epidural catheter with SOFT tip (comes only with ZZOM18G).

While ordering, add the letter "S" to the reference number – for example: ZZOM18GS.

Special order can contain system for catheter fixation.

While ordering, add the letters "EM" to the reference number – for example: ZZOM18GEM.

**Caution:**

The content of this document, including provided data, information, trademarks and logos are the sole property of BALTON. Any reproduction, whether in whole or in part, is strictly prohibited and will be considered a violation of BALTON's copyrights and other intellectual proprietary rights. ©2022 Balton Sp. z o.o. All rights reserved. This document and any images contained therein is intended solely for the use by healthcare professionals and BALTON's distributors.

This publication is for information purposes only and should not be construed as a commercial offer in the meaning of art. 66(1) of the Polish Civil Code. All medical devices made commercially available by BALTON must be used in accordance with the enclosed instructions for use. This publication presents only a small sample of our complete offering of medical devices. For further information please contact your Balton representative or email us at [balton@balton.pl](mailto:balton@balton.pl).

