



VACUSERA SERUM CLOT ACTIVATOR TUBES

Vacusera Serum Clot Activator Tubes are used for clinical biochemistry and immunology. The inner wall of the Serum Clot Activator Tube is specially coated with clot activator, which activate coagulation process.

Reference Number	Description	Draw vol. (ml)	Size (mm)	Additive	Shelf life	Label Type	Pack/Carton	Cap Colour
234201	VACUSERA 1 ml CAT Serum Clot Activator	1	13 x 75	Clot Activator	18 months	Paper	100/1200	●
234202	VACUSERA 2 ml CAT Serum Clot Activator	2	13 x 75	Clot Activator	18 months	Paper	100/1200	●
234203	VACUSERA 3 ml CAT Serum Clot Activator	3	13 x 75	Clot Activator	18 months	Paper	100/1200	●
234204	VACUSERA 4 ml CAT Serum Clot Activator	4	13 x 75	Clot Activator	18 months	Paper	100/1200	●
235206	VACUSERA 6 ml CAT Serum Clot Activator	6	13 x 100	Clot Activator	18 months	Paper	100/1200	●
236209	VACUSERA 9 ml CAT Serum Clot Activator	9	16 x 100	Clot Activator	18 months	Paper	100/1200	●
236200	VACUSERA 10 ml CAT Serum Clot Activator	10	16 x 100	Clot Activator	18 months	Paper	100/1200	●
434201	VACUSERA 1 ml CAT Serum Clot Activator	1	13 x 75	Clot Activator	18 months	Transparent	100/1200	●
434202	VACUSERA 2 ml CAT Serum Clot Activator	2	13 x 75	Clot Activator	18 months	Transparent	100/1200	●
434203	VACUSERA 3 ml CAT Serum Clot Activator	3	13 x 75	Clot Activator	18 months	Transparent	100/1200	●
434204	VACUSERA 4 ml CAT Serum Clot Activator	4	13 x 75	Clot Activator	18 months	Transparent	100/1200	●
435206	VACUSERA 6 ml CAT Serum Clot Activator	6	13 x 100	Clot Activator	18 months	Transparent	100/1200	●
436209	VACUSERA 9 ml CAT Serum Clot Activator	9	16 x 100	Clot Activator	18 months	Transparent	100/1200	●
436200	VACUSERA 10 ml CAT Serum Clot Activator	10	16 x 100	Clot Activator	18 months	Transparent	100/1200	●