

Mirror and endoscope heater Intra-red



An easy to clean (disinfection wipe) heater for warming laryngeal mirrors and endoscopes up to body temperature to avoid fogging.

Via pushing a button, the infra-red lamp will be switched on and last for approx. 10 sec.. During this time, mirrors and endoscopes can be warmed to body temperature without touching any surface. For pre- and re-heating purposes.

The heater will switch off automatically – so you have all the time for your patient.

Technical data:

- Ø 86 x 140 mm
- 35 W
- temperature of focus point approximately 300°C
- Overheating protection