


MUSTANG™ 0.035" Balloon Dilatation Catheter

EXCELLENCE ON EVERY DIMENSION



Non-Compliant



24 ATM



7x200 thru 5 F

Powerful dilatation combined with superior deliverability and a complete size matrix.

Deliverable High Pressure

- + 24 atmosphere rated burst with a low profile
- + NyBax™ balloon material enables great trackability and high pressure in the same balloon


Broadest Size Matrix

- + Balloon lengths 20–200 mm
- + Diameters 3–12 mm
- + Shafts 40 cm, 75 cm, 135 cm

Low Profile

- + Lesion entry only 15% larger than the wire
- + 7x200 through 5 F
- + 10x80 through 6 F






Only one balloon provides total performance.

With superior cross and track, powerful dilatation, longer lengths, and smaller sheath sizes, Mustang™ is the 0.035" workhorse balloon that offers excellence across the board.

Mustang™ Balloon Dilatation Catheter



Data on file PDM 90721880.

MUSTANG™ BALLOON DILATATION CATHETER

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The Mustang Balloon Dilatation Catheter is indicated for Percutaneous Transluminal Angioplasty (PTA) in the peripheral vasculature, including iliac, femoral, popliteal, tibial, peroneal, subclavian, and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The Mustang Balloon Dilatation Catheter is also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature. Mustang Balloon Dilatation Catheters with balloons up to 120 mm in length are indicated for the treatment of biliary strictures.

CONTRAINDICATIONS: None Known.

WARNINGS: The Mustang PTA balloon catheter is not for use in the coronary arteries. Any use for procedures other than those indicated in these instructions is not recommended.

PRECAUTIONS: The Mustang PTA Balloon Dilatation Catheter shall only be used by physicians trained in the performance of percutaneous transluminal angioplasty. The Mustang PTA Balloon Dilatation Catheter should be used with caution for procedures involving calcified lesions or synthetic vascular grafts due to the abrasive nature of these inflation sites. The Mustang PTA Balloon Dilatation Catheters are not intended for injection of contrast medium. Precautions to prevent or reduce clotting should be taken when any catheter is used: • Consider systemic anticoagulation • Flush or rinse all products with sterile saline or a similar solution prior to use.

ADVERSE EVENTS: The complications that may result from a balloon dilatation procedure include, but are not limited to: • Allergic reaction (device, contrast medium and medications) • Arteriovenous fistula • Embolization (air, device, plaque, etc.) • Hematoma • Hemorrhage, including bleeding at puncture site • Pseudoaneurysm • Sepsis/infection • Thromboembolic episodes • Vessel injury, e.g. dissection, perforation, rupture • Vessel occlusion • Vessel spasm

Additionally, the following complications may result from a balloon dilatation procedure within the biliary system: • Cholangitis • Hemobilia • Pancreatitis • Duct injury, e.g. dissection, perforation, rupture • Duct occlusion • Duct spasm

Mustang is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates.

Boston Scientific
Advancing science for life™

Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752 USA
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2014 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-225202-AA NOV2014

Mustang™ Balloon Dilatation Catheter

40 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171030240	39171-03024	08714729794639	3	20	5	24	C1725
H74939171030340	39171-03034	08714729795025	3	30	5	24	C1725
H74939171030440	39171-03044	08714729794646	3	40	5	24	C1725
H74939171030640	39171-03064	08714729795032	3	60	5	24	C1725
H74939171030840	39171-03084	08714729795049	3	80	5	24	C1725
H74939171031040	39171-03104	08714729795308	3	100	5	24	C1725
H74939171040240	39171-04024	08714729794172	4	20	5	24	C1725
H74939171040340	39171-04034	08714729795063	4	30	5	24	C1725
H74939171040440	39171-04044	08714729794103	4	40	5	24	C1725
H74939171040640	39171-04064	08714729795070	4	60	5	24	C1725
H74939171040840	39171-04084	08714729795087	4	80	5	24	C1725
H74939171041040	39171-04104	08714729795315	4	100	5	24	C1725
H74939171050240	39171-05024	08714729794066	5	20	5	24	C1725
H74939171050340	39171-05034	08714729795094	5	30	5	24	C1725
H74939171050440	39171-05044	08714729793731	5	40	5	24	C1725
H74939171050640	39171-05064	08714729795100	5	60	5	24	C1725
H74939171050840	39171-05084	08714729795117	5	80	5	24	C1725
H74939171051040	39171-05104	08714729795124	5	100	5	24	C1725
H74939171051240	39171-05124	08714729794882	5	120	5	24	C1725
H74939171060240	39171-06024	08714729793908	6	20	5	24	C1725
H74939171060340	39171-06034	08714729795131	6	30	5	24	C1725
H74939171060440	39171-06044	08714729793564	6	40	5	24	C1725
H74939171060640	39171-06064	08714729794783	6	60	5	24	C1725
H74939171060840	39171-06084	08714729795148	6	80	5	24	C1725
H74939171061040	39171-06104	08714729794899	6	100	5	24	C1725
H74939171061240	39171-06124	08714729794905	6	120	5	22	C1725
H74939171070240	39171-07024	08714729794370	7	20	5	20	C1725
H74939171070340	39171-07034	08714729795155	7	30	5	20	C1725
H74939171070440	39171-07044	08714729793588	7	40	5	20	C1725
H74939171070640	39171-07064	08714729794790	7	60	5	20	C1725
H74939171070840	39171-07084	08714729794806	7	80	5	20	C1725
H74939171071040	39171-07104	08714729794813	7	100	5	20	C1725
H74939171080240	39171-08024	08714729794660	8	20	6	20	C1725
H74939171080340	39171-08034	08714729794400	8	30	6	20	C1725
H74939171080440	39171-08044	08714729793632	8	40	6	20	C1725
H74939171080640	39171-08064	08714729794820	8	60	6	20	C1725
H74939171080840	39171-08084	08714729794837	8	80	6	20	C1725
H74939171081040	39171-08104	08714729794844	8	100	6	20	C1725
H74939171090240	39171-09024	08714729794912	9	20	6	18	C1725
H74939171090340	39171-09034	08714729795179	9	30	6	18	C1725
H74939171090440	39171-09044	08714729794622	9	40	6	18	C1725
H74939171090640	39171-09064	08714729795193	9	60	6	18	C1725
H74939171090840	39171-09084	08714729795216	9	80	6	18	C1725
H74939171100240	39171-10024	08714729795322	10	20	6	14	C1725

Mustang™ Balloon Dilatation Catheter (cont.)

40 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171100340	39171-10034	08714729795339	10	30	6	14	C1725
H74939171100440	39171-10044	08714729794226	10	40	6	14	C1725
H74939171100640	39171-10064	08714729795247	10	60	6	14	C1725
H74939171100840	39171-10084	08714729794721	10	80	6	14	C1725
H74939171120240	39171-12024	08714729794929	12	20	7	14	C1725
H74939171120340	39171-12034	08714729795346	12	30	7	14	C1725
H74939171120440	39171-12044	08714729794592	12	40	7	14	C1725
H74939171120640	39171-12064	08714729795278	12	60	7	14	C1725
H74939171120840	39171-12084	08714729795292	12	80	7	14	C1725

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171030270	39171-03027	08714729793991	3	20	5	24	C1725
H74939171030370	39171-03037	08714729794981	3	30	5	24	C1725
H74939171030470	39171-03047	08714729793946	3	40	5	24	C1725
H74939171030670	39171-03067	08714729794943	3	60	5	24	C1725
H74939171030870	39171-03087	08714729794950	3	80	5	24	C1725
H74939171031070	39171-03107	08714729794967	3	100	5	24	C1725
H74939171031270	39171-03127	08714729794554	3	120	5	24	C1725
H74939171031570	39171-03157	08714729794561	3	150	5	24	C1725
H74939171031870	39171-03187	08714729794547	3	180	5	24	C1725
H74939171032070	39171-03207	08714729794530	3	200	5	24	C1725
H74939171040270	39171-04027	08714729793496	4	20	5	24	C1725
H74939171040370	39171-04037	08714729794752	4	30	5	24	C1725
H74939171040470	39171-04047	08714729793403	4	40	5	24	C1725
H74939171040670	39171-04067	08714729794349	4	60	5	24	C1725
H74939171040870	39171-04087	08714729794363	4	80	5	24	C1725
H74939171041070	39171-04107	08714729794240	4	100	5	24	C1725
H74939171041270	39171-04127	08714729794158	4	120	5	24	C1725
H74939171041570	39171-04157	08714729794462	4	150	5	24	C1725
H74939171041870	39171-04187	08714729794479	4	180	5	24	C1725
H74939171042070	39171-04207	08714729794455	4	200	5	24	C1725
H74939171050270	39171-05027	08714729793427	5	20	5	24	C1725
H74939171050370	39171-05037	08714729794707	5	30	5	24	C1725
H74939171050470	39171-05047	08714729793366	5	40	5	24	C1725
H74939171050670	39171-05067	08714729793977	5	60	5	24	C1725
H74939171050870	39171-05087	08714729793960	5	80	5	24	C1725
H74939171051070	39171-05107	08714729793892	5	100	5	24	C1725
H74939171051270	39171-05127	08714729793915	5	120	5	24	C1725
H74939171051570	39171-05157	08714729794042	5	150	5	24	C1725
H74939171051870	39171-05187	08714729794486	5	180	5	24	C1725

Mustang™ Balloon Dilatation Catheter (cont)
75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171052070	39171-05207	08714729794202	5	200	5	24	C1725
H74939171060270	39171-06027	08714729793397	6	20	5	24	C1725
H74939171060370	39171-06037	08714729794318	6	30	5	24	C1725
H74939171060470	39171-06047	08714729793335	6	40	5	24	C1725
H74939171060670	39171-06067	08714729793762	6	60	5	24	C1725
H74939171060870	39171-06087	08714729793854	6	80	5	24	C1725
H74939171061070	39171-06107	08714729793625	6	100	5	24	C1725
H74939171061270	39171-06127	08714729794127	6	120	5	22	C1725
H74939171061570	39171-06157	08714729793984	6	150	5	22	C1725
H74939171061870	39171-06187	08714729794493	6	180	5	22	C1725
H74939171062070	39171-06207	08714729794509	6	200	5	22	C1725
H74939171070270	39171-07027	08714729793502	7	20	5	20	C1725
H74939171070370	39171-07037	08714729794271	7	30	5	20	C1725
H74939171070470	39171-07047	08714729793342	7	40	5	20	C1725
H74939171070670	39171-07067	08714729793861	7	60	5	20	C1725
H74939171070870	39171-07087	08714729794264	7	80	5	20	C1725
H74939171071070	39171-07107	08714729793779	7	100	5	20	C1725
H74939171071270	39171-07127	08714729794424	7	120	5	18	C1725
H74939171071570	39171-07157	08714729794578	7	150	5	18	C1725
H74939171071870	39171-07187	08714729794608	7	180	5	18	C1725
H74939171072070	39171-07207	08714729794189	7	200	5	18	C1725
H74939171080270	39171-08027	08714729793519	8	20	6	20	C1725
H74939171080370	39171-08037	08714729793809	8	30	6	20	C1725
H74939171080470	39171-08047	08714729793328	8	40	6	20	C1725
H74939171080670	39171-08067	08714729793885	8	60	6	20	C1725
H74939171080870	39171-08087	08714729793656	8	80	6	20	C1725
H74939171081070	39171-08107	08714729793700	8	100	6	20	C1725
H74939171081270	39171-08127	08714729794431	8	120	6	18	C1725
H74939171081570	39171-08157	08714729794585	8	150	6	18	C1725
H74939171081870	39171-08187	08714729794615	8	180	6	18	C1725
H74939171082070	39171-08207	08714729794196	8	200	6	18	C1725
H74939171090270	39171-09027	08714729793755	9	20	6	18	C1725
H74939171090370	39171-09037	08714729794851	9	30	6	18	C1725
H74939171090470	39171-09047	08714729793489	9	40	6	18	C1725
H74939171090670	39171-09067	08714729794691	9	60	6	18	C1725
H74939171090870	39171-09087	08714729794868	9	80	6	18	C1725
H74939171100270	39171-10027	08714729793601	10	20	6	14	C1725
H74939171100370	39171-10037	08714729794394	10	30	6	14	C1725
H74939171100470	39171-10047	08714729793359	10	40	6	14	C1725
H74939171100670	39171-10067	08714729794004	10	60	6	14	C1725
H74939171100870	39171-10087	08714729794332	10	80	6	14	C1725
H74939171120270	39171-12027	08714729793922	12	20	7	14	C1725

Mustang™ Balloon Dilatation Catheter (cont)
75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171120370	39171-12037	08714729794875	12	30	7	14	C1725
H74939171120470	39171-12047	08714729793441	12	40	7	14	C1725
H74939171120670	39171-12067	08714729794738	12	60	7	14	C1725
H74939171120870	39171-12087	08714729794745	12	80	7	14	C1725

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171030210	39171-03021	08714729794141	3	20	5	24	C1725
H74939171030310	39171-03031	08714729794974	3	30	5	24	C1725
H74939171030410	39171-03041	08714729794035	3	40	5	24	C1725
H74939171030610	39171-03061	08714729794998	3	60	5	24	C1725
H74939171030810	39171-03081	08714729794936	3	80	5	24	C1725
H74939171031010	39171-03101	08714729794295	3	100	5	24	C1725
H74939171031210	39171-03121	08714729793816	3	120	5	24	C1725
H74939171031510	39171-03151	08714729794516	3	150	5	24	C1725
H74939171031810	39171-03181	08714729794523	3	180	5	24	C1725
H74939171032010	39171-03201	08714729794387	3	200	5	24	C1725
H74939171040210	39171-04021	08714729793724	4	20	5	24	C1725
H74939171040310	39171-04031	08714729795056	4	30	5	24	C1725
H74939171040410	39171-04041	08714729793557	4	40	5	24	C1725
H74939171040610	39171-04061	08714729794257	4	60	5	24	C1725
H74939171040810	39171-04081	08714729794028	4	80	5	24	C1725
H74939171041010	39171-04101	08714729793717	4	100	5	24	C1725
H74939171041210	39171-04121	08714729793670	4	120	5	24	C1725
H74939171041510	39171-04151	08714729793694	4	150	5	24	C1725
H74939171041810	39171-04181	08714729794653	4	180	5	24	C1725
H74939171042010	39171-04201	08714729793595	4	200	5	24	C1725
H74939171050210	39171-05021	08714729793618	5	20	5	24	C1725
H74939171050310	39171-05031	08714729794769	5	30	5	24	C1725
H74939171050410	39171-05041	08714729793465	5	40	5	24	C1725
H74939171050610	39171-05061	08714729793878	5	60	5	24	C1725
H74939171050810	39171-05081	08714729793687	5	80	5	24	C1725
H74939171051010	39171-05101	08714729793533	5	100	5	24	C1725
H74939171051210	39171-05121	08714729793410	5	120	5	24	C1725
H74939171051510	39171-05151	08714729793434	5	150	5	24	C1725
H74939171051810	39171-05181	08714729794219	5	180	5	24	C1725
H74939171052010	39171-05201	08714729793373	5	200	5	24	C1725
H74939171060210	39171-06021	08714729793571	6	20	5	24	C1725
H74939171060310	39171-06031	08714729794776	6	30	5	24	C1725
H74939171060410	39171-06041	08714729793472	6	40	5	24	C1725
H74939171060610	39171-06061	08714729793939	6	60	5	24	C1725
H74939171060810	39171-06081	08714729793793	6	80	5	24	C1725
H74939171061010	39171-06101	08714729793540	6	100	5	24	C1725
H74939171061210	39171-06121	08714729793526	6	120	5	22	C1725

Mustang™ Balloon Dilatation Catheter (cont)
135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939171061510	39171-06151	08714729793458	6	150	5	22	C1725
H74939171061810	39171-06181	08714729794059	6	180	5	22	C1725
H74939171062010	39171-06201	08714729793380	6	200	5	22	C1725
H74939171070210	39171-07021	08714729793830	7	20	5	20	C1725
H74939171070310	39171-07031	08714729795001	7	30	5	20	C1725
H74939171070410	39171-07041	08714729793649	7	40	5	20	C1725
H74939171070610	39171-07061	08714729794325	7	60	5	20	C1725
H74939171070810	39171-07081	08714729794288	7	80	5	20	C1725
H74939171071010	39171-07101	08714729794011	7	100	5	20	C1725
H74939171071210	39171-07121	08714729794097	7	120	5	18	C1725
H74939171071510	39171-07151	08714729793786	7	150	5	18	C1725
H74939171071810	39171-07181	08714729794677	7	180	5	18	C1725
H74939171072010	39171-07201	08714729793748	7	200	5	18	C1725
H74939171080210	39171-08021	08714729793847	8	20	6	20	C1725
H74939171080310	39171-08031	08714729794417	8	30	6	20	C1725
H74939171080410	39171-08041	08714729793663	8	40	6	20	C1725
H74939171080610	39171-08061	08714729794714	8	60	6	20	C1725
H74939171080810	39171-08081	08714729794301	8	80	6	20	C1725
H74939171081010	39171-08101	08714729794233	8	100	6	20	C1725
H74939171081210	39171-08121	08714729794448	8	120	6	18	C1725
H74939171081510	39171-08151	08714729794110	8	150	6	18	C1725
H74939171081810	39171-08181	08714729794684	8	180	6	18	C1725
H74939171082010	39171-08201	08714729794080	8	200	6	18	C1725
H74939171090210	39171-09021	08714729794165	9	20	6	18	C1725
H74939171090310	39171-09031	08714729795162	9	30	6	18	C1725
H74939171090410	39171-09041	08714729794134	9	40	6	18	C1725
H74939171090610	39171-09061	08714729795186	9	60	6	18	C1725
H74939171090810	39171-09081	08714729795209	9	80	6	18	C1725
H74939171100210	39171-10021	08714729794073	10	20	6	14	C1725
H74939171100310	39171-10031	08714729795223	10	30	6	14	C1725
H74939171100410	39171-10041	08714729793823	10	40	6	14	C1725
H74939171100610	39171-10061	08714729795230	10	60	6	14	C1725
H74939171100810	39171-10081	08714729795254	10	80	6	14	C1725
H74939171120210	39171-12021	08714729794356	12	20	7	14	C1725
H74939171120310	39171-12031	08714729795261	12	30	7	14	C1725
H74939171120410	39171-12041	08714729793953	12	40	7	14	C1725
H74939171120610	39171-12061	08714729795018	12	60	7	14	C1725
H74939171120810	39171-12081	08714729795285	12	80	7	14	C1725

Charger™ Balloon Dilatation Catheter
75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code
H74939206030270	39206-03027	08714729803942	3	20	5	20	C1725
H74939206030370	39206-03037	08714729803966	3	30	5	20	C1725
H74939206030470	39206-03047	08714729803980	3	40	5	20	C1725