



# AESCULAP® POWER SYSTEMS

CATALOG 2022

ELAN 4 and Acculan 4

# AESCULAP® ELAN 4

PRECISION. PERFORMANCE. PLEASURE.

With its technical concept ELAN 4 addresses important requirements of the surgeon and all parties involved in handling power systems in the hospital.

- Convenient and reliable grip
- Powerful, precise and reliable work
- Easy to set up and simple to use



# AESCULAP® ELAN 4

ELAN 4 electro

ELECTRIC HIGHSPEED SYSTEM



## ELAN 4 electro

### BASIC COMPONENTS


#### GA800

##### ELAN 4 electro CONTROL UNIT

with irrigation pump; touch screen;  
power cord not included  
pump flow 0–65 ml/min

100–240 V AC; 50/60 Hz;  
380 x 330 x 200 mm; 9,500 g

consisting parts:  
GD412804 bottle holder



#### GA808

##### ELAN 4 electro FOOT CONTROL

protection type IPX8; anaesthetic  
proof class AP

5.0 m;  
245 x 144 x 76 mm; 2,240 g



#### GA810

##### ELAN 4 electro WIRELESS FOOT CONTROL

protection type IPX8; anaesthetic  
proof class AP

271 x 266 x 85 mm; 3,800 g

Requires control unit with software  
version V3.00.



#### GA806

#### GB073R

##### ELAN 4 electro MOTOR CABLE

for use with ELAN 4 electro applied  
parts

4.0 m; 270 g

consisting parts:  
GB600860 Sterilit® spray adapter  
ELAN 4 electro



# AESCULAP® ELAN 4

## ELAN 4 electro

### BASIC COMPONENTS, PERFORATOR DRIVER & LOWSPEED MOTOR

#### GA804

GB073R

##### ELAN 4 electro MOTOR CABLE WITH HAND SWITCH

for use with ELAN 4 electro applied  
parts

4.0 m; 374 g

consisting parts:  
GB600860 Sterilit® spray adapter  
ELAN 4 electro



#### GA805

GB073R

##### ELAN 4 electro MOTOR CABLE WITH HAND CONTROL

for use with ELAN 4 electro applied  
parts

4.0 m; 369 g

consisting parts:  
GB600860 Sterilit® spray adapter  
ELAN 4 electro



#### GA822

GB076R

##### ELAN 4 electro PERFORATOR DRIVER

for use with cranial perforators with  
Hudson shank

180 W; 0 - 1,200 rpm; 250 Ncm;  
166 x 35 mm; 410 g



#### GA824


GB077R

##### ELAN 4 electro LOWSPEED MOTOR

with intra coupling

180 W; 0 - 40,000\* rpm;  
10 Ncm; 123 x 23 mm; 150 g

\* Requires control unit with software  
version V3.00. With lower versions  
the motor runs max. 20,000 rpm.



## ELAN 4 electro

### CRANIOTOME & ATTACHMENTS

**GA849**

**GB085R**

#### ELAN 4 electro CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with dura guards and holding sleeves GB941R - GB947R; spray nozzle GA259SU (not included)

140 W;  0 - 80,000 rpm; 2.2 Ncm; 156 x 17 mm;  121 g



**GB941R**

**GB719R**

#### ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding; for craniotome and multifunction handpieces GA749, GA849



**GB942R**

**GB719R**

#### ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



**GB943R**

**GB719R**

#### ELAN 4 DURA GUARD LONG

fixed; 3 dot coding; for craniotome and multifunction handpieces GA749, GA849



**GB945R**

**GB719R**

#### ELAN 4 HOLDING SLEEVE

for craniotome and multifunction handpieces GA749, GA849



**GB947R**

**GB719R**

#### ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



# AESCULAP® ELAN 4

## ELAN 4 electro

### 1-RING HIGHSPEED HANDPIECES

#### GA861

GB083R

##### ELAN 4 electro 1-RING HANDPIECE L4

1-ring coding;  
spray nozzle GA261SU (not included)  
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;  
152 x 16 mm;  80 g



#### GA862

GB083R

##### ELAN 4 electro 1-RING HANDPIECE L7

1-ring coding;  
spray nozzle GA262SU (not included)  
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;  
177 x 16 mm;  83 g



#### GA863

GB083R

##### ELAN 4 electro 1-RING HANDPIECE L10

1-ring coding;  
spray nozzle GA263SU (not included)  
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;  
206 x 16 mm;  87 g



#### GA864

GB083R

##### ELAN 4 electro 1-RING HANDPIECE L13

1-ring coding;  
spray nozzle GA264SU (not included)  
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;  
237 x 16 mm;  91 g



## ELAN 4 electro

### 2-RING HIGHSPEED HANDPIECES

#### GA865

GB083R

##### ELAN 4 electro 2-RING HANDPIECE L7

2-ring coding;  
spray nozzle GA265SU (not included)  
140 W;  ; 0–80,000 rpm; 2.2 Ncm;  
176 x 16 mm;  89 g



#### GA866

GB083R

##### ELAN 4 electro 2-RING HANDPIECE L10

2-ring coding;  
spray nozzle GA266SU (not included)  
140 W;  ; 0–80,000 rpm; 2.2 Ncm;  
206 x 16 mm;  95 g




#### GA867

GB083R

##### ELAN 4 electro 2-RING HANDPIECE L13

2-ring coding;  
spray nozzle GA267SU (not included)  
140 W;  ; 0–80,000 rpm; 2.2 Ncm;  
236 x 16 mm;  100 g



#### GA868

GB083R

##### ELAN 4 electro 2-RING HANDPIECE L17

2-ring coding;  
spray nozzle GA268SU (not included)  
140 W;  ; 0–80,000 rpm; 2.2 Ncm;  
276 x 16 mm;  108 g



#### GA869

GB083R

##### ELAN 4 electro 2-RING HANDPIECE L22

2-ring coding;  
spray nozzle GA269SU (not included)  
140 W;  ; 0–80,000 rpm; 2.2 Ncm;  
326 x 16 mm;  118 g



# AESCULAP® ELAN 4

ELAN 4 electro

MIS HANDPIECE & HANDPIECE SHAFTS



## GA860

GB084R

### ELAN 4 electro MIS HANDPIECE

for use with MIS handpiece shafts  
GB920R - GB927R

140 W;  ; 0 - 80,000 rpm;  
2.2 Ncm; 118 x 17 mm;  100 g



## GB920R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L10 STRAIGHT

spray nozzle GB796SU (not included)



## GB921R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L10 CURVED

spray nozzle GB796SU (not included)



## GB922R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L10 HEAVY CURVED

spray nozzle GB796SU (not included)



## GB925R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L13 STRAIGHT

spray nozzle GB797SU (not included)



## GB926R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L13 CURVED

spray nozzle GB797SU (not included)



## GB927R

GB723R

### ELAN 4 MIS HANDPIECE SHAFT L13 HEAVY CURVED

spray nozzle GB797SU (not included)



# AESCLAP<sup>®</sup> ELAN 4


ELAN 4 electro

SAWS



**GA836****GB086R****ELAN 4 electro MICRO SAGITTAL SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 161 x 23 mm;  
∠ 7°; 190 g**GA837****GB087R****ELAN 4 electro  
MICRO RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 184 x 23 mm;  
1.3 mm; 198 g**GA831****GB088R****ELAN 4 electro SAGITTAL SAW**


spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 129 x 23 mm;  
4.5°; 220 g**GA832****GB087R****ELAN 4 electro RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 194 x 23 mm;  
2.5 mm; 238 g**GA833****GB088R****ELAN 4 electro TRANSVERSAL SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 10 Ncm;  
130 x 23 mm; 10°; 185 gto be used with saw blade shafts  
GP580R, GP582R, GP586R and GP590R

**GA844****GB482R****ELAN 4 electro DRILL**

with fixed motor cable, oscillating and tapping mode

200 W;  3.2 mm;  
122 x 130 x 27 mm;  800 g

consisting parts:

GA344244 K-wire protection guide

TA011944 cleaning brush

GA344211 rinsing adapter

GB600880 Sterilit® spray adapter

Requires control unit with software version V3.00.

**GB870R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with Jacobs chuck 0.5 – 7.4 mm, not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;  
105 x 29 mm;  230 g

consisting parts:

key GA031R

**GB871R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with large AO chuck, not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;  
88 x 29 mm;  170 g

**GB872R****GB716R****INTRAMEDULLARY REAMER ATTACHMENT**

with Hudson/Zimmer chuck, not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;  
86 x 25 mm;  150 g



# AESCULAP® ELAN 4

ELAN 4 electro

ATTACHMENTS

**GB881R**

**GB716R**

## DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;  
89 x 29 mm;  163 g

consisting parts:  
key GA031R



**GB882R**

**GB716R**

## DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;  
76 x 22 mm;  100 g

consisting parts:  
key GA062R



**GB884R**

**GB716R**

## DRILL ATTACHMENT

with small AO chuck

0 – 1,250 rpm; 2 Nm;  2.8 mm;  
60 x 22 mm;  66 g



**GB886R**

**GB716R**

## DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,250 rpm; 2 Nm;  3.2 mm;  
57 x 24 mm;  67 g



**GB887R**

**GB716R**

## DRILL ATTACHMENT

with Dental shank

0 – 1,250 rpm; 2 Nm;  1.8 mm;  
56 x 22 mm;  50 g



**GB878R**

**GB716R**

**DRILL ATTACHMENT**

for Synthes radiolucent drive 511.300

0 - 1,250 rpm; 2 Nm;  
111 x 27 mm; 181 g



**GB894R**

**GB716R**

**K-WIRE ATTACHMENT**

for K-wires from 0.6 - 3.2 mm

0 - 1,250 rpm; 3.2 mm;  
92 x 27 x 124 mm; 179 g



**GB891R**

**GB716R**

**SAGITTAL SAW ATTACHMENT**

for "L-sagittal" sawblades

0 - 15,750 rpm; 74 x 33.5 mm;  
4.5°; 179 g



**GB892R**

**GB716R**

**RECIPROCATING SAW ATTACHMENT**

for reciprocating sawblades

0 - 15,750 rpm; 135 x 24.5 mm;  
2.5 mm; 166 g



# AESCULAP® ELAN 4

ELAN 4 electro

ATTACHMENTS

**GB896R**

**GB716R**

## SCREW ATTACHMENT

with small AO chuck

0 - 165 rpm; 12 Nm;  1.7 mm;  
95 x 28 x 129.5 mm;  235 g



**GB897R**


**GB716R**

## SCREW ATTACHMENT

with AESCULAP® hexagonal chuck

0 - 165 rpm; 12 Nm;  1.7 mm;  
92 x 28 x 129.5 mm;  231 g





# AESCLAP<sup>®</sup> ELAN 4

## ELAN 4 ACCESSORIES

### GB861P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE  
for use with 1-ring handpieces



### GB863P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE  
for use with 2-ring handpieces



### GB862P

GB727R

ELAN 4 MIS SLEEVE  
for use with 1-ring handpieces



### GB866R

ELAN 4 LATERAL BURR GUARD  
for use with 1-ring handpieces and  
burrs up to 4.5 mm diameter



### TE780

POWER CORD  
Europe, 1.5 m

### TE730

POWER CORD  
Europe, 5.0 m

### TE734

POWER CORD  
UK, 5.0 m

### TE735


POWER CORD  
USA, Canada, Japan, 3.5 m



## GA395SU

sterile   10

ELAN 4 electro DISPOSABLE TUBE SET



## ELAN 4 SPRAY NOZZELS

sterile 

### GA259SU

for craniotome attachments

### GA261SU

for 1-ring handpiece L4

### GA262SU

for 1-ring handpiece L7

### GA263SU

for 1-ring handpiece L10

### GA264SU

for 1-ring handpiece L13

### GA265SU

for 2-ring handpiece L7

### GA266SU

for 2-ring handpiece L10

### GA267SU

for 2-ring handpiece L13

### GA268SU

for 2-ring handpiece L17

### GA269SU

for 2-ring handpiece L22

### GB796SU

for MIS handpiece shafts L10

### GB797SU

for MIS handpiece shafts L13

### GA258SU

for saws





# AESCULAP® MECHANICAL PROCESSING

- The solution for storage and mechanical reprocessing of our AESCULAP® Power Systems.
- For detailed information please refer to the instructions for use.



# AESCU LAP<sup>®</sup> MECHANICAL PROCESSING

## ECCOS<sup>®</sup> SETS

### ELAN 4

#### GB261R

##### ELAN 4 electro ECCOS<sup>®</sup> SET

 56 mm

consisting parts:

GB073R

GB076R


GB083R

GB085R

GB719R

GB072R

JF212R



#### GB262R

##### ELAN 4 electro ECCOS<sup>®</sup> SET FOR DRILL GA844

 106 mm

consisting parts:

GB072R

GB073R

GB482R


GB715R

GB716R

2x GB690R

JG342

JG344



ECCOS®  
ELAN 4 electro

**GB073R**

ECCOS® HOLDER

for motor cable GA804, GA805, GA806  
(set consisting of 2 pieces)

 56 mm



**GB076R**

ECCOS® HOLDER

for perforator driver GA822

 56 mm



**GB077R**

ECCOS® HOLDER

for intra motor GA824

 56 mm




**GB082R**

ECCOS® HOLDER

for 1 handpiece

 56 mm




**GB083R**

ECCOS® HOLDER

for 3 handpieces

 56 mm



**GB085R**

ECCOS® HOLDER

for craniotome GA849

 56 mm



**GB086R**

ECCOS® HOLDER

for micro sagittal saw GA836

 56 mm



**GB087R**

ECCOS® HOLDER

for reciprocating saws GA832, GA837

 56 mm



**GB088R**

ECCOS® HOLDER

for saws GA831, GA833



 56 mm



# AESULAP<sup>®</sup> MECHANICAL PROCESSING

ECCOS<sup>®</sup>

ELAN 4 electro

| GB482R                                                                                                                                                                | GB715R                                                                                                                                                                                                                         | GB716R                                                                                                                                                                                                                            |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>ECCOS<sup>®</sup> HOLDER</p> <p>for ELAN 4 electro drill GA844</p> <p> 106 mm</p> | <p>ECCOS<sup>®</sup> HOLDER</p> <p>for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p> | <p>ECCOS<sup>®</sup> HOLDER</p> <p>for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p> |
|                                                                                      |                                                                                                                                               |                                                                                                                                                |
|                                                                                     |                                                                                                                                              |                                                                                                                                               |

## ECCOS® & RINSING DEVICES


### ELAN 4 electro

#### GB084R

##### ECCOS® HOLDER

for MIS handpiece GA860

 56 mm



#### GB723R

##### ECCOS® HOLDER

for MIS handpiece shafts GB920R - GB927R


 56 mm



#### GB692R

##### RINSING DEVICE


for applied parts, incl. ceramic filter GB693  
(not suitable for autoclaving)



#### GB679R

##### RINSING DEVICE

for 4 MIS shafts, incl. ceramic filter GB693  
(not suitable for autoclaving)



#### GB698R

##### RINSING ADAPTER

for applied parts  
(for manual cleaning)



#### GB699R

##### RINSING ADAPTER

for MIS shafts  
(for manual cleaning)



# AESCU LAP<sup>®</sup> MECHANICAL PROCESSING

ECCOS<sup>®</sup>


ELAN 4 ACCESSORIES

## GB718R

ECCOS<sup>®</sup> HOLDER

for 12 1- and 2-ring burrs

 56 mm



## GB719R

ECCOS<sup>®</sup> HOLDER

for 3 craniotome attachments

 56 mm



## GB720R

ECCOS<sup>®</sup> UNIVERSAL HOLDER

for saw blades

 56 mm



## GB729R

ECCOS<sup>®</sup> HOLDER

for 2 reciprocating saw shafts

 56 mm



## GB730R

ECCOS<sup>®</sup> HOLDER

for 2 transversal saw shafts

 56 mm



## GB072R

ECCOS<sup>®</sup> CABLE/HOSE POS. GUIDE

for cable and hoses (set of 4 pieces)

 56 mm



## GB727R

ECCOS<sup>®</sup> HOLDER

for handpiece holding-sleeves  
GB861P, GB862P, GB863P

 56 mm






## GB572R

ECCOS<sup>®</sup> FIXATION

for GB300R - GB304R


 56 mm






1 sterile

Craniotome cutter – spiral type  
 Kraniotomfräser – spiralig  
 Fraise craniotome – en spirale  
 Fresa cranéotomo – helicoidal  
 Fresa di craniotomo – a spirale




1 sterile

Craniotome cutter – straight type  
 Kraniotomfräser – geradverzahnt  
 Fraise craniotome – dents droites  
 Fresa cranéotomo – dientes rectos  
 Fresa di craniotomo – denti dritti




1 sterile

Lindemann




1 sterile

Twist drill  
 Spiralbohrer  
 Fraise spirale  
 Broca espiral  
 Punta elicoidale




1 sterile

Twist drill  
 Spiralbohrer  
 Fraise spirale  
 Broca espiral  
 Punta elicoidale



1 sterile

Pin cutter  
 Pinfräser  
 Fraise pointue  
 Fresa punta  
 Fresa puntuale



1 sterile

Twin-Cut burr  
 Twin-Cut Fräser  
 Fraise Twin-Cut  
 Fresa Twin-Cut  
 Fresa Twin-Cut

# AESCU LAP<sup>®</sup> BURRS & BLADES


## ELAN 4 BURRS – 2-RING

ELAN 4 Fräser – 2-Ring

Fraises ELAN 4 – 2 anneaux

Fresa ELAN 4 – 2 anillos

Fresa ELAN 4 – 2 anelli



### Tungsten Carbide Burrs

| Ø                                                                                 | 2.3 mm                                                                            | 3.0 mm                                                                            | 4.0 mm                                                                            | 5.0 mm                                                                            | 6.0 mm                                                                            | 7.0 mm                                                                            | 8.0 mm                                                                            |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
|  |  |  |  |  |  |  |  |
| II                                                                                | GP301TC                                                                           | GP302TC                                                                           | GP304TC                                                                           | GP306TC                                                                           | GP307TC                                                                           | GP308TC                                                                           | GP309TC                                                                           |



Rosen

| Ø                                                                                   | 1.2 mm                                                                             | 1.6 mm                                                                             |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|  |  |  |
| II                                                                                  | GP395TC                                                                            | GP396TC                                                                            |



Taper burr  
Konusfräser  
Fraise conique  
Fresa cónica  
Fresa conica

| Ø                                                                                   | 3.0 mm                                                                             |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|  |  |
| II                                                                                  | GP367TC                                                                            |




Neuro cutter  
Olivenfräser  
Fraise à olive  
Fresa oliva  
Fresa ad oliva

### Burrs for metal dissection

| Ø                                                                                   | 2.3 mm                                                                              |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |
| II ②                                                                                | GP399TC-SU                                                                          |



Side cutter  
Zylinderfräser  
Fraise cylindrique  
Fresa cilíndrica  
Fresa cilíndrica

| Ø                                                                                   | 25 mm                                                                               |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |
| II ②                                                                                | GP398SU                                                                             |



Diamond disc  
Diamant-Trennscheibe  
Disque cou pant diamant  
Disco cortante de diamante  
Disco tagliente diamante

## CARE AND MAINTENANCE


Pflege und Aufbereitung  
Entretien et maintenance  
Conservación y mantenimiento  
Cura e manutenzione



**GB600860**  
ELAN 4 electro



**GB600870**  
ELAN 4 MIS shafts,  
microspeed uni XS



**GB600880**  
Acculan 4 small drill




**GB600840**  
Acculan 4, Acculan 3Ti



**GB600820**  
HiLAN, HiLAN XS, Hi-Line  
(XS, XXS), Microspeed EC,  
microspeed uni, microLAN,  
miniLAN



**GB600830**  
mini-Line



**GB600810**  
micro-Line



**GB600**



Sterilit® Power Systems  
consisting parts:  
3x GB600890 spray tube

etc.

# Basis - Size 3/4

## Container bodies

**without** base perforation = JK



**with** base perforation = JN



## Lids for individual color combination

Inner lid



Outer lid



|  |  | Outside dimensions<br>Length x Width x Height (mm) |
|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|----------------------------------------------------|
| JK740                                                                              | JN740                                                                              | 470 x 274 x 90                                     |
| JK741                                                                              | JN741                                                                              | 470 x 274 x 120                                    |
| JK742                                                                              | JN742                                                                              | 470 x 274 x 135                                    |

|  | Color  |  | Color  |
|------------------------------------------------------------------------------------|--------|--------------------------------------------------------------------------------------|--------|
| JK785                                                                              | red    | JK775A                                                                               | red    |
| JK786                                                                              | blue   | JK776A                                                                               | blue   |
| JK787                                                                              | green  | JK777A                                                                               | green  |
| JK788                                                                              | gold   | JK778A                                                                               | gold   |
| JK789                                                                              | silver | JK779A                                                                               | silver |

## Dimensions complete - Length x Width x Height (mm)





|  | Outside dimensions<br>Inside dimensions | Outside dimensions<br>Inside dimensions |  | Outside dimensions<br>Inside dimensions | Outside dimensions<br>Inside dimensions |
|-------------------------------------------------------------------------------------|-----------------------------------------|-----------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------|-----------------------------------------|
| Body 3/4                                                                            | with inner lid                          | with outer lid                          | Body 3/4                                                                            | with inner lid                          | with outer lid                          |
| JK740                                                                               | 470 x 285 x 108<br>421 x 258 x 75       | 470 x 291 x 115<br>421 x 258 x 75       | JN740                                                                               | 470 x 285 x 108<br>421 x 258 x 63       | 470 x 291 x 115<br>421 x 258 x 63       |
| JK741                                                                               | 470 x 285 x 138<br>421 x 258 x 105      | 470 x 291 x 145<br>421 x 258 x 105      | JN741                                                                               | 470 x 285 x 138<br>421 x 258 x 93       | 470 x 291 x 145<br>421 x 258 x 93       |
| JK742                                                                               | 470 x 285 x 153<br>421 x 258 x 120      | 470 x 291 x 160<br>421 x 258 x 120      | JN742                                                                               | 470 x 285 x 153<br>421 x 258 x 108      | 470 x 291 x 160<br>421 x 258 x 108      |

All dimensions: Length x Width x Height (mm)




## Basket **with** feet







|                                                                                        |  | Outside dim.<br>Inside dim.       |  | Outside dim.<br>Inside dim.       |  | Outside dim.<br>Inside dim.       |  | Outside dim.<br>Inside dim.       |
|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------|-----------------------------------------------------------------------------------|-----------------------------------|-----------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------------------------------------------------|-----------------------------------|
|  36    | JF111R                                                                            | 243 x 253 x 36<br>241 x 251 x 22  | JF251R                                                                            | 406 x 253 x 36<br>404 x 251 x 22  | JF211R                                                                            | 485 x 253 x 36<br>483 x 251 x 22  | JF221R                                                                              | 540 x 253 x 36<br>538 x 251 x 22  |
|  56    | JF112R                                                                            | 243 x 253 x 56<br>241 x 251 x 42  | JF252R                                                                            | 406 x 253 x 56<br>404 x 251 x 42  | JF212R                                                                            | 485 x 253 x 56<br>483 x 251 x 42  | JF222R                                                                              | 540 x 253 x 56<br>538 x 251 x 42  |
|  76  | JF113R                                                                            | 243 x 253 x 76<br>241 x 251 x 62  | JF253R                                                                            | 406 x 253 x 76<br>404 x 251 x 62  | JF213R                                                                            | 485 x 253 x 76<br>483 x 251 x 62  | JF223R                                                                              | 540 x 253 x 76<br>538 x 251 x 62  |
|  106 | JF114R                                                                            | 243 x 253 x 106<br>241 x 251 x 92 | JF254R                                                                            | 406 x 253 x 106<br>404 x 251 x 92 | JF214R                                                                            | 485 x 253 x 106<br>483 x 251 x 92 | JF224R                                                                              | 540 x 253 x 106<br>538 x 251 x 92 |

Instrument baskets

## Basket **without** feet




|                                                                                       |  | Outside dim.<br>Inside dim.      |  | Outside dim.<br>Inside dim.      |  | Outside dim.<br>Inside dim.      |  | Outside dim.<br>Inside dim.      |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------|----------------------------------|-------------------------------------------------------------------------------------|----------------------------------|---------------------------------------------------------------------------------------|----------------------------------|
|  24 | JG111R                                                                              | 243 x 253 x 24<br>241 x 251 x 22 | JG251R                                                                              | 406 x 253 x 24<br>404 x 251 x 22 | JG211R                                                                              | 485 x 253 x 24<br>483 x 251 x 22 | JG221R                                                                                | 540 x 253 x 24<br>538 x 251 x 22 |
|  44 | JG112R                                                                              | 243 x 253 x 44<br>241 x 251 x 42 | JG252R                                                                              | 406 x 253 x 44<br>404 x 251 x 42 | JG212R                                                                              | 485 x 253 x 44<br>483 x 251 x 42 | JG222R                                                                                | 540 x 253 x 44<br>538 x 251 x 42 |
|  64 | JG113R                                                                              | 243 x 253 x 64<br>241 x 251 x 62 | JG253R                                                                              | 406 x 253 x 64<br>404 x 251 x 62 | JG213R                                                                              | 485 x 253 x 64<br>483 x 251 x 62 | JG223R                                                                                | 540 x 253 x 64<br>538 x 251 x 62 |
|  94 | JG114R                                                                              | 243 x 253 x 94<br>241 x 251 x 92 | JG254R                                                                              | 406 x 253 x 94<br>404 x 251 x 92 | JG214R                                                                              | 485 x 253 x 94<br>483 x 251 x 92 | JG224R                                                                                | 540 x 253 x 94<br>538 x 251 x 92 |

All dimensions: Length x Width x Height (mm)




## Stenter frames



|  Standard<br>1/2 | Height |
|---------------------------------------------------------------------------------------------------|--------|
| JF707R                                                                                            | 110    |

Stenter frame with division  
the swabs and compresses.



|  Standard<br>1/2 | Height |
|---------------------------------------------------------------------------------------------------|--------|
| JF704R                                                                                            | 220    |

Stenter frame with division  
the swabs and compresses.

## Identification labels



|        |             |
|--------|-------------|
| JG785B | red         |
| JG786B | blue        |
| JG787B | green       |
| JG788B | gold        |
| JG790B | black       |
| JG792B | olive-brown |
| JG793B | violet      |
| JG794B | orange      |

Identification labels with lettering one line  
(character height 6 mm) or two lines (character  
height 5 mm), up to 13 characters per line.  
Dimensions: 56 x 18 mm  
for all Standard containers except JK299, JK309 and JK339

*Note: Process aluminum identification labels with  
pH-neutral cleaning agent.*



|        |             |
|--------|-------------|
| JG785U | red         |
| JG786U | blue        |
| JG787U | green       |
| JG788U | gold        |
| JG790U | black       |
| JG792U | olive-brown |
| JG793U | violet      |
| JG794U | orange      |

Identification labels without lettering

Specifications as above

*Note: Process aluminum identification labels with  
pH-neutral cleaning agent.*



Reusable filters



|          |              |
|----------|--------------|
| Standard | Outside dim. |
| JK090    | Ø 190        |

Aesculap reusable filter, round  
for use in up to 1,000 cleaning and  
sterilization cycles.  
*Order unit (PAK): Pack of 10*



|          |              |
|----------|--------------|
| Standard | Outside dim. |
| JK089    | 225 x 225    |

Aesculap reusable filter, square  
for legacy Standard containers, service life  
1,000 cleaning and sterilization cycles.  
*Order unit (PAK): Pack of 10*



|        |              |
|--------|--------------|
| Mini   | Outside dim. |
| Dental |              |
| JK091  | 95 x 215     |

Aesculap reusable filter  
for Mini and Dental containers with filter  
holder JK098, service life 1,000 cleaning and  
sterilization cycles.  
*Order unit (PAK): Pack of 2*



|        |              |
|--------|--------------|
| Dental | Outside dim. |
| JK092  | 171 x 231    |

Aesculap reusable filter  
for Dental containers without filter  
holder JK098, service life 1,000 cleaning  
and sterilization cycles.  
*Order unit (PAK): Pack of 2*

Marking pens



|       |  |
|-------|--|
| JK086 |  |
|-------|--|

Marking pen for Aesculap  
reusable filters  
*Order unit (PAK): Pack of 3*

All dimensions: Length x Width x Height (mm)