

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

## Bacterial/Viral Filter, with HME

Bacterial/Viral Filter with HME and CO<sub>2</sub>-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

**Shape (Ventilation Filter):** Round

**ISO style (patient end):** Male/Female

**Inner diameter (patient end):** 20/15 mm

[read more](#)

Article no.	Application area	Particle retention	ISO style (machine end)	Inner diameter (machine end)	Max. tidal volume (Breathing Filter)	PZN	minima quantity
4300906120	Child	99,99 %	Male	13 mm	900 ml	03656356	50
4300906020	Adult	99,999 %	Male/Female	22/13 mm	1.500 ml	03650000	25
4300905708	Adult	99,999 %	Male/Female	22/13 mm	1.500 ml	06783747	50

