










MHKMEDİKAL TEKSTİL

We, MHK Medikal Tekstil Sanayi Ve Ticaret Limited Sirketi, based in **Korkun Mah. Hidayet Cad. No 23/1 Oguzeli Gaziantep, Turkey**, confirm the technical specification of our products:

Lot nr.	Product name	Description	Picture
50	Sterile surgical gown L	1.size: L 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 132 cm 9. Standard EN 13795.	
51	Sterile surgical gown M	1.size: M 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 125 cm 9. Standard EN 13795.	
52	Reinforced sterile surgical gown L	1.size: L 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping	





MHK MEDİKAL TEKSTİL

		closure at the back by ties 8. Length not less than 132 cm 9. Standard EN 13795.	
53	Reinforced sterile surgical gown M	1.size: M 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 125 cm 9. Standard EN 13795.	
54	Reinforced sterile surgical gown S	1.size: S 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 117 cm 9. Standard EN 13795.	
55	Reinforced sterile surgical gown XL	1.size: XL 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 140 cm 9. Standard EN 13795.	
56	Reinforced sterile surgical gown XXL	1.size: XXL 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 150 cm 9. Standard EN 13795.	



MHK MEDİKAL


























TEKSTİL

57	Sterile surgical gown S	1.size: S 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 117 cm 9. Standard EN 13795.	
58	Sterile surgical gown XL	1.size: XL 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H2O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 140 cm 9. Standard EN 13795.	



MHK MEDİKAL

TEKSTİL

59	Sterile surgical gown XXL	1.size: XXL 2.Non-woven material 3 layers SMS resistance to liquid penetration within 50-80 cm H ₂ O according to EN test 20811: 1981 (test method of the European certification body) 3.with elastic cuffs, knitted 4.sterile 5.individually packaged 6.adjustment at the neck 7. overlapping closure at the back by ties 8. Length not less than 150 cm 9. T Standard EN 13795.																																														
82	Surgical set for labor (long-term) interventions, sterile	<p>Field for instrument table with Reinforced Area from SMS 150x190cm (+/- 10cm) - 1 pc; Field for Mayo table with reinforcement SMS 80x140 cm (+/- 10 cm) - 1 pc;</p> <p>Absorbent wipes - should not leave lint- 30x35 cm (+/- 10 cm) – 2 pcs; U-shaped field 200x310 cm (+/- 10 cm) - 1 pc</p> <p>Incision area - 10x110 cm (+/- 5 cm);</p> <p>Positioning tape - on the entire surface of the incision area, impregnated in the fabric, width > 5 cm;</p> <p>Absorbent area - across the entire incision surface; liquid penetration resistance min. 150 cm H₂O / cm²</p> <p>Tube fixators - 4 pcs. of SMS;</p> <p>SMS material over the entire surface of the field with liquid resistance > 50 cm H₂O / cm²; Upper part of field 150x270 cm (+/- 10 cm) - 1 pc</p> <p>Positioning tape - on the entire surface of the incision area, impregnated in the fabric, width> 5 cm;</p> <p>Absorbent area - across the entire incision surface; liquid penetration resistance min. 150 cm H₂O / cm²</p> <p>SMS material over the entire surface of the field with liquid resistance min. 150 cm H₂O/cm²</p> <p>Side section of field 95</p>	<table border="1"><thead><tr><th>DRAW</th><th>ITEM</th><th>SIZE</th><th>MATERIAL</th><th>PIECE</th></tr></thead><tbody><tr><td></td><td>Field for instrument table with Reinforced Area</td><td>150x190cm</td><td>PE+NW</td><td>1</td></tr><tr><td></td><td>Field for Mayo table with reinforcement</td><td>80x140 cm</td><td>PE+NW</td><td>1</td></tr><tr><td></td><td>Absorbent wipes</td><td>30x35 cm</td><td>Cellulose</td><td>2</td></tr><tr><td></td><td>U-shaped field</td><td>200x310 cm</td><td>SMS+Reinforcement</td><td>1</td></tr><tr><td></td><td>Tube fixators</td><td>st</td><td>110 GR</td><td>4</td></tr><tr><td></td><td>Upper part of field</td><td>150x270 cm</td><td>SMS + reinforcement</td><td>1</td></tr><tr><td></td><td>Side section of field</td><td>95 x 105 cm</td><td>SMS + reinforcement</td><td>2</td></tr><tr><td></td><td>Positioning tape</td><td>10x50 cm</td><td></td><td>1</td></tr></tbody></table>	DRAW	ITEM	SIZE	MATERIAL	PIECE		Field for instrument table with Reinforced Area	150x190cm	PE+NW	1		Field for Mayo table with reinforcement	80x140 cm	PE+NW	1		Absorbent wipes	30x35 cm	Cellulose	2		U-shaped field	200x310 cm	SMS+Reinforcement	1		Tube fixators	st	110 GR	4		Upper part of field	150x270 cm	SMS + reinforcement	1		Side section of field	95 x 105 cm	SMS + reinforcement	2		Positioning tape	10x50 cm		1
DRAW	ITEM	SIZE	MATERIAL	PIECE																																												
	Field for instrument table with Reinforced Area	150x190cm	PE+NW	1																																												
	Field for Mayo table with reinforcement	80x140 cm	PE+NW	1																																												
	Absorbent wipes	30x35 cm	Cellulose	2																																												
	U-shaped field	200x310 cm	SMS+Reinforcement	1																																												
	Tube fixators	st	110 GR	4																																												
	Upper part of field	150x270 cm	SMS + reinforcement	1																																												
	Side section of field	95 x 105 cm	SMS + reinforcement	2																																												
	Positioning tape	10x50 cm		1																																												



MHK MEDİKAL

TEKSTİL

		<p>* 105 cm (+/- 10 cm) - 2 pcs</p> <p>Positioning tape - on the entire surface of the incision area, impregnated in the fabric, width > 5 cm;</p> <p>Absorbent area - across the entire incision surface; liquid penetration resistance min. 150 cm H₂O / cm²</p> <p>SMS material over the entire surface of the field with liquid resistance min. 150 cm H₂O/cm² Positioning tape 10*50 (+/-10 cm) - 1 pc</p>																																				
83	Sterile standard surgical set	<p>Field for instrument table with Reinforced Area from SMS 150x190cm (+/- 10cm) - 1 pc; Field for Mayo table with reinforcement SMS 80x140 cm (+/- 10 cm) - 1 pc</p> <p>Absorbent wipes - should not leave lint- 30x35 cm (+/- 10 cm) – 2 pcs</p> <p>Bottom side field 190 x 195 cm (+/- 10 cm) - 1 pc; Positioning tape - on the entire surface of the incision area, impregnated in the fabric, width > 5 cm;</p> <p>Absorbent area - across the entire incision surface; liquid penetration resistance min. 150 cm H₂O / cm²</p> <p>SMS material over the entire surface of the field with liquid resistance > 50 cm H₂O / cm² Upper part of field 150x270 cm (+/- 10 cm) - 1 pc Positioning tape - on the entire surface of the incision area, impregnated in the fabric, width > 5 cm;</p> <p>Absorbent area - across the entire incision surface; liquid penetration resistance min. 150 cm H₂O / cm²</p> <p>SMMS material over the entire surface of the field with liquid resistance min. 150 cm H₂O/cm²</p> <p>Positioning tape 10*50 (+/-5cm) - 1 pc</p>	<table border="1"><thead><tr><th>DRAW</th><th>ITEM</th><th>SIZE</th><th>MATERIAL</th><th>PIECE</th></tr></thead><tbody><tr><td></td><td>Field for instrument table with Reinforced Area</td><td>150x190cm</td><td>PE+NW</td><td>1</td></tr><tr><td></td><td>Field for Mayo table with reinforcement</td><td>80x140 cm</td><td>PE+NW</td><td>1</td></tr><tr><td></td><td>Absorbent wipes</td><td>30x35 cm</td><td>Celulose</td><td>2</td></tr><tr><td></td><td>Bottom side field</td><td>190 x 195 cm</td><td>SMS+Reinforcement</td><td>1</td></tr><tr><td></td><td>Upper part of field</td><td>150x270 cm</td><td>SMS+Reinforcement</td><td>1</td></tr><tr><td></td><td>Positioning tape</td><td>10x50 cm</td><td></td><td>1</td></tr></tbody></table>	DRAW	ITEM	SIZE	MATERIAL	PIECE		Field for instrument table with Reinforced Area	150x190cm	PE+NW	1		Field for Mayo table with reinforcement	80x140 cm	PE+NW	1		Absorbent wipes	30x35 cm	Celulose	2		Bottom side field	190 x 195 cm	SMS+Reinforcement	1		Upper part of field	150x270 cm	SMS+Reinforcement	1		Positioning tape	10x50 cm		1
DRAW	ITEM	SIZE	MATERIAL	PIECE																																		
	Field for instrument table with Reinforced Area	150x190cm	PE+NW	1																																		
	Field for Mayo table with reinforcement	80x140 cm	PE+NW	1																																		
	Absorbent wipes	30x35 cm	Celulose	2																																		
	Bottom side field	190 x 195 cm	SMS+Reinforcement	1																																		
	Upper part of field	150x270 cm	SMS+Reinforcement	1																																		
	Positioning tape	10x50 cm		1																																		



MHK MEDİKAL TEKSTİL

87	Set for caesarean	material SMMS	<table border="1"><thead><tr><th>DRAW</th><th>ITEM</th><th>SIZE</th><th>MATERIAL</th><th>PIECE</th></tr></thead><tbody><tr><td></td><td>Sheet for the instrument table</td><td>150 x 200cm</td><td>FE-NOV</td><td>1</td></tr><tr><td></td><td>Medical sheet for Caesarean with collector pocket</td><td>200 x 300cm</td><td>Blended</td><td>1</td></tr><tr><td></td><td>Medical sheet for newborn</td><td>75 x 90cm</td><td>Spunlax</td><td>1</td></tr><tr><td></td><td>Umbilical clip</td><td>It</td><td></td><td>1</td></tr><tr><td></td><td>Surgical gown reinforced</td><td>SMS IS GR</td><td></td><td>2</td></tr><tr><td></td><td>Towel for hands</td><td>40 x 40cm</td><td></td><td>2</td></tr><tr><td></td><td>Medical sheet</td><td>100 x 100cm</td><td>SMS</td><td>1</td></tr></tbody></table>	DRAW	ITEM	SIZE	MATERIAL	PIECE		Sheet for the instrument table	150 x 200cm	FE-NOV	1		Medical sheet for Caesarean with collector pocket	200 x 300cm	Blended	1		Medical sheet for newborn	75 x 90cm	Spunlax	1		Umbilical clip	It		1		Surgical gown reinforced	SMS IS GR		2		Towel for hands	40 x 40cm		2		Medical sheet	100 x 100cm	SMS	1
		DRAW		ITEM	SIZE	MATERIAL	PIECE																																				
				Sheet for the instrument table	150 x 200cm	FE-NOV	1																																				
				Medical sheet for Caesarean with collector pocket	200 x 300cm	Blended	1																																				
				Medical sheet for newborn	75 x 90cm	Spunlax	1																																				
				Umbilical clip	It		1																																				
				Surgical gown reinforced	SMS IS GR		2																																				
	Towel for hands	40 x 40cm		2																																							
	Medical sheet	100 x 100cm	SMS	1																																							
1. Sheet for the instrument table ~ 150 x 200cm. - 1 piece																																											
2. Medical sheet for Caesarean with collector pocket ~ 200 x 300cm - 1 pc.																																											
3. Medical sheet for newborn ~ 75 x 90cm. 1 pc.																																											
4. Umbilical clip - 1 pc.																																											
5. Surgical gown (SMS material) reinforced size L (universal size) - 2pcs.																																											
6. Towel for hands 40 x 40cm - 2 pcs.																																											
7. Medical sheet 100 x 100cm - 1 pc.																																											

Place: MHK MEDİKAL TEKSTİL SAN. VE TIC. LTD. ŞTİ.
Date: 14.09.2021

Signed:



MHK MEDİKAL TEKSTİL
SAN. VE TIC. LTD. ŞTİ.

Körkün Mahallesi Hidayet Caddesi No:23 Oğuzeli
Tel.: 0342 341 93 66 - 477 10 03 / 342 341 93 67 GAZİANTEP
Tic.Sicil No: 40345 SUBURCU V.D. 620 056 1920