

# Dispogel®

Urethral Lubricant Gel  
with 2 % Lidocaine HCL

Dispogel® is a sterile lubricat gel containing local anesthetic (Lidocaine 2%) and antiseptic (Chlorhexidine Digluconate 0.05%). Dispogel® is available in 12.5 gr and 8,5 gr single use accordion bottle presented in a sterile package.

Dispogel® is used for instillation into the urethra prior to insertion of a catheter or other medical instruments. The insertion of urinary catheters by the helps of Dispogel® Urethral Lubricant Gel with Lidocaine HCL (12.5 gr and 8,5 gr) health professionals minimise associated trauma, pain and catheter-associated infections.

After instillation of gel - WAIT 5-10 minutes prior to insertion of a catheter or other medical instruments- time needed for achieving the lidocaine anesthetic effect.

Dispogel® is a Class III medical device.



12.5 gr



8,5 gr

8.5 gr	12.5 gr	STERIL R	CE 2292
GMDN Code	Product Description	Brand	Pack Size
37717	8.5 gr Accordion Bottle	Dispogel®	Pack of 100
37717	12.5 gr Accordion Bottle	Dispogel®	Pack of 100

Lotul 98

Method of Use  
For Male



Step 1



Step 2



Step 3



Step 4




Step 5



Step 6



Step 7



Step 8

ANAESTHETIC

with

ANTISEPTIC

with

STERILE LUBRICANT

Method of Use  
For Male



Step 1



Step 2



Step 3



Step 4



Step 5



Step 6



Step 7



Step 8

For Female



Step 1



Step 2



Step 3



Step 4



Step 5

Dispogel®



**Dispogel®**  
Urethral Lubricant Gel with  
Lidocaine HCL  
11 ml & 6 ml

**Dispogel® is a sterile lubricat gel containing local anesthetic (Lidocaine 2%) and antiseptic (Chlorhexidine Digluconate 0.05%). Dispogel is available in both 6m and 11 ml single use pre-filled syringes each presented in a sterile package.**

**Dispogel® is used for instillation into the urethra prior to insertion of a catheter or other medical instruments. The insertion of urinary catheters by the helps of Dispogel® Urethral Lubricant Gel with Lidocaine HCL ( 11 ml & 6 ml ) health professionals minimise associated trauma, pain and catheter-associated infections.**

**After instillation of gel - WAIT 5-10 minutes prior to insertion of a catheter or other medical instruments- time needed for achieving the lidocaine anesthetic effect.**

**Dispogel® is a Class III medical device.**



6ml	11ml	STERILE R	CE 2292
GMDN Code	Product Description	Brand	Pack Size
37717	6 ml Syringe	Dispogel®	Pack of 50
37717	11 ml Syringe	Dispogel®	Pack of 25

**Dispofarma İlaç Sanayi ve Ticaret Ltd. Şti.**  
İvedik Organize Sanayi Bölgesi, Ağaç İşleri Sitesi, 1333. Cad. No: 51  
Yenimahalle - ANKARA / TURKEY  
Tel: +90 (0) 312 394 45 52 Fax: +90 (0) 312 394 45 58  
info@dispofarma.com  
www.dispofarma.com.tr



ANAESTHETIC  
with  
ANTISEPTIC  
with  
STERILE LUBRICANT