

VOMIT BAG



Vomit Bag

A disposable and non-sterile device with a drawstring at the opening facilitates the collection and sampling of a patient's vomit for diagnostic purposes.

Available Sizes

Vomit Bag	DLM-8001
-----------	----------

Features

- There is an absorbent pad that has the capacity to absorb up to 450ml of liquid.
- A simple and easy to use drawstring for closure.

Intended Use


Vomit Bags are disposable and non-sterile devices that help prevent environmental contamination in patients while collecting and examining their stomach contents for diagnostic purposes. It is disposable and non-sterile.



Instructions for Use

1. Provide a vomit bag when needed.
2. Open the bag and place your mouth in the middle.
3. After vomiting, pull the strings on both sides of the bag to close it.
4. The vomit bag absorbs up to 450ml of liquid in a short time.

Technical Data



	Part	Qty	Material
1-	Bag	1	PE
2-	Absorbent pad	1	PP

Manufacturing Methods

- Assembly
- Packaging-Boxing

Specifications

- Manufactured under ISO 13485 Quality Management System standard.
- Latex free.
- Phthalate free.

User Information

Storage:

Storage condition only valid for storing in original packaging. Do not use the product if the original packaging is damaged.

Storage conditions: Keep away from direct sunlight
Keep dry

Packaging:

Packaging Unit: Packaging with nylon bag (20pcs)

Box Qty: 1000pcs

Inner Bag Qty: 20pcs

Regulatory Compliance:



Dilamed Tıbbi Ürünler A.Ş.

Yazıbaşı Mah. Fatma Tonkuş Cad. No:24 Torbalı/İzmir

+90 232 281 11 25

info@dilamedtibbi.com www.dilamedtibbi.com

sales@dilamedtibbi.com

