

SEJOY®





Product description

Blood glucose monitoring system	
Model Number	BG-201
Response Time:	5 Seconds
sample:	Fresh capillary whole blood
Blood size:	Minimum 1 microlitre
Measuring Range	40 to 600 mg/dl (2.2-33.3 mmol/L)
Average calculation:	last 7 14 & 21 days result averaging
Memory	360
LCD Size:	28 x 26 mm (L X W)
Power source:	one 3.0 V alkaline battery .size CR2032. included
Battery life:	About 1000 tests
Dimension:	8.02 x 5.5 x 1.63 cm (L x W x H)
Weight:	Approx. 35.8 grams

Test Strips

BS-101 / BS-102



1Microliter



5 seconds
measuring
time



Accurate



More sensitive



High cost
performance

- 25/50pcs in one bottle
- Required volume of blood : 1ul
- Code free optional
- 24 month expiry time
- Can still be used within 90 days after open box



Parts



Lancets

- 10/25/50pcs in one box
- size: 28g



Lancing Device

- The adjustable tip offer 5 depth levels of skin penetration



Control Solution

- 1*2.0ml (bottle)
- Store between 10-30°C (50-86°F)
- Must be used within 90 days after the first opening



FY0902 Stitch Cutters

All sizes are available in both a carbon or stainless steel sterile
Designed with a rounded bull nose tip, the Index family of Stitch Cutter Blades are ideal for cutting the wide range of modern sutures available today whilst minimising patient discomfort
The range consists of a short standard stitch cutter with a larger grip. The long handle enables a more aseptic technique allowing good visibility whilst working further from the wound site
Size: 65mm, 110mm
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag, 100pcs/box, 5000pcs/ctn



FY0903 Surgical Handle

Index stainless steel metal handles are supplied non-sterile and individually wrapped
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag, 10pcs/box, 500pcs/ctn



FY0904 Surgical Scalpels

Standard sheath guarded and blade protected sterile disposable scalpels offer convenience in a number of user environments: on the ward, in emergency or trauma units or even out in the field with paramedics
Material: Carbon Steel or Stainless Steel
Sterilization: Gamma
Packing: 10pcs/Box, 500pcs/ctn



FY0905 Steel Lancet

Low cost: simple, flat, steel lancet
Sample volume: excellent blood flow
High quality: hardened surgical steel
Minimum discomfort: modern asymmetric design
Sterile: sterilized by gamma radiation
Size: Adult 55mm
Packing: 5pcs of packaged, 200pcs/Box, 20000pcs/ctn



FY0906 Twist Lancets

Conventional twist-off cap design
Sterilized by gamma-radiation
Smooth tri-bevel needle tip for comfortable sampling experience
Versatile with most lancing devices
Size: 21G, 23G, 26G, 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



FY0907 Twist Lancets (Exclusive to Roche Diagnostics)

5 twist lancets (28G) are used with Roche Accu-Chek Softclix Lancing Device and the needle tip is siliconized for easy and smooth motion during lancing
Size: 28G, 30G
Packing: 100pcs/Box, 20000pcs/ctn



FY0910 Safety Lancet

Ergonomic design offers a more secure hold when to activate the device
Tri-bevel needle tip ensures the sharpness and accuracy
The spring buffer the puncture making the procedure a more comfortable experience
The complete range of needle sizes from 18G to 32G ensure the accurate
Blood samples for all popular meters
The needle retracts immediately after use
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0911 Pressure Activated Safety Lancets

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0912 Button Activated Safety Lancets

Needle tip is sealed with cap through injection moulding
Pre-loaded, button operated design
Single-use only, sterilised by gamma irradiation, cannot be reused
High quality, ultra sharp tri-bevel needle, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



FY0908 Lancing Device (Releasable) FY0909 Lancing Device

The user can install the lancet with the device easily. And the device will work with the lancet smoothly without any blockage. The lancet won't fall on its own when to fire the device. Retracting automatically after firing, the launcher is safe and reliable.