



PENTAX Medical Accessories Catalogue

01 Index

02 Advanced Therapeutics

- 02-1 ESD Accessories
- 02-2 TOP Endoscopic Hoods
- 02-3 C2 CryoBalloon™ Ablation System

03 Standard Therapeutic Accessories

- 03-1 Stella Line Single-Use Accessories

04 Endoscope Consumables

- 04-1 Distal Caps / Hoods
- 04-2 Valves
- 04-3 Inlet Seals
- 04-4 Balloons
- 04-5 O-Rings
- 04-6 Others
- 04-7 VERSA Consumables

05 Electrosurgery and Argon Plasma Coagulation (APC)

06 Hygiene Consumables

- 06-1 PROFILE Single-Use Cleaning Brushes
- 06-2 Cleaning Brushes
- 06-3 Reprocessing Accessories
- 06-4 VERSA Consumables

Distributors

EMEA Headquarter

Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

United Kingdom

PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical.uk@pentaxmedical.com

France

PENTAX France Life Care S.A.S.
116 quai de Bezons
B.P. 204
95106 Argenteuil Cedex
Tel: +33 1 / 3025 7575
Fax: +33 1 / 3025 7445
E-mail: clients.fr@pentaxmedical.com

Netherlands

PENTAX Nederland B.V.
Edisonring 6
6669NB Dodewaard
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: info.nl@pentaxmedical.com

Italy

PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.it@pentaxmedical.com

Spain

SIMMEDICA – Sistemas Integrales
de Medicina, S.A.
28830 San Fernando de Henares - Madrid
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

Russia

Moscow Representative office
of PENTAX Europe GmbH (Germany)
13, 4-y Lesnoy pereulok
125047 Moscow
Tel.: +7 495 114 52 31
Fax: +7 495 114 52 31

Turkey

PENTAX Turkey
Veko Giz Plaza, Meydan Sokak No:3/43
343396 Maslak – Istanbul
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL

Advanced Therapeutics

02

Advanced Therapeutics

table of contents

02-1 ESD Accessories

02-2 TOP Endoscopic Hoods

02-3 C2 CryoBalloon™ Ablation System



For time-efficient and safe procedures.

PENTAX Medical ESD accessories are innovative, easy to handle and highly controllable for any kind of dissection of the entire GI tract.

A full range of ESD accessories for a complete dissection.

PENTAX Medical accessories cover all steps of the ESD procedure and can be perfectly combined with the PENTAX Medical RetroView scope, with its unique 210° deflection ability and reduced bending radius of its distal end, for a safe and time-efficient dissection.

Marking, marginal incision, dissection and hemostasis with Splash M-Knife

Splash M-Knife is a unique multi-functional accessory allowing complete ESD procedures with a single device.

It has been developed and designed to achieve clear marking, better hemostasis and smoother and easier operation with fixed knife length and hooking mechanism. Splash M-Knife is designed to improve hemostasis ability with better bleeding control during a procedure without replacing the knife.




Marking




Resection



Dissection



Hemostasis



Hemostasis with the bipolar HemoStat WideCup

The HemoStat WideCup is used to accomplish effective initial hemostasis for gastrointestinal bleeding.

The bipolar forceps device is designed to support fast and efficient hemostasis without causing damage to the deeper tissue through thermal degeneration. The HemoStat WideCup is expected to aid prevention of re-bleeding, especially during ESD.*¹

^{*1} M Kataoka et al. "Clinical evaluation of emergency endoscopic hemostasis with bipolar forceps in non-variceal upper gastrointestinal bleeding" Digestive Endoscopy (2010) 22, 151-155



Dissection with Mucosectom

With its innovative design, Mucosectom overcomes the challenges of dissection procedures safely and efficiently.

The insulated tip can reduce the risk of perforation, while the extended knife length can reduce the amount of time the procedure takes.*² The knife on the side of the 360° rotatable tip and the optical marker are designed to assist the physician in holding the accessory in the direction of the lesion and dissecting it precisely.

^{*2} Hiromitsu Kanzaki MD1, Ryu Ishihara MD1, Takashi Ohta MD1, Kengo Nagai MD1, Fumi Matsui MD1, Takeshi Yamashina MD1, Masao Hanafusa MD1, Sachiko Yamamoto MD1, Noboru Hanaoka MD1, Yoji Takeuchi MD1, Koji Higashino MD1, Noriya Uedo MD1, Hiroyasu Iishi MD1 and Masaharu Tatsuta MD1 "Randomized Study of Two Endo-knives for Endoscopic Submucosal Dissection of Esophageal Cancer" The American Journal of Gastroenterology 2013 Aug;108(8):1293-8. doi: 10.1038/ajg.2013.161. Epub 2013 Jun 4

PENTAX Medical ESD – precise, efficient and safe.

Splash M-Knife: Multi-functional device for a complete ESD



- Disk in the knife helps to hook up the tissue for improved safety and ease of cutting
- Metal plate that closes the knife for clear marking and better hemostasis
- Water jet function that enables additional injection without changing the device

*3 Major bleeding needs to be stopped by a clip or hemostasis forceps such as HemoStat WideCup

INNOVATIVE
DESIGN




Mucosectom: For a safe and fast therapeutic approach



- Insulated tip and single-sided knife designed for a safer dissection
- Rotatable tip with optical marker supports a precisely directed dissection
- Extended knife length designed for a faster procedure*4

*4 Kawahara Y, Imagawa A, Fujiki S, Shiratori Y. "Novel procedure of endoscopic submucosal dissection (ESD) using a new device (mucosectome) for early esophageal cancer" *Gastrointest Endosc.* 2006;63:AB222

INNOVATIVE
DESIGN



Ksnare: All-in-one



- Enables hybrid techniques EMR/ESD especially in the colon
- Combines an adjustable knife on the tip designed for precise marking and incision, and a snare for resection
- The all-in-one accessory for supporting a complete lesion resection



Bipolar HemoStat WideCup: Easy control of bleeding

BIPOLAR

INNOVATIVE
DESIGN

Marking

Marginal
Incision

Submucosal
Dissection

Hemostasis

- Bipolar HemoStat WideCup reduces the impact of deeper tissue thermal degeneration^{*5}
- Supports prevention of re-bleeding^{*6}
- Based on its flexibility, it can be used for various applications such as EMR/ESD, gastrointestinal bleeding and papilla resection

^{*5} Compared with monopolar forceps

^{*6} M Kataoka et al. "Clinical evaluation of emergency endoscopic hemostasis with bipolar forceps in non-variceal upper gastrointestinal bleeding" *Digestive Endoscopy* (2010) 22, 151-155



	Model name	Item number	Working Length [mm]	Min. Channel [mm]	Max. Insertion Portion Width [mm]	Knife Length [mm]	Adjustable Knife Length [mm]	Snare/Cup Opening [mm]	Packaging	
									Single unit	Boxes of 5
Splash M-Knife	DN-D2718B	80387	1,800	2.8	2.7	Extended 2.0			X	
						Retracted 0.5			X	
	DN-D2722B	80396	2,200	2.8	2.7	Extended 2.0			X	
						Retracted 0.5			X	
Mucosectom	DP-D2518	47662	1,800	2.8	2.5	5.0			X	
	DP-D2622	47902	2,200	2.8	2.6	2.5			X	
Ksnare	DO-D2618	47905	1,800	2.8	2.6	6.0	2.0-6.0	22	X	
	DO-D2622	47908	2,200	2.8	2.6	6.0	2.0-6.0	22	X	
HemoStat WideCup ^{*7}	HS-D2618	80942	1,800	2.8	2.6			5.8		X
	HS-D2622	80940	2,200	2.8	2.6			5.8		X

^{*7} Equipped with active cord and plug

Consult your local representative for information on the products currently available in your region.

TOP ENDOSCOPIC HOOD

The slit and hole allows quick fluid outlet from the hood.



©2012 TOP Corporation


Slit and Hole

Fluid in the hood is discharged by capillarity and pressure difference.

How to use

1. Attachment direction


Attach the hood to the scope in the direction shown by the arrow. Avoid attaching upside down. The attachment side is thinner, with matching endoscope diameter indicated.



©2012 TOP Corporation

2. Securing with tape

Secure the loose attachment (e.g. unstable hood) with tape. Avoid covering the side hole with the tape.



©2012 TOP Corporation

Lubricant jelly may clog the side hole, preventing fluid outlet.

Straight type



Model name	Item number	Compliant endoscope O.D.	Drainage function	Packaging unit
SS	306030	8.6~9.8 mm	Off	5pcs./box
SL	306031	10.6~12.6 mm		

Slit and Hole type



Model name	Item number	Compliant endoscope O.D.	Drainage function	Packaging unit
SHLL-B	306032	12.6~14.0 mm	On	5pcs./box

Slit and Hole F(Focus) type



Model name	Item number	Compliant endoscope O.D.	Hood protrusion length	Drainage function	Packaging unit
SHF-010	306025	9.0~9.2 mm	3.5~7.0 mm	On	5pcs./box
SHF-025	306026	9.8~9.9 mm			
SHF-030	306027	9.9~10.5 mm			
SHF-035	306028	10.8 mm			
SHF-050	306033	11.7~12.2 mm			
SHF-055	306034	12.8 mm			
SHF-060	306029	13.2 mm			

Manufactured by



TOP Corporation

Ask your local sales representative about compatibility with the PENTAX Medical endoscopes.

In order to select the compatible TOP Hood for your endoscope, please consider that variation of the nominal scope size may occur after service and repair.

Some of the above products may not be available in all countries. Please consult your local representative for further information.

The C2 CryoBalloon™ Ablation System


For the Endoscopic Cryoablation of Barrett's Esophagus*

*Cleared in US, Europe and Canada



C2 CryoBalloon™ Ablation System has not been commercially released in all European countries. Please contact your local sales rep to learn more.




The C2 CryoBalloon™ Ablation System is used to destroy diseased/ unwanted tissue by applying extreme cold. A balloon catheter is positioned in the area of the diseased tissue and the compliant balloon is then inflated to contact the esophageal wall circumferentially. The diffuser, within the expanded balloon, is then guided for ablation by pointing the diffuser outlet to deliver nitrous oxide directly at the targeted tissue. The extreme cold eradicates the targeted precancerous or abnormal lesions by a cryogenic heat transfer process between the balloon and the esophageal lining delivering a controlled therapy ablating only the targeted lesions. The system allows physicians complete control of the procedure with direct visualization of the treatment area during the procedure, which lasts between 15 to 30 minutes.[1],[2]



Circumferential ablation of the gastroesophageal junction. Courtesy of Dr. Haidry, University College Hospital, London

- [1] Canto MI, Shaheen NJ, Almario JA, Voltaggio L, Montgomery E, Lightdale CJ. Multifocal Nitrous Oxide Cryoballoon Ablation with or without EMR for Treatment of Neoplastic Barrett's Esophagus. *Gastrointest Endosc*. 2018 Apr 4; pii: S0016-5107(18)30259-1. doi: 10.1016/j.gie.2018.03.024. PMID: 29626424
- [2] Künzli HT, Schölvinc DW, Meijer SL, Seldenrijk KA, Bergman JGHM, Weusten BLAM. Efficacy of the CryoBalloon Focal Ablation System for the Eradication of Dysplastic Barrett's Esophagus Islands. *Endoscopy*. 2017 Feb;49(2): 169-175. doi: 10.1055/s-0042-120117. Epub 2016 Dec 20. PMID: 27997963


How does the cryoablation work?



- [3] Baust, J.G; Gage, A.A; Bieklund Johansen, T.E and Baust, J.M. Mechanisms of Cryoablation: Clinical Consequences on Malignant Tumors; *Cryobiology*. 2014 Feb. 68(1): 1-11

The new C2 CryoBalloon™ Ablation System offers several key advantages including the ability to ablate a wide range of patients, enhanced visualization and control as the ablation catheter is deployed through the working channel of the endoscope. “With the C2 CryoBalloon™ Ablation System, you can promote your practice as a more comprehensive Barrett’s center with the latest and best treatment options,” said Harshit S. Khara, MD, FACG, FASGE, Clinical Associate Professor of Medicine, Geisinger Medical Center.

Ablation extending to the submucosa¹



Combines deep ablation with a favorable safety profile²

1. Louie B, et al, Evaluation of a Novel Swipe Cryoballoon Ablation System in Bench, Porcine, and Human Esophagus Models. *Gastrointestinal Endoscopy*, Volume 85, No. 5S, AB513-AB514: 2017.
2. Schölvinck DW, et al, Balloon-Based Esophageal Cryoablation with a Novel Focal Ablation Device: Dose-Finding and Safety in Porcine and Human Models. *Diseases of the Esophagus*, 2017 Nov 1;30(11):1-8.

Effective Endoscopic Treatment For BE-Associated Neoplasia³

95%

Complete Eradication of Dysplasia


88%

Complete Eradication of Intestinal Metaplasia

3. Canto M, et al. Multifocal Nitrous Oxide Cryoballoon Ablation With or Without EMR for Treatment of Neoplastic Barrett’s Esophagus. *Gastrointest Endosc*. 2018 2018; 88: 438-446. doi: 10.1016/j.gie.2018.03.024

Less Pain⁴

PATIENTS REPORTING SEVERE PAIN AT DAY 2



4. Van Munster, et al, Post-Procedural Pain After Endoscopic Ablation Therapy of Barrett’s Esophagus: Is Cryo Ablation Less Painful Compared with Radiofrequency Ablation?, *Gastrointestinal Endoscopy*. 2017 May; Volume 85, Issue 5, Supplement, Page AB568.

All in a Smart and Compact Design

C2 CryoBalloon™ Ablation System Provides:

- Targeted access delivered through the working channel of a therapeutic endoscope
- An intuitive interface
- Smart, compact design that eliminates the need for large tanks or capital equipment



C2 CryoBalloon™ Ablation System portfolio of products

Multiple CryoBalloon shapes and targeting options designed for various anatomic structures

Pear shaped catheter



Model name	Article number	Pack size	Description
FG-1024	300617	1	Pear shaped catheter designed for Ablation of GEJ

Standard shaped catheter



Model name	Article number	Pack size	Description
FG-1028	300613	1	Standard shaped catheter for tubular esophagus

C2 CryoBalloon™ Ablation System controller and cartridges



Model name	Article number	Pack size	Description
FG-1017	300610	1	C2 CryoBalloon™ Ablation System Controller
FG-1013	300612	5	C2 CryoBalloon™ Ablation System Cartridge

C2 CryoBalloon™ Ablation System Foot Pedal



Model name	Article number	Pack size	Description
FG-1018	300611	1	C2 CryoBalloon™ Ablation System Foot Pedal, enables rotational & axial targeting

Manufacturers

Splash M-Knife
Mucosectom
Ksnare
HemoStat WideCup

 **HOYA Corporation**
6-10-1, Nishi-shinjuku
Shinjuku-ku, Tokyo
160-0023 Japan

TOP Endoscopic Hood

 **TOP Corporation**
19-10, Senju Nakai-cho
Adachi-ku, Tokyo
120-0035 Japan

C2 CryoBaloon™
Ablation
System

 **PENTAX Medical**
A Division of PENTAX of America, Inc.
303, Convention Way, Suite 1
Redwood City, CA 94063
USA

 PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany



0123

TÜV SÜD CE 0123 Medical Device Class: IIb

 MPS Medical Product Service GmbH
Borngasse 20
35619 Braunfels, Germany



0123

TÜV SÜD CE 0123 Medical Device Class: IIa

 PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany



0197

TÜV Rheinland CE 0197 Medical Device Class: IIa

Distributor

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL

Standard Therapeutic Accessories

03

Standard Therapeutic Accessories

table of contents

03-1 Stella Line Single-Use Accessories

PENTAX Medical

Product line of state-of-the-art accessories for diagnosis and treatment in endoscopy.



Our commitment to quality

Safe and successful endoscopy depends on the use of a finely tuned endoscopic system of components. This is the underlying principle behind PENTAX Medical's commitment to a product line of therapeutic accessories with the highest possible level of quality. We invite you to examine our carefully selected and meticulously crafted therapeutic accessories for state-of-the-art diagnostic and treatment applications. All PENTAX Medical accessories have been tested in conjunction with PENTAX Medical endoscopes for functionality and compatibility.

How to select the correct accessories.

Before you select new accessories, please consider the following:

- The compatibility of your endoscope with the desired diagnostic/treatment application.
- The diameter of the instrument channel and working length of your endoscope.
- The compatibility of the accessories with your endoscope to ensure the highest standards of endoscopic therapy.

Single-use

StellaClip Hemoclips

Usage

For endoscopic hemostasis and clamping the wound in the GI tract.

Features

- Synchronous rotation supports accurate positioning.
- Repositionable clips.
- The handle's security lock system helps prevent unexpected pain for the patient.
- MR conditional.
- Packaging: boxes of 5.



Description	Model name	Item number	Jaw outer diameter [mm]	Working length [mm]	Maximum opening of clip [mm]	Compatible endoscopes channel diameter* [mm]
Uncoated	VDK-CD-11-230-O	500063	2.5 mm	2300	11	≥ Ø 2.8
	VDK-CD-11-230-P	500064	2.5 mm	2300	11	≥ Ø 2.8
Coated	VDK-CD-13-230-P	500012	2.5 mm	2300	13	≥ Ø 2.8
	VDK-CD-16-230-P	500065	2.5 mm	2300	16	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use

StellaNeedle Injection needles

Usage

Endoscopic injection treatment of esophageal and gastric varices.
Endoscopic injection of submucosa in the GI tract.

Features

- Classic handle designed to easily control the needle movement.
- 3-bevel tip designed for precise injection and prevents lumen from clogging.
- Needle shield designed to prevent damage to endoscope working channel.
- Packaging: boxes of 10.



Needle size	Model name	Item number	Outer diameter [mm]	Working length [mm]	Needle length [mm]	Compatible endoscopes channel diameter* [mm]
19G	VDK-IN-23-230-1904-A	500016	2.3	2300	4	≥ Ø 2.8
	VDK-IN-18-120-2205-A	500013	1.8	1200	5	≥ Ø 2.0
22G	VDK-IN-23-180-2204-A	500068	2.3	1800	4	≥ Ø 2.8
	VDK-IN-23-230-2205-A	500069	2.3	2300	5	≥ Ø 2.8
25G	VDK-IN-23-230-2504-A	500017	2.3	2300	4	≥ Ø 2.8
	VDK-IN-23-230-2506-A	500019	2.3	2300	6	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use



StellaLoop Polyp snares

Usage

For excising polyps and other redundant tissues in the GI tract by high-frequency electricity in combination with an endoscope.

Features

- Synchronous rotation supports accurate positioning (rotatable snares).
- Made of medical stainless steel with good conductivity.
- Designed for safe and precise removal with effective coagulation.
- Braided snare loop.
- Packaging: boxes of 10.



Description	Model name	Item number	Outer diameter [mm]	Working length [mm]	Width of snare loop [mm]	Compatible endoscopes channel diameter* [mm]
Oval	VDK-SD-23-230-15-A1	500062	2.3	2300	15	≥ Ø 2.8
	VDK-SD-23-230-25-A1	500022	2.3	2300	25	≥ Ø 2.8
	VDK-SD-23-230-35-A1	500024	2.3	2300	35	≥ Ø 2.8
	VDK-SD-23-230-50-A1	500025	2.3	2300	50	≥ Ø 2.8
Oval rotatable	VDK-SD-23-230-20-A1R	500026	2.3	2300	20	≥ Ø 2.8
	VDK-SD-23-230-30-A1R	500027	2.3	2300	30	≥ Ø 2.8
Hexagonal	VDK-SD-18-160-20-C1	500028	1.8	1600	20	≥ Ø 2.0


*Figures are based on the validation tests done by PENTAX Medical

Single-use

StellaCold Cold snares

Features

- To be used without electricity.
- Packaging: boxes of 10.



Description	Model name	Item number	Outer diameter [mm]	Working length [mm]	Width of snare loop [mm]	Compatible endoscopes channel diameter* [mm]
Rhombic	VDK-CS-23-230-10-C3	500030	2.3	2300	10	≥ Ø 2.8
	VDK-CS-23-230-15-C3	500031	2.3	2300	15	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Usage

Tissue sampling in GI and respiratory tracts.

Features

- Special cup sharpness designed for easy sampling.
- Wide cup opening enables massive tissue sampling.
- 4-bar linkage connecting structure of the jaws designed for good rigidity and accurate biopsy.*
- Packaging: boxes of 10



Model name	Item number	Description				Jaw diameter [mm]	Working length [mm]	Compatible endoscopes channel diameter** [mm]
		Alli-gator-tooth jaws	Spike	Plastic screen printing scale	Coated			
VDK-FB-18-105-O-P-A1	500048	X			X	1.8	1050	≥ Ø 2.0
VDK-FB-23-230-W-P1-A1	500053	X	X	X	X	2.3	2300	≥ Ø 2.8
VDK-FB-23-230-O-P1-A1	500052	X		X	X	2.3	2300	≥ Ø 2.8
VDK-FB-18-105-O-P-B1	500039				X	1.8	1050	≥ Ø 2.0
VDK-FB-18-105-W-P-B1	500040		X		X	1.8	1050	≥ Ø 2.0
VDK-FB-18-160-O-P-B1	500041				X	1.8	1600	≥ Ø 2.0
VDK-FB-23-180-O-P1-B1	500072			X	X	2.3	1800	≥ Ø 2.8
VDK-FB-23-230-O-P1-B1	500045			X	X	2.3	2300	≥ Ø 2.8
VDK-FB-23-230-W-P1-B1	500046		X	X	X	2.3	2300	≥ Ø 2.8

*Applicable for the models with 2.3 mm Jaw diameter
 **Figures are based on the validation tests done by PENTAX Medical

Single-use

Spray catheters

Usage

Used in conjunction with an endoscope to spray the gastrointestinal mucosa.

Features

- Two types of spray: straight (A type) and dispersive (B type: >120°) to meet various clinical requirements.
- Uniform spray application.
- Packaging: boxes of 10.



Model name	Item number	Spray type	Outer diameter [mm]	Working length [mm]	Compatible endoscopes channel diameter* [mm]
VDK-ST-23-230-A	500038	straight spray	2.3	2300	≥ Ø 2.8
VDK-ST-23-230-B	500070	dispersive spray	2.3	2300	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use

StellaNet

Grasping forceps (net type)

Usage

Facilitates the removal of foreign bodies from the GI tract.

Features

- Packaging: boxes of 10.



Model name	Item number	Net size [mm]	Outer diameter [mm]	Working length [mm]	Compatible endoscopes channel diameter* [mm]
VDK-FG-25-230-C3	500056	30*60 mm	2.5	2300	≥ Ø 2.8
VDK-FG-25-230-C4	500057	40*75 mm	2.5	2300	≥ Ø 2.8
VDK-FG-25-230-C5	500058	40*65 mm	2.5	2300	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use**Grasping forceps (jaw type)****Usage**

Facilitates the removal of foreign bodies from the GI tract.

Features

- Variation of tooth design to facilitate the removal of a variety of foreign bodies.
- Packaging: boxes of 10.



Description	Model name	Item number	Jaw diameter [mm]	Working length [mm]	Compatible endoscopes channel diameter* [mm]
Alligator-teeth shape	VDK-FG-23-230-A1	500066	2.3	2300	≥ Ø 2.8
Rat-teeth shape	VDK-FG-23-230-A2	500067	2.3	2300	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use**Polyp traps****Usage**

Allows for the quick and accurate retrieval of endoscopic tissue samples.

Features

- 5 numbered chambers designed to facilitate accurate identification of specimens.
- Packaging: boxes of 10.



Model name	Item number	Description
VDK-PT-05	500061	5 chambers

Manufacturer

**Jiangsu Vedkang Medical
Science & Technology Co., Ltd**

No. 52, Guoxiang Road, Wujin
Economic Development Zone,
Changzhou, Jiangsu, China



Prolinx GmbH

Brehmstr. 56
40239, Düsseldorf, Germany



0197

TÜV Rheinland CE 0197 Medical Devices Class: IIa, IIb & Is

Distributor

PENTAX Europe GmbH

Julius-Vosseler-Straße 104
22527 Hamburg, Germany

Tel.: +49 40 / 5 61 92 - 0

Fax: +49 40 / 5 60 42 13

E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL

Endoscope Consumables

04

Endoscope Consumables

table of contents

04-1	Distal Caps/Hoods
04-2	Valves
04-3	Inlet Seals
04-4	Balloons
04-5	O-Rings
04-6	Others
04-7	VERSA Consumables

PENTAX Medical

Accessories Programme for cleaning and caring for your PENTAX Medical endoscopes

Original is better

Finely tuned system components are essential for safe and successful endoscopy. That is why PENTAX Medical original accessories aim to provide the highest possible quality in every detail. The following pages describe our carefully compiled accessory programme to ensure safe and effective cleaning of your PENTAX Medical endoscope.



How to choose the correct accessories:

Before you decide on new accessories, please ensure that your PENTAX Medical endoscope is compatible with the desired accessory.

- The illustration below is intended to help you to identify the relevant areas of your endoscope quickly and simply. Our PENTAX Medical specialists will be pleased to advise you if you have any additional questions.
- All PENTAX Medical original accessories in this catalogue have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

Orientation

In order to make your choice of PENTAX Medical accessories as simple and reliable as possible, we have identified fundamental product features with simple-to-understand symbols displayed next to the relevant technical data.



Autoclavable



For single
use only



OE-A63, OE-A65 Single-Use, Sterile Disposable Elevator Cap (DEC™)



Model name OE-A63	Description Single-Use sterile disposable elevator cap
-----------------------------	--

Item number 80924	Pack size 10
-----------------------------	------------------------

Compatible Endoscopes	ED34-i10T2
--------------------------	------------

Model name OE-A65	Description Single-Use sterile disposable elevator cap
-----------------------------	--

Item number 90031	Pack size 10
-----------------------------	------------------------

Compatible Endoscopes	ED32-i10
--------------------------	----------

OE-A55 Distal Rubber Tip (C.A.P.™)



Model name OE-A55	Description Distal rubber tip
-----------------------------	---

Item number 83352	Pack size 20
-----------------------------	------------------------

Compatible Endoscopes	ED34-i10T
--------------------------	-----------



OE-A58 Distal Rubber Hood



Model name OE-A58	Description Distal rubber hood
-----------------------------	--

Item number 80314	
-----------------------------	--

Compatible Endoscopes	EG-2990Zi
--------------------------	-----------



OE-A59 Distal Rubber Hood



Model name OE-A59	Description Distal rubber hood
-----------------------------	--

Item number 80315	
-----------------------------	--

Compatible Endoscopes	EC-3890LZi	EC-3890FZi	EC-3890MZi
--------------------------	------------	------------	------------





OF-B205 Single-Use Suction Control Valve



Model name
OF-B205

Description

Suction control valve for the J10 series

Item number
84587

Pack size
10

Compatible
Endoscopes

EB11-J10 EB19-J10 VNL15-J10

EB15-J10 EB19-J10U

OF-B179 Suction Control Valve



Model name
OF-B179

Description

Suction control valve (red)

Item number
86805

Compatible
Endoscopes

EB-1170K EB-1575K EB-1972K ECY1570K

EB-1570AK EB-1970AK EB-1975K ECY1575K

EB-1570K EB-1970K EB-1990i VNL1570STK

EB-1572K EB-1970TK EB-1970UK

OF-B163 Irrigation Cessation Valve



Model name
OF-B163

Description

Irrigation cessation valve (blue)

Item number
86488

Compatible
Endoscopes

ECY1570K ECY1575K

OF-B140 Suction Control Valve



Model name
OF-B140

Description

Suction control valve (red)

Item number
82688

Compatible
Endoscopes

FB-15BS FCY-15RBS FI-10BS FI-16RBS

FB-15RBS FHY-10RBS FI-10RBS FNL-15RP3

FB-18BS FHY-15RBS FI-13BS FUR-9RBS

FB-18RBS FI-9BS FI-13RBS

FCP-9P FI-9RBS FI-16BS

OF-B53 Suction Control Valve



Model name OF-B53	Description Suction control valve			
Item number 83630				
Compatible Endoscopes	FB-15BS	FCY-15RBS	FI-10BS	FI-16RBS
	FB-15RBS	FHY-10RBS	FI-10RBS	FNL-15RP3
	FB-18BS	FHY-15RBS	FI-13BS	FUR-9RBS
	FB-18RBS	FI-9BS	FI-13RBS	
	FCP-9P	FI-9RBS	FI-16BS	

OF-B136 Suction Control Valve



Model name OF-B136	Description Suction control valve (red)		
Item number 82673			
Compatible Endoscopes	FB-8V	FB-15V	FB-19TV
	FB-10V	FB-18V	FCN-15X




OF-B137 Irrigation Cessation Valve



Model name OF-B137	Description Irrigation cessation valve (blue)
Item number 82672	
Compatible Endoscopes	FCN-15X



OF-B147 Irrigation activation valve



Model name OF-B147	Description Irrigation activation valve
Item number 85515	
Compatible Endoscopes	FCN-15X





OF-B120 Suction Control Valve



Model name
OF-B120

Description
Suction valve (red)

Item number
81724

Compatible
Endoscopes

EC-2990Li	EC-3890Mi2	EC38-i10L	EG-290Kp
EC-2990Fi	EC-3890LK	EC38-i10F	EG-2990i
EC-2990Mi	EC-3890FK	EC38-i10M	EG-2990K
EC-3490Li	EC-3890MK	EC38-i10F2	EG-2990Zi
EC-3490Fi	EC-3890FK2	EC38-i10M2	EG-3490K
EC-3490Mi	EC-3890MK2	EC38-i10NL	EG-3890TK
EC-3490LK	EC-3890LZi	EC38-i10NF	EG16-K10
EC-3490FK	EC-3890FZi	EC38-i10cL	EG17-J10
EC-3490MK	EC-3890MZi	EC38-i10cF	EG27-i10
EC-3490TLi	EC-3890TLK	EC38-i10cM	EG29-i10
EC-3490TFi	EC-3890TFK	EC38-i10cF2	EG29-i10c
EC-3490TMi	EC34-i10L	ED-3490TK	EG34-i10
EC-380LKp	EC34-i10F	ED32-i10	FC-38LV
EC-380FKp	EC34-i10M	ED34-i10T	FC-38FV
EC-380FK2p	EC34-i10TL	ED34-i10T2	FD-34V2
EC-380MK2p	EC34-i10TF	EE17-J10	FG-16V
EC-3890Li	EC34-i10TM	EG-1690K	FG-24V
EC-3890Fi	EC34-i10cL	EG-2490K	FG-29V
EC-3890Mi	EC34-i10cF	EG-2790i	FS-34V
EC-3890Fi2	EC34-i10cM	EG-2790K	VSB-2990i

OF-B188 Air/Water Feeding Valve



Model name
OF-B188

Description
Air/Water valve (purple)

Item number
87143

Compatible
Endoscopes

EC-2990Li	EC-3890LK	EC34-i10cF	EE17-J10
EC-2990Fi	EC-3890FK	EC34-i10cM	EG-1690K
EC-2990Mi	EC-3890MK	EC38-i10L	EG-2490K
EC-3490Li	EC-3890FK2	EC38-i10F	EG-2790i
EC-3490Fi	EC-3890MK2	EC38-i10M	EG-2790K
EC-3490Mi	EC-3890LZi	EC38-i10F2	EG-2990i
EC-3490LK	EC-3890FZi	EC38-i10M2	EG-2990K
EC-3490FK	EC-3890MZi	EC38-i10NL	EG-2990Zi
EC-3490MK	EC-3890TLK	EC38-i10NF	EG-3490K
EC-3490TLi	EC-3890TFK	EC38-i10cL	EG-3890TK
EC-3490TFi	EC34-i10L	EC38-i10cF	EG16-K10
EC-3490TMi	EC34-i10F	EC38-i10cM	EG17-J10
EC-3890Li	EC34-i10M	EC38-i10cF2	EG27-i10
EC-3890Fi	EC34-i10TL	ED-3490TK	EG29-i10
EC-3890Mi	EC34-i10TF	ED32-i10	EG29-i10c
EC-3890Fi2	EC34-i10TM	ED34-i10T	EG34-i10
EC-3890Mi2	EC34-i10cL	ED34-i10T2	VSB-2990i



OF-B194 Gas/Water Feeding Valve



Model name
OF-B194

Description
Gas/water feeding valve (purple)

Item number
85133

Compatible
Endoscopes

EC-2990Li	EC-3890LK	EC34-i10cF	EE17-J10
EC-2990Fi	EC-3890FK	EC34-i10cM	EG-1690K
EC-2990Mi	EC-3890MK	EC38-i10L	EG-2490K
EC-3490Li	EC-3890FK2	EC38-i10F	EG-2790i
EC-3490Fi	EC-3890MK2	EC38-i10M	EG-2790K
EC-3490Mi	EC-3890LZi	EC38-i10F2	EG-2990i
EC-3490LK	EC-3890FZi	EC38-i10M2	EG-2990K
EC-3490FK	EC-3890MZi	EC38-i10NL	EG-2990Zi
EC-3490MK	EC-3890TLK	EC38-i10NF	EG-3490K
EC-3490TLi	EC-3890TFK	EC38-i10cL	EG-3890TK
EC-3490TFi	EC34-i10L	EC38-i10cF	EG16-K10
EC-3490TMi	EC34-i10F	EC38-i10cM	EG17-J10
EC-3890Li	EC34-i10M	EC38-i10cF2	EG27-i10
EC-3890Fi	EC34-i10TL	ED-3490TK	EG29-i10
EC-3890Mi	EC34-i10TF	ED32-i10	EG29-i10c
EC-3890Fi2	EC34-i10TM	ED34-i10T	EG34-i10
EC-3890Mi2	EC34-i10cL	ED34-i10T2	VSF-2990i

OF-B121 Air/Water Feeding Valve



Model name
OF-B121

Description
Air/water valve (blue)

Item number
81725

Compatible
Endoscopes

EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
EC-380FKp	EG-290Kp	FD-34V2	FG-29V
EC-380FK2p	FC-38LV	FG-16V	FG-34V

OF-B183 Air/Water Feeding Valve



Model name
OF-B183

Description
Air/water valve (green)

Item number
87999

Compatible
Endoscopes

ES-3870K



OF-B177 Suction Control Valve



Model name OF-B177	Description Suction valve (yellow)	
Item number 86803		
Compatible Endoscopes	EE-1580K	ES-3870K



OF-B105 Irrigation Cessation Valve



Model name OF-B105	Description Irrigation control valve		
Item number 81564			
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	




OF-B151 Irrigation Cessation Valve



Model name OF-B151	Description Irrigation (deactivation) valve		
Item number 86218			
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	



OF-B152 Water Feeding Valve



Model name OF-B152	Description Irrigation (activation) valve (blue)		
Item number 86217			
Compatible Endoscopes	FCP-9P	FHY-10RBS	FUR-9RBS
	FCY-15RBS	FHY-15RBS	

OF-B171 Suction Control Valve



Model name
OF-B171

Description
Suction control valve (yellow)

Item number
87068

Compatible
Endoscopes

EG-3270UK

EG-3670URK

EG-3870UTK

OF-B172 Air/Water/Balloon Feeding Valve



Model name
OF-B172

Description
Air/water/balloon feeding valve (green)

Item number
87069

Compatible
Endoscopes

EG-3270UK

EG-3670URK

EG-3870UTK

OF-B206 Suction/Balloon Suction Control Valve



Model name
OF-B206

Description
Suction/balloon suction control valve (red)

Item number
84833

Compatible
Endoscopes

EG34-J10U

EG36-J10UR

EG38-J10UT

OF-B207 Air/Water/Balloon Feeding Valve



Model name
OF-B207

Description
Air/water/balloon feeding valve (purple)

Item number
84834

Compatible
Endoscopes

EG34-J10U

EG36-J10UR

EG38-J10UT



OF-B215 PENTAX Medical PROfILE Single-Use Endoscope Inlet Seal



Model name
OF-B215

Description
Single-Use rubber inlet seal

Item number
300430

Pack size
50

Compatible
Endoscopes

EC-2990Li	EC34-i10TL	EG27-i10	EG-3870UTK
EC-2990Fi	EC34-i10TF	EG29-i10	EG34-J10U
EC-2990Mi	EC34-i10TM	EG29-i10c	EG36-J10UR
EC-3490Li	EC34-i10cL	EG34-i10	EG38-J10UT
EC-3490Fi	EC34-i10cF	FC-38LV	ES-3870K
EC-3490Mi	EC34-i10cM	FC-38FV	FB-8V
EC-3490LK	EC38-i10L	FD-34V2	FB-10V
EC-3490FK	EC38-i10F	FG-16V	FB-15V
EC-3490MK	EC38-i10M	FG-24V	FB-18V
EC-3490TLi	EC38-i10F2	FG-29V	FB-19TV
EC-3490TFi	EC38-i10M2	FS-34V	FB-15BS
EC-3490TMi	EC38-i10NL	VSB-2990i	FB-15RBS
EC-380LKp	EC38-i10NF	EB-1170K	FB-18BS
EC-380FKp	EC38-i10cL	EB-1570AK	FB-18RBS
EC-380FK2p	EC38-i10cF	EB-1570K	FCN-15X
EC-380MK2p	EC38-i10cM	EB-1572K	FCP-9P
EC-3890Li	EC38-i10cF2	EB-1575K	FCY-15RBS
EC-3890Fi	ED-3490TK	EB-1970AK	FHY-10RBS
EC-3890Mi	ED32-i10	EB-1970K	FHY-15RBS
EC-3890Fi2	ED34-i10T	EB-1970TK	FI-9BS
EC-3890Mi2	ED34-i10T2	EB-1972K	FI-9RBS
EC-3890LK	EE17-J10	EB-1975K	FI-10BS
EC-3890FK	EG-1690K	EB-1990i	FI-10RBS
EC-3890MK	EG-2490K	EB11-J10	FI-13BS
EC-3890FK2	EG-2790i	EB15-J10	FI-13RBS
EC-3890MK2	EG-2790K	EB19-J10	FI-16BS
EC-3890LZi	EG-290Kp	EB-1970UK	FI-16RBS
EC-3890FZi	EG-2990i	EB19-J10U	FNL-15RP3
EC-3890MZi	EG-2990K	ECY-1570K	FUR-9RBS
EC-3890TLK	EG-2990Zi	ECY-1575K	VNL-1570STK
EC-3890TFK	EG-3490K	EE-1580K	VNL15-J10
EC34-i10L	EG-3890TK	EG-3270UK	
EC34-i10F	EG16-K10	EG-3670URK	
EC34-i10M	EG17-J10		

OF-B190 Rubber Inlet Seal



Model name
OF-B190

Description
Rubber inlet seal

Item number
86665

Pack size
10

Compatible
Endoscopes

EC-2990Li	EC34-i10TL	EG27-i10	EG-3870UTK
EC-2990Fi	EC34-i10TF	EG29-i10	EG34-J10U
EC-2990Mi	EC34-i10TM	EG29-i10c	EG36-J10UR
EC-3490Li	EC34-i10cL	EG34-i10	EG38-J10UT
EC-3490Fi	EC34-i10cF	FC-38LV	ES-3870K
EC-3490Mi	EC34-i10cM	FC-38FV	FB-8V
EC-3490LK	EC38-i10L	FD-34V2	FB-10V
EC-3490FK	EC38-i10F	FG-16V	FB-15V
EC-3490MK	EC38-i10M	FG-24V	FB-18V
EC-3490TLi	EC38-i10F2	FG-29V	FB-19TV
EC-3490TFi	EC38-i10M2	FS-34V	FB-15BS
EC-3490TMi	EC38-i10NL	VSB-2990i	FB-15RBS
EC-380LKp	EC38-i10NF	EB-1170K	FB-18BS
EC-380FKp	EC38-i10cL	EB-1570AK	FB-18RBS
EC-380FK2p	EC38-i10cF	EB-1570K	FCN-15X
EC-380MK2p	EC38-i10cM	EB-1572K	FCP-9P
EC-3890Li	EC38-i10cF2	EB-1575K	FCY15RBS
EC-3890Fi	ED-3490TK	EB-1970AK	FHY10RBS
EC-3890Mi	ED32-i10	EB-1970K	FHY15RBS
EC-3890Fi2	ED34-i10T	EB-1970TK	FI-9BS
EC-3890Mi2	ED34-i10T2	EB-1972K	FI-9RBS
EC-3890LK	EE17-J10	EB-1975K	FI-10BS
EC-3890FK	EG-1690K	EB-1990i	FI-10RBS
EC-3890MK	EG-2490K	EB11-J10	FI-13BS
EC-3890FK2	EG-2790i	EB15-J10	FI-13RBS
EC-3890MK2	EG-2790K	EB19-J10	FI-16BS
EC-3890LZi	EG-290Kp	EB-1970UK	FI-16RBS
EC-3890FZi	EG-2990i	EB19-J10U	FNL15RP3
EC-3890MZi	EG-2990K	ECY1570K	FUR-9RBS
EC-3890TLK	EG-2990Zi	ECY1575K	VNL1570STK
EC-3890TFK	EG-3490K	EE-1580K	VNL15-J10
EC34-i10L	EG-3890TK	EG-3270UK	
EC34-i10F	EG16-K10	EG-3670URK	
EC34-i10M	EG17-J10		



OE-A51 Balloons for PENTAX Ultrasound Endoscopes



Model name OE-A51	Description Balloons for endoscopic ultrasound application	
Item number 87070	Pack size 20	
Compatible Endoscopes	EG-3670URK	EG36-J10UR

OE-A52 Balloons for PENTAX Ultrasound Endoscopes



Model name OE-A52	Description Balloons for endoscopic ultrasound application	
Item number 86758	Pack size 20	
Compatible Endoscopes	EB-1970UK	

OE-A56 Balloons for PENTAX Ultrasound Endoscopes



Model name OE-A56	Description Balloons for endoscopic ultrasound application	
Item number 83385	Pack size 20	
Compatible Endoscopes	EG-3270UK	

OE-A60 Single-Use Balloon



Model name OE-A60	Description Balloons for endoscopic ultrasound application	
Item number 84875	Pack size 20	
Compatible Endoscopes	EG34-J10U	EG38-J10UT

OE-A61 PENTAX Medical Single-Use Balloon



Model name OE-A61	Description Single-Use balloon
Item number 84629	
Compatible Endoscopes	EB19-J10U

OF-A67 Balloon-O-Ring Kit



Model name OF-A67	Description Balloon and o-ring kit for EUS Gastroscope
Item number 87763	Pack size 20
Compatible Endoscopes	EG-3870UTK

OE-A53 Balloon Installation Device



Model name OE-A53	Description Balloon installation device
Item number 83348	
Compatible Endoscopes	EB-1970UK EB19-J10U

OF-B219 O-Ring Set



Model name
OF-B219

Description

O-ring set and installation device

Item number
84835

Compatible
Valves

OF-B207

OF-B214 O-Ring Set



Model name
OF-B214

Description


A set of 3 pcs. of O-rings for the cleaning adapter

Item number
84589

Compatible
Valves

OF-B212

OF-B204 O-Ring Set



Model name
OF-B204

Description

O-ring set


Item number
80248

Pack size
2

Compatible
Light Source

BS-LL1

OF-B195 O-Ring Set



Model name
OF-B195

Description

O-ring set

Item number
85138

Pack size
5

Compatible
Valves

OF-B194

OE-B10

OF-B192 O-Ring Set For Air/Water Feeding Valves



Model name
OF-B192

Description

O-ring set for air/water feeding valves

Item number
87144

Compatible
Valves

OF-B188

OF-B191 O-Ring Set for Air/Water/Balloon Feeding Valves



Model name OF-B191	Description O-ring set for air/water/balloon feeding valves
Item number 87066	
Compatible Valves	OF-B172

OF-B187 O-Ring Set for Suction Control Valves




Model name OF-B187	Description To be attached to the suction/balloon suction control valve
Item number 87065	
Compatible Valves	OF-B171 OF-B206

OF-B184 O-Ring Set For Air/Water Feeding Valves



Model name OF-B184	Description O-Ring-set for A/W valve
Item number 87949	
Compatible Valves	OF-B183 OF-B178

OF-B182 O-Ring Set




Model name OF-B182	Description O-ring set, 3 pcs
Item number 86816	
Compatible Valves	OF-B113

OF-B181 O-Ring Set for Suction Control Valves



Model name OF-B181	Description O-Ring-Set for suction valve, 3pcs
Item number 86807	
Compatible Valves	OF-B179

OF-B143 O-Ring Set



Model name OF-B143	Description O-ring set for A/W valve, 6 pcs
Item number 82534	
Compatible Valves	OF-B121

OF-B127 O-Ring Set for Suction Control Valves



Model name OF-B127	Description O-ring set for suction valve, 3 pcs
Item number 81782	
Compatible Valves	OF-B120 OF-B177

OF-B100 O-Ring Set For Suction Control/Irrigation Cessation Valves



Model name OF-B100	Description O-ring set for suction valve, 2 X 3 pcs
Item number 81376	
Compatible Valves	OF-B95 OF-B140



OE-C12 Water Jet Check Valve Adapter



Model name
OE-C12

Description
Forward water jet check valve adapter

Item number
86894

Compatible
Endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10TL	EC38-i10cM
EC-2990Fi	EC-3890Mi2	EC34-i10TF	EC38-i10cF2
EC-2990Mi	EC-3890LK	EC34-i10TM	EG-2990i
EC-3490Li	EC-3890FK	EC34-i10cL	EG-2990K
EC-3490Fi	EC-3890MK	EC34-i10cF	EG-2990Zi
EC-3490Mi	EC-3890FK2	EC34-i10cM	EG-3490K
EC-3490LK	EC-3890MK2	EC38-i10L	EG-3890TK
EC-3490FK	EC-3890LZi	EC38-i10F	EG29-i10
EC-3490MK	EC-3890FZi	EC38-i10M	EG29-i10c
EC-3490TLi	EC-3890MZi	EC38-i10F2	EG34-i10
EC-3490TFi	EC-3890TLK	EC38-i10M2	VSB-2990i
EC-3490TMi	EC-3890TFK	EC38-i10NL	ES-3870K
EC-3890Li	EC34-i10L	EC38-i10NF	
EC-3890Fi	EC34-i10F	EC38-i10cL	
EC-3890Mi	EC34-i10M	EC38-i10cF	

OE-C15 Check Valve Set



Model name
OE-C15

Description
Rubber check valve set (10 pcs) (black)

Item number
86917

Pack size
10

Compatible
Valves

OE-C12 OF-B121





OF-B118 Water Jet Connector Cap



Model name
OF-B118

Description
Water jet connector cap

Item number
81722

Compatible
Endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10TL	EC38-i10cM
EC-2990Fi	EC-3890Mi2	EC34-i10TF	EC38-i10cF2
EC-2990Mi	EC-3890LK	EC34-i10TM	EG-2990i
EC-3490Li	EC-3890FK	EC34-i10cL	EG-2990K
EC-3490Fi	EC-3890MK	EC34-i10cF	EG-2990Zi
EC-3490Mi	EC-3890FK2	EC34-i10cM	EG-3490K
EC-3490LK	EC-3890MK2	EC38-i10L	EG-3890TK
EC-3490FK	EC-3890LZi	EC38-i10F	EG29-i10
EC-3490MK	EC-3890FZi	EC38-i10M	EG29-i10c
EC-3490TLi	EC-3890MZi	EC38-i10F2	EG34-i10
EC-3490TFi	EC-3890TLK	EC38-i10M2	VSB-2990i
EC-3490TMi	EC-3890TFK	EC38-i10NL	ES-3870K
EC-3890Li	EC34-i10L	EC38-i10NF	
EC-3890Fi	EC34-i10F	EC38-i10cL	
EC-3890Mi	EC34-i10M	EC38-i10cF	

OF-B113 Irrigation Tube



Model name
OF-B113

Description
Irrigation tube for WJ channel


Item number
81762

Compatible
Endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10TL	EC38-i10cM
EC-2990Fi	EC-3890Mi2	EC34-i10TF	EC38-i10cF2
EC-2990Mi	EC-3890LK	EC34-i10TM	EG-2990i
EC-3490Li	EC-3890FK	EC34-i10cL	EG-2990K
EC-3490Fi	EC-3890MK	EC34-i10cF	EG-2990Zi
EC-3490Mi	EC-3890FK2	EC34-i10cM	EG-3490K
EC-3490LK	EC-3890MK2	EC38-i10L	EG-3890TK
EC-3490FK	EC-3890LZi	EC38-i10F	EG29-i10
EC-3490MK	EC-3890FZi	EC38-i10M	EG29-i10c
EC-3490TLi	EC-3890MZi	EC38-i10F2	EG34-i10
EC-3490TFi	EC-3890TLK	EC38-i10M2	VSB-2990i
EC-3490TMi	EC-3890TFK	EC38-i10NL	ES-3870K
EC-3890Li	EC34-i10L	EC38-i10NF	
EC-3890Fi	EC34-i10F	EC38-i10cL	
EC-3890Mi	EC34-i10M	EC38-i10cF	



OF-G11 Gas Adapter



Model name	Description			
OF-G11	CO ₂ gas adapter			
Item number				
83647				
Compatible Endoscopes	EC-2990Li	EC-3890LK	EC38-i10M	EG-2990Zi
	EC-2990Fi	EC-3890FK	EC38-i10F2	EG-3490K
	EC-2990Mi	EC-3890MK	EC38-i10M2	EG-3890TK
	EC-3490Li	EC-3890FK2	EC38-i10NL	EG16-K10
	EC-3490Fi	EC-3890MK2	EC38-i10NF	EG17-J10
	EC-3490Mi	EC-3890LZi	EC38-i10cL	EG27-i10
	EC-3490LK	EC-3890FZi	EC38-i10cF	EG29-i10
	EC-3490FK	EC-3890MZi	EC38-i10cM	EG29-i10c
	EC-3490MK	EC-3890TLK	EC38-i10cF2	EG34-i10
	EC-3490TLi	EC-3890TFK	ED-3490TK	VSB-2990i
	EC-3490TFi	EC34-i10L	ED32-i10	EE-1580K
	EC-3490TMi	EC34-i10F	ED34-i10T	EG-3270UK
	EC-380LKp	EC34-i10M	ED34-i10T2	EG-3670URK
	EC-380FKp	EC34-i10TL	EE17-J10	EG-3870UTK
	EC-380FK2p	EC34-i10TF	EG-1690K	EG34-J10U
	EC-380MK2p	EC34-i10TM	EG-2490K	EG36-J10UR
	EC-3890Li	EC34-i10cL	EG-2790i	EG38-J10UT
	EC-3890Fi	EC34-i10cF	EG-2790K	ES-3870K
	EC-3890Mi	EC34-i10cM	EG-290Kp	
	EC-3890Fi2	EC38-i10L	EG-2990i	
	EC-3890Mi2	EC38-i10F	EG-2990K	

OE-B5 Injection Tube



Model name	Description			
OE-B5	Injection tube for ultrasound endoscopes			
Item number				
81160				
Compatible Endoscopes	EB-1970UK	EB19-J10U	EG-3270UK	EG-3870UTK



OE-B3 Elevator Wire Cleaning Inlet Cap



Model name	Description			
OE-B3	Elevator wire cleaning inlet cap			
Item number				
87089				
Compatible Endoscopes	EB-1970UK	EB19-J10U	EG-3270UK	EG-3870UTK




OF-G1 PENTAX Standard Adapter Sleeve



Model name	Description			
OF-G1	Adapter sleeve			
Item number				
83181				
Compatible Endoscopes	FC-38LV	FB-8V	FNL7RP3	VNL1590STi
	FC-38FV	FB-10V	FNL10RP3	VNL8-J10
	FD-34V2	FB-15V	FNL15RP3	VNL11-J10
	FG-16V	FB-18V	FUR-9RBS	VNL15-J10
	FG-24V	FB-19TV	VNL1070STK	BS-LC1
	FG-29V	FCN-15X	VNL1190STK	
	FS-34V	FCP-9P	VNL1570STK	

OF-B55 T-adapter



Model name	Description			
OF-B55	T-adapter for luer lock			
Item number				
83633				
Compatible Endoscopes	EC-2990Li	EC34-i10TL	EG27-i10	EG-3870UTK
	EC-2990Fi	EC34-i10TF	EG29-i10	EG34-J10U
	EC-2990Mi	EC34-i10TM	EG29-i10c	EG36-J10UR
	EC-3490Li	EC34-i10cL	EG34-i10	EG38-J10UT
	EC-3490Fi	EC34-i10cF	FC-38LV	ES-3870K
	EC-3490Mi	EC34-i10cM	FC-38FV	FB-8V
	EC-3490LK	EC38-i10L	FD-34V2	FB-10V
	EC-3490FK	EC38-i10F	FG-16V	FB-15V
	EC-3490MK	EC38-i10M	FG-24V	FB-18V
	EC-3490TLi	EC38-i10F2	FG-29V	FB-19TV
	EC-3490TFi	EC38-i10M2	FS-34V	FB-15BS
	EC-3490TMi	EC38-i10NL	VSB-2990i	FB-15RBS
	EC-380LKp	EC38-i10NF	EB-1170K	FB-18BS
	EC-380FKp	EC38-i10cL	EB-1570AK	FB-18RBS
	EC-380FK2p	EC38-i10cF	EB-1570K	FCN-15X
	EC-380MK2p	EC38-i10cM	EB-1572K	FCP-9P
	EC-3890Li	EC38-i10cF2	EB-1575K	FCY15RBS
	EC-3890Fi	ED-3490TK	EB-1970AK	FHY10RBS
	EC-3890Mi	ED32-i10	EB-1970K	FHY15RBS
	EC-3890Fi2	ED34-i10T	EB-1970TK	FI-9BS
	EC-3890Mi2	ED34-i10T2	EB-1972K	FI-9RBS
	EC-3890LK	EE17-J10	EB-1975K	FI-10BS
	EC-3890FK	EG-1690K	EB-1990i	FI-10RBS
	EC-3890MK	EG-2490K	EB11-J10	FI-13BS
	EC-3890FK2	EG-2790i	EB15-J10	FI-13RBS
	EC-3890MK2	EG-2790K	EB19-J10	FI-16BS
	EC-3890LZi	EG-290Kp	EB-1970UK	FI-16RBS
	EC-3890FZi	EG-2990i	EB19-J10U	FNL15RP3
	EC-3890MZi	EG-2990K	ECY1570K	FUR-9RBS
	EC-3890TLK	EG-2990Zi	ECY1575K	VNL1570STK
	EC-3890TFK	EG-3490K	EE-1580K	VNL15-J10
	EC34-i10L	EG-3890TK	EG-3270UK	
	EC34-i10F	EG16-K10	EG-3670URK	
	EC34-i10M	EG17-J10		

OF-B161 Instrument Channel Selector



Model name
OF-B161

Description

Instrument channel selector for dual channel scopes

Item number
82002

Compatible
Endoscopes

EC-3890TLK

EC-3890TFK

EG-3890TK

OE-C28 Ventilation Cap



Model name
OE-C28

Description

Ventilation cap for IMAGINA endoscopes

Item number
84235

Compatible
Endoscopes

EC34-i10cL

EC34-i10cM

EC38-i10cF

EC38-i10cF2

EC34-i10cF

EC38-i10cL

EC38-i10cM

EG29-i10c

OE-C29 Venting Connector Adapter



Model name
OE-C29

Description

Venting connector adapter for IMAGINA endoscopes

Item number
84236

Compatible
Endoscopes

EC34-i10cL

EC34-i10cM

EC38-i10cF

EC38-i10cF2

EC34-i10cF

EC38-i10cL

EC38-i10cM

EG29-i10c

OE-C21 Venting Cap-English



Model name
OE-C21

Description

Venting cap-English

Item number
82133

Compatible
Endoscopes

EC34-i10TL

EC38-i10NL

ED34-i10T

EG17-J10

EC34-i10TF

EC38-i10NF

ED34-i10T2

EG34-i10

EC34-i10TM

ED32-i10

EE17-J10

OF-C5 Ventilation Cap



Model name	Description			
OF-C5	Ventilation cap			
Item number				
81522				
Compatible Endoscopes	EC-2990Li	EC34-i10F	EB-1572K	FB-18RBS
	EC-2990Fi	EC34-i10M	EB-1575K	FCN-15X
	EC-2990Mi	EC38-i10L	EB-1970AK	FCP-9P
	EC-3490Li	EC38-i10F	EB-1970K	FCY-15RBS
	EC-3490Fi	EC38-i10M	EB-1970TK	FHY-10RBS
	EC-3490Mi	EC38-i10F2	EB-1972K	FHY-15RBS
	EC-3490LK	EC38-i10M2	EB-1975K	FI-7BS
	EC-3490FK	ED-3490TK	EB-1990i	FI-7RBS
	EC-3490MK	EG-1690K	EB11-J10	FI-9BS
	EC-3490TLi	EG-2490K	EB15-J10	FI-9RBS
	EC-3490TFi	EG-2790i	EB19-J10	FI-10BS
	EC-3490TMi	EG-2790K	EB-1970UK	FI-10RBS
	EC-380LKp	EG-290Kp	EB19-J10U	FI-13BS
	EC-380FKp	EG-2990i	ECY-1570K	FI-13RBS
	EC-380FK2p	EG-2990K	ECY-1575K	FI-16BS
	EC-380MK2p	EG-2990Zi	EE-1580K	FI-16RBS
	EC-3890Li	EG-3490K	EG-3270UK	FNL-7RP3
	EC-3890Fi	EG-3890TK	EG-3670URK	FNL-10RBS
	EC-3890Mi	EG16-K10	EG-3870UTK	FNL-10RP3
	EC-3890Fi2	EG27-i10	EG34-J10U	FNL-15RP3
	EC-3890Mi2	EG29-i10	EG36-J10UR	FUR-9RBS
	EC-3890LK	FC-38LV	EG38-J10UT	VNL-1070STK
	EC-3890FK	FC-38FV	ES-3870K	VNL-1190STK
	EC-3890MK	FD-34V2	FB-8V	VNL-1570STK
	EC-3890FK2	FG-16V	FB-10V	VNL-1590STi
	EC-3890MK2	FG-24V	FB-15V	VNL8-J10
	EC-3890LZi	FG-29V	FB-18V	VNL11-J10
	EC-3890FZi	FS-34V	FB-19TV	VNL15-J10
	EC-3890MZi	VSB-2990i	FB-15BS	PVK-1070Z
	EC-3890TLK	EB-1170K	FB-15RBS	
	EC-3890TFK	EB-1570AK	FB-18BS	
	EC34-i10L	EB-1570K		

BS-LC1 Light Guide Cable



Model name	Description			
BS-LC1	Light guide cable for BS-endoscopes			
Item number				
85736				
Compatible Endoscopes	FB-15BS	FHY-15RBS	FI-10RBS	FUR-9RBS
	FB-15RBS	FI-7BS	FI-13BS	
	FB-18BS	FI-7RBS	FI-13RBS	
	FB-18RBS	FI-9BS	FI-16BS	
	FCY-15RBS	FI-9RBS	FI-16RBS	
	FHY-10RBS	FI-10BS	FNL-10RBS	

BS-LL1 LED Light Source Batt-Type




Model name	Description			
BS-LL1	LED light source battery type			
Item number				
85136				
Compatible Endoscopes	FB-15BS	FHY-15RBS	FI-10RBS	FUR-9RBS
	FB-15RBS	FI-7BS	FI-13BS	
	FB-18BS	FI-7RBS	FI-13RBS	
	FB-18RBS	FI-9BS	FI-16BS	
	FCY-15RBS	FI-9RBS	FI-16RBS	
	FHY-10RBS	FI-10BS	FNL-10RBS	

OF-Z5 Bite Block



Model name	Description
OF-Z5	Bite block, autoclavable
Item number	
83036	

OF-Z11 Silicone Oil



Model name	Description
OF-Z11	Silicone oil
Item number	
83033	

AP-PV1 Adapter for Video Module for PVK-1070Z



Model name	Description			
AP-PV1	Adapter lens for PENTAX Medical video module			
Item number				
81362				
Compatible Endoscopes	FC-38LV	FG-29V	FB-19TV	FI-7BS
	FC-38FV	FB-8V	FB-15BS	FI-9BS
	FD-34V2	FB-10V	FB-18BS	FI-10BS
	FG-16V	FB-15V	FCN-15X	FI-13BS
	FG-24V	FB-18V	FCP-9P	

AP-RV1 Scope Mount Adapter for PVK-1070Z

Picture not available	Model name AP-RV1	Description Adapter lens for PENTAX Medical video module	
	Item number 81363		
	Compatible Endoscopes	FB-15RBS	FB-18RBS

AL-SL Storz Light Source Adapter

Picture not available	Model name AL-SL	Description Storz light source adapter			
	Item number 83241				
	Compatible Endoscopes	FC-38LV	FB-8V	FNL-7RP3	VNL-1590STi
		FC-38FV	FB-10V	FNL-10RP3	VNL8-J10
		FD-34V2	FB-15V	FNL-15RP3	VNL11-J10
		FG-16V	FB-18V	FUR-9RBS	VNL15-J10
		FG-24V	FB-19TV	VNL-1070STK	BS-LC1
		FG-29V	FCN-15X	VNL-1190STK	
		FS-34V	FCP-9P	VNL-1570STK	

AT-ST2 Ocular Adapter to Rigid Type

Picture not available	Model name AT-ST2	Description Ocular adapter to rigid type			
	Item number 83383				
	Compatible Endoscopes	FC-38LV	FS-34V	FB-15BS	FI-10BS
		FC-38FV	FB-8V	FB-18BS	FI-13BS
		FD-34V2	FB-10V	FCN-15X	FI-16BS
		FG-16V	FB-15V	FCP-9P	
		FG-24V	FB-18V	FI-7BS	
		FG-29V	FB-19TV	FI-9BS	

OF-D1 Eyepiece Cap



Model name OF-D1	Description Eyepiece cap			
Item number 83043				
Compatible Endoscopes	FD-34V2	FB-8V	FB-19TV	FI-13BS
	FG-16V	FB-10V	FCN-15X	FI-16BS
	FG-24V	FB-15V	FI-9BS	
	FG-29V	FB-18V	FI-10BS	

OF-D3 Eyeshield



Model name OF-D3	Description Eyeshield			
Item number 83544				
Compatible Endoscopes	FD-34V2	FG-24V	FB-15V	FB-19TV
	FG-16V	FG-29V	FB-18V	FCN-15X

OF-C3 Soaking Cap



Model name OF-C3	Description Soaking cap			
Item number 83665				
Compatible Endoscopes	FCN-15X			

GM00261 Air/Water Valve

Picture not available

Model name
GM00261**Description**

Air/water valve for VERSA endoscopes

Item number
300957Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GM00259 Suction Valve

Picture not available

Model name
GM00259**Description**

Suction valve for VERSA endoscopes

Item number
300958Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

MF00004 Biopsy Valve

Picture not available

Model name
MF00004**Description**

Biopsy valve for VERSA endoscopes

Item number
300959Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GS00006 Water Bottle

Picture not available

Model name
GS00006**Description**

Water bottle for VERSA endoscopes

Item number
3000956Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

Manufacturers

VERSA Consumables

 **PENTAX-AOHUA Medical Technologies Co., Ltd.**
East of 3rd Floor, Block C, Building 1,
No.5 Shenwang Road,
Minhang District, Shanghai 201108,
P. R. China.

All other products

 **HOYA Corporation**
6-10-1, Nishi-shinjuku
Shinjuku-ku, Tokyo
160-0023 Japan



EMEA Headquarter Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg, Germany



Medical Device Class: I



EMEA Headquarter Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg, Germany



0123
TÜV SÜD CE 0123 Medical Device Classes: IIa, IIb, Is and I

Distributor

EMEA Headquarter Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL

Electrosurgery and Argon Plasma Coagulation (APC)

05

Electrosurgery systems



Manufacturer reference	Item number	Item description	Notes
900-250	306050	Electrosurgical unit ENDO ARC	



Manufacturer reference	Item number	Item description	Notes
900-000	306051	ENDO PLUS unit for ENDO ARC	



Manufacturer reference	Item number	Item description	Notes
900-400	306090	Electrosurgical unit ARC 400	Only available in specific markets, please check with your local sales representative.



Manufacturer reference	Item number	Item description	Notes
900-001	306091	ARC Plus unit for ARC 400	Only available in specific markets, please check with your local sales representative.

Argon management



Manufacturer reference	Item number	Item description	Notes
900-806	306052	Gas cylinder Argon, 4.8 l, connection DIN477 no. 6, empty	eg. CEE, Spain
900-807	306112	Gas cylinder Argon, 4.8 l, connection DIN477 no. 10, empty	eg. NL, Nordics
900-808	306113	Gas cylinder Argon, 4.8 l, connection CGA no. 580, empty	eg. Italy, UAE
900-809	306114	Gas cylinder Argon, 4.8 l, connection BS 341 no. 3, empty	eg. UK, Turkey



Manufacturer reference	Item number	Item description	Notes
900-901	306053	Pressure reducer, for ENDO PLUS (REF 900-000), DIN 477 no. 6	eg. CEE, Spain (for ENDO PLUS)
900-902	306054	Pressure reducer, for ENDO PLUS (REF 900-000), DIN 477 no. 10	eg. NL, Nordics (for ENDO PLUS)
900-903	306055	Pressure reducer, for ENDO PLUS (REF 900-000), CGA 580	eg. Italy, UAE (for ENDO PLUS)
900-904	306056	Pressure reducer, for ENDO PLUS (REF 900-000), BS 341 no. 3	eg. UK, Turkey (for ENDO PLUS)



Manufacturer reference	Item number	Item description	Notes
900-906	306057	Pressure reducer, for ARC PLUS (REF 900-001), DIN 477 no. 6	eg. CEE, Spain (for ARC Plus)
900-907	306058	Pressure reducer, for ARC PLUS (REF 900-001), DIN 477 no. 10	eg. NL, Nordics (for ARC Plus)
900-908	306059	Pressure reducer, for ARC PLUS (REF 900-001), CGA 580	eg. Italy, UAE (for ARC Plus)
900-909	306060	Pressure reducer, for ARC PLUS (REF 900-001), BS 341 no. 3	eg. UK, Turkey (for ARC Plus)



Manufacturer reference	Item number	Item description	Notes
830-050	306061	ARC PLUS filter set, single-use, sterile (50 pcs.)	

Footswitch



Manufacturer reference	Item number	Item description	Notes
901-022	306064	Double-pedal footswitch, cable 4 m for ENDO ARC, PENTAX branded	
901-031	306065	Double-pedal footswitch with switch, cable 4 m for ARC 400, BOWA branded	

Equipotential bonding




Manufacturer reference	Item number	Item description	Notes
900-030	306074	Equipotential bonding, 0.5 m	
900-031	306075	Equipotential bonding, 5 m	
900-035	306076	Equipotential bonding, 1.5 m	

Cables




Manufacturer reference	Item number	Item description	Notes
932-146	306077	Cable for flexible Argon probe COMFORT, 2.5 m	
131-145	306081	Monopolar cable, endoscopy 2.8 mm socket, COMFORT, 4.5 m	for ARC 400
431-045	306085	Monopolar cable, endoscopy, 2.8 mm socket, for 4 mm, 4.5 m	for ENDO ARC

Disposable probes



Manufacturer reference	Item number	Item description	Notes
932-248	306086	Flexible Argon probe, Ø 1,5 mm, 1,5 m, single-use, sterile (10 pcs.)	
932-249	306078	Flexible Argon probe, Ø 2.3 mm, 2.2 m, single-use, sterile (10 pcs.)	
932-250	306087	Flexible Argon probe, Ø 3,2 mm, 2,2 m, single-use, sterile (10 pcs.)	
932-251	306088	Flexible Argon probe, Ø 1,5 mm, 3 m, single-use, sterile (10 pcs.)	
932-252	306079	Flexible Argon probe, Ø 2.3 mm, 3 m, single-use, sterile (10 pcs.)	
932-253	306080	Flexible Argon probe, 90°, Ø 2.3 mm, 2.2 m, single-use, ster. (10 pcs.)	

Reusable probes



Manufacturer reference	Item number	Item description	Notes
932-148	306092	Flexible Argon probe, Ø 1,5 mm, 1,5 m, reusable, non-sterile (1 pc.)	
932-151	306093	Flexible Argon probe, Ø 1,5 mm, 3 m, reusable, non-sterile (1 pc.)	
932-154	306094	Flexible Argon probe, Ø 2.3 mm, 1 m, reusable, non-sterile (1 pc.)	
932-149	306095	Flexible Argon probe, Ø 2.3 mm, 2,2 m, reusable, non-sterile (1 pc.)	
932-152	306096	Flexible Argon probe, Ø 2.3 mm, 3 m, reusable, non-sterile (1 pc.)	
932-150	306097	Flexible Argon probe, Ø 3.2 mm, 2,2 m, reusable, non-sterile (1 pc.)	
932-153	306098	Flexible Argon probe, Ø 2.3 mm, 2,2 m, reusable, non-sterile (1 pc.)	

Rinsing adaptor




Manufacturer reference	Item number	Item description	Notes
932-910	306089	Rinsing adaptor for flexible argon probes to connect to washing machines (5 pieces)	

Bipolar adaptor



Manufacturer reference	Item number	Item description	Notes
901-215	306099	Bipolar adaptor	Adaptor for HemoStat WideCup

Disposable Return plate



Manufacturer reference	Item number	Item description	Notes
380-050	306083	Cable, single-use return plate, international, 4.5 m	Cable to be used with single use return plate
816-092	306110	Return plate, 90 cm ² , split, EASY (100 pcs.)	
818-092	306111	Return plate, 90 cm ² , split, with cable, 3 m, EASY (15 pcs.)	
818-112	306124	Return plate, 110 cm ² , split, with cable, 3 m, EASY (15 pcs.)	

Reusable return plate



Manufacturer reference	Item number	Item description	Notes
385-050	306120	Cable, reusable rubber return plate, international, 4.5 m	Cable to be used with reusable rubber return plate
242-003	306118	Rubber reusable return plate, adults, 250 x 150 mm, international	
232-003	306119	Rubber reusable return plate, pediatric, 80mm x 150 mm, international	

Rubber band




Manufacturer reference	Item number	Item description	Notes
194-075	306125	Rubber band with fixation button, 75 cm, latex free	

Fibre optics



Manufacturer reference	Item number	Item description	Notes
901-051	306115	Fibre optics, for ARC PLUS (REF 900-001)	
901-052	306116	Fibre optics, for ENDO ARC (REF 900-000)	

Cart



Manufacturer reference	Item number	Item description	Notes
10825340	306500	Fuego Basic cart II	basic cart
10825350	306501	Fuego push back handle	to be additional to the cart
10825460	306502	Wire basket	to be additional to the cart
10000285	306503	Tray	to be additional to the cart
10826010	306504	Drawer	to be additional to the cart
10825740	306505	Fuego gas bottle holder left, closed version, max bottle diameter 145 mm	to be additional to the cart
10825710	306506	Fuego gas bottle holder left, 5l, open version, max bottle diameter 145 mm	to be additional to the cart
10825730	306507	Fuego gas bottle holder right/left, closed version, max bottle diameter 145 mm	to be additional to the cart
10825720	306508	Fuego gas bottle holder right/left, open version, max bottle diameter 145 mm	to be additional to the cart
30309	306509	Front handle for Fuego cart	to be additional to the cart
30310	306510	POAG-PIN	to be additional to the cart
30313	306506	Fuego holder for gas bottle 5/10 ltrs., left side, max. diam. bottle = 145mm	to be additional to the cart
30315	306511	Fuego upper holder for gas bottle 5/10 ltrs., right side	to be additional to the cart
30316	306512	Fuego upper holder for gas bottle 5/10 ltrs., left side	to be additional to the cart

Manufacturers

Carts

HAEBERLE GmbH & Co. KG

Breitwiesenstr. 13
70565 Stuttgart, Germany

ARC PLUS filter set

Pall Medical - A division of Pall International Sàrl

Avenue de Tivoli 3
CH-1700 Fribourg, Switzerland

All other products

BOWA-electronic GmbH & Co.KG

Heinrich-Hertz-Strasse 4-10
72810 Gomaringen, Germany



Pall Genedisc Technologies

1 Rue du Courtil
Centre d'Affaires Cicea - Bat 1
35170 Bruz, France



1282

ECM CE 1282 Medical Device Class: IIa



0123

TÜV SÜD CE 0123 Medical Device Class: I, IIb

Distributors

EMEA Headquarter

Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

United Kingdom

PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical.uk@pentaxmedical.com

France

PENTAX France Life Care S.A.S.
116 quai de Bezons
B.P. 204
95106 Argenteuil Cedex
Tel: +33 1 / 3025 7575
Fax: +33 1 / 3025 7445
E-mail: clients.fr@pentaxmedical.com

Netherlands

PENTAX Nederland B.V.
Edisonring 6
6669NB Dodewaard
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: info.nl@pentaxmedical.com

Italy

PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.it@pentaxmedical.com

Spain

SIMMEDICA – Sistemas Integrales
de Medicina, S.A.
28830 San Fernando de Henares - Madrid
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

Russia

Moscow Representative office
of PENTAX Europe GmbH (Germany)
13, 4-y Lesnoy pereulok
125047 Moscow
Tel.: +7 495 114 52 31
Fax: +7 495 114 52 31

Turkey

PENTAX Turkey
Veko Giz Plaza, Meydan Sokak No:3/43
343396 Maslak – Istanbul
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL

Hygiene Consumables

06

Hygiene Consumables

table of contents

06-1 PROfILE Single-Use Cleaning Brushes

06-2 Cleaning Brushes

06-3 Reprocessing Accessories

06-4 VERSA Consumables



DESIGNED FOR PENTAX Medical – One product range that fits every PENTAX Medical endoscope


Manually pre-cleaning flexible endoscopes is key to ensuring reliable reprocessing as an end result. Even with magnification, visual assessment of the manual pre-cleaning results is insufficient to detect microorganisms. As such, PENTAX Medical provides a verified and validated solution to improve confidence in manual pre-cleaning, with the new PROfILE single-use cleaning brushes.

1**Proven cleaning efficiency****2****Designed and validated to meet PENTAX Medical quality standards****3****Compliant with the industry's highest reprocessing standards**


Overview PROfILE Single-Use Cleaning Brushes

Long instrument/suction channel brushes

CS5522A



CS5530A



CS2010A




CS2037A



CS4010A



CS3022A



CS5503A



CS2022A




Short instrument channel inlet brushes


CS-C13A



CS-C11A



CS-C12A




Find the right brush with the brush finder on www.the-hygiene-solution-that-fits.com

Proven cleaning efficiency

In combination with PENTAX Medical endoscopes, the PROfile single-use cleaning brushes were tested against residual soil and compared against the performance of legacy brushes.

Through an authorized statistical test method, the internal validation proved both brushes to have a similar performance.



Designed and validated to meet PENTAX Medical quality standards

PROfile single-use cleaning brushes were made to meet the high-quality standards of PENTAX Medical engineering. Extensive internal tests have proven that they are the best possible fit to ensure proper hygiene results, with no impairment to the endoscope, across all types.

- Tailored for PENTAX Medical endoscopes' precise dimensions, structures, and materials
- Improved bonding and pull-off capability at optimal insertion resistance of $\leq 10N$
- No damage to endoscope working channels and reliable even at increased thermal conditions

Compliant with the industry's highest reprocessing standards

How clean is clean enough? This was the leading question of PENTAX Medical engineers when deciding upon the validation criteria for hygienic efficacy of the PROfile single-use cleaning brushes. In order to follow established standards, PENTAX Medical validated against the requirements of the Association for the Advancement of Medical Instrumentation (AAMI)¹.

Article	Hemoglobin [< 2.2µg/cm²]			Protein [< 6.4µg/cm²]			Carbohydrate [< 1.8µg/cm²]		
	µg/article	µg/cm²	Result	µg/article	µg/cm²	Result	µg/article	µg/cm²	Result
Endoscope 1	14.88	1.06	PASS	31.05	2.20	PASS	10.82	0.77	PASS
Endoscope 2	9.92	0.70	PASS	22.68	1.61	PASS	8.05	0.57	PASS
Endoscope 3	19.84	1.41	PASS	59.40	4.21	PASS	6.82	0.48	PASS



Test example: Residual of three markers at scope's distal end after third cleaning cycle, complying with AAMI requirements.

¹ AAMI TIR30:2011

PROfile Single-Use Cleaning Brushes

One product range that matches every channel

Brush				Compatible endoscope							
Model	Length [mm]	Outer Diameter [ø mm]	Color [Shaft / Handle]	Compatible Channel I.D. [mm]	Working Length [mm]	Gastro-enterology	Pulmo-nology	Uro-logy	ENT & Speech	Intu-bation	EUS
CS2010A	1,020	2	Blue	1.1 - 1.5	Up to 720		•	•	•	•	
CS2022A	2,200	2		1.1 - 1.5	Up to 1,900						•
CS2037A	3,720	2		1.1 - 1.5	Up to 3,420	•					
CS3022A	2,220	3	White	1.7 - 2.4	Up to 1,920	•					
CS4010A	1,020	4	Orange	1.8 - 3.2	Up to 720		•	•	•	•	
CS5503A	310	5.5	Yellow	2.0 - 5.0	N/A		•	•	•	•	
CS5522A	2,420	5.5		2.0 - 5.0	Up to 2,120	•					
CS5530A	3,020	5.5		2.0 - 5.0	Up to 2,720	•					
CS-C11A	N/A	6-14	Cobalt Blue	N/A	N/A		•	•	•	•	
CS-C12A	N/A	5	Yellow	N/A	N/A						•
CS-C13A	N/A	6-15	Green	N/A	N/A	•					



PROfile Single-Use Cleaning Brushes

One product range to match every channel



PENTAX Medical PROfile single-use endoscope cleaning brush kits with proven usability are available in a wide range of dimensions to cover all channel sizes.

Compliant with reprocessing guidelines

- Single-use cleaning accessories are internationally recommended by reprocessing guidelines.

Meeting infection control requirements

- PENTAX Medical single-use cleaning brush kits help to ensure safe reprocessing and simplify processes.

Proven cleaning results

- Streamlined product range to match all PENTAX Medical endoscopes in all product segments
- Tailored for PENTAX Medical endoscopes' precise dimensions, structures and materials
- Validated cleaning efficacy

PROfile Shaft Color	Model Name	Working Length [mm]	Endoscope Channel Diam- eter [mm]	Endoscope type (both video and fiberoptic)					
				Gastro- entero- logy	Pulmo- nology	Urology	ENT & Speech	Intuba- tion	EUS
Blue	CS2010A	Up to 720	1.1 to 1.5		•	•	•	•	
	CS2022A	Up to 1,900	1.1 to 1.5						•
	CS2037A	Up to 3,420	1.1 to 1.5	•					
White	CS3022A	Up to 1,920	1.7 to 2.4	•					
Orange	CS4010A	Up to 720	1.8 to 3.2		•	•	•	•	
Yellow	CS5503A	N/A	2.0 to 5.0		•	•	•	•	
	CS5522A	Up to 2,120	2.0 to 5.0	•					
	CS5530A	Up to 2,720	2.0 to 5.0	•					
Cobalt Blue	CS-C11A	N/A	N/A		•	•	•	•	
Yellow	CS-C12A	N/A	N/A						•
Green	CS-C13A	N/A	N/A	•					

Consult your local representative for information on the products currently available in your region.

PENTAX Medical

Accessories Programme for cleaning and caring for your PENTAX Medical endoscopes

Original is better

Finely tuned system components are essential for safe and successful endoscopy. That is why PENTAX Medical original accessories aim to provide the highest possible quality in every detail. The following pages describe our carefully compiled accessory programme to ensure safe and effective cleaning of your PENTAX Medical endoscope.



How to choose the correct accessories:

Before you decide on new accessories, please ensure that your PENTAX Medical endoscope is compatible with the desired accessory.

- The illustration below is intended to help you to identify the relevant areas of your endoscope quickly and simply. Our PENTAX Medical specialists will be pleased to advise you if you have any additional questions.
- All PENTAX Medical original accessories in this catalogue have been tested in conjunction with PENTAX Medical endoscopes for function and compatibility.

Orientation

In order to make your choice of PENTAX Medical accessories as simple and reliable as possible, we have identified fundamental product features with simple-to-understand symbols displayed next to the relevant technical data.



Autoclavable



**For single
use only**



CS6021T Tri-Bristled Cleaning Brushes



Model name
CS6021T

Description

Single-use cleaning brushes for channel

Item number
86147

Pack size
10

Compatible
endoscopes

EC-2990Li	EC-3890Mi2	EC34-i10TM	EG-2990Zi
EC-2990Fi	EC-3890LK	EC38-i10L	EG-3490K
EC-2990Mi	EC-3890FK	EC38-i10F	EG-3890TK
EC-3490Li	EC-3890MK	EC38-i10M	EG16-K10
EC-3490Fi	EC-3890FK2	EC38-i10F2	EG27-i10
EC-3490Mi	EC-3890MK2	EC38-i10M2	EG29-i10
EC-3490LK	EC-3890LZi	EC38-i10NL	EG34-i10
EC-3490FK	EC-3890FZi	EC38-i10NF	EE-1580K
EC-3490MK	EC-3890MZi	ED-3490TK	EG-1870K
EC-3490TLi	EC-3890TLK	ED34-i10T	EG-3270UK
EC-3490TFi	EC-3890TFK	EG-1690K	EG-3670URK
EC-3490TMi	EC34-i10L	EG-2490K	EG-3870UTK
EC-3890Li	EC34-i10F	EG-2790i	ES-3870K
EC-3890Fi	EC34-i10M	EG-2790K	
EC-3890Mi	EC34-i10TL	EG-2990i	
EC-3890Fi2	EC34-i10TF	EG-2990K	

CS3010S Cleaning Brush for Endoscope Channel



Model name
CS3010S

Description

Cleaning brush for channel

Item number
82288

Compatible
endoscopes

EB-1170K	FB-10V	FI-9RBS	FI-13BS
EB-1990i	FHY-10RBS	FI-10BS	FI-13RBS
FB-8V	FI-9BS	FI-10RBS	FUR-9RBS

CS3025S Cleaning Brush for Endoscope Channel



Model name
CS3025S

Description

Cleaning brush for A/W and WJ channel

Item number
82287

Compatible
endoscopes

EE-1580K	FCP-9P
----------	--------

CS6002SN Cleaning Brush for Endoscope Channel



Model name
CS6002SN

Description
Cleaning brush for channel

Item number
82292

Compatible
endoscopes

EB-1170K	EB-1990i	FB-15RBS	FI-10BS
EB-1570AK	EB-1970UK	FB-18BS	FI-10RBS
EB-1570K	ECY-1570K	FB-18RBS	FI-13BS
EB-1572K	ECY-1575K	FCN-15X	FI-13RBS
EB-1575K	FB-8V	FCP-9P	FI-16BS
EB-1970AK	FB-10V	FCY-15RBS	FI-16RBS
EB-1970K	FB-15V	FHY-10RBS	FNL-15RP3
EB-1970TK	FB-18V	FHY-15RBS	FUR-9RBS
EB-1972K	FB-19TV	FI-9BS	VNL-1570STK
EB-1975K	FB-15BS	FI-9RBS	

CS6015ST Cleaning Brush for Endoscope Channel



Model name
CS6015ST

Description
Cleaning brush for channel

Item number
81454

Compatible
endoscopes

FG-16V	EB-1970K	FB-15V	FCN-15X
FG-24V	EB-1970TK	FB-18V	FCY-15RBS
EB-1570AK	EB-1972K	FB-19TV	FHY-15RBS
EB-1570K	EB-1975K	FB-15BS	FI-16BS
EB-1572K	EB-1970UK	FB-15RBS	FI-16RBS
EB-1575K	ECY-1570K	FB-18BS	FNL-15RP3
EB-1970AK	ECY-1575K	FB-18RBS	VNL-1570STK

CS6021SN Cleaning Brush for Endoscope Channel



Model name
CS6021SN

Description
Cleaning brush for working channel

Item number
82290

Compatible
endoscopes

EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
EC-380FKp	FG-290Kp	FD-34V2	FG-29V
EC-380FK2p	FC-38LV	FG-16V	FS-34V



CS6030SN Cleaning Brush for Endoscope Channel



Model name CS6030SN	Description Cleaning brush for working channel, extra long
Item number 82289	
Compatible endoscopes	VSB-2990i



CS-C3S Cleaning Brush for Endoscope Cylinder



Model name CS-C3S	Description Cleaning brush for cylinders or certain valves			
Item number 82293				
Compatible endoscopes	EB-1170K	EB-1990i	FB-15RBS	FI-10BS
	EB-1570AK	EB-1970UK	FB-18BS	FI-10RBS
	EB-1570K	ECY-1570K	FB-18RBS	FI-13BS
	EB-1572K	ECY-1575K	FCN-15X	FI-13RBS
	EB-1575K	FB-8V	FCP-9P	FI-16BS
	EB-1970AK	FB-10V	FCY-15RBS	FI-16RBS
	EB-1970K	FB-15V	FHY-10RBS	FNL-15RP3
	EB-1970TK	FB-18V	FHY-15RBS	FUR-9RBS
	EB-1972K	FB-19TV	FI-9BS	VNL-1570STK
	EB-1975K	FB-15BS	FI-9RBS	



CS-C5S Cleaning Brush for Endoscope Cylinder



Model name CS-C5S	Description Cleaning brush for cylinder			
Item number 82294				
Compatible endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	FG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V

CS-C9S Cleaning Brush for Endoscope Cylinder



Model name
CS-C9S

Description
Cleaning brush for valve cylinders, inlet/channel ports, scope tips/
elevator mechanisms and valves

Item number
87145

Compatible
endoscopes

EC-2990Li	EC-3890Mi2	EC34-i10TM	EG-2990Zi
EC-2990Fi	EC-3890LK	EC38-i10L	EG-3490K
EC-2990Mi	EC-3890FK	EC38-i10F	EG-3890TK
EC-3490Li	EC-3890MK	EC38-i10M	EG16-K10
EC-3490Fi	EC-3890FK2	EC38-i10F2	EG27-i10
EC-3490Mi	EC-3890MK2	EC38-i10M2	EG29-i10
EC-3490LK	EC-3890LZi	EC38-i10NL	EG34-i10
EC-3490FK	EC-3890FZi	EC38-i10NF	VSB-2990i
EC-3490MK	EC-3890MZi	ED-3490TK	EE-1580K
EC-3490TLi	EC-3890TLK	ED34-i10T	EG-1870K
EC-3490TFi	EC-3890TFK	EG-1690K	EG-3270UK
EC-3490TMi	EC34-i10L	EG-2490K	EG-3670URK
EC-3890Li	EC34-i10F	EG-2790i	EG-3870UTK
EC-3890Fi	EC34-i10M	EG-2790K	ES-3870K
EC-3890Mi	EC34-i10TL	EG-2990i	
EC-3890Fi2	EC34-i10TF	EG-2990K	



OF-B212 Cleaning Adapter



Model name
OF-B212

Description
Cleaning adapter

Item number
84588

Compatible
endoscopes

EB11-J10	EB19-J10	VNL15-J10
EB15-J10	EB19-J10U	

Fitting O-ring set is OF-B214

OF-G17 Air-Water Channel Cleaning Adapter



Model name
OF-G17

Description
A/W cleaning adapter

Item number
83189

Compatible
endoscopes

EC-2990Li	EC-3890LK	EC38-i10M	EG-2990Zi
EC-2990Fi	EC-3890FK	EC38-i10F2	EG-3490K
EC-2990Mi	EC-3890MK	EC38-i10M2	EG-3890TK
EC-3490Li	EC-3890FK2	EC38-i10NL	EG16-K10
EC-3490Fi	EC-3890MK2	EC38-i10NF	EG17-J10
EC-3490Mi	EC-3890LZi	EC38-i10cL	EG27-i10
EC-3490LK	EC-3890FZi	EC38-i10cF	EG29-i10
EC-3490FK	EC-3890MZi	EC38-i10cM	EG29-i10c
EC-3490MK	EC-3890TLK	EC38-i10cF2	EG34-i10
EC-3490TLi	EC-3890TFK	ED-3490TK	FC-38LV
EC-3490TFi	EC34-i10L	ED32-i10	FC-38FV
EC-3490TMi	EC34-i10F	ED34-i10T	FD-34V2
EC-380LKp	EC34-i10M	ED34-i10T2	FG-16V
EC-380FKp	EC34-i10TL	EE17-J10	FG-24V
EC-380FK2p	EC34-i10TF	EG-1690K	FG-29V
EC-380MK2p	EC34-i10TM	EG-2490K	FS-34V
EC-3890Li	EC34-i10cL	EG-2790i	VSB-2990i
EC-3890Fi	EC34-i10cF	EG-2790K	EG-1870K
EC-3890Mi	EC34-i10cM	EG-290Kp	EG-3270UK
EC-3890Fi2	EC38-i10L	EG-2990i	EG-3870UTK
EC-3890Mi2	EC38-i10F	EG-2990K	ES-3870K

OE-C20 Water Jet Channel Cleaning Adapter



Model name
OE-C20

Description
Water jet channel cleaning adapter

Item number
82054

Compatible
endoscopes

EC-2990Li	EC-3890Fi2	EC34-i10TL	EC38-i10cM
EC-2990Fi	EC-3890Mi2	EC34-i10TF	EC38-i10cF2
EC-2990Mi	EC-3890LK	EC34-i10TM	EG-2990i
EC-3490Li	EC-3890FK	EC34-i10cL	EG-2990K
EC-3490Fi	EC-3890MK	EC34-i10cF	EG-2990Zi
EC-3490Mi	EC-3890FK2	EC34-i10cM	EG-3490K
EC-3490LK	EC-3890MK2	EC38-i10L	EG-3890TK
EC-3490FK	EC-3890LZi	EC38-i10F	EG29-i10
EC-3490MK	EC-3890FZi	EC38-i10M	EG29-i10c
EC-3490TLi	EC-3890MZi	EC38-i10F2	EG34-i10
EC-3490TFi	EC-3890TLK	EC38-i10M2	VSB-2990i
EC-3490TMi	EC-3890TFK	EC38-i10NL	ES-3870K
EC-3890Li	EC34-i10L	EC38-i10NF	
EC-3890Fi	EC34-i10F	EC38-i10cL	
EC-3890Mi	EC34-i10M	EC38-i10cF	

OF-B153 Air-Water Suction Channel Cleaning Adapter



Model name
OF-B153

Description
Cleaning adapter

Item number
86153

Compatible
endoscopes

EC-2990Li	EC-3890FK	EC38-i10L	EG-2790i
EC-2990Fi	EC-3890MK	EC38-i10F	EG-2790K
EC-2990Mi	EC-3890FK2	EC38-i10M	EG-2990i
EC-3490Li	EC-3890MK2	EC38-i10F2	EG-2990K
EC-3490Fi	EC-3890LZi	EC38-i10M2	EG-2990Zi
EC-3490Mi	EC-3890FZi	EC38-i10NL	EG-3490K
EC-3490LK	EC-3890MZi	EC38-i10NF	EG-3890TK
EC-3490FK	EC-3890TLK	EC38-i10cL	EG16-K10
EC-3490MK	EC-3890TFK	EC38-i10cF	EG17-J10
EC-3490TLi	EC34-i10L	EC38-i10cM	EG27-i10
EC-3490TFi	EC34-i10F	EC38-i10cF2	EG29-i10
EC-3490TMi	EC34-i10M	ED-3490TK	EG29-i10c
EC-3890Li	EC34-i10TL	ED32-i10	EG34-i10
EC-3890Fi	EC34-i10TF	ED34-i10T	VSB-2990i
EC-3890Mi	EC34-i10TM	ED34-i10T2	EG-3670URK
EC-3890Fi2	EC34-i10cL	EE17-J10	
EC-3890Mi2	EC34-i10cF	EG-1690K	
EC-3890LK	EC34-i10cM	EG-2490K	



OF-B218 Cleaning Cap for Instrument Channel Inlet



Model name OF-B218	Description Cleaning cap for instrument channel inlet			
Item number 8000				
Compatible endoscopes	EB11-J10	EB19-J10	EG34-J10U	EG38-J10UT
	EB15-J10	EB19-J10U	EG36-J10UR	VNL15-J10



OE-B11 PENTAX Medical Cleaning Adapter



Model name OE-B11	Description Cleaning adapter			
Item number 84630				
Compatible endoscopes	EB19-J10U			



OE-B10 Cleaning Adapter



Model name OE-B10	Description Cleaning adapter			
Item number 84832				
Compatible endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT	

Fitting O-ring set is OF-B195



OE-G1 Cleaning Adapter




Model name OE-G1	Description Cleaning adapter			
Item number 84831				
Compatible endoscopes	EG34-J10U	EG36-J10UR	EG38-J10UT	

OE-C9 PVE Soaking Cap



Model name	Description			
OE-C9	Soaking cap			
Item number				
86643				
Compatible endoscopes	EC-2990Li	EC-3890LZi	EE17-J10	EB-1975K
	EC-2990Fi	EC-3890FZi	EG-1690K	EB-1990i
	EC-2990Mi	EC-3890MZi	EG-2490K	EB11-J10
	EC-3490Li	EC-3890TLK	EG-2790i	EB15-J10
	EC-3490Fi	EC-3890TFK	EG-2790K	EB19-J10
	EC-3490Mi	EC34-i10L	EG-290Kp	EB-1970UK
	EC-3490LK	EC34-i10F	EG-2990i	EB19-J10U
	EC-3490FK	EC34-i10M	EG-2990K	ECY-1575K
	EC-3490MK	EC34-i10TL	EG-2990Zi	EE-1580K
	EC-3490TLi	EC34-i10TF	EG-3490K	EG-1870K
	EC-3490TFi	EC34-i10TM	EG-3890TK	EG-3270UK
	EC-3490TMi	EC38-i10L	EG16-K10	EG-3670URK
	EC-380LKp	EC38-i10F	EG17-J10	EG-3870UTK
	EC-380FKp	EC38-i10M	EG27-i10	EG34-J10U
	EC-380FK2p	EC38-i10F2	EG29-i10	EG36-J10UR
	EC-380MK2p	EC38-i10M2	EG34-i10	EG38-J10UT
	EC-3890Li	EC38-i10NL	VSB-2990i	ES-3870K
	EC-3890Fi	EC38-i10NF	EB-1170K	VNL-1070STK
	EC-3890Mi	EC38-i10cL	EB-1570AK	VNL-1190STK
	EC-3890Fi2	EC38-i10cF	EB-1570K	VNL-1570STK
	EC-3890Mi2	EC38-i10cM	EB-1572K	VNL-1590STi
	EC-3890LK	EC38-i10cF2	EB-1575K	VNL8-J10
	EC-3890FK	ED-3490TK	EB-1970AK	VNL11-J10
	EC-3890MK	ED32-i10	EB-1970K	VNL15-J10
	EC-3890FK2	ED34-i10T	EB-1970TK	PVK-1070Z
	EC-3890MK2	ED34-i10T2	EB-1972K	

OE-C2 PVE Soaking Cap



Model name	Description
OE-C2	Soaking cap
Item number	
81461	

OE-C27 Soaking Cap for Navigation Connector



Model name	Description
OE-C27	Soaking cap for navigation connector
Item number	
84650	
Compatible endoscopes	EC38-i10NL EC-38-i10NF

OE-U1 Soaking Cap for Scanning Unit Connector



Model name OE-U1	Description Soaking cap for scanning unit connector			
Item number 87080				
Compatible endoscopes	EB-1970UK	EG-3270UK	EG-3670URK	EG-3870UTK

OE-U4 Soaking Cap for PENTAX Medical Ultrasound Endoscopes



Model name OE-U4	Description Soaking cap for PENTAX Medical ultrasound endoscopes			
Item number 84631				
Compatible endoscopes	EB19-J10U	EG34-J10U	EG36-J10UR	EG38-J10UT

OF-B162 Cleaning Cap for 2ch



Model name OF-B162	Description Cleaning cap for 2ch			
Item number 82003				
Compatible endoscopes	EC-3890TLK	EC-3890TFK	EG-3870TK	



OF-B115 Suction & Air-Water Channel Cleaning Adapter



Model name OF-B115	Description Cleaning cap for A/W/suction channel			
Item number 81781				
Compatible endoscopes	EC-380LKp	EC-380MK2p	FC-38FV	FG-24V
	EC-380FKp	EG-290Kp	FD-34V2	FG-29V
	EC-380FK2p	FC-38LV	FG-16V	FS-34V



OF-B166 Air-Water-Suction Channel Cleaning Adapter



Model name OF-B166	Description Cleaning adapter		
Item number 86463			
Compatible endoscopes	EE-1580K	EG-1870K	ES-3870K

OE-C8 Suction Nipple Cleaning Adapter



Model name OE-C8	Description Cleaning adapter suction nipple		
Item number 86150			
Compatible endoscopes	EG-1870K	ES-3870K	

OF-B170 Air Cleaning Brush Adapter



Model name OF-B170	Description Cleaning brush adapter for A/W channel EG-1580K		
Item number 86546			
Compatible endoscopes	EE-1580K		

OE-C5 Air-Water & Suction Channel Adapter



Model name OE-C5	Description Cleaning adapter for A/W channel		
Item number 82148			
Compatible endoscopes	EE-1580K	EG-3670URK	



OF-B103 Suction Channel Cleaning & Disinfection Adapter



Model name	Description			
OF-B103	Suction channel cleaning & disinfection adapter			
Item number	81411			
Compatible endoscopes	FB-15BS	FCY-15RBS	FI-10BS	FI-16RBS
	FB-15RBS	FHY-10RBS	FI-10RBS	FNL-15RP3
	FB-18BS	FHY-15RBS	FI-13BS	FUR-9RBS
	FB-18RBS	FI-9BS	FI-13RBS	
	FCP-9P	FI-9RBS	FI-16BS	



BS-B1 Soaking Cap for BS



Model name	Description			
BS-B1	Soaking cap for BS series			
Item number	81953			
Compatible endoscopes	FB-15BS	FHY-10RBS	FI-9RBS	FI-16BS
	FB-15RBS	FHY-15RBS	FI-10BS	FI-16RBS
	FB-18BS	FI-7BS	FI-10RBS	FNL-10RBS
	FB-18RBS	FI-7RBS	FI-13BS	FUR-9RBS
	FCY-15RBS	FI-9BS	FI-13RBS	



OF-B155 Cleaning Adapter



Model name	Description			
OF-B155	Cleaning cap for suction channel			
Item number	86148			
Compatible endoscopes	EB-1170K	EB-1575K	EB-1972K	ECY-1575K
	EB-1570AK	EB-1970AK	EB-1975K	VNL-1570STK
	EB-1570K	EB-1970K	EB-1990i	
	EB-1572K	EB-1970TK	EB-1970UK	



OF-B75 Suction Channel Cleaning Adapter



Model name	Description		
OF-B75	Cleaning cap for suction channel		
Item number	83127		
Compatible endoscopes	FB-8V	FB-15V	FB-19TV
	FB-10V	FB-18V	FCN-15X

SHA-P6 Leakage Tester



Model name SHA-P6	Description Leakage tester
Item number 90041	

GS00003 Cleaning Kit

Picture not available

Model name
GS00003**Description**
Cleaning kit for VERSA endoscopes**Item number**
300955Compatible
endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GF00089 Cleaning Plug

Picture not available

Model name
GF00089**Description**
Cleaning plug for VERSA endoscopes**Item number**
300953Compatible
endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GF00027 Leakage Tester

Picture not available

Model name
GF00027**Description**
Leakage tester for VERSA endoscopes**Item number**
300950Compatible
endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GM00906 Water Resistant Cap

Picture not available

Model Name
GM00906**Description**
Water resistant cap for VERSA endoscopes**Item number**
300952Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

GS00017 Plug Label

Picture not available

Model Name
GS00017**Description**
Plug label for VERSA endoscopes**Item number**
300954Compatible
Endoscopes

EG27-V10c

EG29-V10c

EC38-V10cM

EC38-V10cL

Manufacturers

VERSA Consumables

 **PENTAX-AOHUA Medical Technologies Co., Ltd.**
East of 3rd Floor, Block C, Building 1,
No.5 Shenwang Road,
Minhang District, Shanghai 201108,
P. R. China.

All other products

 **HOYA Corporation**
6-10-1, Nishi-shinjuku
Shinjuku-ku, Tokyo
160-0023 Japan

 **EMEA Headquarter
Germany**

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg, Germany



Medical Device Class: I

 **EMEA Headquarter
Germany**

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg, Germany



0123
TÜV SÜD CE 0123 Medical Device Classes: I, Is, IIa and IIb

Distributor

EMEA Headquarter Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

This product must be used only by healthcare professionals.

Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

Product availability cannot be guaranteed in all countries. For further information please contact your local sales representative.

PENTAX
MEDICAL