


Test Report

– Eppendorf PCR clean –

Safe-Lock Tubes 2.0 mL



1000 EPPENDORF TUBES®
PCR clean
centrifugation up to 25 000 x g

LOT O2231840

 **2030-09-28**

Order no.: 0030 123.344
Eppendorf SE
22331 Hamburg · Germany
Email: eppendorf@eppendorf.com

• **Human DNA**

Method: *real-time PCR acc. to SOP 99880/V8*
Limit: *< 2 pg; less than one human cell*

• **DNase**

Method: *DNA Digestion acc. to SOP 100570/V4*
Limit: *Not detectable (LOD: 1.0×10^{-7} Kunitz units)*

• **PCR inhibition**

Method: *real-time PCR acc. to SOP 100559/V4*
Limit: *Less than 10 targets amplifiable*

• **RNase**

Method: *RNA Digestion acc. to SOP 100571/V4*
Limit: *Not detectable (LOD: 1.0×10^{-9} Kunitz units)*

Customer: Eppendorf SE, Barkhausenweg 1, D-22339 Hamburg

Date of sample receipt / examination: *22.09.25*

Testing of the rinse solution of the above described lot showed conformity within the limits of detection.

22.09.25
Date of release

Epp 2509-6470
Registration number

S. V. P. Kuhnert
Dr. rer. nat. Sybille Kuhnert (Head)

The test results refer exclusively to the items tested as arrived. The test report must not be copied partially without the approval of LADR GmbH.

page 1 of 1