



RACING FOR LIFE | PRODUCT CATALOG

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City, Huangpu District,
510700, Guangzhou, P.R. China
Tel: (+86) 400-830-8768
E-mail: sales@wondfo.com.cn
Website: en.wondfo.com



© 2024 Guangzhou Wondfo Biotech Co., Ltd.
All rights reserved Rev. A2 Rel: 2024/07/31

Wondfo

Stock Code: 300482.SZ

WONDFO PROFILE

Founded in 1992, headquartered in Guangzhou science city, Wondfo focuses on the R&D, manufacturing, sales and service of POCT (Point-of-care testing) industry, providing professional rapid diagnosis and chronic disease management solutions.

Wondfo has 9 technology platforms. Our rapid diagnostic devices and reagents cover a wide range, including cardiovascular diseases, diabetes, renal functions, respiratory, fertility, DOA, infections, tumor etc. Our products have been distributed to over 150 countries and regions.

In 2015, Wondfo was listed on Shenzhen Stock Exchange, becoming the first public POCT company in China.



5
R&D Centers



30+
Years in POCT



3000+
Employees



WONDFO PRODUCT LINES

Dry Chemistry Platform

- Organophosphates Pesticides Poisoning ;
- Trauma / Anemia ;
- Blood Lipid ;
- Liver Function ;
- Pancreatitis ;
- Blood Glucose ;
- Myocardial Enzyme ;
- Kidney Function ;
- Alcoholism

Immunofluorescence Platform

- Diabetes & Kidney Disease ;
- Craniocerebral Injury ;
- Thyroid Function ;
- Cardiovascular ;
- Coagulation ;
- Inflammation ;
- Fertility ;
- Tumor ;
- Others

Molecular Diagnostics Platform


- Pathogenic Microorganism Nucleic Acid Detection ;
- Tumor Gene Detection ;
- Circulating Tumor Cell Detection

Electrochemistry Platform

- Coagulation Function ;
- Metabolites ;
- Blood Gas ;
- Electrolytes ;
- Hematology

Chemiluminescence Platform

- Inflammation ;
- Infectious Disease ;
- Anemia ;
- Bone metabolism ;
- Allergy ;
- Tumor ;
- Fertility ;
- Cardiovascular ;
- Thyroid Function



Pathology Platform

- Immunohistochemistry ;
- In Situ Hybridization ;
- Automatic Pathological Staining

Colloidal Gold Platform

- Respiratory Tract ;
- Infectious Disease Screening ;
- Tropical Infectious Diseases ;
- Fertility ;
- STDs ;
- Digestive Tract ;
- Cardiac Markers

Trademark

610+
Registered Trademarks

Patent

420+
Valid Patents

Certificate

280+ NMPA	240+ CE
70+ CLIA WAIVED	13 MDALL
2 WHO PQ	



Cardiovascular



Diabetes



Respiratory



Renal Injury

8

MAJOR FIELDS



Fertility



DOA



Infection



Tumor

PRODUCT PORTFOLIO



Immunofluorescence




Colloidal Gold



Coagulation



Chemiluminescence



Blood Gas



Molecular



Dry Chemistry



Immunohistochemistry



CONTENTS

Qualitative Rapid Diagnostic Family 09

Fertility Tests 10

HCG Pregnancy Test
LH Ovulation Test
FSH Urine Test
Sperm Concentration Test

Infectious Disease Tests 11

HBsAg Hepatitis B Surface Antigen Test
HBsAb Hepatitis B Surface Antibody Test
HBV Hepatitis B Virus Test
HCV Hepatitis C Virus Test
HAV Hepatitis A Virus Test
HIV 1/2 Test (2 lines)
HIV 1/2 Test (3 Lines)
HIV 1/2 Oral Mucosal Transudate Test
HIV-HCV-HBsAg-Syphilis Test
HIV-HCV-HBsAg Test
Syphilis Test
Chlamydia Test
Gonorrhea Test
Malaria P.f Test
Malaria P.f/Pan Test
Malaria P.f/P.v Test
H. Pylori Antigen Test
H. Pylori Antibody Test
Adenovirus Antigen Test
Rotavirus Antigen Test
Combined Adenovirus Antigen and Rotavirus Antigen Test
EV71 IgM Antibody Test
Flu A&B Test
Strep A Test
MP IgM Antibody Test
CP IgM Antibody Test
Dengue NS1 Antigen Test
Dengue IgG/IgM Antibody Test
Combined Dengue NS1 Antigen and Dengue IgG/IgM Antibody Test

Drug of Abuse Tests 12

Multi-Drug Tests
DOA Tests (Urine)
DOA Tests (Oral Fluid)

Tumor Marker Tests 14

PSA Prostate Specific Antigen Test
FOB Fecal Occult Blood Test

Other Rapid Tests 14

Troponin I Test

Quantitative Detection Family 15

Electrochemistry Platform 15

Optical Coagulation Analyzer (OCG-102)
Thrombelastography Analyzer (TLT-101)
Blood Gas Analyzer (BGA-102)
Automatic Blood Gas Analyzer (Ucare-6000)

Immunofluorescence Platform 20

Finecare™ platform 21

Finecare™ FIA Meter Plus (FS-113)
Finecare™ FIA Meter II Plus SE (FS-114)
Finecare™ FIA Meter III Plus (FS-205)
Finecare™ FIA Meter Automatic (FS-301)

New Finecare™ platform 25

Finecare™ FIA Meter X1 (FS-1000)
Finecare™ FIA Meter X2 (FS-2000)
Finecare™ FIA Meter X2 PLUS (FS-2100)

Mono-test Chemiluminescence 28

ACCRE 8
ACCRE 90
ACCRE 160

Fully-auto Chemiluminescence 33

SmarLumi™ FC-2000/2100

Molecular Diagnostics Platform 35

Molecular Diagnostic Solutions
Real-Time PCR Assay Test Menu
U-Card Dx™ Molecular Diagnostic System

Immunohistochemistry Platform 40

PA-3600 Fully Automated IHC Staining Instrument

Dry Biochemistry Platform 45

Dry Chemistry Analyzer (DC-101)

Wondfo COVID-19 diagnostic solution 48

Wondfo 9 Technology Platforms

Colloidal Gold

Immunofluorescence

Electrochemistry

Chemiluminescence

Molecular Diagnostics

Dry Biochemistry

Instrument

Biological Raw Material

Pathology



Qualitative Rapid Diagnostic Family

Qualitative Rapid Diagnostic Family

Fertility Tests



Sperm Concentration Test Kit



Pregnancy & Ovulation tests

- Other colors of Midstream caps are available.
- Reliable original equipment manufacturer, high-quality service is provided.

Product	Cat. No.	Specimen	Width	Format	Sensitivity	Qualification
	W1-S	Urine	3.5mm	Strip	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W1-S(20mIU)	Urine	3.5mm	Strip	20mIU/mL	CE ₀₁₂₃
	W1-S(10mIU)	Urine	3.5mm	Strip	10mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W001-S1HJ(20mIU)1	Urine	2.5mm	Strip	20mIU/mL	CE ₀₁₂₃
	W1-C	Urine	4.0mm	Cassette	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W1-CI	Urine	3.5mm	Cassette		
	W1-C(20mIU)	Urine	4.0mm	Cassette	20mIU/mL	CE ₀₁₂₃
	W1-CI(20mIU)	Urine	3.5mm	Cassette		
	W1-C(10mIU)	Urine	4.0mm	Cassette		
	W1-CI(10mIU)	Urine	3.5mm	Cassette	10mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W1-M	Urine	6.0mm	Midstream	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W1-MII	Urine	5.0mm	Midstream		
	W1-MIII	Urine	6.0mm	Midstream		
	W1-M4	Urine	3.5mm	Midstream		
HCG Test	W1-M6	Urine	6.0mm	Midstream	20mIU/mL	CE ₀₁₂₃
	W1-M(20mIU)	Urine	6.0mm	Midstream		
	W1-MII(20mIU)	Urine	5.0mm	Midstream		
	W1-MIII(20mIU)	Urine	6.0mm	Midstream		
	W1-M4(20mIU)	Urine	3.5mm	Midstream		
	W1-M6(20mIU)	Urine	6.0mm	Midstream		
	W1-M(10mIU)	Urine	6.0mm	Midstream	10mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W1-MII(10mIU)	Urine	5.0mm	Midstream		
	W1-MIII(10mIU)	Urine	6.0mm	Midstream		
	W1-M4(10mIU)	Urine	3.5mm	Midstream		
W1-M6(10mIU)	Urine	6.0mm	Midstream	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED	
W1-S3P	Urine/ Serum	3.5mm	Strip			
W001P0720	Urine/ Serum	3.5mm	Strip			
W1-C3P	Urine/ Serum	4.0mm	Cassette	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED	
W001P0721	Urine/ Serum	4.0mm	Cassette			
LH Test	W2-S	Urine	3.5mm	Strip	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W2-C	Urine	4.0mm	Cassette	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W2-CI	Urine	3.5mm	Cassette		
	W2-M	Urine	6.0mm	Midstream	25mIU/mL	CE ₀₁₂₃ , CLIA WAIVED, NMPA
	W2-MII	Urine	5.0mm	Midstream		
	W2-M4	Urine	3.5mm	Midstream		
W2-M6	Urine	6.0mm	Midstream			
FSH Test	W52-S	Urine	3.5mm	Strip	25mIU/mL	CE ₀₁₂₃ , NMPA
	W52-C	Urine	4.0mm	Cassette	25mIU/mL	CE ₀₁₂₃ , NMPA
	W52-M	Urine	6.0mm	Midstream	25mIU/mL	CE ₀₁₂₃ , NMPA
Sperm Concentration Test	W121-C	Semen	4.0mm	Cassette	15x10 ⁶ /mL	CE ₀₁₂₃ , NMPA

Infectious Disease Tests



Category	Cat. No.	Test Item	Specimen	Format	Sensitivity	Qualification
Hepatitis	W003	HBsAg	WB/S/P	Cassette/Strip	5.0 IU/mL	NMPA
	W004	HBsAb	WB/S/P	Cassette/Strip	30 mIU/mL	NMPA
	W040-P	HBV(HBsAg, HBsAb, HBeAg, HBeAb, HBcAb)	WB/S/P	Panel	/	NMPA
	W005	HCV	WB/S/P	Cassette/Strip	2 NCU/mL	NMPA
	W071	HAV-IgM	S/P	Cassette	/	CE
HIV	W752	HIV 1/2 (2 lines)	WB/S/P	Cassette	/	CE; NMPA
	W006	HIV 1/2 (2 lines)	WB/S/P	Cassette/Strip	/	NMPA
	W006-P	HIV 1/2 (2 lines)	WB/S/P	Cassette	/	WHO Pre-qualified; NMPA
	W006P00XX (58~60)	HIV 1/2 (2 lines) (Self-Test)	Fingerstick WB	Cassette	/	WHO Pre-qualified; NMPA
	W050-C	HIV 1/2 (3 lines)	WB/S/P	Cassette	/	
	W063-C	HIV 1/2 Oral Mucosal Transudate	Oral Mucosal Transudate	Cassette	/	NMPA
	W026-C	HIV-HCV-HBsAg-Syphilis	WB/S/P	Cassette	/	NMPA
	W094-C	HIV-HCV-HBsAg	WB/S/P	Cassette	/	NMPA
	W034	Syphilis	WB/S/P	Cassette/Strip	2 NCU/mL	CE; NMPA
	W35-C	Chlamydia	Female Endocervical Swab/ Male Penile Swab/Male Urine	Cassette	5x10 ⁴ IFU/mL	CE; NMPA
W036	Gonorrhea	Female Endocervical Swab/ Male Penile Swab	Cassette	1x10 ⁵ CFU/mL	CE; NMPA	
Malaria	W037	Malaria HRP2 (P.f)	WB	Cassette/Strip	100 p/μL	CE; NMPA
	W054	Malaria HRP2/pLDH (P.f/Pan)	WB	Cassette	100 p/μL	CE; NMPA
	W056	Malaria HRP2/pLDH (P.f/P.v)	WB	Cassette	100 p/μL	CE
Respiratory Infection	W059	Flu A&B	Nasal Swab	Cassette	Flu A: 5.1x10 ⁵ TCID50/mL Flu B: 1.5x10 ⁶ TCID50/mL	CE; NMPA
	W059	Flu A&B	Nasal Swab/ Throat Swab	Cassette	Flu A: 5.1x10 ⁵ TCID50/mL Flu B: 1.5x10 ⁶ TCID50/mL	CE; NMPA
	W71-CH	Flu A&B (Self-Test)	Nasal Swab/ Throat Swab	Cassette	Flu A: 5.1x10 ⁵ TCID50/mL Flu B: 1.5x10 ⁶ TCID50/mL	CE; NMPA
	W039	Strep A	Swab	Cassette	1.5x10 ⁵ CFU/mL	CE; CLIA WAIVED
	W039	Strep A (Self-Test)	Swab	Cassette	1.5x10 ⁵ CFU/mL	CE
	W113	MP IgM	WB/S/P	Cassette	/	CE; NMPA
	W135	CP IgM	WB/S/P	Cassette	/	CE; NMPA
	W771	Flu A&B/2019-nCoV/RSV/ADV/MP	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W770	Flu A&B/2019-nCoV/RSV/ADV	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W769	Flu A&B/RSV/ADV/MP	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W766	RSV/ADV/MP	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W773	RSV Ag	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W774	ADV Ag	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	W772	MP Ag	Nasopharyngeal Swab/ Oropharyngeal Swab	Cassette	/	CE
	Enteric Infection	W028	H. Pylori Antigen	Feces	Cassette/Strip	/
W038		H. Pylori Antibody	S/P	Cassette/Strip	/	CE
W133		ADV	Feces	Cassette	/	CE; NMPA
W112		RV	Feces	Cassette	/	CE; NMPA
W134		ADV-RV	Feces	Cassette	/	CE; NMPA
W141		EV71	WB/S/P	Cassette	/	CE; NMPA

Category	Cat. No.	Test Item	Specimen	Format	Sensitivity	Qualification
Enteric Infection	W158	EV 71-CA16 IgM	WB/S/P	Cassette	/	CE; NMPA
	W187	Norovirus Ag	Feces	Cassette	/	CE
Dengue	W041	Dengue NS1 Antigen	WB/S/P	Cassette	/	CE; NMPA
	W109	Dengue IgG/IgM Antibody	WB/S/P	Cassette	/	CE; NMPA
	W111	Dengue NS1/IgG/IgM	WB/S/P	Cassette	/	CE; NMPA

DOA Urine Tests



Cassette



Strip



Panel



T-Cup



Q-Cup



T-Cube

Multi-Drug Tests

Product	Cat. No.	Format	Test items	Specimen	Sensitivity	Qualification
Multi-Drug Urine Test	W20XX(2-16)-P	Panel	2-16 Items	Urine	Variable	CE; CLIA WAIVED
	W5XX(2-18)-CU2	T-Cup	2-18 Items	Urine	Variable	CE; CLIA WAIVED
	W5XX(2-18)-CU3	Q-Cup	2-18 Items	Urine	Variable	CE
Multi-Drug Oral Fluid	W6XX(2-10)-CU6	T-Cube	2-10 Items	Oral Fluid	Variable	CE
	W605P-W606P	T-Swab	5-6 Items	Oral Fluid	Variable	CE
Multi-Drug Hair Test	W602P-W612P	T-Square	2-12 Items	Oral Fluid	Variable	CE
	W5XX(62-72)P0001	T-Cube	2-12 Items	Hair	Variable	CE
Marijuana (THC) Multi-Level Urine Drug Test	W775P0001	Panel	THC15/50/100/200/300	Urine	Variable	CE
Marijuana(THC) Hair Test	W628P0001	T-Cube	THC100/200	Hair	Variable	CE

Drug of Abuse Tests (Urine)

Cat. No.	Product	Test Item	Specimen	Format	Sensitivity	Qualification
W007	AMP	Amphetamine	Urine	Cassette/Panel/Cup	1000/500 ng/mL	CE
W008	BAR	Barbiturates	Urine	Cassette/Panel/Cup	300 ng/mL	CE
W009	BZO	Benzodiazepines	Urine	Cassette/Panel/Cup	300 ng/mL	CE; CLIA WAIVED
W010	COC	Cocaine	Urine	Cassette/Panel/Cup	300/150 ng/mL	CE; CLIA WAIVED; NMPA
W011	mAMP/MET	Methamphetamine	Urine	Cassette/Panel/Cup	1000/500 ng/mL	CE; CLIA WAIVED; NMPA
W012	MDMA	Methylenedioxymethamphetamine; Ecstasy	Urine	Cassette/Panel/Cup	500 ng/mL	CE; NMPA
W013	BUP	Buprenorphine	Urine	Cassette/Panel/Cup	10 ng/mL	CE; CLIA WAIVED
W014	MOP	Morphine	Urine	Cassette/Panel/Cup	300 ng/mL	CE; NMPA
W015	MTD	Methadone	Urine	Cassette/Panel/Cup	300 ng/mL	CE; CLIA WAIVED
W016	OPI	Opiate	Urine	Cassette/Panel/Cup	2000 ng/mL	CE; CLIA WAIVED
W017	PCP	Phencyclidine	Urine	Cassette/Panel/Cup	25 ng/mL	CE; CLIA WAIVED
W018	TCA	Tricyclic Antidepressants	Urine	Panel/Cup	1000 ng/mL	CE
W019	THC	Marijuana	Urine	Cassette/Panel/Cup	50 ng/mL	CE; CLIA WAIVED
W021	OXY	Oxycodone	Urine	Cassette/Panel/Cup	100 ng/mL	CE; CLIA WAIVED
W055	KET	Ketamine	Urine	Cassette/Panel/Cup	1000 ng/mL	CE; NMPA
W057	PPX	Propoxyphene	Urine	Cassette/Panel/Cup	300 ng/mL	CE
W074	EDDP	Ethylenediamine-dimethylphosphinic	Urine	Cassette/Panel/Cup	100 ng/mL	CE
W086	TRA	Tramadol	Urine	Cassette/Panel/Cup	200 ng/mL	CE
W087	FTY	Fentanyl	Urine	Cassette/Panel/Cup	20 ng/mL	CE
W089	K2	Synthetic Cannabis	Urine	Panel/Cup	50 ng/mL	CE
W139	COT	Cotinine	Urine	Cassette/Panel/Cup	100 ng/mL	CE
W140	ETG	Ethyl Glucuronide	Urine	Panel/Cup	500 ng/mL	CE

DOA Tests (Oral Fluid)

Product	Test Item	Specimen	Format	Sensitivity	Qualification
AMP	Amphetamine	Oral Fluid	T-Cube	50 ng/mL	CE
BAR	Barbiturates	Oral Fluid	T-Cube	60 ng/mL	CE
BUP	Buprenorphine	Oral Fluid	T-Cube	5 ng/mL	CE
BZO	Benzodiazepines	Oral Fluid	T-Cube	30 ng/mL	CE
COC	Cocaine	Oral Fluid	T-Cube	20 ng/mL	CE
EDDP	Ethylenediamine-dimethylphosphinic	Oral Fluid	T-Cube	20 ng/mL	CE
FTY	Fentanyl	Oral Fluid	T-Cube	30 ng/mL	CE
MDMA	Methylenedioxymethamphetamine; Ecstasy	Oral Fluid	T-Cube	100 ng/mL	CE
MET	Methamphetamine	Oral Fluid	T-Cube	50 ng/mL	CE
MTD	Methadone	Oral Fluid	T-Cube	30 ng/mL	CE
OPI	Opiate	Oral Fluid	T-Cube	40 ng/mL	CE
PCP	Phencyclidine	Oral Fluid	T-Cube	10 ng/mL	CE
THC	Marijuana	Oral Fluid	T-Cube	25 ng/mL	CE
OXY	Oxycodone	Oral Fluid	T-Cube	20 ng/mL	CE

Tumor Marker Tests

Cat. No.	Test Item	Specimen	Format	Sensitivity	Qualification
W030	FOB	Feces	Cassette/Strip	200 ng/mL	CE; NMPA
W030	FOB (Self-Test)	Feces	Cassette/Strip	50ng/mL	CE; CLIA WAIVED 0123
W31-C	PSA	S/P	Cassette	/	CE 0123

Other Rapid Tests

Cat. No.	Product	Format	Specimen	Sensitivity	Qualification
W046	Troponin I	Cassette	WB/S/P	1 ng/mL	CE; NMPA

Electrochemistry Platform

Optical Coagulation Analyzer

OCG-102

PT,APTT,TT,FIB,ACT,PT/INR

- **Sample Type:** Venous/Capillary whole blood
- **Specimen volume:** 20µL
- **Long shelf life:** 18 months
- **Precision:** CV≤5%
- **Room-temperature storage:** 2-30°C
- **Portable:** Rechargeable lithium battery



Cat. No.	Item	Category	Specification	Qualification
W955	Optical Coagulation Analyzer	Instrument	/	CE NMPA
W480-484	PT,APTT,TT,FIB,ACT,PT/INR	Test strip	24 T/Kit	CE NMPA

Test Items

Cat. No.	Item	Specimen Type	Test strip Volume	Measuring Range	Clinical Reference
W480	PT	Citrated WB/Capillary WB	20 µL	7~90 s	10~14 s
W481	APTT	Citrated WB	20 µL	18~120 s	22~38 s
W482	FIB	Citrated WB	20 µL	1.5~6 g/L	2~4 g/L
W483	ACT	Citrated WB	20 µL	50~600 s	80~140 s
W483P0007	ACT (Non-anticoagulant)	Fresh WB	20 µL	50~600s	80~140s
W484	TT	Citrated WB	20µL	10~90 s	14~20 s
W486	PT/INR	Citrated WB/Capillary WB	20µL	7~90s/0.5-9	10~14s/0.7-1.3

* Note: INR is calculated by PT.

Thrombelastography Analyzer

TLT-101

Fast, Accurate

- Reliable results can be obtained within 20mins~2hours

Easy to use

- Integrated large touch screen
- Linux operating system
- Immediate, intuitive reading of results

High throughput

- Four independent channels
- Compatible with various cup types

Intelligent operation


- Automatic loading
- Automatic lifting
- Intelligent cover removal

Patented design

- Small size and light weight
- Six technical modules



Perioperative Bleeding Management



Guidance for Blood Transfusion



Monitoring of antiplatelet therapy



Detection of heparin therapy

Test Items	Reportable Parameters
Thrombelastography Kaolin Test Kit (Coagulation Method)	R, K, α , MA, LY30, EPL
Thrombelastography Rapid Test Kit (Coagulation Method)	ACT, K, α , MA, LY30, EPL
Thrombelastography Platelet Function-AA pathway Test Kit (Coagulation Method)	AA% inhibition
Thrombelastography Platelet Function-ADP pathway Test Kit (Coagulation Method)	ADP% inhibition
Thrombelastography Platelet Function-AA&ADP pathways Test Kit (Coagulation Method)	AA&ADP% inhibition
Thrombelastography Heparinase Test Kit (Coagulation Method)	R'=R(Kaolin Cup)-R (Heparinase Cup)
Thrombelastography Quality Control	R, K, α , MA

Blood Gas Analyzer

BGA-102

- Specimen type:** Aeterial blood/Capillary blood
- Specimen volume:** 80 μ L
- Time to result:** 180s
- Portable:** Built-in rechargeable battery

Suitable for Syringe and Arterial/Capillary Blood Sampler
 Good Correlation to Reference Method: $R^2 > 0.95$
 Test cartridge is stable and effective up to 15 months
 Reagent pack can be used for up to 60 days after installation



Testing Parameters

pH, pCO₂, pO₂, K⁺, Na⁺, Ca⁺⁺, Cl⁻, Hct, Glu, Lac

Calculating Parameters

cH⁺, cH⁺(T), pH(T), pCO₂(T), pO₂(T), HCO₃-atd, HCO₃-std, BB(B), BE(B), BE(ecf), ct CO₂, Ca⁺⁺(7.4), AnGap, tHb(est), sO₂(est), pO₂(A-a), pO₂(A-a)(T), pO₂(a/A), pO₂(a/A)(T), RI, RI(T), pO₂/FI O₂, pO₂(T)/FI O₂, POP

Cat. No.	Test Item	Specification	Qualification
W459	pH, pCO ₂ , pO ₂	25T/Kit	CE NMPA
W461	K ⁺ , Na ⁺ , Ca ⁺⁺ , Cl ⁻	25T/Kit	CE NMPA
W465	pH, pCO ₂ , pO ₂ , K ⁺ , Na ⁺ , Cl ⁻ , Ca ⁺⁺ , Hct	25T/Kit	CE NMPA
W466	pH, pCO ₂ , pO ₂ , K ⁺ , Na ⁺ , Cl ⁻ , Ca ⁺⁺ , Hct, Glu, Lac	25T/Kit	CE NMPA

Automatic Blood Gas Analyzer

Ucare-6000



Various specifications
75/150/300-test pack



Room temperature storage
Store at 2-30°C, NO refrigeration required



Long shelf life
7 months storage or 30 days on-board life



Built-in battery
Long standby for 6 hours WITHOUT power supply



One-step blood gas analysis



- All-in-one design of solution pack
- One-step sample loading and code scanning
- Hands-free automatic sample aspiration

10 Testing Parameters :

pH, pCO₂, pO₂, K⁺, Na⁺, Cl⁻, Ca⁺⁺, Glu, Lac, Hct

24 Calculating Parameters :

cH⁺(T), pH(T), pCO₂(T), pO₂(T), pO₂(A-a)(T), pO₂(a/A)(T), RI(T), pO₂(T)/FiO₂, cH⁺, Ca⁺⁺(7.4), HCO₃⁻act, HCO₃⁻std, BE(B), BB(B), BE(ecf), AG, sO₂(est), tHb(est), ctCO₂, pO₂(A-a), pO₂(a/A), RI, POP, pO₂/FiO₂



Ucare-6000 Solution Pack



Immunofluorescence Platform

Immunofluorescence Quantitative Test System

Finecare™ Platform



Finecare™ FIA Meter series is a fluorescence immunochromatographic analysing system which can help diagnose conditions such as infection, diabetes, cardiovascular diseases, renal injury and cancers, etc.

- Reliable**
 - CE/NMPA/NGSP/IFCC Certified
 - Fast and Accurate
 - Full Traceability
 - QC Available
- Room Temperature**
 - Store & Operate at 4~30°C
 - 24 Months Shelf Life
- Competitive**
 - Full Test Line
 - Result in 3~15 mins
 - Cost-efficient
 - Built-in Temperature Control
- Date Management and Transmission**
 - HIS/LIS
- Easy-to-use**
 - LCD Touch Screen
 - Multiple Language Available
 - Internal Thermal Printer
- Application**
 - Laboratory, ER, Cardiology, ICU, Respiratory, Pediatrics, etc.

TEST MENU

COVID-19

- Antibody & Antigen Test

Diabetes

- HbA1c

Kidney

- MAU
- CysC
- NGAL
- β2-MG

Inflammation

- PCT
- CRP
- SAA
- IL-6

RA

- Anti-CCP
- ASO
- RF

Anemia

- Ferritin
- Vitamin B12

Sex Hormone

- β-HCG
- LH
- FSH
- Prog
- PRL
- E2
- Testosterone
- AMH

Thyroid

- T3
- T4
- TSH
- fT3
- fT4

Allergy

- Total IgE

Hypertension

- Cortisol

Gastroenterology

- iFOB
- H.pylori Ag

Tumor

- AFP
- PSA
- CEA
- fPSA

Cardiac

- cTn I
- cTn T
- Myo
- CK-MB
- H-FABP
- NT-proBNP
- BNP
- hsCRP
- Lp-PLA2
- ST2

Bone Metabolism

- Vitamin D

Thrombus

- D-Dimer


Brain Injury

- S100β

Infectious Disease

- Dengue NS1-Ag

Instrument




Reagent



Quality Control



DISPOSABLE REAGENT CARTRIDGES PARAMETER



Cat. No.	Category	Test Item	Diseases	Specimen Type	Specimen Volume	Reaction Time	Measuring Range	Qualification
W241P0002	Bone Metabolism	Vitamin D	Vitamin D deficiency	S/P	75 µL	15 mins	5~100 ng/mL	CE NMPA
W241P0006/ W241PP007		Vitamin D	Vitamin D deficiency	WB/S/P	10 µL	8 mins	5~100 ng/mL	CE
W207		Diabetes	HbA1c	Diabetes Mellitus	WB	10 µL	5 mins	4.0~14.5%
W206	Kidney Function Markers	One Step MAU	Renal Injury	Urine	75 µL	3 mins	5~300 mg/L	CE NMPA
W219		CysC	Renal Injury	WB/S/P	10 µL	5 mins	0.2~10 mg/L	CE NMPA
W228		NGAL	Renal Injury	Urine	75 µL	15 mins	10~1500 ng/mL	CE NMPA
W229		β2-MG	Proximal Renal Tubule Dysfunction	WB/S/P	10 µL	10 mins	0.3~20 mg/L	CE NMPA
W231	Thyroid Function	T3	Thyroid Malfunction	WB/S/P	75 µL	15 mins	0.61~9.22 nmol/L	CE NMPA
W232		T4	Thyroid Malfunction	WB/S/P	75 µL	15 mins	12.87~300 nmol/L	CE NMPA
W220		TSH	Thyroid Malfunction	WB/S/P	75 µL	15 mins	0.1~100 mIU/L	CE NMPA
W200		ft3	Thyroid Malfunction	WB/S/P	75 µL	15 mins	0.40~50.00 pmol/L or 0.26~32.55 pg/mL	CE NMPA
W259		ft4	Thyroid Malfunction	WB/S/P	75 µL	15 mins	1.00~100 pmol/L or 0.078~7.77 ng/dL	CE NMPA
W225	Sex Hormone Markers	β-HCG	Fertility	WB/S/P	20 µL	15 mins	2~200,000 mIU/mL	CE NMPA
W244		LH	Fertility	WB/S/P	75 µL	15 mins	1~100 mIU/mL	CE NMPA
W245		FSH	Fertility	WB/S/P	75 µL	15 mins	1~100 mIU/mL	CE NMPA
W233		Prog	Fertility	WB/S/P	75 µL	15 mins	1.5~60 ng/mL	CE NMPA
W246P0001		PRL	Fertility	WB/S/P	75 µL	15 mins	1~200 ng/mL	CE NMPA
W247-C4P1-E		E2	Fertility	WB/S/P	75 µL	15 mins	9~3000 pg/mL	CE NMPA
W248-C4P1-E		Testosterone	Fertility	WB/S/P	75 µL	15 mins	0.2~15 ng/mL	CE NMPA
W257P0001		AMH	Fertility	WB/S/P	75 µL	15 mins	0.1~16 ng/mL	CE NMPA
W202	Cardiac Markers	NT-proBNP	Heart Failure	WB/S/P	75 µL	15 mins	18~35000 pg/mL	CE NMPA
W203		cTn I	Myocardial Infarction	WB/S/P	75 µL	15 mins	0.1~50 ng/mL	CE NMPA
W204		Myo	Myocardial Infarction	WB/S/P	75 µL	15 mins	2.0~400 ng/mL	CE NMPA
W205		CK-MB	Myocardial Infarction	WB/S/P	75 µL	15 mins	0.3~100 ng/mL	CE NMPA
W217		H-FABP	Myocardial Infarction	WB/S/P	75 µL	15 mins	1~120 ng/mL	CE NMPA
W216		3 in 1 (cTn I/Myo/CK-MB)	Myocardial Infarction	WB/S/P	75 µL	15 mins	Same with single items	CE NMPA
W222		BNP	Heart Failure	WB/P	75 µL	15 mins	5~5000 pg/mL	CE NMPA
W223		2 in 1 (cTn I/ NT-proBNP)	Myocardial Infarction	WB/S/P	75 µL	15 mins	Same with single items	CE NMPA
W240		hsCRP	CVD	WB/S/P	10 µL	3 mins	0.1~10mg/L	CE NMPA
W254P0001		One Step cTnT	Myocardial Infarction	WB/S/P	100µL-WB 75µL-S/P	15 mins	0.03~10 ng/mL	CE NMPA
W242	Lp-PLA2	Atherosclerosis	WB/S/P	75 µL	15 mins	10~1000 ng/mL	NMPA	
W243	ST2	Heart Failure	WB/S/P	75µL	15 mins	10~400ng/mL	CE	
W239	Brain Injury Markers	S100β	Cerebral Ischemia Cerebral Hemorrhage	WB/S/P	75 µL	15 mins	0.05~10 ng/mL	NMPA
W211		Thrombus Markers	D-Dimer	Thrombosis	WB/P	15µL-WB 10µL-P	5 mins	0.1~10 mg/L
W227	One step D-Dimer		Thrombosis	WB/P	100µL-WB 75µL-S/P	15 mins	0.1~10 mg/L	CE NMPA
W210/ W210P0004	Inflammation Markers	PCT	Sepsis	WB/S/P	75µL-WB 50µL-S/P	15 mins	0.1~100 ng/mL	CE NMPA
W201/ W201P0008		CRP (hsCRP+CRP)	Infection	Fingerstick or Venous WB/S/P	8.5µL-WB 5µL-S/P	3 mins	0.5~200 mg/L	CE NMPA
W218		2 in 1 (CRP+PCT)	Infection	WB/S/P	75 µL	15 mins	CRP:0.5~150 mg/L PCT:0.1~100 ng/mL	CE NMPA
W221/ W221P0004		SAA	Infection	Fingerstick or Venous WB/S/P	8.5µL-WB 5µL-S/P	3 mins	1 ~ 500 mg/L	CE NMPA
W251P0001		IL-6	Infection	WB/S/P	75 µL	15 mins	3~4000 pg/mL	CE NMPA

Cat. No.	Category	Test Item	Diseases	Specimen Type	Specimen Volume	Reaction Time	Measuring Range	Qualification
W208	Tumor Markers	AFP	Liver Cancer	WB/S/P	75 µL	15 mins	5~400 ng/mL	CE NMPA
W209		PSA	Prostate Cancer	WB/S/P	75 µL	15 mins	2~100 ng/mL	CE ₀₁₂₃ NMPA
W226		CEA	Cancer	WB/S/P	75 µL	15 mins	1~500 ng/mL	CE NMPA
W230		fPSA	Prostate Cancer	WB/S/P	75 µL	15 mins	0.2~30 ng/mL	CE ₀₁₂₃ NMPA
W279P0001	Gastroenterology	iFOB	Gastrointestinal bleeding	Feces	6 different parts of the feces	5 mins	25~1000 ng/mL	CE
W270P0001		H.pylori Ag	Infectious disease	Feces	10 mins	Quantitative	CE	
W282P0001	Infections Disease	Dengue NS1 Ag	Infectious disease	WB/S/P	75 µL	12 mins	Quantitative	CE
W298P0001	Allergy	Total IgE	Allergic diseases	WB/S/P	75 µL	10 mins	1~1000 IU/mL	CE
W250	Hypertension	Cortisol	Cortisol disorders	WB/S/P	75 µL	15 mins	50~1000 nmol/L	CE NMPA
W269P0001	Anemia Markers	Ferritin	Iron metabolism	WB/S/P	75 µL	10 mins	5 ~ 1000 ng/mL	CE NMPA
WX01P0001		Vitamin B12	VB12 Deficiency	S/P	40µL	10 mins	80~2000pg/mL	CE
W297P0001	Rheumatoid Arthritis	ASO	Rheumatic fever and glomerulonephritis	WB/S/P	10 µL	5 mins	20~1000 IU/mL	CE
W262P0001		RF	Rheumatoid Arthritis	WB/S/P	10 µL	5 mins	10 ~ 650 IU/mL	CE
WW01P0001		Anti-CCP	Rheumatoid Arthritis	WB/S/P	10 µL	10 mins	3.5-300 U/mL	CE
W290	COVID-19	2019-nCoV RBD Antibody Test	COVID-19	Fingerstick or Venous WB/ S/P	20µL	15 mins	Quantitative	CE
W286		2019-nCoV Antigen Test		NP/ OP	/	15 mins	Quantitative	CE
W278		2019-nCoV IgM/IgG Test		WB/S/P	10 µL	10 mins	Quantitative	CE

AUTOMATIC REAGENT KITS PARAMETER



Category	Cat. No.	Test Item	Diseases	Specimen Type	Reaction Time	Measuring Range	Qualification
Inflammation Markers	W210P0017	PCT	Sepsis	WB/S/P	15 mins	0.1-100 ng/mL	CE, NMPA
	W201P0015	CRP (hsCRP+CRP)	Inflammation	WB/S/P	3 mins	0.5-200 mg/L	CE, NMPA
Diabetes	W207P0005	HbA1c	Diabetes Mellitus	WB	5 mins	4-14.5%	CE, NMPA
Cardiac Markers	W203P0012	cTnI	Myocardial Infarction	WB/S/P	15 mins	0.1-50 ng/mL	CE, NMPA
	W205P0007	CK-MB	Myocardial Infarction	WB/S/P	15 mins	0.3-100 ng/mL	CE, NMPA
	W202P0013	NT-proBNP	Heart Failure	WB/S/P	15 mins	18~35000 pg/mL	CE, NMPA
Thrombus Markers	W211P0012	D-Dimer	Thrombosis	P/WB	5 mins	0.1-10 mg/L	CE, NMPA

Immunofluorescence Quantitative Test System

Finecare™ New Platform



Finecare™ FIA Meter X1
(FS-1000)



Finecare™ FIA Meter X2
(FS-2000)



Finecare™ FIA Meter X2 PLUS
(FS-2100)



Stability

TRFIA Methodology
Lower CV higher stability



Convenient

Simple operation procedure
Integrated handle design



Delicate design

Automatic reagent door
Charging indicator light



Instant results

Result in 6s under quick mode
Result in 8 min under standard mode

Specification

- Dimension: 209*290*187 mm
- Weight: 3.5 kg

Built-in Code Scanner

- No need for ID chip
- Ample batch information

Built-in Thermal Printer

- Real-time auto printing
- External printer supported

Rechargeable Internal Battery

- 10400 mAh battery
- 24-hour standby
- Low battery alert

Large Display

- 8-inch capacity touch screen
- Intuitive user interface
- Smooth operation experience like smart-phone

Connectivity

- Wifi connection to Bi-directional LIS/HIS
- RS232 and ethernet port
- 4 USB ports

Linux System

- Smooth operation
- Security system
- Custom intelligent management module
- Multi-user management

1 Faster to Get the Result:


Only 8mins

2 New Items:

Support for hs-cTnl

3 Excellent CV and Accuracy with Optimized FIA Methodology

Satisfying and stable results delivered in short time



Operation Procedure >



Entry of test information



Entry of patient information



Sample loading



Insert the cartridge

Cat. No.	Category	Test Item	Diseases	Specimen Type	Specimen Volume	Reaction Time	Measuring Range	Qualification
WS02P0001	Cardiac Markers	One Step hs-cTnl/CK-MB/Myo	Myocardial Infarction	S/P/WB	100uL whole blood; 75ul serum/plasma	8 min	hs-cTnl: 0.02 ~ 50 ng/mL CK-MB: 1 ~ 80 ng/mL Myo: 10 ~ 600 ng/mL	CE NMPA
WS01P0001		One Step hs-cTnl		S/P/WB	100uL whole blood; 75ul serum/plasma	8 min	0.02 - 50 ng/mL	CE NMPA
WS04P0001		One Step NT-proBNP	Heart Failure	S/P/WB	100uL whole blood; 75ul serum/plasma	8 min	18pg/mL-35000pg/mL	CE NMPA
WS03P0001	One Step BNP	P/WB		100uL whole blood; 75ul plasma	8 min	5-5000pg/mL	CE NMPA	
WS06P0001	Thrombus Markers	One Step D-Dimer	Thrombosis	P/WB	100uL whole blood; 75ul plasma	5 min	0.10 - 10.00mg/L FEU	CE
WS08P0001	Inflammation Markers	CRP	Infection	S/P/WB	8.5uL whole blood; 5ul serum/plasma	3 min	0.3 ~ 200 mg/L	CE
WS14P0001		One Step IL-6		S/P/WB	100uL whole blood; 75ul serum/plasma	8 min	3 ~ 4000 pg/mL	CE
WS09P1001 WS09P1002		SAA		S/P/WB	8.5uL whole blood; 5ul serum/plasma	8 min	1~ 200 mg/L	CE
WS12P0001		One Step PCT	Sepsis	S/P/WB	100uL whole blood; 75ul serum/plasma	8 min	0.020 ~ 100 ng/mL	CE



Chemiluminescence Platform

Mono-test Chemiluminescence Quantitative Test System

ACCRES CLIA Platform



ACCRES 90



ACCRES 8



*ACCRES 120

ACCRES is a series of automatic bench-top chemiluminescence immunoassay (CLIA) systems using mono-test reagent cartridge. It provides reliable test results with less maintenance cost.

Flexible & Efficient

Simultaneous processing of difference parameters

Maintenance-free

No liquid path inside
600 days MTBF

User-friendly

Easy operation requiring only 30-minute training

Fast

10 mins for the first result,
Up to 120 tests/hour



Reliable Accuracy

High sensitivity
Wide measuring range

Compact Design

70% smaller than traditional
CLIA analyzer

* Coming soon



EXCELLENT AS CLIA

- ALP-APS-5 chemiluminescence
- Excellent reagent performance
- Nano magnetic particles for separation
- Broad menu with featured items
- Reliable calibration and quality control system

EASY AS POCT

- Mono-test cartridge, ready-to-use
- Avoid cross-contamination
- Short turnaround time (TAT)
- All-in-one design, no consumables and extra cost
- Built-in scanner reads reagent information automatically

TEST MENU

- Available
- In Development

Tumor

- CEA
- AFP
- CA125
- tPSA
- CA19-9
- CA15-3
- PGI
- PGII
- G-17
- CA72-4
- HE4
- fPSA
- CYFRA 21-1
- CA50
- CA242
- SCCA
- NSE
- ProGRP

Infectious Disease

- Dengue NS1 Ag
- Dengue IgG
- Dengue IgM
- HIV
- HBsAg
- Anti-HCV
- Anti-TP
- H.Pylori (HP)

Inflammation

- PCT
- IL-6
- CRP
- SAA

Bone Metabolism

- 25-OH Vitamin D Total
- PTH

Fertility

- AMH
- Progesterone
- Testosterone
- Prolactin
- LH
- E2
- FSH
- β -HCG
- Inhibin B
- fTesto
- DHEA-S

TORCH

- Toxo IgG
- Toxo IgM
- Rubella IgG
- Rubella IgM
- CMV IgG
- CMV IgM
- HSV-1 IgG
- HSV-1 IgM
- HSV-2 IgG
- HSV-2 IgM

Thrombus

- D-dimer
- TM
- TAT
- FDP
- PIC
- t-PAIC

Anemia

- VB12
- Folate
- Ferritin

Cardiac

- cTnI
- hs cTnI
- CK-MB
- MYO
- BNP
- NT-proBNP
- H-FABP
- ST2
- hs cTnT
- Lp-PLA2
- hs-CRP

Thyroid

- TSH
- TT3
- TT4
- FT3
- FT4
- Tg
- Anti-Tg
- Anti-TPO

Allergy

- IgE

DIABETES

- Insulin
- C-Peptide
- INS-Ab

Others

- S100 β
- Cortisol
- ACTH

ACCRES

PRODUCT CATALOGUE NUMBER



- Specification: 36 tests/kit
- Shelf life: 18 months
- Storage: 2-8°C
- Calibrator and Quality Control included in the kit.

Category	Parameter	Tisenc Cat. No.	Wondfo Cat. No.	Sample Volume (µL)	Sample Type*	Test Time (mins)	Measuring Range	Qualification
Cardiology	cTnl	W121C	F121C	75	S/P/WB	15	0.05-50 ng/mL	CE, NMPA
	hs-cTnl	W122C	F122C	75	S/P/WB	15	0.005-50 ng/mL	CE, NMPA
	MYO	W123C	F123C	75	S/P/WB	15	1-1000 ng/mL	CE, NMPA
	CK-MB	W124C	F124C	75	S/P/WB	15	0.1-300 ng/mL	CE, NMPA
	BNP	W125C	F125C	75	P/WB	15	15-5000 pg/mL	CE, NMPA
	NT-proBNP	W126C	F126C	75	S/P/WB	15	15-35000 pg/mL	CE, NMPA
	H-FABP	W127C	F127C	75	S/P/WB	15	0.5-160 ng/mL	CE, NMPA
	hs cTnT	W86C	F86C	75	S/P/WB	15	5.0-10000.0 ng/L	CE
	ST2	W87C	F87C	75	S/P/WB	15	4.00-400.00 ng/m	CE
Thyroid	TSH	W30C	F30C	75	P/S	25	0.01-100 µIU /mL	CE, NMPA
	TT3	W31C	F31C	75	P/S	25	0.3-10 nmol/L	CE, NMPA
	TT4	W32C	F32C	75	P/S	25	6.44-320 nmol/L	CE
	FT3	W33C	F33C	70	P/S	25	1.0-50 pmol/L	CE, NMPA
	FT4	W34C	F34C	70	P/S	25	2-77.2 pmol/L	CE, NMPA
	Tg	W110C	F110C	55	P/S	25	0.10-500.00 ng/mL	CE
	Anti-TPO	W150C	F150C	55	P/S	25	2.00-1000.00 IU/mL	CE
	Anti-Tg	W151C	F151C	55	P/S	28	1.50~2500.00 IU/mL	CE
	Fertility	FSH	W41C	F41C	70	P/S	15	0.2-200 mIU/mL
LH		W42C	F42C	70	P/S	15	0.1-200 mIU/mL	CE, NMPA
PRL		W43C	F43C	70	P/S	15	4-5000 µIU/ml	CE, NMPA
AMH		W117C	F117C	75	P/S	15	0.01-23 ng/mL	CE, NMPA
Progesterone		W44C	F44C	100	P/S	23	0.5-191 nmol/L	CE, NMPA
Testosterone		W45C	F45C	100	P/S	23	0.1-15 ng/mL	CE, NMPA
E2		W40C	F40C	100	P/S	23	15-3000 pg/mL	CE, NMPA
β-HCG		W46C	F46C	60	P/S	15	1-5000 mIU/mL	CE, NMPA

* S: Serum, P: Plasma, WB: Whole Blood

Category	Parameter	Tisenc Cat. No.	Wondfo Cat. No.	Sample Volume (µL)	Sample Type*	Test Time (mins)	Measuring Range	Qualification
Allergy	IgE	W107C	F107C	65	P/S	15	0.3-2500 IU/mL	CE
	AFP	W60C	F60C	70	P/S	15	0.5-1000 IU/mL	CE
	CEA	W85C	F85C	65	P/S	15	0.2-1000 ng/mL	CE
	CA125	W104C	F104C	65	P/S	15	0.6-5000 U/mL	CE
	CA19-9	W105C	F105C	65	P/S	15	0.6-1000 U/mL	CE
Tumor	tPSA	W84C	F84C	65	P/S	15	0.1-100 ng/mL	
	fPSA	W103C	F103C	55	P/S	15	0.010-50.000 ng/mL	
	CA15-3	W106C	F106C	65	P/S	15	1-300 U/mL	CE
	CA72-4	W63C	F63C	55	P/S	15	0.20-300.00 U/mL	CE
	HE4	W64C	F64C	55	P/S	15	15.00-1500pmol/L	CE
	CYFRA 21-1	W65C	F65C	55	P/S	15	0.10-500.00 ng/mL	CE
	G-17	W93C	F93C	55	P/S	15	1.00-320.00pmol/L	CE
	PG I	W91C	F91C	55	P/S	15	1.00-200.00 ng/mL	CE, NMPA
	PG II	W92C	F92C	55	P/S	15	0.5~100 ng/mL	CE, NMPA
	Anemia	Ferritin	W06C	F06C	60	P/S	15	1-2000 ng/mL
Folate		W119C	F119C	85	P/S	28	1.0-20 ng/mL	CE, NMPA
VB12		W118C	F118C	75	P/S	35	125-2000 pg/mL	CE, NMPA
Bone Metabolism	25-OH VD	W120C	F120C	75	P/S	40	3.0-100 ng/mL	CE, NMPA
Thrombus	D-Dimer	W128C	F128C	75	P/WB	15	0.045-10 µg/mL (FEU)	CE, NMPA
	PCT	W116C	F116C	75	S/P/WB	15	0.04-100 ng/mL	CE, NMPA
Inflammation	IL-6	W90C	F90C	75	S/P/WB	15	1.5-5000 pg/mL	CE
	S100β	W61C	F61C	55	S	21	0.020-39.000 µg/L	CE
Brain Damage	S100β	W61C	F61C	55	S	21	0.020-39.000 µg/L	CE
Hypertention	Cortisol	W180C	F180C	75	S/P	21	0.054-63.400 µg/dL	CE
	Insulin	W186C	F186C	55	S/P	15	0.2-1000.00 µU/mL	CE
Diabetes	C-Peptide	W187C	F187C	55	S/P	15	0.01-40.00 ng/mL	CE

* S: Serum, P: Plasma, WB: Whole Blood

Fully-auto Chemiluminescence Quantitative Test System

Product Features

● Flexibility & Efficiency

- STAT mode supported
- Only 0.64 m² footprint with 60 sample positions
- <12 mins for the first result
- High throughput with 200T per hour
- RFID management, easy to load reagents

● Better User Experience

- Built-in touch screen operation
- Bidirectional LIS communication
- Visual management of consumables and reagents
- Cuvettes loading by dumping

● Smart & Accurate Machine System

- The carry-over is ≤ 0.5ppm
- Constant 2-8°C reagent refrigeration
- Needle collision protection and block detection
- Liquid level detection
- Automatic switching of double substrate bottles during the operation

● Fully Automatic Analyzer

- Primary tube sampling
- Automatic barcode identification
- Automatic dilution function
- Automatic rerun

Specification

Model *	FC-2000	FC-2100
Methodology	Chemiluminescence Immunoassay (CLIA)	
Dimension	900 mm * 710 mm * 590 mm	
Weight	140 kg	
Throughput	200 T/h	
Time for first result	< 12 min	
Reagent position	15	25
Reagent specification	50 or 100 T/Bottle	25 or 50 T/Bottle
Reagent refrigeration	2-8°C on board refrigeration	
Sample position	60	
Sample volume	5-150 µL	
Sample type	Whole blood / Serum / Plasma / Urine/ Cerebrospinal fluid	
Test mode	Batch / STAT	
Probe system	Liquid level detection / Block detection / Collision protection	
Storage capacity	≥ 60000 results	
Operation system	Windows 10 / bidirectional LIS compatible	
Dilution function	Support automatic dilution Fixed multiple dilution Exceeded linearity range remeasurement dilution	

* FC-2000 and FC-2100 share the same analyzer outlook.



FC-2000

- More reagent volume
- Larger test volume

FC-2100

- More reagent positions
- More test combinations

TEST MENU

- Available
- In Development



📄 Cardiac

- cTnT
- BNP
- MYO
- hs-cTnI
- Lp-PLA2
- CK-MB
- NT-proBNP
- MPO
- H-FABP
- ST2
- hs-CRP
- HCY*

📄 Inflammation

- IL-6
- PCT
- SAA

📄 Thyroid

- TSH
- FT3
- FT4
- TT3
- TT4
- Anti-TPO*
- Anti-TG*
- Tg*
- CT*
- T-Uptake*
- TBG*
- TRAb*
- TSAb*
- rT3*

📄 TORCH

- TOXO-IgG*
- TOXO-IgM*
- Rubella IgG*
- Rubella IgM*
- CMV-IgG*
- CMV-IgM*,
- HSV-1-IgG*
- HSV-1-IgM*
- HSV-2-IgG*
- HSV-2-IgM*

📄 Thrombus

- TM
- TAT
- FDP
- PIC
- t-PAIC
- D-Dimer

📄 Hypertension

- REN
- Cortisol
- A I
- A II
- ACTH
- ALD

📄 Gastric function

- G-17
- PGI
- PGI2
- HP-IgG*

📄 Tumor Markers

- CA50
- CA242
- CYFRA21-1
- SCC
- HE4
- NSE
- CA72-4
- CEA
- CA15-3
- CA125
- CA19-9
- ProGRP
- FPSA
- TPSA
- AFP*
- PIVKA-II*

📄 Liver fibrosis

- HA*
- CIV*
- LN*
- PIIINP*

📄 Craniocerebral injury

- S100β

📄 Bone metabolism

- 25-OH VD
- PTH

📄 Pre-eclampsia

- PLGF
- sFlt-1

📄 Kidney function

- NGAL

📄 Anemia

- Ferritin
- Folate
- VB12

📄 Fertility

- AMH
- β-HCG
- TESTO*
- PROG*
- E2*
- PRL*
- LH*
- FSH*
- InhB*

📄 Infectious disease

- HBeAg*
- HBeAb*
- HBcAb*
- HBsAg*
- HBsAb*
- anti-TP*
- HIV Ag/Ab*
- HCV*
- M. Pneumonia IgG*
- M. Pneumonia IgM
- Chlamydia pn. IgG*
- Chlamydia pn. IgM*



Molecular Diagnostics Platform

MOLECULAR DIAGNOSTIC SOLUTIONS

1

Collection

Flocked Swabs



- Our sterile flocked swabs provide optimized specimen collection at a low cost.
- The swabs feature a molded breakpoint that can be easily snapped.

Specimen Collection Buffer



- Optional inactivated and non-inactivated preservation solutions.

2

Extraction

Extraction Instrument



- High throughput automated system with magnetic bead-based purification.
- Compact size with large screen.
- Highly efficient and stable design provides optimized solution.

Extraction Kit



3

Amplification

Portable PCR Instrument



- Flexibility - Open platform with high capacity.
- Flexibility - High precision, Real-Time temperature control.
- Automated - Automatic analysis of results, no external computer required.
- User-Friendly - Visual animation guides operation; quick mode saves time.

Real-Time PCR Assay



Real-Time PCR Assay Test Menu

Category	Product	Targets
2019-nCoV	2019-nCoV PCR Assay	2019-nCoV (N gene) 2019-nCoV (ORF1ab gene)
	2019-nCoV/FluA/FluB/RSV PCR Assay	2019-nCoV (N gene) RSV A/B 2019-nCoV (ORF1ab gene) Flu A Flu B
	2019-nCoV/Flu/RSV PCR Assay	2019-nCoV (ORF1ab gene) Flu A/B RSV A/B
	2019-nCoV PCR Assay for B.1.1.529	2019-nCoV S371L 2019-nCoV E484A
	2019-nCoV Variants PCR Assay	2019-nCoV (ORF1ab gene) HV69-70del K417T T478K 246-252del K417N N501Y
	Mpox Virus	Mpox Virus PCR Assay
Respiratory Infection	FluA/FluB/RSV PCR Assay	RSV A/B Flu A Flu B
	MP/AdV/PIV PCR Assay	<i>M. pneumoniae</i> Adenovirus Parainfluenza virus
	Respiratory-Virus PCR Assay	RSV A/B <i>M. pneumoniae</i> Flu A Adenovirus Flu B Parainfluenza virus
	Respiratory-Bacteria PCR Assay	<i>S. pneumoniae</i> <i>M. pneumoniae</i> <i>H. influenzae</i> <i>C. pneumoniae</i> <i>B. pertussis</i> <i>L. pneumophila</i>
	Tuberculosis	MTB/RIF <i>M. tuberculosis</i> rpoB
	Gastrointestinal Infection	GI-Virus PCR Assay
GI-Bacteria PCR Assay		Shigella/EIEC <i>Vibrio parahaemolyticus</i> <i>Y. enterocolitica</i> Salmonella <i>C. jejuni</i> <i>E. coli</i> O157 <i>C. difficile</i>
C. difficile PCR Assay		<i>C. trachomatis</i> <i>M. hominis</i> <i>N. gonorrhoeae</i> <i>U. urealyticum</i> <i>M. genitalium</i> <i>U. parvum</i> <i>T. vaginalis</i>
Sexually Transmitted Diseases		STD PCR Assay
HPV	HPV16/18 PCR Assay	HPV16 HPV18

U-Card Dx™ Molecular Diagnostic System

MULTIPLEXED

RAPID

ACCURATE

FULLY AUTOMATED


Instrument features

- Designed for use at the point of care
- Portable and fully automated instrument
- Comprehensive results with PCR curves and Ct values
- Simple workflow without sample preparation
- No maintenance required
- Multiple data connection interface options
- Robust memory capability



Reagent features

- Microfluidic cartridge design
- Lyophilized reagents for ambient temperature stability
- Closed-system design avoids contamination
- Multiplex detection for up to 8 targets/controls
- Rapidly expanding test menu




3-STEP OPERATION

- Simple-to-use operation with on-screen guide
- Positive samples identified in as soon as 28 minutes

28~38 mins



Scan Cartridge



Load Sample



Insert Cartridge

Test Menu

• CE available • RUO

Upper Respiratory Tract

- 2019-nCoV
- 2019-nCoV/Flu A/Flu B
- Mpox virus
- 2019-nCoV/Flu A/Flu B/RSV
- Flu A/Flu B/RSV
- Strep A

Sexual Transmitted Infection (STI)

- CT/NG
- MG/TV
- UU
- CT/NG/UU
- CT/NG/MG/TV
- HPV 16/18 +12 highrisk
- GBS
- HSV 1/2



Immunohistochemistry Platform

FULLY AUTOMATED IHC STAINING SYSTEM

KEY TECHNOLOGY



PA-3600 Fully Automated IHC Staining Instrument

Efficient, flexible and stable immunohistochemistry staining platform

—01—

Flexibility

supports customization of incubation time and temperature for each slide

—02—

Optimal quality

Staining cover plate & wave-like mixing technology

—03—

High Accuracy

RFID chip technology, Automatic scanning of reagent information to reduce manual errors and real-time monitoring of reagent capacity

—04—


Eco-friendly design

V-shape design bottom of bottle enables small amount of reagent residue

—05—

High workflow efficiency

Continuous and random slide processing to optimize laboratory workflow




Independent temperature control technology for each slide

- 36 staining modules are equipped with independent heating modules.
- Temperature control range from room temperature to 101°C, accurate and uniform temperature control, temperature fluctuation <0.5°C.



Wave-like mixing technology

- The liquid on the surface of the slide moves in waves to improve the antigen-antibody reaction, shorten the staining time, and ensure thorough cleaning.
- The motor drives the cover plate through the worm gear, and the cover plate can be opened to a larger angle, which is convenient for maintenance and cleaning.



Independent liquid cover film technology

- Works long with IHC staining
- Improve reaction efficiency and ensure staining quality



RFID chip technology

- Reduce the pre-treatment time of the experiment, and place the reagent bottle without restriction.
- 40 reagent bottles are equipped with independent scanning antennas, and the time to identify 40 reagent bottles is less than 10s. It is recognized when inserted, efficient and accurate.



QR code identification technology

- Using QR code technology to identify slides, slides can be placed at will
- Identify 36 slides <30s



Economical reagent design

- Only needs 130 μL to cover the whole slide
- Reduce reagent residue



Uninterrupted operation


- During operation, slides can be loaded continuously without interfering with existing staining procedure.
- In the same way as stain module.

Software Operation Interface

Clear Partition & High Visualization

Easy Operation

Situation Display & Reminding Timely



WORK PROCESS



Secondary Antibody Reagent System

- PA-3600 Immunostaining Solution
- PA-3600 Dewax Solution
- PA-3600 Wash Solution
- PA-3600 Antigen Retrieval Solution (PH9.0)
- PA-3600 Antigen Retrieval Solution (PH6.0)




Dry Chemistry Platform

Design For POCT-Racing For Life





MAIN SPECIFICATIONS

Test channel	1
Dimension	331.5 (L) * 220.5 (W) * 153 (H) mm
Weight	2.3 kg
Screen	8 inches
Record storage	50,000
Incubation zone temperature	37 °C ±0.3°C
Thermal Printer	Auto-print, Support for external printing
Power adapter	AC 220 V ±22 V, 50 Hz ±1 Hz
Cable Interface	USB 2.0*2, RS232 Interface*1, LIS/HIS Interface*1



ONE-STEP OPERATION



APPLICATION SCENARIOS



COMPARED WITH OTHER BIOCHEMICAL ANALYZERS

	Accurate	Less Sample Volume	Easy-to-use	Free Of Maintenance
Fully Automatic Wet Chemistry Analyzer	√	√	×	×
Semi-automatic Wet Chemistry Analyzer	√	√	√	×
Disc-microfluidic Chemistry Analyzer	√	×	√	×
Dry Chemistry Analyzer (DC-101)	√	√	√	√

Test Parameters


Catalog No.	Category	Test Item	Sample Type	Reaction Time	Sample Volume	Measuring Range	Qualification
W416	Liver Function	ALT	Whole blood/Serum/Plasma	2.5 min	15 µL	5~800 U/L	CE NMPA
W420		AST	Serum	2 min	20 µL	10~800 U/L	CE NMPA
W414		ALB	Serum	1 min	7 µL	10~60 g/L	CE NMPA
W441		TP	Serum/Plasma	1.5 min	15 µL	20 g/L~110 g/L	CE NMPA
W444		DBIL	Serum/Plasma	3 min	15 µL	5 µmol/L~342 µmol/L	CE NMPA
W443		TBIL	Serum/Plasma	3 min	15 µL	8 µmol/L~650 µmol/L	CE NMPA
W785		GGT	Serum/Plasma	3 min	15 µL	7~1000U/L	CE
W784		ALP	Serum/Plasma	3 min	15 µL	20~1250U/L	CE
W442		TBA	Serum/Plasma	3 min	15 µL	8 mol/L~300.0 mol/L	CE NMPA
W412		UREA	Serum/Plasma	2 min	15 µL	2.5~40 mmol/L	CE NMPA
W418	Kidney Function	CREA	Serum	5 min	15 µL	30~1000 µmol/L	CE NMPA
W413		UA	Serum/Plasma	2 min	15 µL	120~1200 µmol/L	CE NMPA
W422		CK	Serum	3 min	20 µL	20~1600 U/L	CE NMPA
W423	Myocardial Enzyme	CK-MB	Serum	4 min	20 µL	5~300 U/L	CE NMPA
W419		LDH	Serum	2 min	15 µL	100~1800 U/L	CE NMPA
W426		a-HBDH	Serum	2 min	15 µL	20~1000 U/L	CE NMPA
W411	Pancreatitis	α-AMY	Serum/Urine	2 min	10 µL	20~1000 U/L	CE NMPA
W424		LPS	Serum	3 min	10 µL	5~300 U/L	CE NMPA
W417	Organophosphates Pesticides Poisoning	CHE	Serum	2 min	15 µL	200~12500 U/L	CE NMPA
W421	Alcoholism	ALC	Serum	2 min	15 µL	10.0~300.0 mg/dL	CE
W415	Trauma / Anemia	HB	Whole Blood	1 min	10 µL	45~200 g/L	CE NMPA
W425	Blood Glucose	GLU	Serum	2 min	7 µL	1.00~25.00 mmol/L	CE NMPA
W447	Blood Lipid	TG	Serum/Plasma	2 min	15 µL	0.5 mmol/L~6.5 mmol/L	CE
W448		HDL-C	Serum/Plasma	2 min	15 µL	0.38 mmol/L~2.6 mmol/L	CE
W446		TC	Serum/Plasma	2 min	15 µL	1.07 mmol/L~8.47 mmol/L	CE
W781	Ion	PHOS	Serum/Plasma	3 min	15 µL	0.2 mmol/L~3.5 mmol/L	CE
W783		Ca	Serum/Plasma	1 min	15 µL	0.50 mmol/L~3.50 mmol/L	CE
W782		Mg	Serum/Plasma	2 min	15 µL	0.15mmol/L~2.70mmol/L	CE
W434	Intestine Barrier	DAO	Serum/Plasma	5 min	15 µL	5-80 U/L	NMPA
W435		D-LAC	Serum/Plasma	2 min	15 µL	2-50 mg/L	NMPA

WONDFO COVID-19 DIAGNOSTIC SOLUTION

Speed Up the COVID-19 Control !

COVID-19 TEST TOTAL SOLUTION

*For illustrative purpose only



Molecular Test



• 2019-nCoV Real-time RT-PCR Assay



U-Card Dx™ System
 • 2019-nCoV/Flu Test
 • 2019-nCoV Test

Antigen Test



Colloidal Gold Platform
 • 2019-nCoV Antigen Test (NP/OP/Nasal/Saliva/Sputum)
 • Influenza/2019-nCoV Antigen Combo Test



Finecare™ FIA Platform
 • 2019-nCoV Antigen Test



ACCRE CLIA Platform
 • 2019-nCoV Antigen (CLIA)

Antibody Test



Finecare™ FIA Platform
 • 2019-nCoV IgM/IgG Test
 • 2019-nCoV Antibody Test
 • 2019-nCoV IgM Test
 • 2019-nCoV IgG Test



ACCRE CLIA Platform
 • 2019-nCoV IgG kit (CLIA)
 • 2019-nCoV IgM kit (CLIA)

Neutralizing Antibody Test



Finecare™ FIA Platform
 • 2019-nCoV RBD Test



ACCRE CLIA Platform
 • 2019-nCoV Neutralizing Antibodies (CLIA)

COVID-19 PRODUCT CATALOGUE

Platform	Product Name	Specimen	Reaction Time	Storage Condition	Qualification	Packing Size	Catalog No.
Colloidal Gold	Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)	Nasopharyngeal swab/ Oropharyngeal swab	15 mins	2~30 °C	CE TGA	20 T	W196
	Wondfo 2019-nCoV Antigen Test (For Self-testing)	Nasal Swab	15 mins	2~30 °C	CE 0123 TGA	1/5 T	W634
	Wondfo Influenza/2019-nCoV Antigen Combo Test	Nasopharyngeal swab	15 mins	2~30 °C	CE	20 T	W630
Finecare™ FIA	Finecare™ 2019-nCoV RBD Antibody Test	Fingerstick or Venous whole blood/ Serum/ Plasma	15 mins	4~30 °C	CE	25 T	W290
	Finecare™ 2019-nCoV Antigen Test	Nasopharyngeal swab/ Oropharyngeal swab	15 mins	4~30 °C	CE	25 T	W286
	Finecare™ 2019-nCoV Antibody Test		10 mins	4~30 °C	CE	25 T	W276
	Finecare™ 2019-nCoV IgM/IgG Test	Venous whole blood/ Serum/Plasma	10 mins	4~30 °C	CE	25 T	W278
	Finecare™ 2019-nCoV IgM Test		10 mins	4~30 °C	CE	25 T	W277
	Finecare™ 2019-nCoV IgG Test		10 mins	4~30 °C	CE	25 T	W280
Molecular	Wondfo 2019-nCoV Real-Time RT-PCR Assay	Nasopharyngeal swab/ Oropharyngeal swab/ Sputum	<60 mins after extraction	10~30 °C	CE	24/48T	W275P001-004
	U-Card Dx™ 2019-nCoV/Flu Test	Nasopharyngeal swab	40 mins	2~30 °C	CE	12/24T	WG01P0001/2
	U-Card Dx™ 2019-nCoV Test	Nasopharyngeal swab	40 mins	2~30 °C	CE	12/24 T	WG01P0001/2
ACCRE CLIA	2019-nCoV IgG kit (CLIA)	Serum/Plasma	22 mins	2~8 °C	CE	36/60 T	W78
	2019-nCoV IgM kit (CLIA)	Serum/Plasma	22 mins	2~8 °C	CE	36/60 T	W77
	2019-nCoV Neutralizing Antibodies (CLIA)	Serum/Plasma	15 mins	2~8 °C	CE	36/60 T	W101
	2019-nCoV Antigen (CLIA)	Nasal swab	20 mins	2~8 °C	CE	36/60 T	W94



✓ Part of major markets worldwide

Vision

Racing for life

Mission

Make biotechnology benefit the public

Core Values

We strive to achieve goals with full passion to meet every customer's needs; we respect every employee and give them opportunities to constantly improve themselves; we dare to dream and act to bring success to every partner. We pursue excellence to safeguard everyone's health.

WONDFO MASCOT – FOFO

Where there is a whale, there is a way.



To raise awareness of protecting environment and blue whales.