

BLOOD BAG WITH LEUCOCYTE FILTER CPDA-1 SOLUTION													
 STERILE • NON PYROGENIC • NON TOXIC For collection of 450ml of Blood. 63ml of Anticoagulant Citrate Phosphate Dextrose Adenine Solution (USP) Each 100ml of CPDA-1 Contains: <table style="width: 100%; border-collapse: collapse;"> <tr> <td>Citric Acid monohydrate</td> <td style="text-align: right;">0,327g</td> </tr> <tr> <td>Sodium Citrate dihydrate</td> <td style="text-align: right;">2,63g</td> </tr> <tr> <td>Sodium Biphosphate monohydrate</td> <td style="text-align: right;">0,222g</td> </tr> <tr> <td>Dextrose monohydrate</td> <td style="text-align: right;">3,19g</td> </tr> <tr> <td>Adenine</td> <td style="text-align: right;">0,0275g</td> </tr> <tr> <td>Water for injection</td> <td style="text-align: right;">q.s</td> </tr> </table>	Citric Acid monohydrate	0,327g	Sodium Citrate dihydrate	2,63g	Sodium Biphosphate monohydrate	0,222g	Dextrose monohydrate	3,19g	Adenine	0,0275g	Water for injection	q.s	Donor Number (Name) Blood Group Rh - Type Collection Date Do not infuse after Serology: non-reactive to Test
Citric Acid monohydrate	0,327g												
Sodium Citrate dihydrate	2,63g												
Sodium Biphosphate monohydrate	0,222g												
Dextrose monohydrate	3,19g												
Adenine	0,0275g												
Water for injection	q.s												
CAUTION 1. Do not use unless the anticoagulant is clear. 2. Store blood between + 2°C and + 6°C. 3. Crossmatch before transfusion. 4. Mix blood thoroughly immediately before use. 5. Do not add medication to this blood. 6. Transfusion set must have a filter.													
DO NOT VENT LOT : : :													
Shandong Weigao Group Medical Polymer Co.,Ltd. Add:No.312 Shichang Road, Weihai, Shandong, China.													

Collection Date	Donation Number	EXPIRES
Place Product Label here after bag is filled.		Store Platelets 3 days at 20°C - 24°C Store Red Blood Cells 35 days at 2°C - 6°C
DO NOT VENT		
STERILE, NON-PYROGENIC FLUID PATH This product may transmit infectious agents. PROPERLY IDENTIFY INTENDED RECIPIENT		LOT NO.: