




CATALOGUE 2021

CARDIAC CATHETERS
AND CANNULAE



Cardiovascular cannulation



Index

- 4 Arterial cannulae
- 10 Aortic catheters
- 12 Flex line venous cannulae
- 14 Flex line two stage venous
cannulae

Arterial cannulae

REINFORCED, CURVED TIP WITH FLANGE



TYPE A



TYPE B

TYPE A - BODY LENGTH: 30 CM

TYPE B - BODY LENGTH: 21 CM

Description	Size	Standard	With vent plug	With vented red cap	Standard	With vent plug	With vented red cap
Without connector	20 Fr	A2016	A2016V	-	B2016	B2016V	-
	22 Fr	A2216	A2216V	-	B2216	B2216V	-
	24 Fr	A2416	A2416V	-	B2416	B2416V	-
With connector	20 Fr	A20161	-	A20161V	B20161	-	B20161V
	22 Fr	A22161	-	A22161V	B22161	-	B22161V
	24 Fr	A24161	-	A24161V	B24161	-	B24161V
With luer connector	20 Fr	A20162	-	A20162V	B20162	-	B20162V
	22 Fr	A22162	-	A22162V	B22162	-	B22162V
	24 Fr	A24162	-	A24162V	B24162	-	B24162V

PRODUCT INFORMATION


Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.

Q Connection site: 3/8"
Qty/box: 20 Pcs

+ Kink resistant reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Arterial cannulae

REINFORCED, STRAIGHT TIP WITH RING



TYPE A



TYPE B

TYPE A - BODY LENGTH: 30 CM

TYPE B - BODY LENGTH: 21 CM

Description	Size	Standard	With vent plug	With vented red cap	Standard	With vent plug	With vented red cap
Without connector	20 Fr	A2013	A2013V	-	B2013	B2013V	-
	22 Fr	A2213	A2213V	-	B2213	B2213V	-
	24 Fr	A2413	A2413V	-	B2413	B2413V	-
With connector	20 Fr	A20131	-	A20131V	B20131	-	B20131V
	22 Fr	A22131	-	A22131V	B22131	-	B22131V
	24 Fr	A24131	-	A24131V	B24131	-	B24131V
With luer connector	20 Fr	A20132	-	A20132V	B20132	-	B20132V
	22 Fr	A22132	-	A22132V	B22132	-	B22132V
	24 Fr	A24132	-	A24132V	B24132	-	B24132V

PRODUCT INFORMATION

AndoCor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.




Connection site: 3/8"
Qty/box: 20 Pcs



Kink resistant reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Flex line venous cannulae

REINFORCED WITH LIGHTHOUSE TIP



PRODUCT INFORMATION


Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.

Q Connection site: 3/8" - 36 Fr: 1/2"
Length: 39 cm
Qty/box: 20 Pcs


+ Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

Size	Straight tip	Right angled 90° tip
20 Fr	01V201L8	01V201L9
22 Fr	01V221L8	01V221L9
24 Fr	01V241L8	01V241L9
26 Fr	01V261L8	01V261L9
28 Fr	01V281L8	01V281L9
30 Fr	01V301L8	01V301L9
32 Fr	01V321L8	01V321L9
34 Fr	01V341L8	01V341L9
36 Fr	01V361L8	01V361L9

VENOUS CANNULA PRESSURE DROP VS FLOW



VENOUS CANNULA PRESSURE DROP VS FLOW



Special configurations

CONNECTION SITE 3/8"



Size	Straight tip	Right angled 90° tip
36 Fr	01V362L8	01V362L9

NON-REINFORCED



Size	Straight tip
32 Fr	VC3208

Flex line two stage venous cannulae

PROXIMAL AND DISTAL REINFORCED WITH LIGHTHOUSE TIP



ROUND BODY



FLAT BODY



ROUND BODY					FLAT BODY		
Size	Length	Without connector	With connector	With luer connector	Without connector	With connector	With luer connector
32/40	36	01V32L40	01V32L401	01V32L402	01V32L40FB	01V32L401FB	01V32L402FB
34/40	36	01V34L40	01V34L401	01V34L402	01V34L40FB	01V34L401FB	01V34L462FB
34/46	40	01V34L46	01V34L461	01V34L462	01V34L46FB	01V34L461FB	01V34L462FB
36/46	40	01V36L46	01V36L461	01V36L462	01V36L46FB	01V36L461FB	01V36L462FB
36/51	40	01V36L51	01V36L511	01V36L512	01V36L51FB	01V36L511FB	01V36L512FB

PRODUCT INFORMATION


Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.

- 


Connection site: 1/2"
Qty/box: 20 Pcs
- 

Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

VENOUS CANNULA PRESSURE DROP VS FLOW




Cardioplegia



Index

- 18 Aortic root cannulae
- 20 Retrograde cardioplegia cannulae
- 21 Air aspiration needles
- 22 Ostial perfusion cannulae
- 24 Cardioplegia adapters

Aortic root cannulae



Description	Tip size	Needle size	Reference
Type A Without vent line	6 Fr	16 G	ARC07
	7 Fr	14 G	ARC04
	9 Fr	12 G	ARC01
Type B Without vent line With suture ring	6 Fr	16 G	ARC07R
	7 Fr	14 G	ARC04R
	9 Fr	12 G	ARC01R
Type C With vent line	6 Fr	16 G	ARC08
	7 Fr	14 G	ARC05
	9 Fr	12 G	ARC02
Type D With vent line	6 Fr	16 G	ARC08S
	7 Fr	14 G	ARC05S
	9 Fr	12 G	ARC02S


PRODUCT INFORMATION

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart. The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.

Q Length: 15 cm
Connection size vent line: 1/4"
Qty/box: 20 Pcs

+ Easy of use
Anti-bleeding valve helps reduce blood loss

AORTIC ROOT CANNULAE PRESSURE DROP VS FLOW



Paediatric cannulae



Index

- 30 Paediatric arterial cannulae
- 31 Paediatric venous cannulae

Paediatric venous cannulae

REINFORCED WITH 90° METAL TIP



PRODUCT INFORMATION

BODY LENGTH: 35 CM

Andocor paediatric venous cannulae are intended for drainage from the superior and inferior vena cava during cardiopulmonary bypass surgery up to six hours or less.

Q Connection: 10-12 Fr: 3/16" - 1/4" | Other: 1/4"
Qty/box: 20 Pcs

+ Kink resistant wire reinforced PVC tubing
Steel tip with thin wall for high flow rates
Additional side holes for improved drainage

Description	Size	Reference
90° metal tip	10 Fr	01V101L7
	12 Fr	01V121L7
	14 Fr	01V141L7
	16 Fr	01V161L7
	18 Fr	01V181L7
	20 Fr	01V201L7

METAL TIP VENOUS CANNULA PRESSURE DROP VS FLOW

