

ERCP Catheters

Single Use



Guide Wire Compatibility

Contrast injection while guide wire is in place

Color-Marked

Macroscopically controlled cannulation depth

ERCP Catheters – Single Use



Standard Tip



Tapered Tip



Metal Cone



Round Metal



Metal Ball

ERCP Catheter – with Luer-Lock Stylet

With colored markers; compatible with guide wire .035"

Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GER-01-18-035	Standard	1.8	200	2.0
GER-02-18-035	5 mm Tapered	1.8	200	2.0
GER-04-18-035	Metal Cone	1.8	200	2.0
GER-05-18-035	Round Metal	1.8	200	2.0
GER-07-18-035	Metal Ball	1.8	200	2.0



Tapered Tip

ERCP Catheter – with Luer-Lock Stylet

With colored markers; compatible with guide wire .021"

Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GER-03-18-021	10mm Tapered	1.8	200	2.0



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
 Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
 Email: sales@medi-globe.de • www.medi-globe.de