

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

## Nebulizer Kits, funnel/funnel

Nebulizer Kit for inhalation of medicines and liquids by mask and safety tubing with funnel. The set is for single patient use only and should be disposed of after 24 hours. Aerosol capacity of nebulizer according to EN 13544-1:2007.

**Length (oxygen tube):** 210 cm

**Set content:** Mask with Lips style, nebulizer chamber 6 ml with lateral oxygen connection, safety connecting tube with funnel/funnel connector

**Sterile:** No

[read more](#)

Article no.	Application area	minimal quantity
0100008852	Child	50
<b>0100008853</b>	Adult	50

